Seraline[®] Anti-Borrelia-8 IgG

Seraline[®] Anti-Borrelia-8 IgM

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Seraline[®] Anti-Borrelia-8 IgG

Seraline[®] Anti-Borrelia-8 IgM

Line Immunoassay for the qualitative detection of IgG or IgM antibodies directed against *Borrelia burgdorferi* sensu lato in serum or plasma of human origin

REF	LIA-006-8 G		20
REF	LIA-006-8 G-12		12 x 20
REF	LIA-006-8 M		20
REF	LIA-006-8 M-12		12 x 20
IVD	In-vitro-diagnostic medical device		

Seramun Diagnostica GmbH • Sprenhagener Str. 1 • 15754 Heidesee • Germany •
 T +49 33767 791-10 • info@seramun.com • www.seramun.com

In-vitro diagnostic medical device	Unique device identifier	Manufacturer
Country of manufacture and date of manufacture	Article number	Serial number
Keep away from sunlight	Humidity limitation	Batch code
Consult instructions for use	Temperature limit	Do not reuse
Sufficient for <i>n</i> tests	Biohazard	Use-by date
		Attention

Intended Use

Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia-8 IgM are IVD tests for the qualitative determination of antibodies of the IgG or IgM isotype against the antigens p100 Bafz, VlsE Bafz, p58 Bgar, p41 Bbur, p39 Bafz (BmpA), OspA Bafz (p31), OspC mix (p23) and DbpA mix (p18) of various *Borrelia* species in serum or plasma (citrate, EDTA, heparin) of human origin by a laboratory professional user.

The tests can be used in combination with accessory Seraline®scan (software).

They are intended to aid in the diagnosis or for the exclusion of a disease caused by *Borrelia* spp. in samples from patients with suspected Lyme disease.

The tests must not be used with specimen materials other than serum or plasma (citrate, EDTA, heparin) of human origin, for screening, monitoring, prediction, diagnosis, prognosis, as companion diagnostic, in the near-patient setting and by lay persons.

Principle of the Test

Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia-8 IgM are line immunoassays (membrane-based immunoassays).

Detection of specific antibodies is performed in 3 steps:

Step 1

Incubation of 1 : 101 diluted samples and test strips for 45 min at room temperature. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with diluted wash and incubation buffer (WIB).

Step 2

Incubation of test strips with HRP-labeled conjugate antibodies for 45 min at room temperature. Removal of unbound conjugate antibodies by aspiration and 3 wash cycles with WIB.

Step 3

Incubation of test strips at room temperature for 10 min with substrate. The reaction is stopped by aspiration and 3 wash cycles with deionized water. Removing of the residual liquid by drying the test strips between filter paper. The coloration of the developed strips is stable when stored in the dark.

Fixation of the developed test strips on the evaluation template. When doing so, placement of the functional control of the test strips exactly on the separating line pre-printed in the evaluation template. Identification of the stained antigen bands on the dried test strips with help of the evaluation template. Determination of the test result according to the evaluation criteria and documentation of the identified bands in the evaluation template.

Alternatively, image acquisition with a conventional flatbed scanner and evaluation with software Seraline®scan is possible.

Test Components (Delivery Scope)

		For 20 determinations	For 12 x 20 determinations
1	TESTSTR	Test strips with evaluation template Specific antigens on nitrocellulose membrane (see antigens used)	20 strips color coding: dark blue in a folding booklet in a pressure seal bag
2	WIB (5x)	Wash and incubation buffer (5x) Seramun® Wash and incubation buffer (5x) Aqueous buffer with salt additives	20 strips color coding: dark blue in a folding booklet in a pressure seal bag
3	CONJ HRP IgG CONJ HRP IgM	Anti-human IgG-HRP-Conjugate anti-human IgG-HRP or Anti-human IgM-HRP-Conjugate anti-human IgM-HRP	20 strips color coding: dark blue in a folding booklet in a pressure seal bag
4	SUBSTR	Substrate SeramunBlau® prec <0.1 % 3,3',5,5'-tetramethylbenzidine	70 mL concentrate for 350 mL diluted buffer, colorless, black cap
5	INCUTRAY	Anti-human IgG-HRP-Conjugate anti-human IgG-HRP or Anti-human IgM-HRP-Conjugate anti-human IgM-HRP	12 x 70 mL concentrate for 350 mL diluted buffer, colorless, black cap
6	TEMPLATE	Anti-human IgG-HRP-Conjugate anti-human IgG-HRP or Anti-human IgM-HRP-Conjugate anti-human IgM-HRP	12 x 35 mL IgG: ready-to-use solution, colored red, red cap 12 x 35 mL IgM: ready-to-use solution, colored green, green cap
7		Substrate SeramunBlau® prec <0.1 % 3,3',5,5'-tetramethylbenzidine	35 mL ready-to-use solution blue cap
8		Incubation tray with cover	12 x 35 mL ready-to-use solution blue cap
		Evaluation template	-
		Certificate of Analysis	1 piece
		Instructions for Use	1 piece

Antigens Used

Name*	Characterization / Specificity		Origin
p100 Bafz	Protein of membrane vesicles on the surface	Highly specific	<i>B. afzelii</i>
VlsE Bafz	Surface lipoprotein	Specific	<i>B. afzelii</i>
p58 Bgar	Not characterized	Specific	<i>B. garinii</i>
p41 Bbur	Structural protein of Endoflagella	Unspecific	<i>B. burgdorferi sensu stricto</i>
p39 Bafz (BmpA)	Flagellin complex Borrelia membrane protein A	Highly specific	<i>B. afzelii</i>
OspA Bafz (p31)	Surface protein A Outer surface protein A	Specific	<i>B. afzelii</i>
OspC mix (p23)	Surface protein C Outer surface protein C	Highly specific	<i>B. burgdorferi sensu stricto</i> , <i>B. garinii</i> , <i>B. bavariensis</i> , <i>B. afzelii</i> , <i>B. spielmanii</i>
DbpA mix (p18)	Surface protein Decorin binding protein A	Specific	<i>B. burgdorferi sensu stricto</i> , <i>B. garinii</i> , <i>B. afzelii</i> , <i>B. spielmanii</i>

* In literature, the nomenclature of the antigens concerning molecular mass and specificity is not consistent in different stages of the disease.

Additional Materials and Aids Required for the Test Procedure

Single-channel micropipette and pipette tips for 15 µL and 1500 µL • plastic tweezers • graduated cylinders and beakers • test tubes for sample dilution • rocking shaker • filter paper • glue stick • aspiration system with collection vessel for infectious solutions • deionized water • stopwatch

Optional:

Flatbed scanner • evaluation software Seraline® scan • automatic processing device

Important Information

This device is for *in vitro* diagnostic use only. The kit may be performed by trained laboratory personnel only.

Follow the instructions for use carefully.

The shelf life specified must be observed. Do not mix components with reagents from other manufacturers. Do not use reagents from damaged packages or bottles. Do not mix components with reagents from other manufacturers.

Mixing of test kit components of different lots is only allowed for wash and incubation buffer (5x) and substrate.

All serious incidents occurring in relation with Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia IgM must be reported to the manufacturer and the competent authority of the EU member state in which the user and/or the patient are located.

Workplace Requirements

A clean and fiber-free workspace is required to perform line immunoassays. Particles or fibers which may adhere to the test strips can lead to erroneous results. The workstation and kit components should not be exposed to direct sunlight.

Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all components to room temperature before use.

Store substrate protected from light!

The order of the pipetting steps and duration of the wash and incubation steps must be observed.

Avoid bubble formation during pipetting, this will lead to evaluation errors.

Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious. Additional information may be taken from the Material Safety Data Sheet. The product contains the following hazard component/-s.

TESTSTR	-	Contains material of animal and microbiological origin.
WIB (5x)	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May cause allergic reactions.
	EUH210	Safety data sheet available on request. Contains material of animal origin.
CONJ HRP	EUH 210	Safety data sheet available on request Contains material of animal and microbiological origin.
SUBSTR	EUH210	Safety data sheet available on request.

Limitations of the Procedure

Result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. In an early stage of infection antibodies may not be present yet or below detection limit.

Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

The use of contaminated samples may increase the background color of the test strips and lead to incorrect results.

In rare cases, e.g., hypergammaglobulinemia, circulating immune complexes, or milk antibodies, the test strips may turn blue very quickly. In such cases the color reaction should be stopped early by rinsing the test strips 3 times with deionized water.

Sample Treatment

Sample Collection

Collect serum or plasma (citrate, EDTA, heparin) of human origin in suitable sampling container.

Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 48 hours. For longer periods samples have to be stored at < -15 °C. Repeated freezing and thawing of samples should be avoided.

It should be noted that diluted antibodies may not be stable, which may lead to a loss of activity. Required follow-up testing at higher dilutions therefore should be performed same day.

Sample Preparation

Serum/Plasma:

Bring samples to RT before use. In general, homogeneity should be ensured by brief shaking. Dilute serum and/or plasma samples 1 : 101 (v/v) (15 µL sample and 1500 µL diluted wash and incubation buffer (prepared from WIB (5x)) in the incubation well.

Reagent Treatment

Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and test strips is stable until the printed expiration date when stored at 2...8 °C. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash and incubation buffer can be stored at 2...8 °C for up to 1 month.

Reagent preparation

Allow all components to reach **room temperature** prior to use. Return unused test strips to the pouch and seal.

Dilute **wash and incubation buffer WIB (5x)** 1 : 5 with deionized water.

Example: 70 mL WIB (5x) + 280 mL deionized water. The prepared wash and incubation buffer must be thoroughly mixed before use. Avoid foam formation!

The **substrate** must be protected from direct light.

Assay Procedure

Test must be performed at room temperature (18...25 °C). Bring **all** ready-to-use test reagents and the diluted wash buffer to **room temperature**.

Place test strips in incubation wells using plastic tweezers. Imprinted numbers should point upwards. All incubation steps should be performed on a rocking shaker with a recommended frequency of 20 to 30 rotations per minute.

Follow the pipetting scheme and time schedules of the protocol

1. Incubate test strip **TESTSTR** with **1.5 mL** diluted wash and incubation buffer for **5 min** on the shaker
2. Add **15 µL** of sample.
3. Cover tray with lid **INCUTRAY** and incubate on rocking shaker for **45 min**.
4. Aspirate solutions and wash test strips **3x 5 min** with **1.5 mL** diluted wash and incubation buffer.
5. **Add 1.5 mL** conjugate **CONJ HRP IgG** or **CONJ HRP IgM**.
6. Cover tray with lid **INCUTRAY** and incubate for **45 min** on rocking shaker.
7. Aspirate solutions and wash test strips **3x 5 min** with **1.5 mL** diluted wash and incubation buffer.
8. Incubate test strips with **1.5 mL** **SUBSTR** for **10 min** on rocking shaker.
9. Aspirate substrate and rinse test strips **3x** with **1.5 mL** deionized water to stop color reaction.
10. **Dry** test strips between filter paper, then evaluate.

Evaluation of Results

Validity Criteria for the Test

Each test strip of Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia IgM contains 5 control bands:

1. The function control must be visible.
2. The conjugate controls IgG and IgM must be visible (positive band with the conjugate IgG or IgM used).

Caution: The respective other conjugate control can also show a pale, non-specific band depending on the sample.

3. The intensity of the negative control must be lower than that of the cut-off control.
4. The cut-off control must be visible (intensity of the band is used to discriminate between positive and negative results).

Test may be evaluated if function control and cut-off control are visible.

If the above-mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash and incubation buffer dilution, wash steps, etc.). In case of repeated failure of the quality criteria contact the manufacturer.

Evaluate dry test strips only. Protect test strips from direct sunlight. Assign protein bands using the enclosed evaluation template.

Alternatively, image acquisition with a conventional flatbed scanner and automatic evaluation with Seraline®scan evaluation software is possible.

Interpretation of Results

Bands with an intensity \geq Cut-off control only are considered positive.

Evaluation	IgG	IgM
Positive	≥ 2 bands except 1 band + p41 Bbur	1 band OspC mix (p23) or DbpA mix (p18) or ≥ 2 other bands
Borderline	1 band VisE Bafz or 1 band + p41 Bbur	1 band except OspC mix (p23) or p41 Bbur or DbpA mix (p18)
Negative	Max. 1 band except VisE Bafz	No band except 1 band p41 Bbur

Performance Characteristics

Interfering Substances

All tested substances do not have a significant impact on the test results even if present in elevated concentrations in the serum of 20 mg/dL bilirubin C and F each (simulation of icteric samples); 500 mg/dL hemoglobin (simulation of hemolytic samples), 1000 mg/dL lipids (triglycerides) (simulation of lipemic samples) and 500 IU/mL rheumatoid factor.

Cross-Reactivity

Sera (n = 195) from patients with potential cross-reactive antibodies were analyzed by Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia-8 IgM.

The evaluation of the IgG and IgM antibody analysis is shown in the table:

n = 195		Seraline® Anti-Borrelia-8		
		positive	borderline	negative
ANA	(n = 51)	7	3	41
EBV	(n = 30)	3	0	27
CMV	(n = 67)	13	2	52
Treponema pallidum	(n = 47)	4	1	42

ANA: antinuclear antibodies, EBV: Epstein-Barr Virus, CMV: Cytomegalovirus

Diagnostic Sensitivity

Pre-characterized sera (n = 78) were analyzed to determine the diagnostic sensitivity of Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia-8 IgM:

n = 78		Reference test	
		positive	negative
Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia-8 IgM	positive	51	10
	negative	1	16
Samples total		52	26

Sensitivity**: 98.1%

** All n = 11 discrepantly determined samples were examined by two different western blots. 8 of the 10 positive Seraline® Anti-Borrelia-8 IgM results could be confirmed by the blot results. One sample determined negative with Seraline® Anti-Borrelia-8 IgM and positive with the comparative test turned out weakly positive in the western blot

Sera from patients with suspected acute Borrelia infection (n = 32) were analyzed by Seraline® Anti-Borrelia-8 IgM and a reference test for Borrelia-specific IgM antibodies:

n = 32		Reference test	
		positive	negative
Seraline® Anti-Borrelia-8 IgM	positive	9	8
	negative	0	15
Samples total		9	23

Sensitivity***: 100%

***All n = 8 discrepantly determined samples were tested in western blot. 6 of the 8 positive Seraline® Anti-Borrelia-8 IgM results could be confirmed by the blot results.

Diagnostic Specificity

A study population of 400 blood donors was investigated to determine the specificity in Seraline® Anti-Borrelia-8 IgG and Seraline® Anti-Borrelia-8 IgM. 47 samples (39 positive and 8 borderline) tested positive with the IgG test strip, 22 (17 positive and 5 borderline) tested positive with the IgM test strip.

n = 400		Reference test	
		positive	negative
Seraline® Anti-Borrelia-8 IgG	positive	0	47
	negative	0	353
Samples total		0	400

Specificity****: 88,2 %

**** Borderline samples were considered positive.

n = 400		Reference test	
		positive	negative
Seraline® Anti-Borrelia-8 IgM	positive	0	22
	negative	0	378
Samples total		0	400

Specificity*****: 94,5 %

***** Borderline samples were considered positive.

Application

Software-based Evaluation

Evaluation is performed visually using the enclosed template. Alternatively, image acquisition with a conventional flatbed scanner and evaluation with Seraline®scan evaluation software is possible.

Change History

Version	Section	Modifications
2022-08_v02_DE_EN	Evaluation of results	Correction of the validity criteria for the test

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