Seraline[®] Vaskulitis-3 IgG

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Literatur:

1. Jenette JC, Falk RJ, Andrassy K. et al, Nomenclature of systemic vasculitides: proposal of an international conference. *Arthritis Rheumatol* 1994, 37, 187-192.
2. Charles LA, Caldas MLR, Falk RJ et al, Antibodies against granule proteins activate neutrophils in vitro. *J Leukocyte Biol* 1991, 50, 539-546.
3. Churg J, Churg A, Idiopathic and secondary vasculitis: a review. *Mod. Pathol.* 1989, 2, 144-160.
4. Cohen Tervaert JW, van der Waude FJ, Fauci, AS et al, Association between active Wegener's Granulomatosis and anticytoplasmic antibodies. *Arch Intern Med* 1989, 149, 2461-2465.
5. Weiss SJ, Tissue destruction by neutrophils. *New Engl. J Med* 1989, 320, 365-376.
6. Ozaki S, ANCA associated vasculitis: Diagnostic and therapeutic strategy. *Allergol. Int.* 2007, 56, 87-96.
7. Chen M, Yu F, Wang SX et al, Renal histology in Chinese patients with anti-myelo-peroxidase autoantibody-positive Wegener's granulomatosis. *Nephrol. Dial. Transplant.* 2007, 22, 139-145.
8. Papiris SA, Manali ED, Kalomenidis I et al, Bench-to-bedside review: Pulmonary-renal syndromes- an update for the intensivist. *Critical Care* 2007, 213.
9. Salama AD, Levy J, Lightstone L et al, Goodature's disease. *Lancet* 2001, 358, 917-920.
10. Hudson BG, Tryggvason K, Sundaramoorthy M et al, Alport's syndrome, Goodpasture's syndrome, and type IV collagen. *N Engl J Med*, 2003, 348, 2543- 2556.

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Seraline[®] Vaskulitis-3 IgG

Line Immunoassay for detection of IgG- antibodies in systemic vasculitis
in human serum or plasma

REF LIA-003-3 G  20  *In-vitro*-Diagnostic device 

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
Telefon +49 (0) 33767 79110 · Fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Primary systemic vasculitis of small and middle large vessels (1) of unknown etiopathogenesis is associated with anti-neutrophil cytoplasmatic antibodies (ANCA). Because ANCA seem to activate directly neutrophil granulocytes which are accumulated in the lesions of small vessels and thereby initiate an inflammatory destructive process no immunoglobulins or complement components are seen in the vascular lesions (so called pauci-immune vasculitis) (2, 3).

Major antigens of ANCA are:

Proteinase 3 (PR3) (synonym: azurophil granule protein) a member of the serine proteases family, is the main target of cytoplasmic ANCA (cANCA) seen by immunofluorescence and is strongly associated with Granulomatosis with polyangiitis (**GPA**, previously known as Wegener's granulomatosis, 85% anti-PR3 antibody positive) (4).

Myeloperoxidase (MPO) is a strongly cationic protein of the azurophil granula and is involved in the generation of reactive oxygen species (5). It is the main target of perinuclear ANCA (pANCA) which are serological markers for microscopic polyangiitis (**MPA**, 45%), Eosinophilic Granulomatosis with Polyangiitis (**EGPA** previously known as Churg-Strauss syndrome, 60%) and idiopathic glomerulonephritis (65%). MPO-ANCA associated vasculitis is more common in Japanese and Chinese and PR3-ANCA associated vasculitis more common in Europeans (6). In PR3-ANCA associated vasculitis glomerulonephritis is more common and more severe than in patients with MPO-associated vasculitis (7). Although MPO is the main target of pANCA a similar immunofluorescence pattern can be produced by antibodies to elastase, cathepsin and lysozyme. Therefore methods (ELISA or LIA) with purified antigens as targets should be used to determine fine specificity of antibodies classified by immunofluorescence.

Glomerular basement membrane (GBM) is the target of autoantibodies which bind to the non-collagenous domain of the $\alpha 3$ chain of collagen type IV (10). This structure is also present in other membranes. Anti-GBM antibodies detectable in about 90% of patients activate complement after binding. Cellular proteases are attracted and cause membrane disruption and destruction of the filtration barrier. Due to knowledge of the target antigen highly purified or recombinant antigen in ELISAs or LIAs enables a rapid diagnosis and an early onset of adequate therapy which dramatically improves prognosis. The good pasture's syndrome (**GPS**) is extremely rare (one case per 1.000.000 inhabitants per year, 3.8% of patients with pulmonary-renal syndrome in intensive care units) (8) and predominantly affects Caucasians. However, due to rapidly progressive glomerulonephritis it is responsible for about 20% of acute renal failure cases (9). Untreated, mortality is between 75 and 90 %. Genetic (> 80% of patients are positive for HLA DR15 or DR4) and environmental factors (smoking, infections) have been implicated in the pathogenesis. Sometimes the good pasture's syndrome can be superimposed to ANCA positive systemic vasculitis (**GPS/ AAV Overlap**).

Literature:

1. Jenette JC, Falk RJ, Andrassy K. et al, Nomenclature of systemic vasculitides: proposal of an international conference. *Arthritis Rheumatol* 1994, 37, 187-192.
2. Charles LA, Caldas MLR, Falk RJ et al, Antibodies against granule proteins activate neutrophils in vitro. *J Leukocyte Biol* 1991, 50, 539-546.
3. Churg J, Churg A, Idiopathic and secondary vasculitis: a review. *Mod. Pathol.* 1989, 2, 144-160.
4. Cohen Tervaert JW, van der Waude FJ, Fauci, AS et al, Association between active Wegener's Granulomatosis and anticytoplasmic antibodies. *Arch Intern Med* 1989, 149, 2461-2465.
5. Weiss SJ, Tissue destruction by neutrophils. *New Engl. J Med* 1989, 320, 365-376.
6. Ozaki S, ANCA associated vasculitis: Diagnostic and therapeutic strategy. *Allergol. Int.* 2007, 56, 87-96.
7. Chen M, Yu F, Wang SX et al, Renal histology in Chinese patients with anti-myelo-peroxidase autoantibody-positive Wegener's granulomatosis. *Nephrol. Dial. Transplant.* 2007, 22, 139-145.
8. Papisir SA, Manali ED, Kalomenidis I et al, Bench-to-bedside review: Pulmonary-renal syndromes- an update for the intensivist. *Critical Care* 2007, 213.
9. Salama AD, Levy J, Lightstone L et al, Goodature's disease. *Lancet* 2001, 358, 917-920.
10. Hudson BG, Tryggvason K, Sundaramoorthy M et al, Alport's syndrome, Goodpasture's syndrome, and type IV collagen. *N Engl J Med*, 2003, 348, 2543- 2556.

Intended use

The *Seraline*[®] Vaskulitis-3 IgG is a Line Immunoassay (LIA) for detection of autoantibodies of the IgG isotype directed against the following antigens: proteinase 3, myeloperoxidase and glomerular basement membrane in human serum or plasma samples.

Principle of the test

The antigens are plotted on nitrocellulose membranes and blocked to avoid unspecific reactions. After this procedure the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

Step 1

The test strips are incubated in the diluted (1 : 101) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates. The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template. Alternatively the strips can be evaluated with the *Seraline*[®] scan software.

Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 101 (15 µl sample and 1500 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

Kit components

1	TESTSTR	test strips incl. evaluation template	20 strips colour coded: red
2	WIB CONC 5x	wash and incubation buffer 5-fold	70 ml concentrate, transparent bottle, black cap
3	CONJ HRP IgG	anti-human IgG-HRP-conjugate	35 ml , ready to use, IgG , coloured red, transparent bottle, red cap
4	SUBSTR TMB	substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	35 ml , ready to use, black bottle, blue cap
5	INCUTRAY	incubation tray with cover	2

Materials required but not provided

- glassware
- pipettes and tips for 15 µl and 1500 µl
- plastic forceps
- rocking platform (vertical)
- filter paper
- collecting devices for infectious material
- distilled or deionized water

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 20 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

Working steps

1. Incubate the test strips with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **15 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **1.5 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **1.5 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **1.5 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.


Test evaluation

Test validity

The *Seraline*[®] Vaskulitis-3 IgG test strips contain 3 consecutive control lines beneath the test strip number:

- a) Function control (positive reaction with every sample).
- b) Negative control (band intensity must be less than the Cut off control).
- c) Cut off control (band intensity is used for result evaluation of the diagnostic bands).

The test run is valid if the function control and the Cut off control are clearly visible.

Result interpretation

Generally the evaluation of the test has to be done with dried strips. To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet. Alternatively the strips can be evaluated with the *Seraline*[®] scan software. Judgement of the bands is performed according to the following classification:

Result interpretation	IgG
negative	Colour intensity of bands < Cut off control
positive	Colour intensity of bands ≥ Cut off control

Inverse bands (bright bands on dark background) should be evaluated as negative. The respective serum should always be examined using other serological methods.

Performance characteristics

Sensitivity

Clinically defined serum samples have been tested with the *Seraline*[®] Vaskulitis-3 IgG and another commercially available Line Immunoassay. The test results of both assays have been interpreted by comparing every single antigen band.

MPO:

		<i>Seraline</i> [®] Vaskulitis-3 IgG		Agreement: 92%
		positive	negative	
competitor	positive	20	0	
	negative	2	2	

PR3:

		<i>Seraline</i> [®] Vaskulitis-3 IgG		Agreement: 90%
		positive	negative	
competitor	positive	17	0	
	negative	2	0	

GBM:

		<i>Seraline</i> [®] Vaskulitis-3 IgG		Agreement: 88%
		positive	negative	
competitor	positive	7	1	
	negative	0	0	

For an external study serum samples from 117 patients have been tested with another commercially available test. The results have been evaluated for every single antigen and specificity and sensitivity of the comparative test were postulated as 100%.

Antigens	Sensitivity	Specificity	Agreement
MPO	97.2% (35/36)	100% (81/81)	99.1% (116/117)
PR3	100% (24/24)	91.4% (85/93)	93.2% (109/117)
GBM	100% (16/16)	100% (101/101)	100% (117/117)

Serum samples with discrepant results have been retested with a second commercially available test. The results of the 8 discrepant samples have been in concordance with the *Seraline*[®] Vaskulitis-3 IgG results.

Specificity

Sera from 120 blood donors have been investigated. From the results a specificity of 95% has been calculated.

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other Immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only. Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

History of changes

Version	Section	Modifications
2018-03-15	Common advices and precautions	Update

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use