# Seraline<sup>®</sup> Anti-Yersinia-6 IgG Seraline<sup>®</sup> Anti-Yersinia-6 IgA

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Line Immunoassay for detection of IgG or IgA antibodies to *Yersinia enterocolitica*  
in human serum or plasma

 LIA-005-6 G	 20	 LIA-005-6 G-12	 12x 20	 <i>In-vitro</i> -Diagnostikum	
 LIA-005-6 A	 20	 LIA-005-6 A-12	 12x 20	 <i>In-vitro</i> -Diagnostikum	

Seramun Diagnostica GmbH · Sprenthagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · [info@seramun.com](mailto:info@seramun.com)

## Introduction

The genus *Yersinia* comprises gram-negative bacteria belonging to the family Enterobacteriaceae. There are at least three *Yersinia* species which are pathogenic to human: *Y. pestis*, *Y. pseudotuberculosis* and *Y. enterocolitica*. More than 60 serovariants are known from *Y. enterocolitica* (6).

The human pathogenic *Yersinia* species are characterized by a 70 kb plasmid encoding for cytosolic, membraneous and releasing proteins. The main virulence markers are the *Yersinia* outer proteins (Yop) which are released by the bacteria (2, 8). These proteins are highly immunogenic and induce the generation of specific antibodies of the IgA, IgG and IgM isotype (3, 4). IgM and IgA antibody levels decrease within months after infection whereas IgG antibodies may persist for years. In cases of immunopathological complications and chronic forms of *Yersinia* infections IgA antibodies are detectable for a longer period of time (1).

Although no precise data are available the assumed incubation period for *Y. enterocolitica* ranges from 2 to 10 days. The pathogens are transferred through contaminated food (meat and sausages) and less frequently through direct contact. Many mammalian species are known to serve as reservoir hosts, but rodents, domestic and farm animals are of major epidemiological importance (6).

The clinical symptoms of an infection with *Y. enterocolitica* and *Y. pseudotuberculosis* are very similar. In infants an infection induces a self limiting gastroenteritis frequently accompanied by fever and vomiting. Juveniles often react with a lymphadenitis characterized by abdominal pain often mistaken as appendicitis. Immunopathological complications often occur as reactive arthritis (HLA-B27 associated) and erythema nodosum (8).

## Literature:

1. Cremer, J., Putzker, M., Faulde, M., Zöller, L., Immunoblotting of *Yersinia* plasmid-encoded released proteins: A tool for serodiagnosis, *Electrophoresis* 1993, 14: 952-959
2. Iriarte, M., Cornelis, G.R, Identification of SycN, YscX, and YscY, Three new elements of the *Yersinia* Yop Virulon, *J. Bacteriol.* 1999, 181/2: 675-680
3. Heesemann, J., Gross, U., Schmidt, N., Laufs, R., Immunochemical analysis of plasmid-encoded proteins released by enteropathogenic *Yersinia* sp. grown in calcium-deficient media *Infect. Immun.* 1986, 54/2: 561-567
4. Heesemann, J., Grüter, L., Genetic evidence that the outer membrane protein Yop1 of *Yersinia enterocolitica* mediates adherence and phagocytosis resistance to human epithelial cells *FEMS Microbiol. Lett.* 1987, 40, 37-41
5. Kendrick, C.J., Baker, B., Morris, A.J., O'Toole, P.W., Identification of *Yersinia*-infected blood donors by anti-Yop IgA immunoassay, *Transfusion*, 2001, 41: 1365-1372
6. Brandis, H., Köhler, W., Eggers, H.J., Pulverer, G., *Lehrbuch der Medizinischen Mikrobiologie* Gustav Fischer Verlag Stuttgart Jena New York, 1994, 7. Auflage
7. Cornelis, G.R., Boland, A., Boyd, A.P., Geuijen, C., Iriarte, M., Neyt, C., Sory, M.P., Stainier, I. The virulence plasmid of *Yersinia*, an antihost genome *Microbiol. Mol. Biol. Rev.* 1998, 62/4, 1315-1352
8. Paerregaard, Interaction between *Yersinia enterocolitica* and the host with special reference to virulence plasmid encoded adhesion and humoral immunity *Dan. Med. Bull.* 1992, 39/2, 155-172

## Intended use

The *Seraline*<sup>®</sup> Anti-Yersinia-6 IgG / *Seraline*<sup>®</sup> Anti-Yersinia-6 IgA is a Line Immunoassay (LIA) for detection of specific IgG or IgA antibodies against Yersinia outer proteins (Yops) of *Yersinia enterocolitica* in human serum or plasma samples.

## Principle of the test

The antigens are plotted on nitrocellulose membranes and blocked to avoid unspecific reactions. After this procedure the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

### Step 1

The test strips are incubated in the diluted (1 : 101) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

### Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

### Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates. The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template. Alternatively the strips can be evaluated with the *Seraline*<sup>®</sup>scan software.

## Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 101 (15 µl sample and 1500 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

## Kit components

			for 20 determinations	for 12 x 20 determinations
1	<b>TESTSTR</b>	test strips incl. evaluation template	20 strips colour coded: light blue	12 x 20 strips colour coded: light blue
2	<b>WIB CONC 5x</b>	wash and incubation buffer 5-fold	70 ml concentrate, transparent bottle, black cap	12 x 70 ml concentrate, transparent bottle, black cap
3	<b>CONJ HRP IgG</b>	anti-human IgG-HRP-conjugate	35 ml, ready to use, IgG, coloured red, transparent bottle, red cap	12 x 35 ml, ready to use, IgG, coloured red, transparent bottle, red cap
	<b>CONJ HRP IgA</b>	anti-human IgA-HRP-conjugate	IgA, coloured purple, transparent bottle, purple cap	IgA, coloured purple, transparent bottle, purple cap
4	<b>SUBSTR TMB</b>	substrate 3,3',5,5'- Tetramethylbenzidine and hydrogen peroxide	35 ml, ready to use, black bottle, blue cap	12 x 35 ml, ready to use, black bottle, blue cap
5	<b>INCUTRAY</b>	incubation tray with cover	2	-

## Materials required but not provided

- glassware
- pipettes and tips for 15 µl and 1500 µl
- plastic forceps
- rocking platform (vertical)
- filter paper
- collecting devices for infectious material
- distilled or deionized water
- incubation trays (▽ 12 x 20)

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 20 or 12 x 20 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

## Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

### Working steps

1. Incubate the test strips with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **15 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **1.5 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **1.5 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **1.5 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

### Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.


## Test evaluation

### Test validity

The *Seraline*<sup>®</sup> Anti-Yersinia-6 IgG / *Seraline*<sup>®</sup> Anti-Yersinia-6 IgA test strips contain 5 consecutive control lines under the strip number:

- Function control (positive reaction with every sample).
- Conjugate controls IgG and IgA (positive reaction with the corresponding conjugate).  
**Note:** The other conjugate control band may develop a weak, non-specific colour.
- The intensity of the negative control must be less than the Cut off control.
- Cut off control (band intensity is used for result evaluation of the diagnostic bands).

The test run is valid if the function control, the corresponding conjugate control and the Cut off control are clearly visible.

### Result interpretation

Generally the evaluation of the test has to be done with dried strips. To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet. Alternatively the strips can be evaluated with the *Seraline*<sup>®</sup> scan software. Judgement of the bands is performed according to the following classification:

Result interpretation	IgG	IgA
<b>negative</b>	all bands < Cut off control	all bands < cut off-control
<b>borderline</b>	Yop D = Cut off control or 1 Yop band $\geq$ Cut off control <b>(but NOT Yop D)</b>	1 Yop band $\geq$ Cut off control <b>(but NOT Yop D)</b>
<b>positive</b>	Yop D > Cut-off control or 2 bands $\geq$ Cut off control <b>(but NOT Yop D=Cut-off)</b>	Yop D $\geq$ Cut off control or 2 bands $\geq$ Cut off control

### Used Yop antigens

Nomenclature*, MW	Specificity
Yop H, 48 kD	<b>all listed proteins (Yops) are highly specific for human pathogenic Yersiniae</b>
Yop M, 45 kD	
Yop B, 42 kD	
Yop D, 36 kD	
Yop N, 34 kD	
Yop E, 27 kD	

\* The nomenclature for some bands is not consistently used in the literature concerning molecular weight (MW) and specificity, respectively (2, 5, 7).

## Performance characteristics

### Sensitivity

Results from the comparative investigation of serum samples from patients suspected of Yersiniosis using *Seraline*<sup>®</sup> Anti-Yersinia-6 IgG / *Seraline*<sup>®</sup> Anti-Yersinia-6 IgA and another commercial assay.

Detection of IgG (n = 150)		<i>Seraline</i> <sup>®</sup> Anti-Yersinia-6 IgG		
		positive	borderline	negative
Commercial-Line-Immunoassay	positive	130	7	3
	borderline	0	0	0
	negative	0	2	8

Detection of IgA (n = 150)		<i>Seraline</i> <sup>®</sup> Anti-Yersinia-6 IgA		
		positive	borderline	negative
Commercial-Line-Immunoassay	positive	71	4	7
	borderline	2	1	3
	negative	5	4	53

Sensitivity IgG: 98.6%    Sensitivity IgA: 88.6%

### Specificity

Results from the comparative investigation of randomly collected blood donor samples using *Seraline*<sup>®</sup> Anti-Yersinia-6 IgG / *Seraline*<sup>®</sup> Anti-Yersinia-6 IgA and another commercial assay.

Detection of IgG (n = 100)		<i>Seraline</i> <sup>®</sup> Anti-Yersinia-6 IgG		
		positive	borderline	negative
Commercial-Line-Immunoassay	positiv	25	3	4
	grenzwertig	0	3	6
	negativ	1	6	52

Detection of IgA (n = 100)		<i>Seraline</i> <sup>®</sup> Anti-Yersinia-6 IgA		
		positive	borderline	negative
Commercial-Line-Immunoassay	positive	9	1	5
	borderline	0	0	2
	negative	0	6	77

Specificity IgG: 88.1%    Specificity IgA: 92.8%

### Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other Immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

### Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

## Common advices and precautions

**This kit is for *in-vitro* use only.** Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only.

Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

**Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

## History of changes

Version	Section	Modifications
2018-03-15	Common advices and precautions	Update

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use