Serazym[®] Adenovirus

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Serazym[®] Adenovirus

Enzyme immunoassay for detection of *Adenovirus* in faecal samples

REF E-017 ▽ 96 REF E-017-A2 ▽ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Spreehagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Adenovirus is the causative agent of respiratory tract infections, conjunctivitis or enteritis. *Adenovirus* infections are spread via faecal-oral transmission or aerosols (1). Most infections are mild to moderate and do not last longer than one week. *Adenovirus* is responsible for 2 - 8% of respiratory tract infections and 7 - 17% of diarrheal diseases in children (2). Thirty percent of viral diarrheas in immunocompromised patients are caused by *Adenovirus* infections (1). The diagnosis of *Adenovirus* infections is preferably done by direct detection in faecal specimens or swabs. Due to the difficulty of *Adenovirus* culture on tissue or cell culture virus detection so far has been performed by electron microscopy. Meanwhile immunological methods like enzyme immunoassay (ELISA) have been developed for antigen detection (3). ELISA techniques are based on poly- and/or monoclonal antibodies to the protein hexon representing the major part of the virus capsid.

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Intended Use

Serazym® Adenovirus is an *in-vitro*-diagnostic medical device for direct detection of *Adenovirus* in faecal samples.

Principle Of The Test

Serazym® Adenovirus is a one-step enzyme immunoassay on the basis of monoclonal antibodies against an epitope of the capsid protein hexon, common to all human pathogenic *Adenovirus* serotypes. Diluted stool specimens and horseradish peroxidase (HRP) labelled monoclonal anti-Adenovirus-antibodies are dispensed simultaneously into the wells of a microtitration plate coated with monoclonal anti-Adenovirus-antibodies. After an incubation time of 60 min at room temperature (RT) unbound components are removed by a washing step. HRP converts the subsequently added colourless substrate solution of 3,3',5,5'-Tetramethylbenzidine (TMB) within a 10 min reaction time at room temperature protected from light into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / ≥ 620 nm is directly proportional to the specifically bound amount of *Adenovirus*. Considering the cut-off value results are interpreted as positive or negative.

Test Components

| | | | For 96 Wells | For 2x 96 Wells |
|---|-------------------------|--|---|--|
| 1 | WELLS | Microtitration plate coated with monoclonal anti-Adenovirus-antibodies (mouse) | 12 single breakable 8-well strips colour coding violet vacuum-sealed with desiccant | 2x 12 single breakable 8-well strips colour coding violet vacuum-sealed with desiccant |
| 2 | WASHBUF CONC 10x | Wash buffer 10-fold | 100 ml concentrate for 1000 ml solution white cap | 2x 100 ml concentrate for 2x 1000 ml solution white cap |
| 3 | DIL | Sample diluent | 100 ml · ready to use coloured yellow black cap | 2x 100 ml · ready to use coloured yellow black cap |
| 4 | CONTROL + | Positive control <i>Adenovirus</i> antigen Ad 41 (inactivated) | 1.5 ml · ready to use coloured blue red cap | 3.0 ml · ready to use coloured blue red cap |
| 5 | CONTROL - | Negative control <i>Adenovirus</i> negative sample | 1.5 ml · ready to use coloured blue green cap | 3.0 ml · ready to use coloured blue green cap |
| 6 | CONJ HRP | HRP-conjugate HRP-labelled, monoclonal anti-Adenovirus-antibodies (mouse) | 12 ml · ready to use coloured green brown cap | 24 ml · ready to use coloured green brown cap |
| 7 | SUBSTR TMB | Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide | 15 ml · ready to use blue cap | 28 ml · ready to use blue cap |
| 8 | STOP | Stop solution 0.25 M sulphuric acid | 15 ml · ready to use yellow cap | 28 ml · ready to use yellow cap |

Preparation And Storage Of Samples

Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8°C before testing in the ELISA.

Preparation

Quickly thaw frozen samples. Warm samples to room temperature and mix well.

The *Serazym*[®] Adenovirus can be performed with 1 : 6 or 1 : 11 diluted specimens. In case of additional testing of the same sample in the *Serazym*[®] Norovirus, the *Serazym*[®] Campylobacter or the *Serazym*[®] Clostridium difficile Toxin A+B the 1 : 6 dilution is recommended.

Preparation of a 1 : 11 sample dilution:

Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 100 mg (diameter about 2 - 3 mm) of faeces if solid or pipette 100 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Preparation of a 1 : 6 sample dilution:

Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or pipette 200 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Materials Required But Not Provided

Micropipettes · multi-channel pipette or multi-pipette · reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation

Preparation And Storage Of Reagents

Kit size and expiry

One kit is designed for 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least 30 days when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

Dilute samples with sample diluent (3) 1 : 11 or 1 : 6, e.g. 100 mg or 100 µl stool + 1.0 ml (1 : 11) sample diluent (3) or 200 mg or 200 µl stool + 1.0 ml (1 : 6) sample diluent (3).

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Dispense 2 drops (or 75 µl) **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 75 µl **CONTROL +** positive control (4)
75 µl **CONTROL -** negative control (5)
50 µl **diluted sample**, mix gently.
4. Cover plate and incubate for 60 min at RT.
5. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
6. Dispense 2 drops (or 75 µl) **SUBSTR TMB** substrate (7) per well.
7. Incubate for 10 min at RT protected from light.
8. Dispense 2 drops (or 75 µl) **STOP** stop solution (8) per well, mix gently.
9. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction Stop.

Result Interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.20

Samples with OD values equal with or higher than the cut-off are considered positive, samples with OD values below the cut-off are considered negative for *Adenovirus* antigen.

Reference Values

| Serazym® Adenovirus | |
|----------------------------|----------------|
| Positive | \geq Cut-off |
| Negative | $<$ Cut-off |

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters, too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

Performance Characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym*[®] Adenovirus from 8-fold determinations of samples:

| sample | mean OD | standard deviation | CV (%) |
|--------|---------|--------------------|--------|
| 1 | 2.792 | 0.160 | 5.7 |
| 2 | 2.059 | 0.167 | 8.1 |
| 3 | 1.368 | 0.094 | 6.9 |
| 4 | 0.718 | 0.068 | 9.4 |

Inter-assay coefficient of variation (CV) in the *Serazym*[®] Adenovirus in 6 different test runs from 8-fold determinations of samples:

| sample | mean OD | standard deviation | CV (%) |
|--------|---------|--------------------|--------|
| 1 | 1.850 | 0.107 | 5.8 |
| 2 | 1.057 | 0.069 | 6.5 |
| 3 | 0.574 | 0.042 | 7.3 |
| 4 | 0.312 | 0.030 | 9.6 |

Lower detection limit

The lower detection limit of *Adenovirus* antigen in the *Serazym*[®] Adenovirus was determined by titration of purified *Adenovirus* antigen (hexon). Lower detection limit: 6 ng / ml

Specificity and sensitivity

A total of 330 stool samples were investigated in parallel in the *Serazym*[®] Adenovirus and in another commercially available ELISA.

| | comparative ELISA positive | comparative ELISA negative |
|--|----------------------------|----------------------------|
| <i>Serazym</i>[®] ELISA positive | 55 | 1 |
| <i>Serazym</i>[®] ELISA negative | 2 | 272 |

Specificity: 99.6% · Sensitivity: 96.5%

Cross reactivity

Stool samples positive for one of the subsequent pathogens have been tested with the *Serazym*[®] Adenovirus and showed no cross reactivity:

Rotavirus (n = 10), *Astrovirus* (n = 8), *Norovirus* (n = 31), *Clostridium difficile* (n = 11), *Campylobacter jejuni* (n = 7), *Campylobacter coli* (n = 1), *Salmonella enteritidis* (n = 18), *Giardia lamblia* (n = 1).

Negative stool specimens have been spiked with $\geq 10^8$ colony forming units of the following microorganisms and tested negative with the *Serazym*[®] ELISA (OD 450 / 620 nm < Cut-Off):

| | |
|-------------------------------|--------------|
| <i>Aeromonas hydrophila</i> | (ATCC 7966) |
| <i>Bacillus cereus</i> | (ATCC 11778) |
| <i>Bacillus subtilis</i> | (ATCC 6633) |
| <i>Bacteroides fragilis</i> | (ATCC 25285) |
| <i>Candida albicans</i> | (ATCC 10231) |
| <i>Campylobacter coli</i> | (ATCC 33559) |
| <i>Campylobacter jejuni</i> | (ATCC 33291) |
| <i>Citrobacter freundii</i> | (ATCC 8090) |
| <i>Clostridium sordellii</i> | (ATCC 9714) |
| <i>Enterobacter aerogenes</i> | (ATCC 13048) |
| <i>Enterobacter cloacae</i> | (ATCC 13047) |
| <i>Enterococcus faecalis</i> | (ATCC 29212) |
| <i>Escherichia coli</i> | (ATCC 25922) |

| | |
|---|-------------------|
| <i>Klebsiella pneumoniae</i> | (ATCC 13883) |
| <i>Peptostreptococcus anaerobius</i> | (ATCC 27337) |
| <i>Proteus vulgaris</i> | (ATCC 8427) |
| <i>Pseudomonas aeruginosa</i> | (ATCC 10145) |
| <i>Salmonella enterica</i> Serovar <i>enteritidis</i> | (ATCC 13076) |
| <i>Salmonella enterica</i> Serovar <i>typhimurium</i> | (ATCC 14028) |
| <i>Shigella flexneri</i> | (ATCC 12022) |
| <i>Shigella sonnei</i> | (ATCC 25931) |
| <i>Staphylococcus aureus</i> | (ATCC 25923) |
| <i>Staphylococcus epidermidis</i> | (ATCC 12228) |
| <i>Vibrio parahaemolyticus</i> | (ATCC 17802) |
| <i>Vibrio cholerae</i> | clinical isolate |
| <i>Yersinia enterocolitica</i> Serotyp <i>O3, O9</i> | clinical isolates |

Interference

None of the following substances added to positive and negative stool samples showed a significant impact on the test result:

Barium sulfate (5%), Buscopan® (2 mg/ml), Cyclamate (5%), Diclofenac (2 mg/ml), Hemoglobine human (5 mg/ml), Blood human (1,25%), Hylak® N (5%), Iberogast® (5%), Immodium® akut duo (0.2/12.5 mg/ml), Loperamide (0.2 mg/ml), Metronidazole (2 mg/ml), Mucin (5 mg/ml), Nexium® (2 mg/ml), Palmitic acid (20%), Pentofuryl® (2 mg/ml), Pepto-Bismol (1 mg/ml), Perenterol (2.5 mg/ml), Rennie® (8 mg/ml), Simagel® (2 mg/ml), Stearic acid (20%), Vancomycin (2 mg/ml).

Common Advices and Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*® stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

History of Changes

| Version | Section | Modifications |
|------------|--------------------------------|--------------------------------------|
| 2020-11-04 | Common Advices and Precautions | Update |
| | History of Changes | Adding section "History of Changes " |
| | Test Components | Update |

Incubation Scheme *Serazym*[®] Adenovirus (E-017)

1.  2 drops (or 75 µl) **CONJ HRP** (6)
 +
 pipette
 75 µl **CONTROL +** (4)
 75 µl **CONTROL -** (5)
 50 µl **diluted stool sample**, mix gently
 60 min incubation (room temperature)

 5 x wash with wash solution

2.  2 drops (or 75 µl) **SUBSTR TMB** (7)
 10 min incubation (room temperature) protected from light

3.  2 drops (or 75 µl) **STOP** (8)

Read OD at 450 / ≥ 620 nm

| | | | | |
|--|---|--|---|---------------------------|
|  Manufacturer |  Date of manufacture |  Use by | LOT Batch code | REF Catalog number |
|  Keep away from sunlight |  Temperature limits |  Biological risks |  Do not reuse | |
|  Consult instructions for use |  Caution | IVD <i>In-vitro</i> -diagnostic medical device |  Contains sufficient for <n> tests | |