Serazym[®] *Cryptosporidium parvum*

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Literatur:

1. Current, William L. and Lynne s. Garcia (1991): Cryptosporidiosis. *Clinical Microbiology Reviews* Vol. 4 (3): p325-358.
2. Clark, Douglas P. (1999): New Insights into Human Cryptosporidiosis. *Clinical Microbiology Reviews* Vol. 12(4): p554-563.
3. Parisi, Mary T. and Philip M. Tierno, JR (1995): Evaluation of new rapid commercial enzyme immunoassay for detection of *Cryptosporidium* Oocysts in untreated stool specimens. *Journal of Clinical Microbiology* Vol. 33(7): p1963-1965.
4. Bialek, R. et al. (2002): Comparison of Fluorescence, Antigen and PCR Assays to Detect *Cryptosporidium parvum* in Faecal Specimens. *Diagnostic Microbiology And Infectious Disease* Vol. 43(4): p283-288.

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Serazym[®] Cryptosporidium parvum

Enzyme immunoassay for detection of *Cryptosporidium parvum* in faecal samples

REF E-039-A ▼ 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostics GmbH · Spreehagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

The genus *Cryptosporidium* belongs to the phylum Apicomplexa, protozoan parasites characterized by secretory cell organelles enabling the invasion of host cells. Among different human pathogenic species *Cryptosporidium parvum* represents the most important cause of diarrhea in man but also in several animal species (Zoonosis) (1, 2). The transmission of *Cryptosporidium* is effected by food or drinking water contaminated with faeces. Person to person transmission has also been described. After ingestion of the infective stage, the oocyst, the vegetative motile sporozoites will be released in the small intestine, where they invade into the enterocytes and stay in a parasitophorous vacuole beneath the microvilli. The sporozoites multiply by bisection (asexually) as well as sexually. In the course of sexual multiplication micro- and macrogamonts emerge and fuse to form the infective oocyst. The oocysts can be divided into thin-wall and thick-wall oocysts: the thin-wall oocysts are responsible for autoinfections within the same host whereas the thick-wall oocysts are excreted with the faeces and are therefore responsible for transmission to other hosts. The clinical picture of cryptosporidiosis is usually watery self-limiting diarrhea. However immunocompromised persons can develop prolonged and life-threatening diarrhea. In rare cases extraintestinal infections of e.g. the respiratory tract or bile-duct have been reported (1, 2). *Cryptosporidium* infections are usually diagnosed by direct pathogen detection from faecal specimens either by microscopy of dyed smears, immune fluorescence microscopy or by immunological methods like ELISA. Nucleic acid amplification techniques (PCR) are also established but still restricted to special laboratories. Enzyme immunoassays used for routine diagnostics are based on mono- or polyclonal antibodies directed to not further characterized antigens derived from the oocyst stage of the parasite (3, 4).

References:

1. Current, William L. and Lynne s. Garcia (1991): Cryptosporidiosis. Clinical Microbiology Reviews Vol. 4 (3): p325-358.
2. Clark, Douglas P. (1999): New Insights into Human Cryptosporidiosis. Clinical Microbiology Reviews Vol. 12(4): p554-563.
3. Parisi, Mary T. and Philip M. Tierno, JR (1995): Evaluation of new rapid commercial enzyme immunoassay for detection of *Cryptosporidium* Oocysts in untreated stool specimens. Journal of Clinical Microbiology Vol. 33(7): p1963-1965.
4. Bialek, R. et al. (2002): Comparison of Fluorescence, Antigen and PCR Assays to Detect *Cryptosporidium parvum* in Faecal Specimens. Diagnostic Microbiology And Infectious Disease Vol. 43(4): p283-288.

Intended Use

***Serazym*[®] *Cryptosporidium parvum* is an *in-vitro*-diagnostic medical device for direct detection of *Cryptosporidium parvum* in faecal samples.**

Principle Of The Test

The *Serazym*[®] *Cryptosporidium parvum* is based on monoclonal antibodies to *Cryptosporidium parvum* specific proteins. Within the first incubation step (60 min, room temperature, RT) diluted stool specimens as well as positive and negative controls react with the solid-phase adsorbed antibodies. Unbound components are removed by a subsequent washing step. Within the second incubation step (30 min, RT) solid-phase bound immune complexes react with the horseradish (HRP)-labelled antibodies of the conjugate. Unbound reagents are separated from the solid-phase antibody-antigen-antibody immune complexes by a further washing step. HRP converts the subsequently added colourless chromogenic substrate solution into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells after 10 min incubation at RT turning the solution from blue to yellow. The optical density (OD) is read at 450 / \geq 620 nm with a microplate photometer.

Test Components

For 96 Wells

1	WELLS	Microtitration plate coated with monoclonal anti-C. parvum antibodies (mouse)	12 single breakable 8 well strips colour coding black vacuum sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control <i>Cryptosporidium parvum</i> Oocysts, inactivated	2.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control <i>Cryptosporidium parvum</i> negative sample	2.0 ml · ready to use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled monoclonal anti-C. parvum antibodies (mouse)	15 ml · ready to use coloured green violet cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap

Preparation And Storage Of Samples

Collection And Storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at -20°C. Repeated freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8°C before testing. Faecal samples that are already diluted in transportation media (PARA-PAK PLUS, PARA-PAK PLUS SAF) are also suitable in the *Serazym*[®] *Cryptosporidium parvum*.

Preparation

Untreated samples: Warm samples to room temperature and mix well. Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or 200 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Conserved samples, prediluted in transportation media: Mix sample thoroughly and apply 100 µl without any further dilution in the *Serazym*[®] *Cryptosporidium parvum*.

Materials required but not provided

Micropipettes · multi-channel pipette or multi-pipette · reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation

Preparation And Storage Of Reagents

Kit size and expiry

One kit is designed for 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least two months, provided proper storage. The ready to use wash solution can be used for at least one month when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

Dilute samples with sample diluent (3) 1 : 6, e.g. 200 mg or 200 μ l faeces + 1.0 ml sample diluent (3).

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that residual fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CONTROL +** positive control (4)
100 µl **CONTROL -** negative control (5)
100 µl **diluted sample**.
3. Cover plate and incubate for 60 min at RT.
4. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
5. Dispense 3 drops (or 100 µl) **CONJ HRP** HRP-conjugate (6) per well.
6. Cover plate and incubate for 30 min at RT.
7. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
8. Dispense 3 drops (or 100 µl) **SUBST TMB** substrate (7) per well.
9. Incubate for 10 min at RT protected from light.
10. Dispense 3 drops (or 100 µl) **STOP** stop solution (8) per well, mix gently.
11. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

Result Interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.10

Samples with OD values equal with or higher than the cut-off are considered positive, samples with OD values below the cut-off are considered negative for *Cryptosporidium parvum* antigen.

Reference Values

Serazym® Cryptosporidium parvum	
Positive	\geq Cut-off
Negative	$<$ Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters, too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid, if:

- the mean OD of the negative control is ≤ 0.20 (manual performance)
 ≤ 0.30 (automatic performance)
- the mean OD of the positive control is: ≥ 0.80

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure


There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate the absorbance of a sample with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the diluted sample can cause false negative as well as false positive results. A negative test result not necessarily excludes a *Cryptosporidium parvum* infection. Inhomogeneous antigen distribution in the sample can cause false negative results. Therefore in case of a negative test result but clinical suspicion of an infection a further sample of the respective person should be investigated. A final interpretation of the test result should always consider clinical findings as well.

Automatic processing

Performing the *Serazym*[®] *Cryptosporidium parvum* on fully automated microplate processors (e.g. DS2, DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x - 8x.

Correlation: manual - automatic processing

A panel of 90 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.976$.

Reference Values

Precision

Intra-Assay coefficient of variation (CV) in the *Serazym*[®] *Cryptosporidium parvum* from 8-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.732	0.081	4.98
2	0.855	0.089	11.14
3	0.697	0.036	5.54
4	0.129	0.005	4.33

Inter-Assay coefficient of variation (CV) in the *Serazym*[®] *Cryptosporidium parvum* in 8 different test runs from twofold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	1.704	0.066	3.89
2	0.969	0.098	10.21
3	0.519	0.023	4.45
4	0.159	0.012	7.67

Lower detection limit

The lower detection limit of *Cryptosporidium parvum* antigen in the *Serazym*[®] *Cryptosporidium parvum* was determined by titration of oocysts diluted in a negative stool sample: 1.6 x 10⁵ Oocysts per ml stool suspension.

Specificity and sensitivity

Specificity

A sample panel of in all n = 327 stool specimens from a routine microbiological laboratory was tested in the *Serazym*[®] *Cryptosporidium parvum*: Negative: n = 325 · Positive: n = 2 · Specificity: 99.4%
The two ELISA positive samples were confirmed true positive by direct immunofluorescence:
Specificity amended: 100% for this sample panel.

Sensitivity

Faecal samples designated *Cryptosporidium* positive were investigated by *Serazym*[®] *Cryptosporidium parvum* in comparison to two commercially available ELISAs.

n = 34	ELISA 1 positive	ELISA 1 negative
<i>Serazym</i> [®] ELISA positive	32	0
<i>Serazym</i> [®] ELISA negative	1*	1

Sensitivity compared to ELISA 1: 96.9%

* The *Serazym*[®] ELISA negative and ELISA 1 positive sample was negative by direct immunofluorescence test.

n = 19	ELISA 2 positive	ELISA 2 negative
<i>Serazym</i> [®] ELISA positive	18	0
<i>Serazym</i> [®] ELISA negative	0	1

Sensitivity compared to ELISA 2: 100%

Cross reactivity

Stool samples positive for one of the subsequent pathogens have been tested with the *Serazym*[®] *Cryptosporidium parvum* and did not show any cross reactivity:

Adenovirus (n = 2), *Rotavirus* (n = 7), *Astrovirus* (n = 3), *Norovirus* (n = 5), *Clostridium difficile* (n = 2), *Campylobacter jejuni* (n = 2), *Helicobacter pylori* (n = 7), *Salmonella ssp.* (n = 9), *Ancylostoma duodenale* (n = 1), *Ascaris lumbricoides* (n = 2), *Blastocystis hominis* (n = 2), *Giardia lamblia* (n = 4), *Entamoeba dispar/histolytica* (n = 5).

Negative stool samples have been spiked with ≥ 10⁸ colony forming units per ml of the following microorganisms and were tested negative with the *Serazym*[®] ELISA (OD 450 / 620nm < Cut-Off):

<i>Aeromonas hydrophila</i>	(ATCC 7966)
<i>Bacillus cereus</i>	(ATCC 11778)
<i>Bacillus subtilis</i>	(ATCC 6633)
<i>Bacteroides fragilis</i>	(ATCC 25285)
<i>Candida albicans</i>	(ATCC 10231)
<i>Campylobacter coli</i>	(ATCC 33559)
<i>Campylobacter jejuni</i>	(ATCC 33291)
<i>Citrobacter freundii</i>	(ATCC 8090)
<i>Clostridium sordellii</i>	(ATCC 9714)
<i>Enterobacter aerogenes</i>	(ATCC 13048)
<i>Enterobacter cloacae</i>	(ATCC 13047)
<i>Enterococcus faecalis</i>	(ATCC 29212)
<i>Escherichia coli</i>	(ATCC 25922)
<i>Klebsiella pneumoniae</i>	(ATCC 13883)

<i>Peptostreptococcus anaerobius</i>	(ATCC 27337)
<i>Proteus vulgaris</i>	(ATCC 8427)
<i>Pseudomonas aeruginosa</i>	(ATCC 10145)
<i>Salmonella enterica Serovar enteritidis</i>	(ATCC 13076)
<i>Salmonella enterica Serovar typhimurium</i>	(ATCC 14028)
<i>Shigella flexneri</i>	(ATCC 12022)
<i>Shigella sonnei</i>	(ATCC 25931)
<i>Staphylococcus aureus</i>	(ATCC 25923)
<i>Staphylococcus epidermidis</i>	(ATCC 12228)
<i>Vibrio parahaemolyticus</i>	(ATCC 17802)
<i>Vibrio cholerae</i>	Clinical isolate
<i>Yersinia enterocolitica Serotyp O3, O9</i>	Clinical isolates

Interference

None of the following substances added to positive and negative stool samples showed a significant impact on the test result:

Barium sulfate (5%), Buscopan® (2 mg/ml), Cyclamate (5%), Diclofenac (2 mg/ml), Hemoglobine human (5 mg/ml), Blood human (5%), Hylak® N (5%), Iberogast® (5%), Immodium® akut duo (0.2/12.5 mg/ml), Loperamide (0.2 mg/ml), Metronidazole (2 mg/ml), Mucin (5 mg/ml), Nexium® (2 mg/ml), Palmitic acid (20%), Pentofuryl® (2 mg/ml), Pepto-Bismol (1 mg/ml), Perenterol (2.5 mg/ml), Rennie® (8 mg/ml), Simigel® (2 mg/ml), Stearic acid (20%), Vancomycin (2 mg/ml).

Common Advices And Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*® stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

History of Changes

Version	Section	Modifications
2021-01-26	Common Advices and Precautions, symbolbox	Update
	History of Changes	Adding section "History of Changes"
	Intended use	Correction

Incubation Scheme *Serazym*[®] *Cryptosporidium parvum* (E-039-A)

1.  pipette
 100 µl **CONTROL +** (4)
 100 µl **CONTROL -** (5)
 100 µl **diluted stool sample**

 60 min incubation (room temperature)

 5 x wash with wash solution

2.  3 drops (or 100 µl) **CONJ HRP** (6)

 30 min incubation (room temperature)

 5 x wash with wash solution

3.  3 drops (or 100 µl) **SUBSTR TMB** (7)

10 min incubation (room temperature, protected from light)

4.  3 drops (or 100 µl) **STOP** (8)

Read OD at 450 / ≥ 620 nm

 Manufacturer	 Date of manufacture	 Use by	 Batch code	 Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	 <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	