Serazym[®] ANA screen

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Lagertemperatur

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Serazym[®] ANA screen

Enzyme immunoassay for detection of IgG antibodies in connective tissue diseases
in human serum

REF E-065 ▽ 96 IVD *In-vitro*-diagnostic device CE

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Connective tissue diseases are systemic inflammatory rheumatic diseases usually characterized by a chronic course (systemic connective tissue diseases) with overlapping symptoms (Overlap-Syndromes). Connective tissue diseases include:

- Systemic lupus erythematosus (SLE) and subsets
- Sjögren's syndrome
- Systemic sclerosis
- Idiopathic (autoimmune) myositis
- Mixed connective tissue disease (MCTD or Sharp syndrome)
- Overlap syndromes

dsDNA-antibodies:

IgG antibodies to dsDNA are marker antibodies and ACR-criterion for systemic lupus erythematoses (SLE). They are regarded as activity and prognostic marker of SLE. The detection frequency varies in dependence of the activity of the disease and organ manifestation: patients suffering from

- active SLE with renal involvement >95%,
- active SLE without renal involvement 50-70%,
- inactive SLE <40%.

In sera from patients with rheumatoid arthritis, juvenile chronic arthritis, Sjögren's syndrome, sclerosis and other diseases antibodies to dsDNA may also be detected temporarily with low titers. Due to their higher specificity in comparison to anti-nucleosome antibodies anti-dsDNA antibodies undoubtedly are of major importance for the diagnosis of SLE.

U1-snRNP antibodies:

Anti-U1-snRNP antibodies are marker antibodies and diagnostic criterion of mixed connective tissue disease (MCTD) with a sensitivity of 100% and a specificity of 98 % for high titers of anti-U1-RNP antibodies and in the absence of antibodies to Sm and dsDNA.

Antibodies to U1-RNP with often low titers are found in 13-32% of patients with SLE and in 10% of patients with systemic sclerosis.

Using the three RNP-proteins A, C and 68 kD *Serazym*[®] ANA screen detects all U1-RNP antibodies.

Sm antibodies:

Sm-antibodies are diagnostic marker and ACR-criterion of SLE with a specificity of 99%, but a sensitivity of only 10-15% in SLE patients of caucasian origin (and 30% to over 40% in asiatic patients).

Since an association to several organ manifestations (kidneys, central nervous system) has been proven, antibodies to Sm are considered as prognostic marker of SLE.

Ro/SS-A antibodies:

Ro/SS-A antibodies are predominantly detected in patients with Sjögren's syndrome and the different forms of lupus erythematoses. They serve as diagnostic marker and are part of the classification criteria of primary and secondary Sjögren's syndrome with a sensitivity of 96% and 80% respectively. Ro/SS-A antibodies are considered as markers for the early diagnosis of Sjögren's syndrome since they may occur years before clinical manifestation.

Anti-Ro/SS-A antibodies are detectable in 40-60% of patients who suffer from SLE. They may occur together with antibodies typical for SLE (anti-dsDNA, anti-Sm antibodies) or as isolated antibodies indicating a relatively mild form of SLE.

Antibodies to Ro/SS-A are a diagnostic marker of subacute cutaneous lupus erythematosus (SCLE) and detectable in 90-100% of the cases.

Ro/SS-A antibodies are detectable in nearly 100% of neonatal lupus erythematosus (NLE) cases. The coincidence of anti- Ro/SS-A and La/SS-B antibodies is associated with congenital heart block (CHB).

Ro/SS-A antibodies are additionally detectable in patients with rheumatoid arthritis (5-8%) and systemic sclerosis (9%).

La/SS-B antibodies:

Antibodies to La/SS-B serve as important diagnostic marker and belong to the classification criteria of Sjögren's syndrome. The diagnostic sensitivity is about 70% for the primary Sjögren's syndrome and about 50% for the secondary Sjögren's syndrome. In case of simultaneous detection of La/SS-B and Ro/SS-A antibodies diagnostic specificity for Sjögren's syndrome is higher in comparison to isolated detection of Ro/SS-A antibodies. Anti-La/SS-B antibodies almost always occur together with antibodies to Ro/SS-A.

La/SS-B antibodies are regarded as early diagnostic markers of Sjögren's syndrome, since they may be detectable years before clinical symptoms become evident. Antibodies to La/SS-B are detectable in 25% of systemic lupus erythematosus (SLE) cases, in 70% of neonatal lupus erythematosus (NLE) cases and in 80% of subacute cutaneous lupus erythematosus (SCLE) cases.

Scl 70 antibodies:

Scl 70- or Topoisomerase I antibodies are diagnostic marker of systemic sclerosis with a sensitivity of

- 18-30% in general
- 10-15% in limited sclerosis (CREST-syndrome)
- 40-65% in diffuse forms of systemic sclerosis

The specificity reaches nearly 100%.

Patients who develop anti-Scl 70 antibodies usually suffer from a more severe disease in comparison to patients with anti-centromer antibodies. Antibodies to Scl 70 may occur already years before specific symptoms of sclerosis are detectable.

CENP-B-antibodies:

CENP-B antibodies are considered as diagnostic marker of the limited form of systemic sclerosis (sensitivity 50-70%). Patients who develop antibodies to CENP-B usually show a rather mild course of sclerosis. Antibodies to CENP-B may occur already years before specific symptoms of sclerosis are detectable.

CENP-B antibodies are detectable in only 10-30% of patients with primary biliary cirrhosis (PBC) and very rarely in patients suffering from other connective tissue diseases.

Jo-1 antibodies:

Antibodies to Jo-1 represent a diagnostic marker for idiopathic myositis with a diagnostic specificity of nearly 100%. Diagnostic sensitivity reaches 18-46% for polymyositis and 25% for dermatomyositis. About 60% of patients with detectable levels of antibodies to Jo-1 suffer from fibrosing alveolitis.

Jo-1 antibodies serve as prognostic marker of myositis, since these patients often develop a severe course of the disease.

The cytoplasmatic localization of the target autoantigen (histidyl-tRNA-synthetase) actually excludes anti-Jo-1 antibodies from anti-nuclear antibodies in the original meaning of the word.

References:

1. Conrad, K., Schößler, W., Hiepe, F., Fritzler, M.J.; Autoantibodies in Systemic Autoimmune Diseases; Pabst Science Publishers Lengerich, 2015
2. Peter, J.B., Shoenfeld Y. (Eds); Autoantibodies; Elsevier Amsterdam 1996
3. Tan, E.M.; Antinuclear antibodies: diagnostic markers for autoimmune disease and probes for cell biology; Adv. Immunol. 44, 1989: 93-151

Intended use

Serazym[®] ANA screen is an *in-vitro* diagnostic device for simultaneous determination of IgG antibodies against the following nuclear and cytoplasmatic antigens: dsDNA, RNP, Sm, Ro/SS-A/60kD, La/SS-B, Scl-70, CENP and Jo-1 in human serum.

Principle of the test

In the first incubation step diluted samples, ready-to-use calibrator and negative control react with the solid-phase adsorbed highly purified antigens. Unbound components are removed after 60 minutes incubation at room temperature (RT, 18...25 °C) by aspirating the samples followed by a washing step.

In the second incubation step horseradish peroxidase (HRP)-labelled anti-human-IgG antibodies bind to sample antibodies. Unbound conjugate is removed after 30 minutes incubation at RT by aspirating followed by a washing step.

In the third step the HRP converts the subsequently added colourless substrate solution of hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB) within a 15 min reaction time at RT into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow.

The absorbance of the solution read at 450/≥620 nm (OD) is directly proportional to the amount of specifically bound antibodies. The cut-off is established by multiplying the absorbance of the calibrator with a corresponding factor. Patient ratio (binding index) are calculated by dividing the absorbance of the sample by the calculated cut-off.

Test components

		For 96 Wells
1	WELLS	Microtitration plate coated with 8 different antigens / well
		12 single breakable 8-well strips vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold
		100 ml concentrate for 1000 ml solution white cap
3	DIL	Sample diluent
		100 ml · ready-to-use coloured red black cap
4	CAL	Calibrator Diluted serum samples factor: see certificate of analysis enclosed
		1.0 ml ready-to-use coloured blue red cap
5	CONTROL -	Negative control Diluted serum sample
		1.0 ml · ready-to-use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled, anti-human-IgG antibodies
		15 ml · ready-to-use coloured red red cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide
		15 ml · ready-to-use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid
		15 ml · ready-to-use yellow cap
9		Adhesive film
		2

Preparation and storage of samples

Collection and storage

Blood is taken by venipuncture. Serum is separated after clotting by centrifugation. Lipaemic, haemolytic and contaminated samples should not be used. Repeated freezing and thawing should be avoided. The samples may be kept at 2...8 °C for up to two days. Long-term storage requires -20 °C.

Preparation before use

Allow samples to reach RT prior to assay. Take care to agitate samples gently in order to ensure homogeneity.

Note: Patient samples have to be diluted 1 : 101, e.g. 10 µl sample + 1000 µl sample diluent (3).

Materials required but not provided

Adjustable one-channel micropipette 0.100 - 1.000 ml and 0.010 - 0.100 ml · Adjustable 8-channel micropipette 0.050 - 0.200 ml · Pipette tips · Reagent containers for dispensing with 8-channel pipette · Graduated measuring flasks 10 ml and 100 ml · Test tubes 2.0 ml for sample dilution · Microtitration plate washer (automatic or hand wash head) · Microtitration plate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · Distilled or deionized water · Stop-watch

Preparation and storage of reagents

Kit size and expiry

The *Serazym*[®] ANA screen has been designed for 12 x 8 determinations.

The complete kit with unopened reagent bottles and microtitration strips is stable until the expiry date printed on the kit box in case of storage at 2...8 °C. After opening all kit components are stable for at least 2 months, provided proper storage. The ready-to-use wash solution can be used for at least one month when stored at 2...8°C.

Reagent preparation

The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For example: 10 ml wash buffer (2) + 90 ml distilled or deionized water.

Besides wash buffer all components of the device are ready-to-use.

Assay procedure

Allow all components to reach room temperature prior to use in the assay.

Dilute patient samples with sample diluent (3) 1 : 101 (v/v), e.g. 10 µl sample + 1000 µl sample diluent (3).

The calibrator and the negative control are ready-to-use.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CAL** calibrator (4)
100 µl **CONTROL -** negative control (5)
100 µl **diluted samples** resp.
into the intended wells.
3. Cover plate (9) and incubate for 60 min at RT.
4. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
5. Dispense 100 µl **CONJ HRP** HRP-conjugate (6) per well.
6. Cover plate (9) and incubate for 30 min at RT.
7. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
8. Dispense 100 µl **SUBSTR TMB** substrate (7) per well.
9. Incubate for 15 min at RT protected from light.
10. Dispense 100 µl **STOP** stop solution (8) per well, mix gently.
11. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

Result interpretation

Results can be interpreted qualitatively by calculating a cut-off or semi-quantitatively by calculating the binding index (BI) for each sample.

Therefore the cut-off value has to be determined.

The cut-off is established by multiplying the absorbance of the calibrator with a corresponding calculation factor.

Absorbance _{Calibrator} **x factor = absorbance** _{cut-off}

The calculation factor is lot-specific and stated on the certificate of analysis of the kit.

For calculation of the binding index (BI) the absorbance of the sample is divided by the calculated cut-off value.

BI sample = absorbance _{sample} / **absorbance** _{cut-off}

Reference values

Serazym® ANA screen	Binding index
positive	≥ 1.0
negative	< 1.0

Test validity

The test run is valid if:

- absorbance of calibrator CAL is ≥ 0.80
- absorbance of negative control CONTROL- is ≤ 0.35

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation volumes, times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the method


Any clinical diagnosis should not be based on the results of *in-vitro* diagnostic methods alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis. False results may be caused by cross contaminations of the kit reagents and samples, bacterial or fungal contaminations of reagents and/or samples, incorrect washing and incorrect incubation times. In case of background problems when using automatic microplate ELISA systems it may be recommendable to wash wells 5 times (instead of 3 times) in every wash cycle.

Performance characteristics

Due to the lack of an international reference preparation results are interpreted as binding index (BI).

Frequency distribution

The frequency distribution of 128 unselected blood donor samples in the *Serazym*[®] ANA screen is illustrated:

Precision

Intra-assay coefficient of variation (CV [%]) in the *Serazym*[®] ANA screen from 6-fold determinations of samples:

sample	BI	CV [%]
1	4.87	3.26
2	3.76	2.01
3	2.00	1.80
4	1.31	2.71

Inter-assay coefficient of variation (CV [%]) in the *Serazym*[®] ANA screen in 16 different test runs from 2-fold determination of samples:

sample	BI	CV [%]
A	4.49	5.21
B	3.66	4.07
C	2.22	3.48
D	1.74	5.16

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for diluted reagents. Do not use or mix reagents from different lots, damaged packages or bottles or reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8 °C before use.

Do not use or mix reagents from different lots except for sample diluent, wash buffer, substrate and stop solution.

Some of the reagents contain small amounts of Thimerosal (< 0.01% w/v) and Kathon (1.0% v/v) as preservative. They may neither be swallowed nor come into contact with skin or mucous membranes. In case of contact rinse with plenty of water. Handle all samples as potentially infectious. Since the kit contains potentially hazardous materials, the following precautions must always be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

Incubation scheme *Serazym*[®] ANA screen

- | | | | |
|----|---|----------|--|
| 1. |  | 100 µl | CAL (4) |
| | | 100 µl | CONTROL - (5) |
| | | 100 µl | diluted samples resp. |
| |  | 60 min | incubation (room temperature) |
| | | 3 x wash | with wash solution |
| 2. |  | 100 µl | CONJ HRP (6) |
| | | 30 min | incubation (room temperature) |
| |  | 3 x wash | with wash solution |
| 3. |  | 100 µl | SUBSTR TMB (7) |
| | | 15 min | incubation (room temperature) protected from light |
| 4. |  | 100 µl | STOP (8) |

Read OD at 450 / \geq 620 nm

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use