Serazym[®] Anti-PR3 (cANCA) IgG

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Serazym[®] Anti-PR3 (cANCA) IgG

Enzyme immunoassay for detection of IgG antibodies to Proteinase 3 (PR3)
in human serum or plasma

r E-088 n 96 v *In-vitro*-diagnostic medical device c

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Anti-neutrophil cytoplasmic antibodies (ANCA) serve as important diagnostic markers for systemic vasculitis and can be detected by indirect immunofluorescence (IIF) on ethanol-fixed neutrophils. Depending from the staining pattern, ANCA are divided into c (cytoplasmic) and p (perinuclear) ANCA.

A positive indirect immunofluorescence test (IIF) should be confirmed by antigen specific ELISA. Proteinase 3 (PR3), a protein of the azurophil granules of neutrophilic granulocytes has been confirmed as the primary target antigen of cANCA. cANCA are detectable in more than 90 % of sera from patients suffering from Wegener's granulomatosis (WG). The anti-PR3 antibody titer correlates with the activity of the disease and may support optimization of therapy. Antibodies to PR3 may be detectable less frequently in sera from patients with other vasculitis forms like microscopic polyangiitis, Churg-Strauss syndrome or panarteritis nodosa.

References:

1. Conrad, K., Schöbler, W., Hiepe, F., Fritzler, M.J.: Autoantibodies in Systemic Autoimmune Diseases. Pabst Science Publishers Lengerich, 3rd edition 2015
2. Peter, J.B., Shoenfeld Y. (Eds): Autoantibodies. Elsevier Amsterdam, 1996
3. Tan, E.M.; Antinuclear antibodies: diagnostic markers for autoimmune disease and probes for cell biology. Adv. Immunol. 44, 1989: 93-151

Intended use

Serazym® Anti-PR3 (cANCA) IgG is an *in-vitro* diagnostic medical device for quantitative determination of IgG antibodies to proteinase 3 in human serum or plasma.

Principle of the test

In the first incubation step diluted samples, ready-to-use calibrators and positive control react with the solid-phase adsorbed antigen. Unbound components are removed after 60 minutes incubation at room temperature (RT, 18...25 °C) by aspirating the samples followed by a washing step.

In the second incubation step horseradish peroxidase (HRP)-labelled anti-human-IgG antibodies bind to sample antibodies. Unbound conjugate is removed after 30 minutes incubation at RT by aspiration followed by a washing step.

In the third step the HRP converts the subsequently added colorless substrate solution of hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB) within a 15 min reaction time at RT into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow.

The absorbance of the solution read at 450/≥620 nm (OD) is directly proportional to the amount of specifically bound antibodies. A calibration curve is established by plotting the concentrations of the antibodies of the calibrators (x-axis) and their corresponding absorbance values (y-axis). The antibody concentrations of the unknown specimens are directly read from the calibration curve.

Test components

| | | | For 96 Wells |
|---|-------------------------|---|--|
| 1 | WELLS | Microtitration plate coated with antigen | 12 single breakable 8-well strips vacuum-sealed with desiccant |
| 2 | WASHBUF CONC 10x | Wash buffer 10-fold | 100 ml concentrate for 1000 ml solution white cap |
| 3 | DIL | Sample diluent | 100 ml · ready-to-use coloured red black cap |
| 4 | CAL 0 - 4 | Calibrator 0 - 4 diluted serum sample CAL 0 = 1 U/ml CAL 1 = 10 U/ml CAL 2 = 30 U/ml CAL 3 = 100 U/ml CAL 4 = 300 U/ml | 1.0 ml per calibrator ready-to-use coloured blue white cap |
| 5 | CONTROL + | Positive control diluted serum sample concentration see certificate of analysis enclosed | 1.0 ml · ready-to-use coloured blue red cap |
| 6 | CONJ HRP | HRP-conjugate HRP-labelled, anti-human-IgG antibodies | 15 ml · ready-to-use coloured red red cap |
| 7 | SUBSTR TMB | Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide | 15 ml · ready-to-use blue cap |
| 8 | STOP | Stop solution 0.25 M sulphuric acid | 15 ml · ready-to-use yellow cap |
| 9 | | Adhesive film | 2 |

Preparation and storage of samples

Collection and storage

Blood is taken by venipuncture. Serum or plasma is separated after clotting by centrifugation. Serum or plasma can be used in the *Serazym*[®] Anti-PR3 (cANCA) IgG. Contaminated samples should not be used. Repeated freezing and thawing should be avoided. The samples may be kept at 2...8 °C for up to two days. Long-term storage requires -20 °C.

Preparation before use

Allow samples to reach RT prior to assay. Take care to agitate samples gently in order to ensure homogeneity.

Note: Patient samples have to be diluted 1 : 101, e.g. 10 µl sample + 1000 µl sample diluent (3).

Materials required but not provided

Adjustable one-channel micropipette 0.100 - 1.000 ml and 0.010 - 0.100 ml · Adjustable 8-channel micropipette 0.050 - 0.200 ml · Pipette tips · Reagent containers for dispensing with 8-channel pipette · Graduated measuring flasks 10 ml and 100 ml · Test tubes 2.0 ml for sample dilution · Microtitration plate washer (automatic or hand wash head) · Microtitration plate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · Distilled or deionized water · Stop-watch

Preparation and storage of reagents

Kit size and expiry

The *Serazym*[®] Anti-PR3 (cANCA) IgG has been designed for 12 x 8 determinations. The complete kit with unopened reagent bottles and microtitration strips is stable until the expiry date printed on the kit box in case of storage at 2...8 °C. After opening all kit components are stable for at least 2 months, provided proper storage. The ready-to-use wash solution can be used for at least one month when stored at 2...8 °C.

Reagent preparation

The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For example: 10 ml wash buffer (2) + 90 ml distilled or deionized water.

Besides wash buffer all components of the device are ready-to-use.

Assay procedure

Allow all components to reach room temperature prior to use in the assay.

Dilute patient samples with sample diluent (3) 1 : 101 (v/v), e.g. 10 μ l sample + 1000 μ l sample diluent (3).

The calibrators and the positive control are ready-to-use.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash solution in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CAL 0** calibrator 0 (4) if desired*,
 100 µl **CAL 1 - 4** calibrators 1, 2, 3, 4 (4)
 100 µl **CONTROL +** positive control (5)
 100 µl **diluted samples** resp.
 into the intended wells.
3. Cover plate (9) and incubate for 60 min at RT.
4. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
5. Dispense 100 µl **CONJ HRP** HRP-conjugate (6) per well.
6. Cover plate (9) and incubate for 30 min at RT.
7. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
8. Dispense 100 µl **SUBSTR TMB** substrate (7) per well.
9. Incubate for 15 min at RT protected from light.
10. Dispense 100 µl **STOP** stop solution (8) per well, mix gently.
11. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction stop.

* Note:

Using CAL 0 (1 U/ml) enables quantification of concentrations below the value of CAL 1 (10 U/ml), otherwise CAL 0 is not needed.

Result interpretation

Create a calibration curve by plotting the mean absorbances of the calibrators CAL 0-4 (Y-axis) to their corresponding antibody-concentrations (X-axis).
 Alternatively create a calibration curve by plotting the mean absorbances of the calibrators CAL 1-4 (Y-axis) to their corresponding antibody-concentrations (X-axis).

For the calibration curve use the 4-Parameter-regression-model.

Determine the antibody-concentrations of the unknown samples by referring their mean absorbances to the calibration curve. In case of using the prescribed sample dilution of 1 : 101, antibody-concentrations are not to be determined by multiplication with the dilution factor, since the calibrators are already prediluted.

If a sample after proper dilution of 1:101 produced a higher absorbance than calibrator CAL 4 please retest the sample at higher dilution and calculate its concentration in the following way:

| Sample volume | Added volume of diluent | Final volume | Dilution of original sample | Multiplication factor to determine analyte concentration of the original sample |
|------------------------|-------------------------|----------------|-----------------------------|---|
| <i>10 µl</i> | <i>1000 µl</i> | <i>1010 µl</i> | <i>1 : 101 = Sample 1</i> | <i>No multiplication required</i> |
| 500 µl Sample 1 | 500 µl | 1000 µl | 1 : 202 = Sample 2 | factor 2 |
| 500 µl Sample 2 | 500 µl | 1000 µl | 1 : 404 = Sample 3 | factor 4 |
| 500 µl Sample 3 | 500 µl | 1000 µl | 1 : 808 = Sample 4 | factor 8 |

Please be aware that antibodies may not be stable in diluted samples resulting in activity loss. By that reason reinvestigation at higher dilution should be done at the same day.

Reference values

| Serazym® Anti-PR3 (cANCA) IgG | |
|-------------------------------|-----------|
| positive | ≥ 10 U/ml |
| negative | < 10 U/ml |

Test validity

The test run is valid if:

- absorbance of calibrator CAL 0 is < CAL 1
- absorbance of calibrator CAL 1 is ≤ 0.70
- absorbance of calibrator CAL 4 is ≥ 1.20

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation volumes, times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the method

Any clinical diagnosis should not be based on the results of *in-vitro* diagnostic methods alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis. False results may be caused by cross contaminations of the kit reagents and samples, bacterial or fungal contaminations of reagents and/or samples, incorrect washing and incorrect incubation times. In case of background problems when using automatic microplate ELISA systems it may be recommendable to wash wells 5 times (instead of 3 times) in every wash cycle.

Interference

Haemolytic, lipaemic and icteric samples do not interfere with the test up to a concentration of 500 mg/dl (haemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors do not interfere up to a concentration of 500 IU/ml.

Performance characteristics

Diagnostic sensitivity and specificity

A total of 82 serum samples were tested both in Serazym® Anti-PR3 (cANCA) IgG and in another commercially available test. The results are shown in the following crosstabulation.

| | comparative test | |
|----------------------------------|------------------|------------|
| | positive | negative |
| Serazym® Anti-PR3 (cANCA) IgG | positive 37 | negative 3 |
| | negative 0 | 42 |


In reference to this panel the diagnostic sensitivity was determined with 100 % and diagnostic specificity with 93.3 %.

Control groups

A total of 247 serum samples were tested in Serazym® Anti-PR3 (cANCA) IgG. The 9 control groups consist of serum samples from the following diseases: Rheumatoid arthritis, SLE, Scleroderma, MCTD (Mixed connective tissue disease), Myositis, Sjögren Syndrome, Colitis ulcerosa, Morbus Crohn and other collagenoses. The Serazym® Anti-PR3 (cANCA) IgG determined 17 serum samples of this collective positive with a clinical specificity of 93.1 %.

Frequency distribution

The frequency distribution of the *Serazym*[®] Anti-PR3 (cANCA) IgG antibody concentrations measured in 97 unselected human serum samples is illustrated:

Precision

Intra-assay coefficient of variation (CV [%]) in the *Serazym*[®] Anti-PR3 (cANCA) IgG from 8-fold determinations of samples:


| sample | concentration [U/ml] | CV [%] |
|--------|----------------------|--------|
| 1 | 195.8 | 5.96 |
| 2 | 88.2 | 4.72 |
| 3 | 25.6 | 5.18 |
| 4 | 16.8 | 2.47 |

Inter-assay coefficient of variation (CV [%]) in the *Serazym*[®] Anti-PR3 (cANCA) IgG from 8 different runs with twofold determinations of samples:

| sample | concentration [U/ml] | CV [%] |
|--------|----------------------|--------|
| A | 177.5 | 7.82 |
| B | 101.2 | 8.81 |
| C | 31.4 | 8.63 |
| D | 19.2 | 8.58 |

Typical reference curve

The typical reference curve in the *Serazym*[®] Anti-PR3 (cANCA) IgG is illustrated:

Common advices and precautions

This kit is for *In-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for diluted reagents. Do not use or mix reagents from different lots, damaged packages or bottles or reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8 °C before use.

Do not use or mix reagents from different lots except for sample diluent, wash buffer, substrate and stop solution.

Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They may neither be swallowed nor come into contact with skin or mucous membranes. Handle all samples as potentially infectious. Since the kit contains potentially hazardous materials, the following precautions must always be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

History of Changes

| Version | Section | Modifications |
|------------|--------------------------------|--------------------------------------|
| 2021-03-11 | Common Advices and Precautions | Update |
| | History of Changes | Adding section "History of Changes " |
| | Test Components | Update |

Incubation scheme *Serazym*[®] Anti-PR3 (cANCA) IgG

- | | | | |
|----|---|----------|--|
| 1. |  | 100 µl | CAL 0 - 4 (4) |
| | | 100 µl | CONTROL + (5) |
| | | 100 µl | diluted samples resp. |
| |  | 60 min | incubation (room temperature) |
| |  | 3 x wash | with wash solution |
| 2. |  | 100 µl | CONJ HRP (6) |
| | | 30 min | incubation (room temperature) |
| |  | 3 x wash | with wash solution |
| 3. |  | 100 µl | SUBSTR TMB (7) |
| | | 15 min | incubation (room temperature) protected from light |
| 4. |  | 100 µl | STOP (8) |

Read OD at 450 / \geq 620 nm

| | | | | |
|--|---|--|---|---------------------------|
|  Manufacturer |  Date of manufacture |  Use by | LOT Batch code | REF Catalog number |
|  Keep away from sunlight |  Temperature limits |  Biological risks |  Do not reuse | |
|  Consult instructions for use |  Caution | IVD <i>In-vitro</i> -diagnostic medical device |  Contains sufficient for <n> tests | |