SeraSpot® ANA-12 IgG

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SeraSpot® ANA-12 IgG

Spotimmunoassay for detection of IgG antibodies to nuclear and cytoplasmatic antigens
in human serum or plasma

 SP-002-12 G-S6  48  SP-002-12 G-S12  96  SP-002-12 G-S24  2x 96
 *In-vitro*-diagnostic medical device 

Seramun Diagnostica GmbH · Spreehagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 33767 791-10 · fax +49 33767 791-99 · info@seramun.com

Introduction

Connective tissue diseases are systemic inflammatory rheumatic diseases usually characterized by a chronic course (systemic connective tissue diseases) with overlapping symptoms (Overlap-Syndromes). Connective tissue diseases include:

- Systemic lupus erythematosus (SLE) and subsets
- Sjögren`s syndrome
- Systemic sclerosis
- Idiopathic (autoimmune) myositis
- Mixed connective tissue disease (MCTD or Sharp syndrome)
- Overlap syndromes

Connective tissue diseases are characterized by typical autoantibody profiles, all detectable with SeraSpot® ANA-12 IgG.

dsDNA antibodies:

IgG antibodies to dsDNA are marker antibodies and ACR-criterion for systemic lupus erythematosus (SLE). They are regarded as activity and prognostic marker of SLE. The detection frequency varies in dependence of the activity of the disease and organ manifestation: patients suffering from

- active SLE with renal involvement > 95 %,
- active SLE without renal involvement 50 - 70 %,
- inactive SLE < 40 %.

In sera from patients with rheumatoid arthritis, juvenile chronic arthritis, Sjögren`s syndrome, sclerosis and other diseases antibodies to dsDNA may also be detected temporarily with low titers. Due to their higher specificity in comparison to anti-nucleosome antibodies anti-dsDNA antibodies undoubtedly are of major importance for the diagnosis of SLE.

Nucleosome antibodies:

Autoantibodies to nucleosomes, a complex consisting of four homodimers of the histones H2A, H2B, H3 and H4 resulting in an octamer surrounded by helical DNA are detectable in 56 to 90 % of SLE patients but also in cases of drug induced Lupus. Sometimes they can be detected earlier than anti-dsDNA antibodies and therefore are useful as early diagnostic indicator of an onset of SLE.

Sm antibodies:

Sm-antibodies are diagnostic marker and ACR-criterion of SLE with a specificity of 99 %, but a sensitivity of only 10 - 15 % in SLE patients of caucasian origin (and 30 % to over 40 % in asiatic patients).

Since an association to several organ manifestations (kidneys, central nervous system) has been proven, antibodies to Sm are considered as prognostic marker of SLE.

P0 antibodies:

Autoantibodies to ribosomal phosphoproteins (P0, P1, P2) are regarded as highly specific for SLE. Autoantibodies to the major target antigen P0 are known as diagnostic marker of SLE reaching a diagnostic sensitivity of 10 - 20 % (and up to 40 % in asiatic patients) and a diagnostic specificity of nearly 100 %. Anti-P0 antibodies are predominantly detectable during the active phase of SLE and associated with renal and liver involvement.

Histone antibodies:

Antibodies to histones are detectable in many, preferably rheumatic diseases and therefore not specific for a disease. Nevertheless high titers of anti-histone antibodies are found almost exclusively in patients with SLE and drug-induced lupus. In the absence of SLE -markers like anti-dsDNA and anti-Sm antibodies high titers of anti-histone antibodies are characteristic for drug-induced lupus.

U1-snRNP antibodies:

Anti-U1-snRNP antibodies are marker antibodies and diagnostic criterion of mixed connective tissue disease (MCTD) with a sensitivity of 100 % and a specificity of 98 % for high titers of anti-U1-snRNP antibodies and in the absence of antibodies to Sm and dsDNA.

Antibodies to U1-RNP with often low titers are found in 13 - 32 % of patients with SLE and in 10 % of patients with systemic sclerosis.

Using the three RNP-proteins A, C and 68 kD *SeraSpot*[®] ANA-12 IgG detects all U1-RNP antibodies.

Ro/SS-A antibodies:

Anti-Ro/SS-A antibodies belong to the group of antinuclear antibodies (ANA), although the proteins Ro60 and Ro52 are localized in both, nucleus and cytoplasm. Ro60 is known as protein of the hY-RNP-complex, whereas Ro52 has been identified as E3-Ubiquitin-Ligase. Antibodies to Ro/SS-A usually occur together with La/SS-B antibodies.

Ro/SS-A antibodies are predominantly detected in patients with Sjögren's syndrome and the different forms of lupus erythematosus with an assumed higher specificity of the antibodies to Ro60 in comparison to anti-Ro52 antibodies. They serve as diagnostic marker and are part of the classification criteria of primary and secondary Sjögren's syndrome with a sensitivity of 96 % and 80 % respectively. Ro/SS-A antibodies are considered as markers for the early diagnosis of Sjögren's syndrome since they may occur years before clinical manifestation.

Anti-Ro/SS-A antibodies are detectable in 40 - 60 % of patients who suffer from SLE. They may occur together with antibodies typical for SLE (anti-dsDNA, anti-Sm antibodies) or as isolated antibodies indicating a relatively mild form of SLE.

Antibodies to Ro/SS-A are a diagnostic marker of subacute cutaneous lupus erythematosus (SCLE) and detectable in 90 - 100 % of the cases. About 10 % of them turn into SLE.

Ro/SS-A antibodies are detectable in nearly 100 % of neonatal lupus erythematosus (NLE) cases. The coincidence of anti- Ro/SS-A-and La/SS-B antibodies is associated with congenital heart block (CHB).

Ro/SS-A antibodies are additionally detectable in patients with rheumatoid arthritis (5 - 8 %) and systemic sclerosis (9 %).

La/SS-B antibodies:

Antibodies to La/SS-B serve as important diagnostic marker and belong to the classification criteria of Sjögren's syndrome. The diagnostic sensitivity is about 70 % for the primary Sjögren's syndrome and about 50 % for the secondary Sjögren's syndrome. In case of simultaneous detection of La/SS-B- and Ro/SS-A antibodies diagnostic specificity for Sjögren's syndrome is higher in comparison to isolated detection of Ro/SS-A antibodies. Anti-La/SS-B antibodies almost always occur together with antibodies to Ro/SS-A.

La/SS-B antibodies are regarded as early diagnostic markers of Sjögren's syndrome, since they may be detectable years before clinical symptoms become evident.

Antibodies to La/SS-B are detectable in 25 % of systemic lupus erythematosus (SLE) cases, in 70 % of neonatal lupus erythematosus (NLE) cases and in 80 % of subacute cutaneous lupus erythematosus (SCLE) cases.

La/SS-B antibodies are only rarely detected in other connective tissue diseases.

Scl-70 antibodies:

Scl-70 or Topoisomerase I antibodies are diagnostic marker of systemic sclerosis with a sensitivity of

- 18 - 30 % in general
- 10 - 15 % in limited sclerosis (CREST-syndrome)
- 40 - 65 % in diffuse forms of systemic sclerosis

The specificity reaches nearly 100 %.

Patients who develop anti-Scl-70 antibodies usually suffer from a more severe disease in comparison to patients with anti-centromer antibodies. Antibodies to Scl-70 may occur already years before specific symptoms of sclerosis are detectable.

CENP-B antibodies:

CENP-B antibodies are considered as diagnostic marker of the limited form of systemic sclerosis (sensitivity 50 - 70 %). Patients who develop antibodies to CENP-B usually show a rather mild course of sclerosis. Antibodies to CENP-B may occur already years before specific symptoms of sclerosis are detectable.

CENP-B antibodies are detectable in only 10 - 30 % of patients with primary biliary cirrhosis (PBC) and very rarely in patients suffering from other connective tissue diseases.

Jo-1 antibodies:

Antibodies to Jo-1 represent a diagnostic marker for idiopathic myositis with a diagnostic specificity of nearly 100 %. Diagnostic sensitivity reaches 18 - 46 % for polymyositis and 25 % for dermatomyositis. About 60 % of patients with detectable levels of antibodies to Jo-1 suffer from fibrosing alveolitis.

Jo-1 antibodies serve as prognostic marker of myositis, since these patients often develop a severe course of the disease.

The cytoplasmatic localization of the target autoantigen (histidyl-tRNA-synthetase) actually excludes anti-Jo-1 antibodies from anti-nuclear antibodies in the original meaning of the word.

Literature:

1. Conrad, K., Schöbner, W., Hiepe, F., Autoantikörper bei systemischen Autoimmunerkrankungen Pabst Science Publishers Lengerich u.a. 2006, 2. Aufl.
2. Peter, J.B., Shoenfeld Y. (Eds), Autoantibodies Elsevier Amsterdam u.a. 1996
3. Tan, E.M., Antinuclear antibodies: diagnostic markers for autoimmune disease and probes for cell biology, Adv. Immunol. 44, 1989: 93-151

Intended use

The *SeraSpot*[®] ANA-12 IgG is an *in vitro* diagnostic medical device (spotimmunoassay, SIA) for the detection of antinuclear (ANA) and cytoplasmic antibodies of the IgG isotype directed against the following antigens: dsDNA, Nucleosome, Sm, P0, Histone, U1-snRNP, Ro/SS-A 60, Ro/SS-A 52, La/SS-B, Scl-70, CENP-B and Jo-1 in human serum or plasma samples.

Principle of the test

SeraSpot[®] ANA-12 IgG test is a solid phase immunoassay based on the use of recombinant or purified native proteins / protein complexes and dsDNA as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96-well microtiter plates. The antigens serve as capture molecules for autoantibodies to nuclear antigens (ANA). Bound antibodies are detected by Horseradish Peroxidase-(HRP)-labeled antibodies against human antibodies of IgG-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-Tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots develop as precipitated product from colorless substrate solution. Pale blue to dark blue spots are visible by eye.

The test procedure for antibody detection is performed in 3 steps:

Step 1

Incubation of 1 : 101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are sucked off and unbound components are removed by 3 wash cycles with wash buffer.

Step 2

Incubation of wells with the anti-human IgG-HRP-conjugate for 30 minutes at room temperature. After incubation the antibody is sucked off and unbound components are removed by 3 wash cycles with wash buffer.

Step 3

Incubation of wells with substrate solution SeramunBlau[®] spot dark for 30 minutes at room temperature. The reaction is stopped by aspirating the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the Seramun *SpotSight*[®] plate or Seramun *SpotSight*[®] strip scanner. The obtained images are interpreted by using the Software Seramun *SpotSight*[®] scan. Taking account of the test-specific evaluation criteria, the test result is determined.

Kit components

		For 48 determinations	For 96 determinations	For 2x 96 determinations	
1	WELLS	Wells with arrays	6 single breakable 8-well strips in frame vacuum-sealed with desiccant	12 single breakable 8-well strips in frame vacuum-sealed with desiccant	2x 12 single breakable 8-well strips in frame separately vacuum-sealed with desiccant
Color coding: purple					
2	WASHBUF CONC 10X	Wash buffer Seramun® Wash buffer A	100 ml concentrate transparent bottle, white cap	100 ml concentrate transparent bottle, white cap	2x 100 ml concentrate transparent bottle, white cap
3	DIL	Sample diluent Seramun® Sample diluent B	55 ml ready-to-use solution, colored red transparent bottle, black cap	2x 55 ml ready-to-use solution, colored red, transparent bottle, black cap	4x 55 ml ready-to-use solution, colored red, transparent bottle, black cap
4	CONJ HRP IgG	Anti-human IgG-HRP-conjugate	8 ml ready-to-use solution, colored red, transparent bottle, red cap	8 ml ready-to-use solution, colored red, transparent bottle, red cap	2x 8 ml ready-to-use solution, colored red, transparent bottle, red cap
5	SUBSTR TMB	Substrate SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	8 ml ready-to-use solution, black bottle, blue cap	8 ml ready-to-use solution, black bottle, blue cap	2x 8 ml ready-to-use solution, black bottle, blue cap
6	COVER	Adhesive film	2x	2x	4x
7	SWAB	Swab with 70 % (v/v) Isopropyl alcohol	3x2	6x2	12x2

Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 7 days. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided. Bring the samples to room temperature and mix well.

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (e.g. 10 µl sample and 1000 µl buffer) with the ready-to-use sample diluent.

Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- stop watch
- Washer for 96-well microtiter plates (optional)
- Collecting devices for infectious material
- Seramun *SpotSight*® plate or Seramun *SpotSight*® strip scanner with evaluation software Seramun *SpotSight*® scan

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used up to 4 weeks when stored at 2...8 °C and up to 2 weeks when stored at room temperature.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtiter plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or deionized water. For Example: 10 ml **WASHBUF CONC 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

Workplace requirements

Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

Assay procedure

Performance at room temperature (18...25 °C) and at the specified incubation times.

For the handling of SeraSpot® tests the use of the starter kit (order number: SP-000-1) is recommended. The kit is designed for multiple use.

Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers, which may adhere to the underside of the bottom of the wells, by using the swab **SWAB** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

Working steps


1. Warm all reagents and required wells **WELLS** to room temperature (RT) before use. Mix gently without causing foam.
2. Add **100 µl** of 1 : 101 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 1000 µl sample diluent **DIL**.
3. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **CONJ HRP IgG** to each well.
6. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard conjugate), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **SUBSTR TMB** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent paper.
11. Clean bottom of wells with swab **SWAB**.
12. Perform image acquisition using the Seramun *SpotSight*® plate or the Seramun *SpotSight*® strip scanner and evaluate the results with the Seramun *SpotSight*® scan software. If image acquisition is performed with the scanner Seramun *SpotSight*® strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

Aspiration and washing steps can be performed manually or with a washer for 96-well microtiter plates.

After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when the plate is stored in the dark

Test evaluation

Array layout

Parameter	Controls
05 U1-snRNP	01 Positive control (PC)
06 CENP-B	02 Negative control (NC), 0 rel. units
07 Nucleosome	03 Cut-off control (CO), 30 rel. units
08 Histone	04 Reference 3 (R3), 60 rel. units
10 dsDNA	09 Reference 2 (R2), 100 rel. units
11 Jo-1	15 Reference 1 (R1), 300 rel. units
12 Scl-70	19 Serum control (SC)
13 Sm	
14 P0	
16 La/SS-B	
17 Ro/SS-A 52	
18 Ro/SS-A 60	P Well position marker

Validity criteria for the test

The *SeraSpot*[®] ANA-12 IgG test includes the following control spots:

1. Positive control (PC). Intensively stained spot, darker than Cut-off control. Always stained.
2. Cut-off control (CO), corresponding to 30 rel. units. Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC), corresponding to 0 rel. units. Pale spot with intensity lower than Cut-off control.
4. Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.
5. Reference 1 (R1), corresponding to 300 rel. units. Intensively stained spot. Always stained.
6. Reference 2 (R2), corresponding to 100 rel. units. Intensity lower than R1. Always stained.
7. Reference 3 (R3), corresponding to 60 rel. units. Intensity lower than R2. Always stained.

The spots NC, CO, R1 ... R3 are used to create a reference curve (four parameter non-linear regression) to calculate the staining intensities of the spots in relative units (rel. units) or in international units (IU) for dsDNA.

The test cannot be evaluated if one of the validity criteria listed in points 1 to 7 is not fulfilled.

Result interpretation

Evaluation of the test has to be performed by using the Seramun *SpotSight*[®] scan software. A visual plausibility check is possible with the aid of the magnifying lens contained in the starter kit (order number: SP-000-1).

Judgment of spot staining is performed in accordance to the following classification:

Result Interpretation	IgG
negative	color intensity of antigen spots \leq Cut-off control
positive	color intensity of one or more antigen spots $>$ Cut-off control

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

Interference

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) do not interfere with the test. Rheumatoid factors also do not interfere with the test.

Assay procedure

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere underside the bottom of the wells or to the array area cause erroneous results.

Performance characteristics

Diagnostic sensitivity

Serologically predetermined samples (reference test 1) of patients with suspected connective tissue diseases were tested in the *SeraSpot*[®] ANA-12 IgG test. Samples with discrepant results were retested in a second assay (reference test 2) for detection of antibodies against autoantigens. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results ("Amended sensitivity").

U1-snRNP	n = 47	Reference test 1		Initial sensitivity:	82.1 %
		positive (n = 39)	negative (n = 8)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	32	0	Amended sensitivity:	91.4 %
	negative	7	8		
CENP-B	n = 112	Reference test 1		sensitivity:	100.0 %
		positive (n = 18)	negative (n = 94)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	18	2		
	negative	0	92		
Nucleosome	n = 103	Reference test 1		Initial sensitivity:	45.3 %
		positive (n = 75)	negative (n = 28)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	34	6	Amended sensitivity:	80.4 %
	negative	41	22		
Histone	n = 103	Reference test 1		Initial sensitivity:	67.3 %
		positive (n = 49)	negative (n = 54)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	33	15	Amended sensitivity:	94.7 %
	negative	16	39		
dsDNA	n = 103	Reference test 1		Initial sensitivity:	83.1 %
		positive (n = 77)	negative (n = 26)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	64	14	Amended sensitivity:	97.2 %
	negative	13	12		
Jo-1	n = 60	Reference test 1		Initial sensitivity:	93.3 %
		positive (n = 15)	negative (n = 45)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	14	3	Amended sensitivity:	93.3 %
	negative	1	42		
Scl-70	n = 112	Reference test 1		Initial sensitivity:	92.2 %
		positive (n = 64)	negative (n = 48)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	59	3	Amended sensitivity:	92.2 %
	negative	5	45		
Sm	n = 103	Reference test 1		Initial sensitivity:	53.3 %
		positive (n = 15)	negative (n = 88)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	8	1	Amended sensitivity:	98.9 %
	negative	7	87		
P0	n = 103	Reference test 1		sensitivity:	100.0 %
		positive (n = 6)	negative (n = 97)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	6	5		
	negative	0	92		
La/SS-B	n = 103	Reference test 1		Initial sensitivity:	57.9 %
		positive (n = 19)	negative (n = 84)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	11	6	Amended sensitivity:	94.3 %
	negative	8	78		
Ro/SS-A 52	n = 103	Reference test 1		Initial sensitivity:	96.4 %
		positive (n = 28)	negative (n = 75)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	27	0	Amended sensitivity:	100.0 %
	negative	1	75		
Ro/SS-A 60	n = 103	Reference test 1		Initial sensitivity:	87.5 %
		positive (n = 48)	negative (n = 55)		
<i>SeraSpot</i> [®] ANA-12 IgG	positive	42	0	Amended sensitivity:	95.5 %
	negative	6	55		

Diagnostic specificity

A total of 200 blood donor samples were tested by the SeraSpot® ANA-12 IgG test ("Initial specificity"). Positive samples were followed up with a reference test (reference test 1).

U1-snRNP		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	18	specificity: 91.0 %
	negative		0	182	
CENP-B		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	2	specificity: 99.0 %
	negative		0	198	
Nucleosome		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	11	specificity: 94.5 %
	negative		0	189	
Histone		n = 200	Reference test 1		
			positive (n = 3)	negative (n = 197)	
SeraSpot® ANA-12 IgG	positive		3	10	specificity: 94.9 %
	negative		0	187	
dsDNA		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	2	specificity: 99.0 %
	negative		0	198	
Jo-1		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	1	specificity: 99.5 %
	negative		0	199	
Scl-70		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	7	specificity: 96.5 %
	negative		0	193	
Sm		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	1	specificity: 99.5 %
	negative		0	199	
PO		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	1	specificity: 99.5 %
	negative		0	199	
La/SS-B		n = 200	Reference test 1		
			positive (n = 2)	negative (n = 198)	
SeraSpot® ANA-12 IgG	positive		2	7	specificity: 96.5 %
	negative		0	191	
Ro/SS-A 52		n = 200	Reference test 1		
			positive (n = 2)	negative (n = 198)	
SeraSpot® ANA-12 IgG	positive		2	2	specificity: 99.0 %
	negative		0	196	
Ro/SS-A 60		n = 200	Reference test 1		
			positive (n = 0)	negative (n = 200)	
SeraSpot® ANA-12 IgG	positive		0	1	specificity: 99.5 %
	negative		0	199	

Cross-reactivity

Samples from patients with rheumatoid arthritis

A collection of 113 samples from patients suffering from rheumatoid arthritis with a positive anti-CCP (n=30) or RF (n=83) serology (reference test 1) were tested in the *SeraSpot*[®] ANA-12 IgG test and followed up with a second reference test ("Amended specificity").

U1-snRNP	n = 113	Reference test 1		Initial specificity: 86.7 % Amended specificity: 86.7 %	
		positive (n = 0)	negative (n = 113)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		0
		negative	0	98	
CENP-B	n = 113	Reference test 1		Initial specificity: 97.3 % Amended specificity: 99.1 %	
		positive (n = 0)	negative (n = 113)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		0
		negative	0	110	
Nucleosome	n = 113	Reference test 1		Initial specificity: 98.2 % Amended specificity: 98.2 %	
		positive (n = 4)	negative (n = 109)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		2
		negative	2	107	
Histone	n = 113	Reference test 1		Initial specificity: 86.5 % Amended specificity: 86.5 %	
		positive (n = 9)	negative (n = 104)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		9
		negative	0	90	
dsDNA	n = 113	Reference test 1		Initial specificity: 96.4 % Amended specificity: 97.3 %	
		positive (n = 1)	negative (n = 112)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		1
		negative	0	108	
Jo-1	n = 113	Reference test 1		Specificity: 100.0 %	
		positive (n = 0)	negative (n = 113)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		0
		negative	0	113	
Scl-70	n = 113	Reference test 1		Initial specificity: 97.2 % Amended specificity: 98.1 %	
		positive (n = 6)	negative (n = 107)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		6
		negative	0	104	
Sm	n = 113	Reference test 1		Initial specificity: 99.1 % Amended specificity: 99.1 %	
		positive (n = 0)	negative (n = 113)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		0
		negative	0	112	
P0	n = 113	Reference test 1		Initial specificity: 98.2 % Amended specificity: 99.1 %	
		positive (n = 0)	negative (n = 113)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		0
		negative	0	111	
La/SS-B	n = 113	Reference test 1		Initial specificity: 93.8 % Amended specificity: 93.8 %	
		positive (n = 0)	negative (n = 113)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		0
		negative	0	106	
Ro/SS-A 52	n = 113	Reference test 1		Initial specificity: 97.3 % Amended specificity: 99.1 %	
		positive (n = 3)	negative (n = 110)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		3
		negative	0	107	
Ro/SS-A 60	n = 113	Reference test 1		Initial specificity: 98.2 % Amended specificity: 99.1 %	
		positive (n = 1)	negative (n = 112)		
		<i>SeraSpot</i> [®] ANA-12 IgG	positive		1
		negative	0	110	

Epstein-Barr-virus (EBV) antibody positive samples

EBV antibody positive samples (reference test 1) were tested in the SeraSpot® ANA-12 IgG test.

U1-snRNP	n = 22	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	22	
CENP-B	n = 22	Reference test 1		Specificity: 100.0 %	
			positive (n = 2)		negative (n = 20)
		SeraSpot® ANA-12 IgG	positive		2
		negative	0	20	
Nucleosome	n = 22	Reference test 1		Specificity: 95.5 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	21	
Histone	n = 22	Reference test 1		Specificity: 77.3 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	17	
dsDNA	n = 22	Reference test 1		Specificity: 95.2 %	
			positive (n = 1)		negative (n = 21)
		SeraSpot® ANA-12 IgG	positive		1
		negative	0	20	
Jo-1	n = 22	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	22	
Scl-70	n = 22	Reference test 1		Specificity: 95.5 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	21	
Sm	n = 22	Reference test 1		Specificity: 95.5 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	21	
P0	n = 22	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	22	
La/SS-B	n = 22	Reference test 1		Specificity: 86.4 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	19	
Ro/SS-A 52	n = 22	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	22	
Ro/SS-A 60	n = 22	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 22)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	22	

Cytomegalovirus (CMV) antibody positive samples

CMV antibody positive samples (reference test 1) were tested in the SeraSpot® ANA-12 IgG test.

U1-snRNP	n = 27	Reference test 1		Specificity: 88.9 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	24	
CENP-B	n = 27	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	27	
Nucleosome	n = 27	Reference test 1		Specificity: 92.6 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	25	
Histone	n = 27	Reference test 1		Specificity: 100.0 %	
			positive (n = 1)		negative (n = 26)
		SeraSpot® ANA-12 IgG	positive		0
		negative	1	26	
dsDNA	n = 27	Reference test 1		Specificity: 96.2 %	
			positive (n = 1)		negative (n = 26)
		SeraSpot® ANA-12 IgG	positive		1
		negative	0	25	
Jo-1	n = 27	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	27	
Scl-70	n = 27	Reference test 1		Specificity: 92.6 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	25	
Sm	n = 27	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	27	
P0	n = 27	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	27	
La/SS-B	n = 27	Reference test 1		Specificity: 88.9 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	24	
Ro/SS-A 52	n = 27	Reference test 1		Specificity: 96.3 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	26	
Ro/SS-A 60	n = 27	Reference test 1		Specificity: 100.0 %	
			positive (n = 0)		negative (n = 27)
		SeraSpot® ANA-12 IgG	positive		0
		negative	0	27	

Precision

Samples of known antibody titer were assayed by *SeraSpot*[®] ANA-17 IgG test. The staining intensities of spots were determined as relative units (rel. units) according to the reference curve and the variation coefficients (CV) were calculated in order to determine the within-run precision (intra-assay-CV), between-run precision (inter-assay-CV) and Lot-to-Lot precision (Interbatch-CV).

	Intra-assay-CV		Inter-assay-CV		Interbatch-CV	
	rel. units \bar{x} n = 40	CV [%]	rel. units \bar{x} n = 80	CV [%]	rel. units \bar{x} n = 240	CV [%]
U1-snRNP	64.9	8.9	55.8	9.0	56.3	15.5
CENP-B	39.0	8.1	37.4	12.3	35.3	12.7
Nucleosome	43.6	11.9	48.5	19.3	48.9	24.2
Histone	67.6	8.8	57.3	14.8	55.6	12.3
dsDNA	63.3	7.8	72.9	11.2	72.4	11.7
Jo-1	48.7	5.4	47.4	10.5	49.1	15.2
Scl-70	37.1	11.6	40.4	15.9	42.3	16.4
Sm	77.4	10.5	93.2	12.8	91.8	13.3
P0	57.9	7.9	48.7	10.1	54.6	21.0
La/SS-B	54.7	11.9	53.1	10.2	52.8	13.9
Ro/SS-A 52	68.0	9.3	84.2	12.9	88.9	19.6
Ro/SS-A 60	96.3	13.3	130.7	15.9	143.6	18.7
Procedure	1 operator, 40x determinations, 1 batch		5 operators, 2x determinations, 2x testings per day 20 days, 1 batch		5 operators, 2x determinations, 2x testings per day, 20 days, 3 batches	

Automation

Samples (n = 96) of known antibody titer were assayed manually and (in parallel) by use of the microplate processor DS2[®] (Dynerx Technologies; manual sample dilution). Color intensity of antigen spots was recorded and used for calculation of the coefficient of determination r^2 for selected antibody reactivities and the different methods of operation.

	Dynex DS2 [®] vs. manual processing
	r^2
U1-snRNP	0.978
CENP-B	0.989
Nucleosome	0.962
Histone	0.758
dsDNA	0.996
Jo-1	0.998
Scl-70	0.995
Sm	0.999
P0	0.963
La/SS-B	0.991
Ro/SS-A 52	0.997
Ro/SS-A 60	0.986

r^2 was not determined if more than 90 % of the total number (n = 96) of the samples showed no antibody reactivity for the particular antigens.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

Do not use or mix reagents from different lots except for sample diluent, wash buffer and TMB substrate solution. These can be used across batches and parameters.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**


History of changes

Version	Section	Modifications
2020-08-12	General	General revision
	Preparation and storage of samples	Stability of samples
	Preparation and storage of reagents	Storage of ready-to-use wash solution
	Kit components	Update
	Limitations of the procedure	Interference
	Performance characteristics	Update of data

Incubation scheme *SeraSpot*[®] ANA-12 IgG

-  100 µl
30 min
 3 x wash
diluted sample (1 : 101)
incubation (room temperature)
with wash solution, each 400 µl per well
-  50 µl
30 min
 3 x wash
conjugate **CONJ HRP IgG**
incubation (room temperature)
with wash solution, each 400 µl per well
-  50 µl
30 min
 Aspiration
substrate **SUBSTR TMB**
incubation (room temperature, protected from light)
- imaging of wells and image analysis
scanner Seramun *SpotSight*[®] plate / Seramun *SpotSight*[®] strip and software Seramun *SpotSight*[®] scan

 Manufacturer	 Date of manufacture	 Use by	 LOT Batch code	 REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	 IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

SeramunBlau®, SeraSpot® und Seramun SpotSight® are registered trademarks of Seramun Diagnostica GmbH.