## SeraSpot<sup>®</sup> Anti-EBV-4 IgG

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## SeraSpot<sup>®</sup> Anti-EBV-4 IgG

Spot immunoassay for the qualitative detection of IgG antibodies against 4 Epstein-Barr-virus (EBV) specific antigens in serum or plasma of human origin

<b>REF</b>	SP-013-4 G-S6		48
<b>REF</b>	SP-013-4 G-S12		96
<b>REF</b>	SP-013-4 G-S24		2 x 96
<b>IVD</b>	In-vitro-diagnostic medical device		

Seramun Diagnostica GmbH • Sprenhagener Str. 1 • 15754 Heidesee • Germany •  
T +49 33767 791-10 • [info@seramun.com](mailto:info@seramun.com) • [www.seramun.com](http://www.seramun.com)

In-vitro diagnostic medical device

Unique device identifier

Manufacturer

Country of manufacture and date of manufacture

Article number

Serial number

Keep away from sunlight

Humidity limitation

Batch code

Consult instructions for use

Temperature limit

Do not reuse

Sufficient for  $n$  tests

Biohazard

Use-by date

Attention

## Intended Use

SeraSpot® Anti-EBV-4 IgG is an IVD test for the qualitative determination of antibodies of the IgG isotype against the antigens EBNA-1 (Epstein-Barr Nuclear Antigen-1), VCA p18 (Viral Capsid Antigen p18), VCA gp125 (Viral Capsid Antigen Glycoprotein 125) and EA-D p54 (Early Antigen Diffuse p54) in serum or plasma (citrate, EDTA, heparin) of human origin by a laboratory professional user.

The test is used in combination with the Seramun SpotSight® plate mono / Seramun SpotSight® strip device for image acquisition and the software Seramun SpotSight® scan for image analysis.

It is intended to aid in the diagnosis of infectious mononucleosis (Pfeiffer's glandular fever) in specimen materials from patients with suspicion of EBV infection.

The test must not be used with specimen materials other than serum or plasma (citrate, EDTA, heparin) of human origin, for diagnosis, monitoring, screening, prediction, prognosis, as companion diagnostic, in the near patient setting and by lay persons and in transplantation medicine for the detection of transmissible agents in blood, blood components, cells, tissues, organs, or any of their derivatives.

## Principle of the Test

SeraSpot® Anti-EBV-4 IgG is a solid phase immunoassay (spot immunoassay) using recombinant antigens printed in array format (spot array) on the bottom of wells of 96-well microtiter plates. Antibodies will bind to the immobilized antigens EBNA-1, VCA p18, VCA gp125 and EA-D p54. After incubation and wash steps bound antibodies are detected by horseradish peroxidase (HRP)-labeled antibodies against human antibodies of IgG-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). Blue spots indicate the formation of immune complexes. Spots range from pale to dark blue and are visible by eye.

### Detection of Specific Antibodies is Performed in 3 steps:

#### Step 1

Incubation of 1 : 101 diluted samples for 30 min at room temperature in selected wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with diluted wash buffer.

#### Step 2

Incubation of wells with the HRP-labeled anti-human IgG for 30 min at room temperature. Removal of unbound conjugate antibodies by aspiration and 3 wash cycles with diluted wash buffer.

#### Step 3

Incubation of wells with substrate SeramunBlau® spot dark for 30 min at room temperature. The reaction is stopped by aspirating the substrate, followed by tapping the plate / strip dry onto lint-free absorbent paper. The developed arrays are to be stored light-protected until image analysis.

## Test Components (Delivery Scope)

			For 48 determinations	For 96 determinations	For 2 x 96 determinations
1	<b>WELLS</b>	<b>Microtiter plate (wells with arrays)</b> EBV-specific antigens and controls immobilized as spots in array layout	6 single breakable 8-well strips in frame color coding: red vacuum sealed with desiccant	12 single breakable 8-well strips in frame color coding: red vacuum sealed with desiccant	2 x 12 single breakable 8-well strips in frame color coding: red vacuum sealed with desiccant
2	<b>WASHBUF (10x)</b>	<b>Wash buffer (10x)</b> Seramun® Wash buffer A 10x TRIS-based buffer	100 mL <b>concentrate</b> for 1000 mL buffer colorless white cap	100 mL <b>concentrate</b> for 1000 mL buffer colorless white cap	2 x 100 mL <b>concentrate</b> for each 1000 mL buffer colorless white cap
3	<b>DIL</b>	<b>Sample diluent</b> Seramun® Sample diluent B	55 mL ready-to-use solution colored red black cap	2 x 55 mL ready-to-use solution colored red black cap	4 x 55 mL ready-to-use solution colored red black cap
4	<b>CONJ HRP IgG</b>	<b>Conjugate</b> anti-Human IgG-HRP conjugate (Sheep)	8.0 mL ready-to-use solution colored red red cap	8.0 mL ready-to-use solution colored red red cap	2 x 8.0 mL ready-to-use solution colored red red cap
5	<b>SUBSTR</b>	<b>Substrate</b> SeramunBlau® spot dark 3,3',5,5'- Tetramethylbenzidine	8.0 mL ready-to-use solution colorless blue cap	8.0 mL ready-to-use solution colorless blue cap	2 x 8.0 mL ready-to-use solution colorless blue cap
6	<b>COVER</b>	<b>Covering film</b>	2 pieces	2 pieces	4 pieces
7	<b>SWAB</b>	<b>Swab</b> 70 % (v/v) Isopropyl alcohol	3 x 2 pieces	6 x 2 pieces	12 x 2 pieces
8		<b>Certificate of Analysis</b>	1 piece	1 piece	1 piece
9		<b>Instructions for Use</b>	1 piece	1 piece	1 piece

## Antigens Used

Name	Description	Clinical relevance
EBNA-1	Epstein-Barr Nuclear Antigen-1	IgG antibodies against EBNA-1 are specific markers for a past EBV infection. In some patients the EBNA-1-IgG antibody titer may decrease over time (secondary EBNA-1 loss).
VCA p18	Virus Capsid Antigen p18	VCA p18- IgM and IgG antibodies are specific markers for an EBV infection. IgM antibodies are detectable early after infection, decline in the course of the infection and are replaced by the IgG antibody response. Typically, VCA p18-specific IgG antibodies are detectable throughout life.
VCA gp125	Virus Capsid Antigen Glycoprotein125	Glycoprotein gp125 is the second of the "virus capsid antigens". Proteins gp125 and p18 are immunodominant. VCA gp125-specific IgM antibodies are often produced in a very early stage of a primary infection. Typically, VCA gp125-specific IgG antibodies are detectable throughout life.
EA-D p54	Early Antigen-Diffuse p54	EA-D has been described as an early antigen. It is involved in the replication of the viral DNA. EA-D-specific antibodies can occur during primary infection.

## Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micro-pipettes • adjustable 8-channel micro-pipette or multi-pipette • reagent container for multi-channel micro-pipettes • measuring cylinder • test tubes for sample dilution • test tube rack • washer for 96-well microtiter plates • laboratory system: Seramun SpotSight® plate mono / strip scanner with evaluation software Seramun SpotSight® scan • deionized water • lint-free absorbent paper • stop watch • collecting devices for infectious material • opaque cover (substrate reaction)

## Important Information

**This device is for *in-vitro* diagnostic use only.** The kit may be performed by laboratory professional user only.

Follow the instructions carefully. The shelf life specified must be observed. Do not use reagents from damaged packages or bottles. Do not use reagents from other manufacturers.

**Mixing of test kit components of different lots is only allowed for sample diluent, wash buffer and substrate.**

All serious incidents occurring in relation with SeraSpot® Anti-EBV-4 IgG must be reported to the manufacturer and the competent authority of the EU member state in which the user and/or the patient are located.

### Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all components to room temperature before use.

For larger sample series, pipetting reagents from liquid reservoirs using a multi-channel pipette is recommended to avoid time delays. Reagents that appear contaminated should not be used.

Avoid time delays when dispensing reagents.

Substrate must be protected from direct light!

Incorrect dilution of samples, incorrect washing of wells and incorrect timing may lead to erroneous results.

**Air bubbles** generated by forceful pipetting of samples and/or reagents may cause uneven signal formation of arrays. These arrays cannot be evaluated.

Damaged arrays, e.g. by **scratching** the bottom of the well with pipet tips or washer needles, are not suitable for evaluation.

Before taking images of the wells make sure to remove particles or fibers which may adhere to the bottom side of the wells using the swab provided in the kit.

Images of dried spots may appear more intense, which can lead to slight deviations of the measured values when scanned repeatedly. Assessment of the individual parameters is not changed.

### Workplace Requirements

Processing spot immunoassays requires a clean workplace. Adhesion of fibers to the plastic surface of the wells should be avoided. Fibers may cause interferences when taking images of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

## Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

<b>WELLS</b>	-	Contains material of animal origin.
<b>WASHBUF (10x)</b>	EUH208 EUH210	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request.
<b>DIL</b>	EUH208 EUH210 -	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request. Contains material of animal origin.
<b>CONJ HRP IgG</b>	Hazard components EUH208 H360D P202 P280 P308+P313 P501 -	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Restricted to professional users. May damage the unborn child. Do not handle until all safety precautions have been read and understood. Wear protective gloves / protective clothing / eye protection / face protection. IF exposed or concerned: Get medical advice / attention. Dispose of contents / container in accordance with local regulations. Contains material of animal origin.
<b>SWAB</b>	Hazard components H225 H319 H336 P210 P280 P304+P340 P305+P351+P338	Propan-2-ol; isopropyl alcohol; isopropanol Highly flammable liquid and vapour. Causes serious eye irritation. May cause drowsiness or dizziness. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Wear protective gloves / protective clothing / eye protection / face protection / hearing protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### Limitations of the Procedure

Result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. In an early stage of infection antibodies may not be present yet or below detection limit.

Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

### Interference

In rare cases, specimens may contain antibodies to BSA (bovine serum albumin) and / or AGE (advanced glycation end products), which may cause nonspecific reactions, so that the test result is read as "not available" ("n.a.") by the Seramun SpotSight® scan software.

## Sample Treatment

### Sample Collection

Collect serum or plasma (citrate, EDTA, heparin) of human origin in suitable sampling container.

### Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 7 days. For longer periods samples have to be stored at < -15 °C. Repeated freezing and thawing of samples should be avoided!

### Sample Preparation

Bring samples to room temperature before use. Homogeneity should be ensured by shaking components briefly.

Samples have to be diluted 1 : 101 (v/v) with sample diluent.

Example: 10 µL sample and 1000 µL sample diluent.

## Reagent Treatment

### Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips is stable until the printed expiration date when stored at 2...8 °C. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

### Reagent Preparation

Allow all components to reach **room temperature** prior to use. The **microtiter plate** is vacuum sealed with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored at 2...8 °C and protected from moisture in the original cover carefully resealed.

Dilute **wash buffer (10x)** 1 : 10 with deionized water.

Example: 10 mL wash buffer (10x) + 90 mL deionized water. The prepared wash buffer must be thoroughly mixed before use.

The **substrate** must be protected from direct light. The test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

## Assay Procedure

Performance at room temperature (RT, 18...25 °C). Bring **all** ready-to-use test reagents, prepared wash buffer and microplate to **room temperature**. Mix all reagents before use by gently shaking, avoiding foaming.

The aspiration and washing steps can be carried out manually or with the aid of a microtiter plate washer.

### Important Notes on Test Procedure:

**Avoid** mechanical contact (**scratching**) on the bottom of the wells with pipette tips or washer needles. This will irreparably damage the array!

All liquid reagents (sample, conjugate and substrate) have to be pipetted **without** causing **air bubbles** into the wells.


1. Provide the required number of wells **WELLS** according to the number of samples.
2. Preparation of working dilution of the sample: 1 : 101 dilution  
e.g. 10 µL sample to 1000 µL sample diluent **DIL**
3. Pipette **100 µL** diluted sample per well.
4. Cover plate with covering film **COVER** and incubate for **30 min** at RT.
5. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,  
Tap plate dry onto lint-free absorbent paper.
6. **Add 50 µL** of the conjugate **CONJ HRP IgG** (anti-Human IgG-HRP) to each well.
7. Cover plate with covering film **COVER** and incubate for **30 min** at RT.
8. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,  
Tap plate dry onto lint-free absorbent paper.
9. **Add 50 µL** of the substrate **SUBSTR** to each well.
10. Cover plate and protect from light. Incubate for **30 min** at RT.
11. Remove liquid. Tap plate dry onto lint-free absorbent paper.
12. Clean bottom of wells with swab **SWAB** shortly before scanning the images.
13. **Take images** using the Seramun SpotSight® plate mono / strip und evaluate the results with the Seramun SpotSight® scan software.

If image acquisition is performed with the scanner Seramun SpotSight® strip, transfer single 8-well microtiter plate strips to the scanner's strip holder.

After aspiration of the substrate the color of developed spots is stable for 24 h when the plate is stored protected from light.

## Evaluation of Results

### Arraylayout

### Antigens

- 04 EBNA-1
- 05 VCA p18
- 06 EA-D p54
- 07 VCA gp125

### Controls

- 01 Positive control (PC)
- 02 Cut-off control (CO)
- 03 Negative control (NC)
- 08 IgG conjugate control (GC)
- 09 Serum control (SC)

P Well position marker

### Qualitative evaluation

#### Validity criteria for the test

SeraSpot® Anti-EBV-4 IgG includes the following control spots:

1. Positive control (PC). Intense spot, stained darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than cut-off control
4. Serum control (SC). Intense spot, always stained in presence of sample. Absence of spot indicates absence of sample.
5. IgG conjugate control (GC). Intensively stained spot. Serves as antibody isotype control.

The test cannot be evaluated if one of the quality criteria listed in points 1. to 5. is not met.

If the above mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the validity criteria contact the manufacturer.

## Interpretation of Results

Evaluation of the test has to be performed by using the Seramun SpotSight® plate mono / strip and in combination with the Seramun SeraSpot® scan software.

The results are interpreted as follows:

Evaluation	Conditions
Positive	Color intensity of at least one antigen spot > cut-off control
Negative	Color intensity of the antigen spot ≤ cut-off control

For Interpretation results of the IgG and the IgM determinations should be considered.

<b>IgG</b>	<b>IgM</b>	<b>Interpretation</b>
EBNA-1 negative	Positive	Indication of an acute or primary EBV infection
EBNA-1 positive	Negative	Indication of a past EBV infection
EBNA-1 negative and VCA p18 + VCA gp125 positive	Negative	Indication of a past EBV infection
Negative	Negative	No antibodies against EBV antigens are detectable
All other reactivities*	All other reactivities*	results not clear (control recommended).

\*During EBV reactivation different antibody reactivities may occur. EBV reactivation cannot be diagnosed by a specific antibody reactivity profile.

Determination of 3 parameters (VCA p18 IgG or VCA gp125 IgG, VCA p18 IgM or VCA gp125 IgM and EBNA-1 IgG) allows for the discrimination between acute and past infection in immunocompetent patients.

## Performance Characteristics

### Precision

Samples with known antibody reactivity were examined in the SeraSpot® Anti-EBV-4 IgG and the color intensities (RAW) were determined. These values were used for calculation of the Ratios. The Ratios were then used to determine the coefficients of variation (CV) as a measure of the precision within a test run (intra-assay CV), between different test runs (inter-assay CV) and between different test lots (lot-to-lot CV).

Antigen	Intra-assay coefficient of variation		Inter-assay coefficient of variation		Lot-to-lot coefficient of variation	
	$\bar{x}$ Ratio n=48	CV [%]	$\bar{x}$ Ratio n=80	CV [%]	$\bar{x}$ Ratio n=240	CV [%]
EBNA-1	1.23	8.94	1.30	9.27	1.39	11.58
VCA p18	1.83	8.16	1.63	12.98	1.63	14.61
EA-D p54	1.56	9.03	1.24	19.18	1.10	24.20
VCA gp125	1.04	8.71	1.34	16.92	1.28	18.86
Procedure	1 operator 48 determinations 1 batch		2 operators 2 determinations 2 testings per day 20 days 1 batch		2 operators 2 determinations 2 testings per day 20 days 3 batches	

### Determination of the Cut-off Value

Data from sensitivity and specificity studies compared to reference tests were used to determine the value for the cut-off control by ROC analysis (Receiver Operating Characteristics). The cut-off range is determined specifically for each test.

### Interfering Substances

All tested substances do not have a significant impact on the test results even if present in elevated concentrations in the serum: 20 mg/dL bilirubin C and F each (simulation of icteric samples), 500 mg/dL hemoglobin (simulation of hemolytic samples), 1000 mg/dL lipids (triglycerides) (simulation of lipemic samples) and 500 IU/mL rheumatoid factor.

### Sensitivity and Specificity

Pre-characterized samples (n=130) were analyzed to determine the sensitivity and specificity of SeraSpot® Anti-EBV-4 IgG in comparison with IgG reference tests. Samples with discrepant results were retested in a second reference test.

IgG		Reference test	
		positive	negative
SeraSpot®	positive	83	1
	negative	2	44
Samples total		85	45

Sensitivity: 97.7 %

Specificity: 97.8 %

## Application

### Automated Workflow

Manually diluted positive samples were assayed side by side in SeraSpot® Anti-EBV-4 IgG, by hand and by use of an automatic microplate processor. Data were generated with the device DS2® (Dynex Technologies). The correlation was calculated with  $R^2 > 0.8$  for the antigens EBNA-1, VCA p18, VCA gp125, and EA-D p54. Assay procedure using an automatic microplate processor such as DS2® is possible.

## Change History

Version	Section	Modifications
2022-10_v02_DE_EN	Important Information	Insertion note on dealing with serious incidents
	Entire document	Updating of wording and correction of spelling mistakes

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