# Serazym<sup>®</sup> H. pylori 2nd Gen.

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#### Literatur:

1. Pepleteiro B, Bastos A, Ferro A et al. (2014): Prevalence of *Helicobacter pylori* infection worldwide: a systemic review of studies with national coverage. *Dig Dis Sci* 59: 1698-1709.
2. Lambert, J R, Lin S K, Sievert W, Nicholson L, Schembri M and Guest C (1995): High prevalence of *Helicobacter pylori* antibodies in an institutionalized population: evidence for person-to-person transmission. *Am J Gastroenterol* 90: 2167-2171.
3. Fischbach W, et al. (2016): S2k-Leitlinie *Helicobacter pylori* und gastroduodenale Ulkuskrankheit. *Z Gastroenterol* 54: 327-363.
4. Marshall B J, and Warren J R (1984): Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. *The Lancet* 16: 1311-1314.
5. Graham D Y (2015): Reviews in basic and clinical gastroenterology and hepatology. *Helicobacter pylori* Update: Gastric Cancer, Reliable Therapy, and Possible Benefits. *Gastroenterology* 148: 719-731.
6. Gisbert J P, de la Morena P and Abaira V (2006): Accuracy of Monoclonal Stool Antigen Test for the Diagnosis of *H. pylori* Infection: A Systematic Review and Meta-Analysis. *Am J Gastroenterol* 101: 1921-1930.

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Lagertemperatur

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# Serazym<sup>®</sup> *H. pylori* 2nd Gen.

Enzyme immunoassay for detection of *Helicobacter pylori* antigen in faecal samples

REF E-114-A  96  *In-vitro*-diagnostic device 

Seramun Diagnostica GmbH · Sprehagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · [info@seramun.com](mailto:info@seramun.com)

## Introduction

The human pathogen *Helicobacter pylori* (*H. pylori*) is a gram-negative microaerophilic bacterium that colonizes the gastric mucosa of humans with lifelong persistence. It is estimated that about 50 % of the world population is infected with *H. pylori* with varying prevalences in dependence of the geographic region and socio-economic factors (1). Infection via oral-oral or fecal-oral route is assumed as the most likely mode of transmission (2). Besides asymptomatic infections, *H. pylori* causes chronic active gastritis in about 10 % of infected persons. Gastroduodenal ulcer, gastric adenocarcinoma or MALT Lymphoma may arise as secondary complications (3). The pathogenic significance of *H. pylori* for the development of gastritis and peptic ulcer has only been realized by the research of Warren and Marshall in 1984, which has been awarded with the Nobel Prize in 2005 (4). The WHO has even classified *H. pylori* as a class I pathogen.

The pathogenic effectiveness of *H. pylori* specific extra cellular products like Urease, Vac A (vacuolating cytotoxin A), and Cag A (cytotoxin coded by the cytotoxin associated gen A and injected into the host cell by the type IV secretion apparatus) is based on direct epithel cell damage. Additionally the infection causes a chronic inflammation with granulocyte infiltration and increased production of inflammatory mediators like Interleucin 8 (5).

Different invasive (culture, histology, urease rapid test, PCR from biopsy material) and non-invasive (urea breath test, stool antigen test, antibody detection in serum) diagnostic methods are available. The "S2k-Leitlinie Helicobacter pylori und gastroduodenale Ulkuskrankheit" recommends to perform tests which detect an actual *H. pylori* infection by direct detection of the pathogen or its metabolites, e.g. histology, culture, PCR, stool antigen test, urease test or urea breath test (3). In the case of gastrointestinal bleedings, gastric cancer or after therapy with proton pump inhibitors, these methods may be impaired and serological antibody detection is the method of choice (3).

For direct antigen detection from stool specimens, ELISA tests based on monoclonal antibodies have been established and are preferable to polyclonal assays (3, 6).

*Serazym*<sup>®</sup> H. pylori 2nd Gen. uses highly specific monoclonal antibodies for rapid and specific detection of an infection with *Helicobacter pylori*.

#### References:

1. Pepleteiro B, Bastos A, Ferro A et al. (2014): Prevalence of *Helicobacter pylori* infection worldwide: a systemic review of studies with national coverage. *Dig Dis Sci* 59: 1698-1709.
2. Lambert, J R, Lin S K, Sievert W, Nicholson L, Schembri M and Guest C (1995): High prevalence of *Helicobacter pylori* antibodies in an institutionalized population: evidence for person-to-person transmission. *Am J Gastroenterol* 90: 2167-2171.
3. Fischbach W, et al. (2016): S2k-Leitlinie *Helicobacter pylori* und gastroduodenale Ulkuskrankheit. *Z Gastroenterol* 54: 327-363.
4. Marshall B J, and Warren J R (1984): Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. *The Lancet* 16: 1311-1314.
5. Graham D Y (2015): Reviews in basic and clinical gastroenterology and hepatology. *Helicobacter pylori* Update: Gastric Cancer, Reliable Therapy, and Possible Benefits. *Gastroenterology* 148: 719-731.
6. Gisbert J P, de la Morena P and Abaira V (2006): Accuracy of Monoclonal Stool Antigen Test for the Diagnosis of H. pylori Infection: A Systematic Review and Meta-Analysis. *Am J Gastroenterol* 101: 1921-1930.

#### Intended Use

The *Serazym*<sup>®</sup> H. pylori 2nd Gen. is an *in-vitro*-diagnostic device for direct detection of *Helicobacter pylori* antigen in faecal samples.

#### Principle Of The Test

*Serazym*<sup>®</sup> H. pylori 2nd Gen. is an indirect two-site-immunoassay for the qualitative determination of *Helicobacter pylori* antigen based on monoclonal antibodies.

Within the first incubation step (60 min, room temperature) diluted stool specimens as well as positive and negative controls react with the solid-phase adsorbed monoclonal antibodies. Unbound components are removed by a subsequent washing step. Within the second incubation step (30 min, room temperature) solid-phase bound immune complexes react with the horseradish (HRP)-labelled monoclonal antibodies of the conjugate. Unbound reagents are separated from the solid-phase antibody-antigen-antibody immune complexes by a further washing step. HRP converts the subsequently added colourless chromogenic substrate solution into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells after 10 min incubation at room temperature turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 /  $\geq$  620 nm is directly proportional to the specifically bound amount of *Helicobacter pylori* antigen. Referring to the cut-off value results are interpreted as positive or negative.

## Test Components

		For 96 Wells	
1	<b>WELLS</b>	<b>Microtitration plate</b> coated with monoclonal anti- <i>H. pylori</i> antibodies (mouse)	12 single breakable 8-well strips colour coding wine red vacuum-sealed with desiccant
2	<b>WASHBUF CONC 10x</b>	<b>Wash buffer</b> 10-fold	100 ml concentrate for 1000 ml solution white cap
3	<b>DIL</b>	<b>Sample diluent</b>	100 ml · ready to use coloured yellow black cap
4	<b>CONTROL +</b>	<b>Positive control</b> <i>H. pylori</i> antigen positive sample (inactivated)	2.0 ml · ready to use coloured blue red cap
5	<b>CONTROL -</b>	<b>Negative control</b> <i>H. pylori</i> antigen negative sample	2.0 ml · ready to use coloured blue green cap
6	<b>CONJ HRP</b>	<b>HRP-conjugate</b> HRP-labelled monoclonal anti- <i>H. pylori</i> antibodies (mouse)	15 ml · ready to use coloured green violet cap
7	<b>SUBSTR TMB</b>	<b>Substrate</b> 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap
8	<b>STOP</b>	<b>Stop solution</b> 0.25 M sulphuric acid	15 ml · ready to use yellow cap

## Preparation And Storage Of Samples

### Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 hours. Longer storage is possible at -20°C. Repeated (> 3x) freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*<sup>®</sup> sample diluent can be stored for up to 72 h at 2...8°C before testing in the ELISA.

### Preparation

Mix samples thoroughly. Pipette 500 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 100 mg (diameter about 2 - 3 mm) of faeces if solid or pipette 100 µl if liquid into the tube and suspend thoroughly. Spin down floating particles in a micro centrifuge at maximum speed for 1 min if necessary.

## Materials Required But Not Provided

Micropipettes · multi-channel pipette or multi-pipette · Reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with 450 nm filter for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation

## Preparation And Storage Of Reagents

### Kit size and expiry

One kit is designed for 96 determinations. The expiry date of each component is reported on its respective label and that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least one month when stored at 2...8°C.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

## Assay Procedure

Dilute samples with sample diluent (3) 1 : 6, e.g. 100 mg or 100 µl faeces + 0.5 ml sample diluent (3).

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash solution in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid direct light exposure of the TMB substrate solution!

### Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CONTROL +** positive control (4)  
100 µl **CONTROL -** negative control (5)  
100 µl **diluted sample**.
3. Cover plate and incubate for 60 min at RT.
4. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
5. 3 drops (or 100 µl) **CONJ HRP** HRP-conjugate (6) per well.
6. Cover plate and incubate for 30 min at RT.
7. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
8. 3 drops (or 100 µl) **SUBSTR TMB** substrate (7) per well.
9. Incubate for 10 min at RT protected from light.
10. 3 drops (or 100 µl) **STOP** stop solution (8) per well, mix gently.
11. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

## Result Interpretation

### Qualitative evaluation

#### Cut-off determination: OD negative control + 0.10

Samples with OD values equal with or higher than the cut-off are considered positive, samples with OD values below the cut-off are considered negative for *Helicobacter pylori* antigen.

## Reference Values

<b>Serazym® H. pylori 2nd Gen.</b>	
Positive	≥ Cut-off
Negative	< Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

### Test validity

The test run is valid if:

- the mean OD of the negative control is ≤ 0.20 (manual performance)  
≤ 0.30 (automatic performance)
- the mean OD of the positive control is ≥ 1.20

If the above mentioned quality criteria are not met, repeat the test and make sure the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps, mixing after addition of controls and samples to the biotin-conjugate in the first incubation step etc.). In case of repeated failure of the quality criteria contact your supplier.

### Limitations of the procedure


There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the suspension can cause false negative as well as false positive results. A negative test result in the *Serazym®* H. pylori 2nd Gen. does not exclude an infection.

### Automatic Processing

Performing the *Serazym®* H. pylori 2nd Gen. on fully automated microplate processors (e.g. DS2 or DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x - 8x.

### Correlation: manual – automatic processing

A panel of 93 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with  $r = 0.928$ .

## Performance Characteristics

### Precision

Intra-assay coefficient of variation (CV) in the *Serazym*<sup>®</sup> H. pylori 2nd Gen. calculated from 24-fold determinations of samples:

sample	mean absorbance	standard deviation	CV (%)
1	2.377	0.126	5.3
2	0.871	0.065	7.4
3	0.124	0.014	11.0

Lot-to-Lot-assay coefficient of variation (CV) in the *Serazym*<sup>®</sup> H. pylori 2nd Gen. in 3 different batches and test runs with 3-fold determinations of samples:

sample	mean absorbance	standard deviation	CV (%)
1	2.626	0.284	10.4
2	1.265	0.138	17.4
3	0.162	0.030	18.5

Inter-assay coefficient of variation (CV) in the *Serazym*<sup>®</sup> H. pylori 2nd Gen. in 10 different test runs with 3-fold determinations of samples:

sample	mean absorbance	standard deviation	CV (%)
1	2.483	0.157	6.3
2	1.044	0.138	13.2
3	0.142	0.015	10.3

### Lower detection limit


The lower detection limit of the *Serazym*<sup>®</sup> H. pylori 2nd Gen. has been determined by titration of faecal samples spiked with purified *H. pylori* antigen.

The lower detection limit is  $< 31.3$  ng/ml.

### Determination of the cut-off value

The frequency distribution of the absorbances of 169 faecal samples has been investigated in order to fix the cut-off value for the *Serazym*<sup>®</sup> H. pylori 2nd Gen.. The cut-off has been determined with absorbance negative control + 0.1.

### Specificity and sensitivity

Specificity and sensitivity of the *Serazym*<sup>®</sup> H. pylori 2nd Gen. have been determined in a retrospective study in comparison to another commercially available ELISA.

	comparative ELISA positive	comparative ELISA negative
<i>Serazym</i> <sup>®</sup> ELISA positive	140	6**
<i>Serazym</i> <sup>®</sup> ELISA negative	2*	85

Specificity: 93.4%    Sensitivity: 98.6%

The samples labeled with \* and \*\* were re-tested in a commercial available Lateral Flow Assay. This results in a corrected specificity of 97.7% and a corrected sensitivity of 99.3%

## Cross reactivity

The following microorganisms were diluted in sample diluent with  $\geq 10^8$  colony forming units and tested negative with the *Serazym*<sup>®</sup> H. pylori 2nd Gen. (OD 450 / 620 nm < Cut-Off):

<i>Aeromonas hydrophila</i>	ATCC (7966)
<i>Bacillus cereus</i>	ATCC (11778)
<i>Bacillus subtilis</i>	ATCC (6633)
<i>Bacteroides fragilis</i>	ATCC (25285)
<i>Citrobacter freundii</i>	ATCC (8090)
<i>Clostridium sordelli</i>	ATCC (9714)
<i>Enterobacter aerogenes</i>	ATCC (13048)
<i>Enterobacter cloacae</i>	ATCC (13047)
<i>Enterococcus faecalis</i>	ATCC (29212)
<i>Escherichia coli</i>	ATCC (25922)
<i>Klebsiella pneumonia</i>	ATCC (13883)
<i>Peptostreptococcus anaerobius</i>	ATCC (27337)
<i>Proteus vulgaris</i>	ATCC (8427)
<i>Pseudomonas aeruginosa</i>	ATCC (10145)
<i>Salmonella enterica</i> Serovar <i>thyphimurium</i>	ATCC (14028)
<i>Salmonella enterica</i> ssp. <i>galolyticus/enteritidis</i>	ATCC (13076)
<i>Shigella flexneri</i>	ATCC (12022)
<i>Shigella sonnei</i>	ATCC (25931)
<i>Staphylococcus aureus</i>	ATCC (25923)
<i>Staphylococcus epidermidis</i>	ATCC (12228)
<i>Vibrio parahaemolyticus</i>	ATCC (17802)
<i>Candida albicans</i>	ATCC (10231)

<i>Campylobacter coli</i>	ATCC (33559)
<i>Campylobacter jejunii</i>	ATCC (32291)
<i>Campylobacter fetus</i>	ATCC (27374)
<i>campylobacter upsaliensis</i>	ATCC (43954)
<i>Campylobacter lari</i>	ATCC (35221)
<i>Vibrio cholerae</i>	clinical isolate
<i>Yersinia enterocolitica</i> O:3	clinical isolate
<i>Yersinia enterocolitica</i> O:9	clinical isolate
<i>Yersinia enterocolitica</i> Y11	clinical isolate
<i>Yersinia enterocolitica</i> RKI 0803733	clinical isolate
<i>Clostridium difficile</i>	VPI 10463
<i>Salmonella infantis</i>	ATCC (51741)
<i>Salmonella anatum</i>	ATCC (9270)
<i>Salmonella paratyphi</i> A	ATCC (11511)
<i>Salmonella paratyphi</i> B	ATCC (8759)
<i>Salmonella paratyphi</i> C	Nr. 2 Pasteur
<i>Lactococcus lactis</i>	DSM (20481)
<i>Proteus mirabilis</i>	ATCC (29906)
<i>Pseudomonas fluorescens</i>	ATCC (13525)
<i>Pseudomonas putida</i>	ATCC (49128)
<i>Streptococcus agalactiae</i>	ATCC (13813)
<i>Morganella morganii</i>	ATCC (25830)

## Interference

None of the following substances added to *H. pylori* positive and negative stool samples showed a significant impact on the test result: The mentioned concentrations refer to 100 mg faeces.

(-)-Scopolamine N-butyl bromid (0.5%, Buscopan<sup>®</sup>), Barium sulfate (5%), Bismuth(III) subsalicylate (0.5%, Pepto-Bismol), Cyclamat (5%), Diclofenac (0.5%), Hemoglobin human (5%), Blood human (5%), Hylak<sup>®</sup> N (5%), Iberogast<sup>®</sup> (5%), Immodium<sup>®</sup> akut duo (0.03/1.9%), Loperamid hydrochloride (5%), Loperamid-CT akut), Metronidazole (0.5%), Mucin (5%), Nexium<sup>®</sup> (0.03%), Nifuroxazide (0.5%, Pentofuryl<sup>®</sup>), Palmitic acid (20%), Perenterol forte (0.5%), Rennie<sup>®</sup> (20%), Simigel<sup>®</sup> (1%), Stearic acid (20%), Vancomycin (0.5%)

## Common Advices and Precautions

**This kit is for *in-vitro* use only.** Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed.

**Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

**The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*<sup>®</sup> stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), Helicobacter pylori (E-081), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).**

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents (2 and 7) contain small amounts of Thimerosal (0.01% w/v) and Kathon (1.0% v/v) as preservative. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

**Do not smoke, eat or drink while handling kit material!**

**Always use protective gloves!**

**Never pipette material by mouth!**

**Note safety precautions of the single test components!**

## Incubation Scheme *Serazym*<sup>®</sup> H. pylori 2nd Gen. (E-114)

1.  pipette  
 100 µl **CONTROL +** (4)  
 100 µl **CONTROL -** (5)  
 100 µl **diluted stool sample**  
 60 min incubation (room temperature)  
 5 x wash with wash solution
2.  3 drops (or 100 µl) **CONJ HRP** (6)  
 30 min incubation (room temperature)  
 5 x wash with wash solution
3.  3 drops (or 100 µl) **SUBSTR TMB** (7)  
 10 min incubation (room temperature) protected from light
4.  3 drops (or 100 µl) **STOP** (8)

**Read OD at 450 / ≥ 620 nm**

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use