# Seraline<sup>®</sup> HepAk-7 IgG

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Lagertemperatur

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# Seraline<sup>®</sup> HepAk-7 IgG

Line Immunoassay for detection of IgG-autoantibodies in autoimmune liver diseases  
in human serum or plasma

REF LIA-004-7 G  20  *In-vitro*-Diagnostic device 

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
Telefon +49 (0) 33767 79110 · Fax +49 (0) 33767 79199 · [info@seramun.com](mailto:info@seramun.com)

## Introduction

Primary autoimmune liver disease (PAL) can be subdivided into autoimmune hepatitis (AIH) type I to III, primary biliary cirrhosis (PBC) and primary sclerosing cholangitis (PSC). Beside these clearly defined diseases, variant forms described as overlap syndrome additionally exist. These forms include characteristics of AIH and PBC or PSC. Clinical symptoms of PAL resemble that of other chronic liver diseases. The progress of chronic liver diseases may result in the development of liver cirrhosis. Only an early diagnosis and treatment may prevent such an outcome. In early stages liver diseases are characterized by increased levels of serum transaminases. Unspecific symptoms as sickness, fatigue and arthralgia can additionally occur. About 15% of chronic liver disease cases underlie an autoimmune pathogenesis. Therefore the determination of different autoantibodies is recommended after exclusion of a virus hepatitis. The early therapy as a result of a correct diagnosis correlates with a positive prognosis.

The *Seraline*<sup>®</sup> HepAk-7 IgG test enables the fast determination of autoantibodies with diagnostic relevance for autoimmune liver diseases.

Autoantigen	Clinical relevance / specificity	Description
AMA/M2	PBC (95%)	antigens in the mitochondrial membrane, mainly the M2 protein of the pyruvatdehydrogenase complex. AMA/M2
LKM1	AIH II (95 - 100%)	epitopes of the Cytochrome p450 protein.
LC1	AIH II (50%)	Cytosolic 62 kDa antigen
SLA	AIH III	50 kDa soluble liver antigen
F-Aktin	AIH I less specific	polymerized filaments of the Actin protein.
gp210	PBC (10 - 25%) Highly specific	membrane protein of nuclear pore complex
Sp100	PBC (30%)	Protein of "nuclear dots"

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## Intended use

The *Seraline*<sup>®</sup> HepAk-7 IgG is a Line Immunoassay (LIA) for detection of autoantibodies of the IgG isotype directed against the following antigens: AMA/M2, LKM1, LC1, SLA, F-Actin, gp210 and Sp100 in human serum or plasma samples.

## Principle of the test

The antigens are plotted on nitrocellulose membranes and blocked to avoid unspecific reactions. After this procedure the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

### Step 1

The test strips are incubated in the diluted (1 : 101) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

### Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

### Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates. The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template. Alternatively the strips can be evaluated with the *Seraline*<sup>®</sup>scan software.

## Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 101 (15 µl sample and 1500 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

## Kit components

1	<b>TESTSTR</b>	<b>test strips incl. evaluation template</b>	<b>20 strips</b> colour coded: green
2	<b>WIB CONC 5x</b>	<b>wash and incubation buffer 5-fold</b>	<b>70 ml</b> concentrate, transparent bottle, black cap
3	<b>CONJ HRP IgG</b>	<b>anti-human IgG-HRP-conjugate</b>	<b>35 ml</b> , ready to use, <b>IgG</b> , coloured red, transparent bottle, red cap
4	<b>SUBSTR TMB</b>	<b>substrate</b> 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	<b>35 ml</b> , ready to use, black bottle, blue cap
5	<b>INCUTRAY</b>	<b>incubation tray with cover</b>	<b>2</b>

## Materials required but not provided

- glassware
- pipettes and tips for 15 µl and 1500 µl
- plastic forceps
- rocking platform (vertical)
- filter paper
- collecting devices for infectious material
- distilled or deionized water

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 20 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

## Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

## Working steps

1. Incubate the test strips with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **15 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **1.5 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **1.5 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **1.5 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

### Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.


## Test evaluation

### Test validity

The *Seraline*<sup>®</sup> HepAk-7 IgG test strips contain 3 consecutive control lines under the test strip number:

- a) Function control (positive reaction with every sample).
- b) The intensity of the negative control must be less than the Cut off control.
- c) Cut off control (band intensity is used for result evaluation of the diagnostic bands).

The test run is valid if the function control and the Cut off control are clearly visible.

## Result interpretation

Generally the evaluation of the test has to be done with dried strips. To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet. Alternatively the strips can be evaluated with the *Seraline*<sup>®</sup> scan software. Judgement of the bands is performed according to the following classification:

Result interpretation	IgG
negative	Colour intensity of bands < Cut off control
positive	Colour intensity of bands ≥ Cut off control

## Performance characteristics

### Sensitivity

73 Sera from patients with clinically defined autoimmune liver diseases have been tested with the *Seraline*<sup>®</sup> HepAk-7 IgG and other commercially available tests. The results for the different antigens are shown in the following tables.

#### AMA/M2

		<i>Seraline</i> <sup>®</sup> HepAk-7 IgG		Agreement: 85%
		positive	negative	
competitor	positive	39	5	
	negative	6	23	

#### LKM1

		<i>Seraline</i> <sup>®</sup> HepAk-7 IgG		Agreement: 96%
		positive	negative	
competitor	positive	3	0	
	negative	3	67	

#### LC1

		<i>Seraline</i> <sup>®</sup> HepAk-7 IgG		Agreement: 100%
		positive	negative	
competitor	positive	6	0	
	negative	0	67	

#### SLA

		<i>Seraline</i> <sup>®</sup> HepAk-7 IgG		Agreement: 97%
		positive	negative	
competitor	positive	8	1	
	negative	1	63	

#### F-Aktin

		<i>Seraline</i> <sup>®</sup> HepAk-7 IgG		Agreement: 92%
		positive	negative	
competitor	positive	4	4	
	negative	2	63	

gp210

		<i>Seraline</i> <sup>®</sup> HepAk-7 IgG	
		positive	negative
competitor	positiv	16	0
	negativ	2	55

Agreement: 97%

SP100

		<i>Seraline</i> <sup>®</sup> HepAk-7 IgG	
		positive	negative
competitor	positiv	6	0
	negativ	0	67

Agreement: 93%

Serum samples from 341 patients with PBC were tested for AMA/M2, Sp100 and gp210 antigens in comparison with 3 different competitor test kits. The estimated sensitivity was 81 % for the *Seraline*<sup>®</sup> HepAk-7 IgG versus 71 %, 81 % and 70 % for the competitor kits.

### Specificity

The screening of 120 serum samples from blood donors and additional 40 serum samples of patients with virus induced hepatitis (HBV, HCV) resulted in an estimated specificity of 92%. In only two serum samples LKM1 and F-Actin specific antibodies were detected, respectively.

### Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other Immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

### Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

## Common advices and precautions

**This kit is for *in-vitro* use only.** Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only. Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

**Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

## History of changes

Version	Section	Modifications
2018-03-15	Common advices and precautions	Update

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use