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Verfallsdatum


Lagertemperatur

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Seraline[®] Anti-Borrelia-8 IgG Seraline[®] Anti-Borrelia-8 IgM

Line Immunoassay for detection of IgG or IgM antibodies to
Borrelia burgdorferi sensu lato in human serum or plasma

 LIA-006-8 G	 20	 LIA-006-8 G-12	 12x 20	 <i>In-vitro</i> -Diagnostikum	
 LIA-006-8 M	 20	 LIA-006-8 M-12	 12x 20	 <i>In-vitro</i> -Diagnostikum	

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Lyme borreliosis is a systemic infectious disease with various clinical manifestations. The disease can be divided into 3 stages involving different organ systems. Cardinal symptoms of the first stage (4 to 8 weeks p. i.) are the erythema chronicum migrans and local or general lymphadenopathies (lymphadenosis cutis benigna). The clinical manifestations of the second stage (1 to 12 months p. i.) range from meningitis, meningopolyneuritis, encephalitis up to hemiparesis, myoneuralgia and arthralgia, especially knee joint arthritis. Cardiac manifestations like life-threatening myocarditis / pancarditis are rarely detected. The third stage (months to years p. i.) is characterized by chronic infection of the nervous system (neuroborreliosis, progressive encephalomyelitis), the skin (acrodermatitis chronica atrophicans) and the joints (chronic erosive arthritis). Especially late manifestations of Lyme borreliosis can reduce quality of life and antibiotic treatment is difficult. Therefore an early diagnosis of borrelia infections is of great importance.

The infectious agent of Lyme disease was isolated from ticks in 1982 by BURGDORFER. It belongs to spirochetes and represents the independent genus *Borrelia*. The tick species *Ixodes ricinus* is the main vector of *Borrelia burgdorferi* in Europe. The occurrence of ticks within densely wooded rural areas correlates with the frequency of Lyme disease cases with a peak in the summer and autumn months. While the tick is taking blood from the host the *Borrelia* are transmitted from the contents of the tick's intestine. Persons at risk for *Borrelia* infections are people who work in forestry or agriculture, hunters and campers, soldiers and all people who occasionally spend time in wooded areas. The incidence of Lyme disease reaches from 2 to 40 / 100 000 in middle Germany to 300 / 100 000 in Austria. The estimated number of unrecognized cases is presumably higher, because a tick bite is not always recognized and the typical erythema chronicum migrans only develops in about 50% of *Borrelia* infections.

The direct detection of the infectious agent via culture from blood, cerebral spinal fluid (CSF) or skin biopsies undoubtedly represents the safest method for diagnosis of Lyme disease. Unfortunately this method is hampered by the long reproduction time of the *Borrelia*, the complexity of the culture medium and the relatively low sensitivity in case of culture from patient material. Thus culture as well as *Borrelia* detection from tissue sections via immunohistological methods are not suitable for routine diagnosis. Methods of choice are detection of IgG and/or IgM antibodies via immunofluorescence or enzyme immunoassay (ELISA). The ELISAs use *Borrelia* sonicates, partially purified or recombinant antigens for plate coating. Because of the regionally differing occurrence of subtypes (genospecies) and the known variability of the cell surface proteins of *Borrelia burgdorferi*, antigen mixtures are often preferred. Confirmatory tests like western blot, dot blot or line assay are essential and the serological test results have to be interpreted together with the clinical picture. Negative test results do not exclude Lyme disease (early stage of infection, seronegative cases).

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Intended use

The *Seraline*[®] Anti-Borrelia-8 IgG / *Seraline*[®] Anti-Borrelia-8 IgM is an in vitro diagnostic device (Line Immunoassay, LIA) for the detection of specific IgG or IgM antibodies against *Borrelia burgdorferi sensu lato* in human serum or plasma samples. It should be used as a confirmatory test after a positive ELISA result.

Principle of the test

The antigens are plotted on nitrocellulose membranes and blocked to avoid unspecific reactions. After this procedure the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

Step 1

The test strips are incubated in the diluted (1 : 101) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates. The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template. Alternatively the strips can be evaluated with the *Seraline*[®] scan software.

Used Borrelia antigens

Nomenclature*	Properties / Specificity	Origin
p100	Protein of surface membrane-vesicles	high specific <i>B. afzelii</i>
VlsE	variable major protein-like sequence Expressed	specific <i>B. afzelii</i>
p58	not characterized	specific <i>B. garinii</i>
p41 (Flagellin)	structural protein of endoflagellin	unspecific <i>B. burgdorferi sensu stricto</i>
p39 (BmpA)	Flagella complex Borrelia membrane protein A	specific <i>B. afzelii</i>
p31 (OspA)	Outer surface protein A	specific <i>B. afzelii</i>
p23 (OspC)	Outer surface protein C, Mix from different Borrelia subtypes	high specific <i>B. burgdorferi sensu stricto</i> , <i>B. garinii</i> , <i>B. garinii II</i> , <i>B. afzelii</i> , <i>B. spielmanii</i>
p18 (DbpA)	Decorin binding protein A	specific <i>B. burgdorferi sensu stricto</i> , <i>B. garinii</i> , <i>B. afzelii</i> , <i>B. spielmanii</i>

* The nomenclature for some bands is not consistently used in the literature concerning molecular mass and specificity, respectively (2, 5, 7)

Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 101 (15 µl sample and 1500 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

Kit components

			for 20 determinations	for 12 x 20 determinations
1	TESTSTR	test strips incl. evaluation template	20 strips colour coded: dark blue	12 x 20 strips colour coded: dark blue
2	WIB CONC 5x	wash and incubation buffer 5-fold	70 ml concentrate, transparent bottle, black cap	12 x 70 ml concentrate, transparent bottle, black cap
3	CONJ HRP IgG CONJ HRP IgM	anti-human IgG-HRP-conjugate or anti-human IgM-HRP-conjugate	35 ml, ready to use, IgG, coloured red, transparent bottle, red cap IgM, coloured green, transparent bottle, green cap	12 x 35 ml, ready to use, IgG, coloured red, transparent bottle, red cap IgM, coloured green, transparent bottle, green cap
4	SUBSTR TMB	substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	35 ml, ready to use, black bottle, blue cap	12 x 35 ml, ready to use, black bottle, blue cap
5	INCUTRAY	incubation tray with cover	2	-

Materials required but not provided

• glassware • pipettes and tips for 15 µl and 1500 µl • plastic forceps • rocking platform (vertical)
 • filter paper • collecting devices for infectious material • distilled or deionized water • incubation trays (▽ 12 x 20)

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 20 or 12 x 20 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

Working steps

1. Incubate the test strips with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **15 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **1.5 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **1.5 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **1.5 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.


Test evaluation

Test validity

The *Seraline*[®] Anti-Borrelia-8 IgG / *Seraline*[®] Anti-Borrelia-8 IgM test strips contain 5 consecutive control lines beneath the strip number:

- Function control (positive reaction with every sample).
- Conjugate controls IgG and IgM (positive reaction with the corresponding conjugate).
Note: The other conjugate control band may develop a weak, non-specific colour.
- Negative control (band intensity must be less than the Cut off control).
- Cut off control (band intensity is used for result evaluation of the diagnostic bands).

The test run is valid if the function control, the corresponding conjugate control and the Cut off control are clearly visible.

Result interpretation

Generally the evaluation of the test has to be done with dried strips. To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet. Alternatively the strips can be evaluated with the *Seraline*[®] scan software. Judgement of the bands is performed according to the following classification:

Result interpretation	IgG	IgM
negative	max. 1 band except for VlsE > Cut off control	no bands except for p41 > Cut off control
borderline	VlsE or 1 band + p41 ≥ Cut off control	1 band except for OspC or p41 or p18 ≥ Cut off control
positive	2 bands except for p41 ≥ Cut off control	OspC or p18 or 2 other bands ≥ Cut off control

Inverse bands (bright bands on dark background) should be evaluated as negative. The respective serum should always be examined using other serological methods.

The cross-reaction between *Borrelia* and other spirochaetes may lead to the occurrence of *Borrelia*-associated bands in the Line Immunoassay, what may finally cause a false positive result. Samples of patients with the following infections may cross-react: Syphilis (*Treponema pallidum*), Framboesie (*Treponema pertenue*), Leptospirosis (*Leptospira spec.*) and relapsing fever (*Borrelia spec.*).

Within the context of EBV primary infections antibody reactivities against *Borrelia burgdorferi* sensu lato antigens due to polyclonal b-cell stimulation, may occur.

Performance characteristics

Sensitivity

Panel 1:

Serum samples (n = 78) from patients with suspect of Lyme Borreliosis were investigated in a comparative study using *Seraline*[®] Anti-Borrelia-8 IgG / *Seraline*[®] Anti-Borrelia-8 IgM and another commercial Line Immunoassay.

n = 78		commercial – Line Immunoassay	
		positive	negative
<i>Seraline</i> [®] Anti-Borrelia-8	positive	51	10
	negative	1	16

A total of 67 samples delivered concordant results in both assays which means an overall agreement of 85.9% and a comparative sensitivity of 98.1% for the *Seraline*[®] Anti-Borrelia-8 IgG/IgM in relation to the commercial assay. Samples with discordant results (n=11) were reinvestigated by Immunoblot (*Serablot*[®] Anti-Borrelia IgG/ IgM and another commercial Line Immunoassay). The *Seraline*[®] results were confirmed by Immunoblot for 8 of 10 positive samples. One discrepant sample (*Seraline*[®] negative, comparative test positive) was determined weak positive by Western Blot (*Serablot*[®]).

Panel 2:

Serum samples from persons with suspect of acute Borrelia infection (n = 32) were investigated for the presence of IgM antibodies to *Borrelia burgdorferi* by use of the *Seraline*[®] Anti-Borrelia-8 IgM and another commercially available Line Immunoassay.

n = 32		commercial – Line Immunoassay IgM	
		positive	negative
<i>Seraline</i> [®] Anti-Borrelia-8 IgM	positive	9	8
	negative	0	15

24 samples show identical results. In comparison to the commercial assay a sensitivity of 100% is achieved. All 8 samples showing discrepant results were reinvestigated by Immunoblot (*Serablot*[®] Anti-Borrelia IgM). Six of these samples were confirmed positive.

Specificity

A population of 400 healthy blood donors was investigated to determine the specificity of *Seraline*[®] Anti-Borrelia-8 IgG and *Seraline*[®] Anti-Borrelia-8 IgM Line Immunoassays.

Out of this population 47 (39 positive and 8 borderline) serum samples were determined positive for IgG and 22 (17 positive and 5 borderline) serum samples were determined positive for IgM antibodies with the *Seraline*[®] Anti-Borrelia-8 IgG and *Seraline*[®] Anti-Borrelia-8 IgM resp. resulting in a specificity of

Specificity IgG: 88.2% Specificity IgM: 94.5%

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other Immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

Cross reactivity

Serum samples (n = 195) from patients with potential cross reactive antibodies were investigated in a study using *Seraline*[®] Anti-Borrelia-8 IgG / *Seraline*[®] Anti-Borrelia-8 IgM.

The result evaluation of IgG and IgM antibody detection is shown in the following table.

n=195		<i>Seraline</i> [®] Anti-Borrelia-8		
		positive	borderline	negative
ANA	(n=51)	7	3	41
EBV	(n=30)	3	0	27
CMV	(n=67)	13	2	52
Syphilis	(n=47)	4	1	42

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only. Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

History of changes

Version	Section	Modifications
2018-03-15	Common advices and precautions	Update

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use