Seraline[®] Anti-Helicobacter-6 IgG Seraline[®] Anti-Helicobacter-6 IgA

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Verfallsdatum

Lagertemperatur

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Seraline[®] Anti-Helicobacter-6 IgG Seraline[®] Anti- Helicobacter-6 IgA

Line Immunoassay for detection of IgG or IgA antibodies to *Helicobacter pylori*
in human serum or plasma

 REF	LIA-007-6 G	 20	 IVD	<i>In-vitro</i> -Diagnostikum	 CE
 REF	LIA-007-6 A	 20	 IVD	<i>In-vitro</i> -Diagnostikum	 CE

Seramun Diagnostica GmbH · Sprenthagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

The human pathogen *Helicobacter pylori* (*H. pylori*) is a gram-negative spiral shaped bacterium that has been found in the stomachs of humans in all parts of the world. In developing countries, 70 to 90% of the population carries *H. pylori* mostly already acquired before the age of 10 years. In developed countries the prevalence of infection ranges from 25 to 50%. *H. pylori* infections are transmitted via the oral-oral, the faecal-oral route or iatrogenic. A substantial reservoir for *H. pylori* aside from the human stomach has not been confirmed so far (1).

H. pylori is characterized by a strong urease activity and some strains (type I) additionally produce toxins like VacA (vacuolating toxin) and CagA (cytotoxin-associated antigen A) (2, 3). Further studies have identified HcpC (Helicobacter cysteine-rich protein C) and GroEL (Chaperonin) as new independent virulence factors (4, 5), whereas GroEL is more important than HcpC.

These extracellular products contribute to the pathogenesis by direct damage of the gastric epithelium accompanied by a chronic inflammation with enhanced levels of inflammation mediators like interleukin 8, interleukin 1 or tumor necrosis factor alpha. Patients infected with *H. pylori* often develop autoantibodies to parietal cells. The key-role of these autoantibodies in the development of atrophic gastritis as precursor of gastric cancer has been reviewed (1, 6).

About 10% of the infected persons develop *H. pylori* associated gastritis and secondary diseases like chronic active gastritis, Ulcus ventriculi, Ulcus duodeni, gastric cancer and MALT lymphoma. Patients suffering from gastric malignomas and Ulcus duodeni are infected with *H. pylori* in nearly 100% of the cases. About 80% of patients with chronic gastritis, 70% of patients with Ulcus ventriculi and 60% of patients with gastric cancer are infected with *H. pylori* (6).

Diagnosis is usually based on gastroendoscopy combined with the detection of the pathogen in biopsy material by culture, histology, PCR or rapid urease test. Culture from biopsy material is difficult and not always successful. The ¹³C breath test which is often performed as follow up test, detects CO₂ released by the bacterial urease from radioactive labeled urea in patients breath. This method is not invasive but the need for special equipment and the uptake of radioactive urea by the patients are disadvantageous.

Immunological tests on the basis of specific anti-*H. pylori* antibodies (mono- or polyclonal) are available. They enable the direct detection of *H. pylori* antigens from faecal specimens and may be used for therapeutic surveillance.

Serology is the third method commonly used as a non-invasive method to diagnose *H. pylori* infection. Furthermore, serology is the only test that is not affected by local changes in the stomach that could lead to false-negative results in the other tests.

Invasive as well as non-invasive tests for detection of *H. pylori* infections contain methodological shortcomings. Therefore the S3 guidelines of the DGVS (8) recommend to include at least two positive test results or a positive culture to ensure a safe diagnosis and before initiating drug therapy.

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Intended use

The *Seraline*[®] Anti-Helicobacter-6 IgG / *Seraline*[®] Anti-Helicobacter-6 IgA is a Line Immunoassay (LIA) for detection of specific IgG or IgA antibodies against *Helicobacter pylori* in human serum or plasma samples.

Principle of the test

The antigens are plotted on nitrocellulose membranes and blocked to avoid unspecific reactions. After this procedure the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

Step 1

The test strips are incubated in the diluted (1 : 101) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates. The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template. Alternatively the strips can be evaluated with the *Seraline*[®]scan software.

Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 101 (15 µl sample and 1500 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

Kit components

1	TESTSTR	test strips incl. evaluation template	20 strips colour coded: brown
2	WIB CONC 5x	wash and incubation buffer 5-fold	70 ml concentrate, transparent bottle, black cap
3	CONJ HRP IgG CONJ HRP IgA	anti-human IgG-HRP-conjugate or anti-human IgA-HRP-conjugate	35 ml , ready to use, IgG , coloured red, transparent bottle, red cap IgA , coloured purple, transparent bottle, purple cap
4	SUBSTR TMB	substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	35 ml , ready to use, black bottle, blue cap
5	INCUTRAY	incubation tray with cover	2

Materials required but not provided

- glassware
- pipettes and tips for 15 µl and 1500 µl
- plastic forceps
- rocking platform (vertical)
- filter paper
- collecting devices for infectious material
- distilled or deionized water

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 20 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

Working steps

1. Incubate the test strips with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **15 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **1.5 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **1.5 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **1.5 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.


Test evaluation

Test validity

The *Seraline*[®] Anti-Helicobacter-6 IgG / *Seraline*[®] Anti-Helicobacter-6 IgA test strips contain 5 consecutive control lines under the strip number:

- Function control (positive reaction with every sample).
- Conjugate controls IgG and IgA (positive reaction with the corresponding conjugate).
Note: The other conjugate control band may develop a weak, non-specific colour.
- The intensity of the negative control must be less than the Cut off control.
- Cut off control (band intensity is used for result evaluation of the diagnostic bands).

The test run is valid if the function control, the corresponding conjugate control and the Cut off control are clearly visible.

Result interpretation

Generally the evaluation of the test has to be done with dried strips. To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet. Alternatively the strips can be evaluated with the *Seraline*[®] scan software. Judgement of the bands is performed according to the following classification:

Result interpretation	IgG / IgA
negative	all bands < Cut off control
borderline	1 band (UreA, HcpC or gGT) ≥ Cut off control
positive	CagA or VacA or GroEL ≥ Cut off control or 2 bands ≥ Cut off control

Used antigens

antigen	molecular weight	Incidence*
CagA	120-145kDa	type I (high virulent)
VacA	87-97kDa	type I (high virulent)
GroEL	54-62kDa	type II
UreA	26kDa	type II
HcpC	28-32kDa	type II
gGT	59-61kDa	type II

*strain typing only with IgG antibody detection

Performance characteristics

Comparative Sensitivity and Specificity

For the determination of the comparative sensitivity and specificity serum samples from patients suspected of an infection with *Helicobacter pylori* and blood donor samples were tested with the *Seraline*[®] Anti-*Helicobacter*-6 IgG / *Seraline*[®] Anti-*Helicobacter*-6 IgA and another commercial assay.

Detection of IgG (n = 174)		<i>Seraline</i> [®] Anti- <i>Helicobacter</i> -6 IgG		
		positive	borderline	negative
Commercial-Line-Immunoassay	positive	96	0	3
	borderline	0	1	0
	negative	4	0	70

Agreement: 96.0%

Comparative Sensitivity: 97.0%

Comparative Specificity: 94.6%

Detection of IgA (n = 167)		<i>Seraline</i> [®] Anti- <i>Helicobacter</i> -6 IgA		
		positive	borderline	negative
Commercial-Line-Immunoassay	positive	50	1	4
	borderline	1	6	0
	negative	13	0	92

Agreement: 89.8%

Comparative Sensitivity: 93.5%

Comparative Specificity: 87.6%

Prevalence

A population of 400 healthy blood donors was investigated to determine the prevalence of Anti-*Helicobacter* IgG and IgA antibodies by using the *Seraline*[®] Anti-*Helicobacter*-6 Line Immunoassays. Out of this population 139 (131 positive and 8 borderline) serum samples were determined positive for IgG and 78 (67 positive and 11 borderline) serum samples were determined positive for IgA antibodies, corresponding to a prevalence of 34.8% for IgG and 19.5% for IgA antibodies.

n=400 healthy blood donors	<i>Seraline</i> [®] Anti- <i>Helicobacter</i> -6 IgG		<i>Seraline</i> [®] Anti- <i>Helicobacter</i> -6 IgA	
	quantity	%	quantity	%
positive	131	32.8	67	16.8
borderline	8	2	11	2.8
negative	261	65.3	322	80.5

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other Immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only. Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

History of changes

Version	Section	Modifications
2018-03-15	Common advices and precautions	Update

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use