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Verfallsdatum

Lagertemperatur

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Seraline[®] Anti-Treponema-4 IgG Seraline[®] Anti-Treponema-4 IgM

Line Immunoassay for detection of IgG or IgM antibodies to *Treponema pallidum*
in human serum or plasma

REF LIA-010-4 G ▼ 20 IVD *In-vitro*-diagnostic device CE
REF LIA-010-4 M ▼ 20 IVD *In-vitro*-diagnostic device CE

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Treponema pallidum ssp. pallidum is the causative agent of Syphilis, a sexually transmitted globally spread disease in humans. The incidence varies depending on the region and a consistent increase is observed since 2000 in western industrial nations with an actual level of 2-4 cases per 100,000 population per year (1).

The clinical course of Syphilis is characterized by different stages with varying symptoms. Infection occurs by sexual contacts. Skin lesions but also intact mucosa in the genital, oral or anal region are the most likely sites of entry for the pathogen and from there it spreads via blood and lymph vessels. During this primary stage a painless ulcer (ulcus durum) develops at the site of entry and a swelling of the regional lymph nodes (lymphadenopathy). Without treatment the infection will enter into the secondary stage after 9-10 weeks. Clinical signs are patchy papular exanthems or other ulcerative skin irritations, accompanied by several less characteristic symptoms like lymphadenopathy, fatigue, weight loss, subfebrile temperatures etc. The papular skin rashes contain a large number of pathogens. In the course of the secondary stage, which lasts between weeks to months, specific antibodies are produced. Nevertheless reinfections may occur. After decline of the symptoms and a latency period of up to several years the tertiary stage with a wide range of symptoms depending on the respective organ manifestations develops. Involvement of the central nervous system results in the clinical picture of Neurosyphilis or other neurological manifestations (1, 2).

Connatal *T. pallidum* infection by transplacental transmission during pregnancy causes spontaneous abortion, premature birth or stillbirth in 40% of the cases. About 50-70% of infected neonates are clinically unsuspectacular. Otherwise symptoms like bloody rhinitis, skin and mucosa irritations, fever, oedema, hepatosplenomegaly, lymphadenopathy, enteritis, hydrocephalus and meningitis may occur. In rare cases clinical symptoms do not develop before early childhood (3, 4).

Antibiotic treatment with Penicillin G is the method of choice for Syphilis therapy. Alternative antibiotics are tetracyclines and cephalosporines (3).

Treponema pallidum infections are diagnosed either by direct (pathogen detection or PCR) or indirect (serology) methods.

Direct pathogen detection from aspirates or smears of primary lesions, lymph nodes or exanthems by dark-field microscopy or direct immune fluorescence is hampered by the necessity to test the samples immediately after collection. Nucleic acid amplification techniques like PCR from aspirates, smears, blood or liquor are available but the diagnostic relevance is currently questionable (2, 5).

Laboratory diagnosis of Syphilis is preferably performed by detection of *T. pallidum* specific antibodies in serum, plasma or liquor. A two-step algorithm comprising a first screening and a subsequent confirmatory test is recommended (1, 2, 6). *Treponema pallidum* hemagglutination test (TPHA) or *Treponema pallidum* particle agglutination test and a polyvalent, IgG, IgA and IgM antibody detecting enzyme immunoassay (ELISA) are suitable screening methods. A positive or questionable screening result must be checked by a confirmatory test.

The Fluorescence *Treponema* Antibody Absorption test (FTA-ABS test) is an established confirmatory test. *Treponema* fixed to the surface of glass slides serve as target antigens for antibody detection from patient samples. In order to reach satisfying test specificity potentially cross reactive antibodies have to be eliminated by absorption with *Treponema phagedenis*. Immunoblots or line immunoassays are also suitable confirmatory tests (2, 4, 7). The major advantage of line immunoassays is the use of defined, highly purified recombinant antigens (Tpp47, Tpp17, Tpp15 and TmpA), that are confirmed as *T. pallidum* specific (8). Therefore cross reactivities appear unlikely and an absorption step is not needed (7). For differentiation between acute and recent infections confirmatory testing should always include the separate detection of IgG and IgM antibodies. Either immunoblot or line immunoassay or 19S-IgM-FTA-ABS test can be used to detect specific IgM antibodies for evaluation of the activity of the infection. Since IgM antibodies may persist for months up to years, a positive IgM result must always consider clinical findings for therapeutic decisions (2).

Quantitative detection methods like IgM-ELISA, 19S-IgM-Test or Lipoid antibody detection (e.g. VDRL-test, Cardiophilin-KBR) are recommended for serological follow-up after therapy (2, 4, 6).

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Intended use

The *Seraline*[®] Anti-Treponema-4 IgG / *Seraline*[®] Anti-Treponema-4 IgM is an in vitro diagnostic device (Line Immunoassay, LIA) for the qualitative detection of specific IgG or IgM antibodies against *Treponema pallidum* in human serum or plasma samples. It should be used as a confirmatory test after a positive ELISA result.

Principle of the test

The antigens are plotted on nitrocellulose membranes and blocked to avoid unspecific reactions. After this procedure the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

Step 1

The test strips are incubated in the diluted (1 : 101) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates. The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template. Alternatively the strips can be evaluated with the *Seraline*[®]scan software.

Used Treponema antigens

Nomenclature*	Properties / Specificity	Origin
Tpp47	Main membrane protein, Lipoprotein, Highly immunogenic, activates endothelial cell, M _r =47 kDa	Highly specific <i>T. pallidum Nichols</i>
TmpA	Membrane associated lipoprotein, signal peptide, M _r =37 kDa	Highly specific <i>T. pallidum Nichols</i>
Tpp17	Main membrane protein, Lipoprotein, M _r =17 kDa	Highly specific <i>T. pallidum Nichols</i>
Tpp15	Main membrane protein, Lipoprotein, immunogenic, M _r =15 kDa	Highly specific <i>T. pallidum Nichols</i>

* The nomenclature for some bands is not consistently used in the literature concerning molecular mass and specificity, respectively.

Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 101 (15 µl sample and 1500 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

Kit components

1	TESTSTR	test strips incl. evaluation template	20 strips colour coded: purple
2	WIB CONC 5x	wash and incubation buffer 5-fold	70 ml concentrate, transparent bottle, black cap
3	CONJ HRP IgG	anti-human IgG-HRP-conjugate	35 ml , ready to use, coloured red, transparent bottle, red cap coloured green, transparent bottle, green cap
	CONJ HRP IgM	anti-human IgM-HRP-conjugate	
4	SUBSTR TMB	substrate 3,3',5,5'- Tetramethylbenzidine and hydrogen peroxide	35 ml , ready to use, black bottle, blue cap
5	INCUTRAY	incubation tray with cover	2

Materials required but not provided

- glassware
- pipettes and tips for 15 µl and 1500 µl
- plastic forceps
- rocking platform (vertical)
- filter paper
- collecting devices for infectious material
- distilled or deionized water

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 20 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

Working steps

1. Incubate the test strips with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **15 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **1.5 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **1.5 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **1.5 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **1.5 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.


Test evaluation

Test validity

The *Seraline*[®] Anti-Treponema-4 IgG / *Seraline*[®] Anti-Treponema-4 IgM test strips contain 5 consecutive control lines beneath the strip number:

- Function control (positive reaction with every sample).
- Conjugate controls IgG and IgM (positive reaction with the corresponding conjugate).
Note: The other conjugate control band may develop a weak, non-specific colour.
- Negative control (band intensity must be less than the Cut off control).
- Cut off control (band intensity is used for result evaluation of the diagnostic bands).

The test run is valid if the function control, the corresponding conjugate control and the Cut off control are clearly visible.

Result interpretation

Generally the evaluation of the test has to be done with dried strips. To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet. Alternatively the strips can be evaluated with the *Seraline*[®] scan software. Judgement of the bands is performed according to the following classification:

Result interpretation	IgG	IgM
negative	max. 1 band ≥ Cut off control	no bands ≥ Cut off control
positive	min. 2 bands ≥ Cut off control	min. 1 band ≥ Cut off control

Inverse bands (bright bands on dark background) should be evaluated as negative. The respective serum should always be examined using other serological methods.

Performance characteristics

Sensitivity

Serum samples (n = 230) from patients with suspect of syphilis were investigated in a study using **Seraline® Anti-Treponema-4 IgG**. The serum samples have been pre-determined with different reference methods (TPHA, TPPA, Immunoblot, VDRL).

n = 230		Reference results	
		positive	negative
Seraline® Anti-Treponema-4 IgG	positive	217	3
	negative	6	4

Sensitivity of the **Seraline® Anti-Treponema-4 IgG** in comparison to the reference results: **97.3%** (217/223).

The following table presents the **Seraline®** sensitivity results compared to the different reference methods.

Panel	Sensitivity			
	Seraline® Anti-Treponema-4 IgG	Immunoblot	TPHA	TPPA
1 n=136	96.3% (131/136)	100% (136/136)	-	-
2 n=45	97.8% (44/45)	97.8% (44/45)	100% (45/45)	-
3 n=42	100% (42/42)	100% (42/42)	-	100% (42/42)

Serum samples (n = 211) from patients with suspect of syphilis were investigated in a study using **Seraline® Anti-Treponema-4 IgM**. The serum samples have been pre-determined with an Immunoblot (reference method).

n = 211		Reference results	
		positive	negative
Seraline® Anti-Treponema-4 IgM	positive	101	15
	negative	4	91

Sensitivity of the **Seraline® Anti-Treponema-4 IgM** in comparison to the reference results: **96.2%** (101/105).

The following table presents the **Seraline®** sensitivity results compared to the different reference methods.

Panel	Sensitivity		
	Seraline® Anti-Treponema-4 IgM	Immunoblot A	Immunoblot B
1 n=48	93.8% (45/48)	100% (48/48)	-
2 n=42	97.6% (41/42)	100% (42/42)	100% (42/42)
3 n=15	100% (15/15)	100% (15/15)	86.6% (13/15)

Specificity

A population of 407 healthy blood donors and 50 samples from pregnant women were investigated to determine the specificity of *Seraline*[®] Anti-Treponema-4 IgG and *Seraline*[®] Anti-Treponema-4 IgM Line Immunoassays.

Out of this population no serum samples were determined positive for IgG and 2 serum samples (blood donor collective) were determined positive for IgM antibodies with the *Seraline*[®] Anti-Treponema-4 IgG and *Seraline*[®] Anti-Treponema-4 IgM resp. resulting in a specificity of

***Seraline*[®] Anti-Treponema-4 IgG: 100%**

***Seraline*[®] Anti-Treponema-4 IgM: 99.6%**

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other Immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only. Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

History of changes

Version	Section	Modifications
2018-03-15	Common advices and precautions	Update

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use