Serablot[®] Anti-Francisella tularensis IgG

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Lagertemperatur

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Serablot[®] Anti-*Francisella tularensis* IgG

Immunoblot for detection of IgG antibodies to *Francisella tularensis*
in human serum or plasma

REF WB-003 G  27  *In-vitro*-Diagnostic device 

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Francisella tularensis is the causative agent of the zoonotic disease tularaemia (deer fly fever) in humans and animals. Additionally to transmission via infected animals especially rodents or their ectoparasites the disease can be transmitted via contaminated water or food and aerosols. After 3 to 5 days lasting incubation period the onset of disease is rather uncharacteristic with symptoms like headache, fatigue and fever. The various clinical manifestations as ulceroglandular, pulmonary or typhoidal form depend on the site of bacterial entry. The pulmonary course of infection is characterized by pneumonia and pleuritis and the abdominal course of infection by severe diarrhea, splenomegaly and intermittent fever.

Francisella tularensis can be identified directly by immunofluorescence microscopy or agglutination assays after isolation from affected tissue (e.g. lymph nodes) and subsequent cultivation. This procedure however is time consuming and detection of specific antibodies induced in the infected individuals e.g. by ELISA techniques gain significance. A positive ELISA result should be confirmed by a different second method which normally is an immunoblot. By that reason an immunoblot procedure was developed to confirm or to rule out antibodies to *Francisella tularensis* as reason for a positive ELISA result.

Intended use

The *Serablot*[®] Anti-*Francisella tularensis* IgG is an Immunoblot for detection of IgG antibodies directed to the lipopolysaccharide of the gram negative bacterium *Francisella tularensis* in human serum or plasma samples.

Principle of the test

Highly purified lipopolysaccharide from *Francisella tularensis* is separated by SDS-PAGE and subsequently transferred to nitrocellulose membrane. After blocking of remaining binding sites the membrane is cut into strips. The reaction for detection of specific antibodies is performed in 3 steps:

Step 1

The test strips are incubated in the diluted (1 : 201) serum or plasma samples for 45 minutes. Specific antibodies present in the samples react with the membrane bound antigens during that time. After incubation unbound proteins are removed by 3 washing cycles.

Step 2

The second incubation starts with the addition of enzyme labelled, isotype specific conjugate antibodies to the nitrocellulose strips. After 45 minutes of incubation unbound conjugate is removed by 3 washing cycles.

Step 3

The addition of the colourless precipitating substrate stains the membrane bound immune complexes by formation of blue precipitates.

The reaction is stopped after 10 minutes by aspirating the substrate followed by rinsing the strips with sufficient quantity of distilled or deionized water.

After drying the strips, interpretation of the resulting pattern is done by comparison with the specific template.

Preparation and storage of samples

Serum or plasma (EDTA, citrate or heparin plasma) can be used. Samples have to be diluted 1 : 201 (10 µl sample and 2000 µl WIB) with the ready to use wash and incubation solution (prepared from concentrated WIB (2)) in the incubation tray. Sample collection should be done in a sterile manner. Samples can be stored at 2...8°C for a maximum of 48 hours. For longer storage times samples have to be stored at -20°C. Frozen samples have to be warmed to room temperature and mixed well before starting the test run. Repeated freezing and thawing of samples should be avoided.

Kit components

1	TESTSTR	test strips incl. evaluation template	27 strips
2	WIB CONC 5x	wash and incubation buffer 5-fold	100 ml concentrate, transparent bottle, black cap
3	CONJ HRP IgG	anti-human IgG-HRP-conjugate	55 ml , ready to use, IgG , coloured red, transparent bottle, red cap
4	SUBSTR TMB	substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	55 ml , ready to use, black bottle, blue cap
5	INCUTRAY	incubation tray with cover	2

Materials required but not provided

- glassware
- pipettes and tips for 10 µl and 2000 µl
- plastic forceps
- rocking platform (vertical)
- filter paper
- collecting devices for infectious material
- distilled or deionized water

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 27 determinations. The complete test kit with closed reagent bottles is stable until the expiry date printed on the labels when stored at 2...8°C. After opening and proper storage as requested the reagents have to be used within 3 months.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. Ready to use wash and incubation solution has to be prepared from the 5-fold concentrated WIB (2) by dilution of 1 : 5 with distilled water. Mix the concentrate before dilution.

For example: 20 ml WIB (5-fold) + 80 ml distilled water

Assay procedure

- The test has to be performed at room temperature (18...25°C).
- The test strips have to be transferred into the incubation tray by use of plastic forceps in a way that the printed strip number on the top is visible.
- All incubation steps should be done on a rocking platform with a recommended frequency of rocking of about 20 - 30 rotations per minute.
- The following procedure has to adhere strictly to the time table.

Working steps

1. Incubate the test strips with **2 ml** wash and incubation solution (ready to use, prepared from (2)) on a rocking platform for **5 min**.
2. Add **10 µl** sample.
3. Cover the incubation tray and incubate for **45 min** on a rocking platform.
4. Aspirate and wash **3 x 5 min** with **2 ml** wash and incubation solution (ready to use, prepared from (2)).
5. Add **2 ml** **CONJ HRP**
6. Cover the incubation tray and incubate **45 min** on a rocking platform.
7. Aspirate and wash **3 x 5 min** with **2 ml** wash and incubation solution (ready to use, prepared from (2)).
8. Incubate the test strips with **2 ml** **SUBSTR TMB** **10 min** on a rocking platform.
9. Aspirate substrate and rinse the test strips **3 x** with **2 ml** distilled water to stop the reaction.
10. The developed test strips have to be dried between filter paper before evaluation.

Please Note:

With so-called „problematic“ samples (e.g. from patients with hypergammaglobulinaemia, circulating immune complexes, antibodies to milk proteins) the background develops rapidly and the test strips discolour blue. In such cases the substrate reaction has to be terminated earlier.

Test evaluation

Test validity

The internal control band on the top of the teststrip must be present with a strong colour. Absence or very low colour intensity indicate an incorrect test procedure.

Result interpretation

Generally the evaluation of the test has to be done with dried strips.

To support identification of stained bands the evaluation template should be used. The identified bands have to be documented on the data sheet.

Positive:

Visible image of a regular band pattern with systematically growing (see evaluation template provided with the kit).

Negative:

No or only a few bands without the characteristic regular band pattern.

Performance characteristics

Sensitivity

Results from the comparative investigation of serum samples (n = 186) that had been characterized by the German reference laboratory using *Serablot*[®] Anti- Francisella tularensis IgG and *Serazym*[®] Anti-Francisella tularensis IgG + IgM + IgA.

n = 186 samples		<i>Serablot</i> [®] Anti-Francisella tularensis IgG	
		positive	negative
<i>Serazym</i> [®] Anti-F. tularensis IgG+IgM+IgA	positive	114	6
	negative	3	63

Sensitivity: 95.2%

Specificity

A population of 100 healthy blood donors was investigated to determine the specificity of *Serablot*[®] Anti-Francisella tularensis IgG. Out of this population 5 serum samples were determined positive.

Specificity: 95.0%

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. As in other immunoassays, impurities and cross contamination of reagents and samples by fungi or bacteria can produce false positive as well as false negative results. Incorrect dilution of samples, incorrect washing of the test strips or insufficient rinsing after substrate reaction and also incorrect timing can produce erroneous results.

Interference

Hemolytic, lipemic and icteric samples did not interfere with the test up to a concentration of 500 mg/dl (hemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere up to a concentration of 500 IU/ml.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The test should be performed by trained technical laboratory staff only. The test instructions have to be followed strictly. Do not use reagents from damaged packages or bottles. Handle the test strips with plastic forceps only. Prior to testing the identification numbers of the samples, the numbers of the test strips and the lot number of the kit should be documented on the data sheet. The expiry dates stated on the respective labels are to be observed.

Do not use or mix reagents from different lots except for wash and incubation buffer and for substrate.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

Since the kit contains potentially hazardous materials, the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

History of changes

Version	Section	Modifications
2018-03-15	Common advices and precautions	Update

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use