



Lyme Borreliose

Seraline[®] Anti-Borrelia IgG **Seraline[®] Anti-Borrelia IgM**

Line - Immunoassay for detection of IgG- or IgM-antibodies against *Borrelia burgdorferi* in human serum or plasma

- ▲ sensitive and specific
- ▲ ready-to-use colour-coded reagents
- ▲ short incubation times
- ▲ automatic or visual evaluation

Indication

The investigation of serum or plasma for the presence of anti-Borrelia antibodies by Line immunoassay is recommended in case of a positive anti-Borrelia ELISA result.

- to control the specificity of the ELISA results
- for discrimination of stages of the disease according to the specific band pattern in the case of low reactivity in anti-Borrelia ELISA
- for detection of a seroconversion
- to control the specificity of the ELISA results in the case of repeated reactivity in IgM-ELISA

Principle

The *Seraline[®] Anti-Borrelia* uses highly purified recombinant Borrelia proteins, which are plotted onto strips of nitrocellulose-membranes in defined distances.

Specific antibodies present in the sample react with the respective antigens during the first incubation step. These immune complexes are detected in the subsequent conjugate (anti-human IgG- or anti-human IgM-HRP) incubation step.

Membrane bound immune complexes are finally detected by formation of blue precipitates with the substrate solution.

One function control, 2 conjugate controls and a Cut off control support the evaluation.

Procedure

- ◆ incubate nitrocellulose strip with 1.5 ml wash- and incubation buffer (WIB)
- ◆ add 15 µl sample to each strip
- ◆ incubate for 45 min on a rocking platform
- ◆ wash 3 times with 1.5 ml WIB (5 min)
- ◆ add 1.5 ml conjugate (IgG or IgM)
- ◆ incubate for 45 min on a rocking platform
- ◆ wash 3 times with 1.5 ml WIB (5 min)
- ◆ add 1.5 ml substrate
- ◆ incubate for 10 min on a rocking platform
- ◆ wash 3 times with 1.5 ml pure water
- ◆ dry strip before evaluation

Evaluation and interpretation

The single bands are assigned to the respective Borrelia antigens by usage of the evaluation template or with the software *Seraline[®] scan*.

Evaluation criteria


IgG

positive	2 bands except for p41 \geq Cut off control
borderline	VlsE, 1 band + p41 \geq Cut off control
negative	max. 1 band except for VlsE \geq Cut off control

IgM

positive	OspC, p18, 2 other bands \geq Cut off control
borderline	1 band except for OspC, p41, p18 \geq Cut off control
negative	no bands except for p41 \geq Cut off control

Band pattern and characterization



antigens	strain
p100	B. afzelii
VlsE	B. afzelii
p58	B. garinii
p41	B. burgdorferi. s.s.
p39	B. afzelii
OspA	B. afzelii
OspC	Mix*
p18	Mix*

Mix*: B. burgdorferi sensu stricto, B. garinii I + II (only Osp C), B. afzelii, B. spielmanii

Sensitivity

IgG and IgM detection

n = 78 sera		comparative – Line Immunoassay	
		positive	negative
Seraline®	positive	51	10
	negative	1	16

A total of 67 samples delivered concordant results in both assays which means an overall agreement of 85.9% and a comparative sensitivity of 98.1% for the Seraline® Anti-Borrelia IgG/IgM in relation to the commercial assay.

Samples with discordant results (n=11) were reinvestigated by immunoblot (Serablott® Anti-Borrelia IgG/ IgM). The Seraline® results were confirmed by immunoblot for 8 of 10 positive samples.


Specificity

The investigation of 400 healthy blood donors revealed a specificity of

88.2 % Seraline® Anti-Borrelia IgG

94.5 % Seraline® Anti-Borrelia IgM

Serology of different stages of Lyme Borreliosis in the Seraline® Anti-Borrelia IgG/IgM



Katalog- Nr. LIA-006-8 G, LIA-006-8 M

Seramun Diagnostica GmbH
Sprenhagener Str. 1
15754 Heidesee, GERMANY
(Stand 01.03.12)



• Certified Management System
• EN ISO 9001
• EN ISO 13485

Tel.: +49-33767-791 10
Fax: +49-33767-791 99
E-Mail: info@seramun.com
www.seramun.com