Serazym[®] Bovine Serum Albumin (BSA) sensitive

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Referenzen:

1. Khomehchian, S.: Development of a sandwich enzyme-linked immunosorbant assay (ELISA) for determining of bovine serum albumin (BSA) in trivalent measles-mump-rubella (MMR) vaccines, Human Vaccines 4:5, 375-378; September/October 2008
2. Granström, M.: A sandwich ELISA for bovine serum in viral vaccines, Journal of Biological Standardization (1987) 15, 193-197
3. Zhang, K.: The establishment of a highly sensitive ELISA for detecting bovine serum albumin (BSA) based on a specific pair of monoclonal antibodies (mAb) and its application in vaccine quality control, Human Vaccines 6:8, 652-658; August 2010

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Serazym[®] Bovine Serum Albumin (BSA) sensitive

Enzyme-linked immunosorbent assay for detection of bovine serum albumin in fluids

REF E-108 ▽ 96

Seramun Diagnostica GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Intended Use

The *Serazym[®]* Bovine Serum Albumin (BSA) sensitive is an enzyme-linked immunosorbent assay for the detection of bovine serum albumin (BSA) in fluids.

Principle Of The Test

The *Serazym[®]* Bovine Serum Albumin (BSA) sensitive is a fast immune enzymometric one-step assay based on affinity purified polyclonal antibodies (rabbit) against bovine serum albumin (BSA).
Specimens

and horseradish peroxidase (HRP) labelled anti-BSA-antibodies are dispensed into the microtitration plate wells, which are coated with anti-BSA-antibodies and incubated at 20...25°C / 68...77°F for 60 minutes. After the incubation unbound components are removed by washing the wells five times with wash buffer. Then substrate solution (tetramethylbenzidine and hydrogen peroxide) is added to each well and incubated for 15 minutes at 20...25°C / 68...77°F protected from light. The presence of specifically bound enzyme labelled antibodies in the wells is indicated by the development of a blue colour. Reaction stop by addition of sulphuric acid to the wells results in a colour change to yellow. Absorbances are read at 450 nm wavelength. For optimum results a reference filter (620 - 690 nm wavelength) should be used. A standard curve is created by plotting the absorbances of the different BSA standards versus their concentrations. The absorbances of the unknown samples can be transformed into their corresponding BSA concentrations by reading from the standard curve.

References:

4. Khomehchian, S. : Development of a sandwich enzyme-linked immunosorbant assay (ELISA) for determining of bovine serum albumin (BSA) in trivalent measles-mump-rubella (MMR) vaccines, Human Vaccines 4:5, 375-378; September/October 2008
5. Granström, M.: A sandwich ELISA for bovine serum in viral vaccines, Journal of Biological Standardization (1987) 15, 193-197
6. Zhang, K.: The establishment of a highly sensitive ELISA for detecting bovine serum albumin (BSA) based on a specific pair of monoclonal antibodies (mAB) and its application in vaccine quality control, Human Vaccines 6:8, 652-658; August 2010

Test Components

			For 96 Wells
1	WELLS	Microtitration plate coated with polyclonal anti-bovine serum albumin- antibodies (rabbit)	12 single breakable 8-well strips vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	50 ml concentrate for 500 ml solution white cap
3	DIL	Sample diluent	70 ml · ready to use coloured red black cap
4	STD 1 - 7	Standard 1-7 S1 = 50.0 ng/ml; S2 = 40.0 ng/ml; S3 = 30.0 ng/ml; S4 = 20.0 ng/ml; S5 = 10.0 ng/ml; S6 = 5.0 ng/ml S7 = 2.5 ng/ml	1.0 ml per standard ready to use coloured red transparent cap
5	CONTROL	Control 25 ng/ml	1.0 ml · ready to use coloured red green cap
6	CONJ HRP	HRP-conjugate HRP-labelled, polyclonal anti- bovine serum albumin - antibodies (rabbit)	15 ml · ready to use coloured red red cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap

Preparation And Storage Of Samples

Collection and storage

The *Serazym*[®] Bovine Serum Albumin sensitive is intended for the determination of the BSA content in diluted fluids. Samples should not be stored longer than 48 hours at 2...8°C / 35.6...46.4°F prior to use. Otherwise a storage temperature of -20°C / -4°F is recommended. Frozen samples should be rapidly warmed to room temperature and mixed thoroughly before testing. Repeated freezing and thawing of samples should be avoided.

Preparation

Dilute samples with sample diluent (3) (usually 1 : 6 to 1 : 21). The dilution varies in dependence from the material and expected BSA concentration. Pipette 1000 µl of sample diluent into a clean tube and add 200 µl (1 : 6) or 50 µl (1 : 21) of the sample. After it mix the solution gently.

Materials Required But Not Provided

Adjustable one-channel micropipette 0.100 - 1.000 ml and 0.010 - 0.100 ml · Adjustable 8-channel micropipette 0.050 - 0.200 ml · Pipette tips · Reagent containers for dispensing with 8-channel pipette · Graduated measuring flasks 10 ml and 100 ml · Eppendorf tubes 2.0 ml · Microplate washer (automatic or hand wash head) · Microplate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · Distilled or deionized water · Stop-watch

Preparation And Storage Of Reagents

Kit size and expiry

One kit (1x 96 wells) enables BSA quantification in a maximum of 40 samples when samples, standards and control are run in duplicate. The complete kit with unopened reagent bottles and microtitration strips is stable until the expiry date printed on the kit box in case of storage at 2...8°C / 35.6...46.4°F. Once opened all kit components are stable for up to 2 months under appropriate storage conditions 2...8°C / 35.6...46.4°F. When stored at 2...8°C / 35.6...46.4°F the ready to use wash solution can be used for at least one month.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Sample diluent

The sample diluent is ready-to-use. During storage at 2...8°C precipitates may occur, that will dissolve when the solution is warmed to room temperature. Make sure the sample diluent has reached room temperature before starting the assay.

Wash solution

Prepare a sufficient amount of wash solution by diluting the 10fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

Dilute samples with sample diluent (3) (usually 1 : 6 to 1 : 21). Avoid any time shift during dispensing of reagents and samples.

Attention: The *Serazym*[®] Bovine Serum Albumin (BSA) sensitive is a very sensitive assay detecting less than 1.0 ng BSA per ml sample material. It is recommended to use new disposable reagent containers for pipetting the reagents. Make sure that the glassware used for buffer preparation and the equipment used for running the assay are absolutely free of BSA. Any possibility of contamination of working equipment with BSA should be strictly avoided! Therefore, special dispensers, pipettes and washers have to be used only for this test and not for other ELISAs.

In case of using a washer make sure that the wells are completely filled (at least 300 µl / well) and drained in every single washing cycle. In case of manual washing avoid foam and air bubbles in the wells. It is recommended to tap the plate onto absorbent paper after each washing cycle.

Make sure that the soak time of the wash solution in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature before use. Mix gently without causing foam.
2. Dispense 100 µl **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 100 µl **diluted samples**
100 µl **STD 1 – 7** standard S1-S7 (4)
100 µl **CONTROL** control (5)
into the intended wells, mix gently.
It is recommended to run samples, standards and control in duplicate.
4. Cover plate and incubate for 60 ± 1 min at $20...25^{\circ}\text{C}$ / $68...77^{\circ}\text{F}$.
5. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
6. Dispense 100 µl **SUBSTR TMB** substrate (7) per well.
7. Incubate for 15 min at $20...25^{\circ}\text{C}$ / $68...77^{\circ}\text{F}$, protected from light.
8. Dispense 100 µl **STOP** stop solution (8) per well, mix gently.
9. Read absorbance at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

Result Interpretation

A standard curve is created by plotting the measured average values of absorbance of the standards 1-7 (y-axis) against the declared BSA concentration (x-axis). Alternatively, the diluent can be included as standard 8 to construct the standard curve. In this case, a concentration of 0.1 ng/ml is to be selected for standard 8 to enable a mathematical model of the standard curve. Plot the mean absorbances of the standards 1 – 8 (y-axis) against the corresponding BSA concentrations of these standards (x-axis). Determine the BSA concentration of the unknown samples by referring their mean absorbances to the reference curve and multiply the values with the dilution factor. Samples with BSA concentrations exceeding those of standard 1 (50 ng/ml) should be retested in higher predilution. In case of sample predilution values have to be multiplied with the dilution factor.

Reference Values

Test validity

A test run is valid if:

- the mean absorbance of standard 1 (50.0 ng/ml) is ≥ 1.50
- the mean absorbance of standard 7 (2.5 ng/ml) is ≤ 0.50
- the mean of the control is determined between 20 ng/ml and 30 ng/ml

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is correct (incubation times and temperatures, sample and wash buffer dilution, washing steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure

Incorrect results may be caused by contamination of reagents or working equipment with BSA. Cross contaminations of the kit reagents and samples, bacterial or fungal contaminations of reagents and/or samples, incorrect washing and incorrect incubation times may also cause erroneous results. Due to the high-dosis-hook-effect samples with BSA concentrations higher than 1 µg/ml may be determined too low. In these cases it is recommended to test several higher dilutions of the respective samples.

Performance Characteristics

Precision

Intra-assay coefficients of variation (CV) in the *Serazym*[®] Bovine Serum Albumin (BSA) sensitive from 8-fold determinations of samples.

sample	mean absorbance	standard deviation	CV [%]
1	2.437	0.033	1.43
2	2.060	0.035	1.79
3	1.670	0.050	3.18
4	1.243	0.051	4.41
5	0.737	0.029	4.24
6	0.473	0.013	2.86
7	0.322	0.002	0.51
8	1.494	0.037	2.61

sample	mean concentration [ng/ml]	standard deviation	CV [%]
1	50.1	0.924	1.97
2	39.8	0.905	2.43
3	30.0	1.197	4.26
4	20.2	1.122	5.93
5	9.9	0.561	6.08
6	5.0	0.221	4.69
7	2.5	0.025	1.06
8	25.9	0.843	3.48

Inter-assay coefficients of variation (CV) in the *Serazym*[®] Bovine Serum Albumin (BSA) sensitive in 8 different test runs from twofold determinations of samples.

sample	mean absorbance	standard deviation	CV [%]
1	2.087	0.099	5.07
2	1.205	0.067	5.93
3	0.715	0.037	5.50
4	0.466	0.027	6.74

sample	mean concentration [ng/ml]	standard deviation	CV [%]
1	50.9	1.701	3.57
2	24.2	0.750	3.31
3	12.0	0.297	2.64
4	6.3	0.299	5.06


Cross reactivity

Samples with different concentrations (200 ng/ml – 3 ng/ml) of Human serum albumin and Ovalbumin were tested in the *Serazym*[®] Bovine Serum Albumin (BSA) sensitive. No cross reactivity was observed.


Lower detection limit

The lower detection limit was determined to be 0.5 ng/ml.

Typical calibration curve in the *Serazym*[®] Bovine Serum Albumin (BSA) sensitive


Dilution linearity and analytical sensitivity

An analytical sensitivity of 5 ng BSA/ml for the reliable detection of a difference in concentration of BSA between samples was determined.

High-Dose-Hook-Effect

To determine the high-dose hook effect that occurs in one-step assays (simultaneous incubation of sample and conjugate), concentrations of 3.9 - 1,000,000 ng / ml were determined. The level of the signal intensity as a function of the BSA concentration shows a high dose hook effect, which begins at a concentration of $\geq 8 \mu\text{g} / \text{ml}$ BSA.

Due to the high-dose hook effect, samples with BSA concentrations $\geq 8 \mu\text{g} / \text{ml}$ can be determined too low. In these cases it is recommended to test several different dilutions of the same sample.

Common Advices and Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiry date on the respective labels is appropriate to note. The same applies for the specified shelf life of the reconstituted reagents.

Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution. Do not use reagents from other manufacturers.

Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They may neither be swallowed nor come into contact with skin or mucous membranes. All reagents should be kept at $2...8^{\circ}\text{C} / 36...46^{\circ}\text{F}$ before use.

Use a new pipette tip for each sample to avoid contaminations and absolutely BSA-free equipment to get correct results. Handle all samples as potentially infectious. Since the kit contains potentially hazardous materials, the following precautions must always be observed:

- Do not smoke, eat or drink while handling kit material!**
- Always use protective gloves!**
- Never pipette material by mouth!**
- Note safety precautions of the single test components!**

History of Changes

Version	Section	Modifications
2020-11-12	Test Components	Correction

Incubation Scheme *Serazym*[®] Bovine Serum Albumin (BSA) sensitive (E-108)

1.  100 µl **CONJ HRP** (6)

+
100 µl **diluted sample**

100 µl **STD S1 – S7** (4) alternatively S1 – S8

100 µl **CONTROL** (5)

 60 min incubation (20...25°C / 68...77°F)

5 x wash with wash solution

2.  100 µl **SUBSTR TMB** (7)

15 min incubation (20...25°C / 68...77°F) protected from light

3.  100 µl **STOP** (8)

Read absorbances at 450 / ≥ 620 nm

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	 Contains sufficient for <n> tests		