# Serazym<sup>®</sup> Bovine Serum Albumin (BSA) sensitive

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## Referenzen:

1. Khomehchian, S. : Development of a sandwich enzyme-linked immunosorbant assay (ELISA) for determining of bovine serum albumin (BSA) in trivalent measles-mump-rubella (MMR) vaccines, Human Vaccines 4:5, 375-378; September/October 2008
2. Granström, M.: A sandwich ELISA for bovine serum in viral vaccines, Journal of Biological Standardization (1987) 15, 193-197
3. Zhang, K.: The establishment of a highly sensitive ELISA for detecting bovine serum albumin (BSA) based on a specific pair of monoclonal antibodies (mAb) and its application in vaccine quality control, Human Vaccines 6:8, 652-658; August 2010

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Lagertemperatur

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## Serazym<sup>®</sup> Bovine Serum Albumin (BSA) sensitive

Enzyme-linked immunosorbent assay for detection of bovine serum albumin in biological fluids

REF E-108 ▼ 96 CE

Seramun Diagnostica GmbH · Sprenthagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · [info@seramun.com](mailto:info@seramun.com)

### Intended Use

The *Serazym<sup>®</sup> Bovine Serum Albumin sensitive* is an *in-vitro*-test for the detection of bovine serum albumin (BSA) in biological fluids.

### Principle Of The Test

The *Serazym<sup>®</sup> Bovine Serum Albumin sensitive* is a fast immune enzymometric one-step assay based on affinity purified polyclonal antibodies (rabbit) against bovine serum albumin (BSA). Specimens and horseradish peroxidase (HRP) labelled anti-BSA-antibodies are dispensed into the microtitration plate wells, which are coated with anti-BSA-antibodies and incubated at 20...25°C / 68...77°F for 60 minutes. After the incubation unbound components are removed by washing the wells five times with wash buffer. Then substrate solution (tetramethylbenzidine and hydrogen peroxide) is added to each well and incubated for 15 minutes at 20...25°C / 68...77°F protected from light. The presence of specifically bound enzyme labelled antibodies in the wells is indicated by the development of a blue colour. Reaction stop by addition of sulphuric acid to the wells results in a colour change to yellow. Absorbances are read at 450 nm wavelength. For optimum results a reference filter (620 - 690 nm wavelength) should be used. A standard curve is created by plotting the absorbances of the different BSA standards versus their concentrations. The absorbances of the unknown samples can be transformed into their corresponding BSA concentrations by reading from the standard curve.

## References:

4. Khamsehchian, S. : Development of a sandwich enzyme-linked immunosorbant assay (ELISA) for determining of bovine serum albumin (BSA) in trivalent measles-mump-rubella (MMR) vaccines, Human Vaccines 4:5, 375-378; September/October 2008
5. Granström, M.: A sandwich ELISA for bovine serum in viral vaccines, Journal of Biological Standardization (1987) 15, 193-197
6. Zhang, K.: The establishment of a highly sensitive ELISA for detecting bovine serum albumin (BSA) based on a specific pair of monoclonal antibodies (mAb) and its application in vaccine quality control, Human Vaccines 6:8, 652-658; August 2010

## Test Components

			For 96 Wells
1	<b>WELLS</b>	<b>Microtitration plate</b> coated with polyclonal anti-bovine serum albumin- antibodies (rabbit)	12 single breakable 8-well strips vacuum-sealed with desiccant
2	<b>WASHBUF CONC 10x</b>	<b>Wash buffer</b> 10-fold	50 ml concentrate for 500 ml solution white cap
3	<b>DIL</b>	<b>Sample diluent</b>	70 ml · ready to use coloured red black cap
4	<b>STD 1 - 7</b>	<b>Standard 1-7</b> S1 = 50.0 ng/ml; S2 = 40.0 ng/ml; S3 = 30.0 ng/ml; S4 = 20.0 ng/ml; S5 = 10.0 ng/ml; S6 = 5.0 ng/ml S7 = 2.5 ng/ml	1.0 ml per standard ready to use coloured red transparent cap
5	<b>CONTROL</b>	<b>Control</b> 25 ng/ml	1.0 ml · ready to use coloured red green cap
6	<b>CONJ HRP</b>	<b>HRP-conjugate</b> HRP-labelled, polyclonal anti- bovine serum albumin - antibodies (rabbit)	15 ml · ready to use coloured green red cap
7	<b>SUBSTR TMB</b>	<b>Substrate</b> 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap
8	<b>STOP</b>	<b>Stop solution</b> 0.25 M sulphuric acid	15 ml · ready to use yellow cap

## Preparation And Storage Of Samples

### Collection and storage

The *Serazym*<sup>®</sup> Bovine Serum Albumin sensitive is intended for the determination of the BSA content in diluted biological fluids. Samples should not be stored longer than 48 hours at 2...8°C / 35.6 ... 46.4°F prior to use. Otherwise a storage temperature of -20°C / -4°F is recommended. Frozen samples should be rapidly warmed to room temperature and mixed thoroughly before testing. Repeated freezing and thawing of samples should be avoided.

### Preparation

Dilute samples with sample diluent (3) (usually 1 : 6 to 1 : 21). The dilution varies in dependence from the material and expected BSA concentration.

## Materials Required But Not Provided

Adjustable one-channel micropipette 0.100 - 1.000 ml and 0.010 - 0.100 ml · Adjustable 8-channel micropipette 0.050 - 0.200 ml · Pipette tips · Reagent containers for dispensing with 8-channel pipette · Graduated measuring flasks 10 ml and 100 ml · Eppendorf tubes 2.0 ml · Microplate washer (automatic or hand wash head) · Microplate reader with optical filters of 450 nm for measurement and  $\geq 620$  nm for reference · Distilled or deionized water · Stop-watch

## Preparation And Storage Of Reagents

### Kit size and expiry

One kit (1x 96 wells) enables BSA quantification in a maximum of 40 samples when samples, standards and control are run in duplicate. The complete kit with unopened reagent bottles and microtitration strips is stable until the expiry date printed on the kit box in case of storage at 2...8°C / 35.6...46.4°F. Once opened all kit components are stable for up to 2 months under appropriate storage conditions 2...8°C / 35.6...46.4°F. When stored at 2...8°C / 35.6...46.4°F the ready to use wash solution can be used for at least 30 days.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

### Sample diluent

The sample diluent is ready-to-use. During storage at 2...8°C precipitates may occur, that will dissolve when the solution is warmed to room temperature. Make sure the sample diluent has reached room temperature before starting the assay.

### Wash solution

Prepare a sufficient amount of wash solution by diluting the 10fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

## Assay Procedure

Dilute samples with sample diluent (3) (usually 1 : 6 to 1 : 21). Avoid any time shift during dispensing of reagents and samples.

**Attention: The *Serazym*<sup>®</sup> Bovine Serum Albumin sensitive is a very sensitive assay detecting less than 1.0 ng BSA per ml sample material. It is recommended to use new disposable reagent containers for pipetting the reagents. Make sure that the glassware used for buffer preparation and the equipment used for running the assay are absolutely free of BSA. Any possibility of contamination of working equipment with BSA should be strictly avoided! Therefore special dispensers, pipettes and washers have to be used only for this test and not for other ELISAs.**

In case of using a washer make sure that the wells are completely filled (at least 300  $\mu$ l / well) and drained in every single washing cycle. In case of manual washing avoid foam and air bubbles in the wells. It is recommended to tap the plate onto absorbent paper after each washing cycle.

Make sure that the soak time of the wash solution in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

## Working steps

1. Warm all reagents to room temperature before use. Mix gently without causing foam.
2. Dispense 100 µl **CONJ HRP** HRP-conjugate (6) per well and
3. Pipette: 100 µl **diluted samples**  
100 µl **STD 1 – 7** standard S1-S7 (4)  
100 µl **CONTROL** control (5)  
into the intended wells, mix gently.  
It is recommended to run samples, standards and control in duplicate.
4. Cover plate and incubate for  $60 \pm 1$  min at  $20...25^{\circ}\text{C}$  /  $68...77^{\circ}\text{F}$ .
5. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
6. Dispense 100 µl **SUBSTR TMB** substrate (7) per well.
7. Incubate for 15 min at  $20...25^{\circ}\text{C}$  /  $68...77^{\circ}\text{F}$ , protected from light.
8. Dispense 100 µl **STOP** stop solution (8) per well, mix gently.
9. Read absorbance at 450 nm /  $\geq 620$  nm with a microplate reader within 30 min after reaction stop.

## Result Interpretation

Calculate the mean absorbance of standards, control and samples. Create a reference curve from the mean absorbances of the standards S1 – S7 (y-axis) and the corresponding BSA concentrations of these standards (x-axis). Alternatively the diluent can be included as standard S8 to construct the standard curve (S1–S8). In this case, a concentration of 0.1 ng/ml is to be selected for S8 to enable a mathematical model of the standard curve. Plot the mean absorbances of the calibrator S1 – S8 (y-axis) against the corresponding BSA concentrations of these standards (x-axis). Determine the BSA concentration of the unknown samples by referring their mean absorbances to the reference curve. Samples with BSA concentrations exceeding those of standard S 1 (50 ng/ml) should be retested in higher predilution. In case of sample predilution values have to be multiplied with the dilution factor.

## Reference Values

### Test validity

A test run is valid if:

- the mean absorbance of standard 1 (50.0 ng/ml) is  $\geq 1.50$
- the mean absorbance of standard 7 (2.5 ng/ml) is  $\leq 0.50$
- the mean of the control is determined between 20 ng/ml and 30 ng/ml

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is correct (incubation times and temperatures, sample and wash buffer dilution, washing steps etc.). In case of repeated failure of the quality criteria contact your supplier.

### Limitations of the procedure

Incorrect results may be caused by contamination of reagents or working equipment with BSA. Cross contaminations of the kit reagents and samples, bacterial or fungal contaminations of reagents and/or samples, incorrect washing and incorrect incubation times may also cause erroneous results. Due to the high-dosis-hook-effect samples with BSA concentrations higher than 1 µg/ml may be determined too low. In these cases it is recommended to test several higher dilutions of the respective samples.

## Performance Characteristics

### Precision

Intra-assay coefficients of variation (CV) in the *Serazym*<sup>®</sup> Bovine Serum Albumin sensitive from 8-fold determinations of samples.

sample	mean absorbance	standard deviation	CV [%]
1	2.437	0.033	1.43
2	2.060	0.035	1.79
3	1.670	0.050	3.18
4	1.243	0.051	4.41
5	0.737	0.029	4.24
6	0.473	0.013	2.86
7	0.322	0.002	0.51
8	1.494	0.037	2.61

sample	mean concentration [ng/ml]	standard deviation	CV [%]
1	50.1	0.924	1.97
2	39.8	0.905	2.43
3	30.0	1.197	4.26
4	20.2	1.122	5.93
5	9.9	0.561	6.08
6	5.0	0.221	4.69
7	2.5	0.025	1.06
8	25.9	0.843	3.48

Inter-assay coefficients of variation (CV) in the *Serazym*<sup>®</sup> Bovine Serum Albumin sensitive in 8 different test runs from twofold determinations of samples.

sample	mean absorbance	standard deviation	CV [%]
1	2.087	0.099	5.07
2	1.205	0.067	5.93
3	0.715	0.037	5.50
4	0.466	0.027	6.74

sample	mean concentration [ng/ml]	standard deviation	CV [%]
1	50.9	1.701	3.57
2	24.2	0.750	3.31
3	12.0	0.297	2.64
4	6.3	0.299	5.06


### Cross reactivity

Samples with different concentrations (200 ng/ml – 3 ng/ml) of Human serum albumin and Ovalbumin were tested in the *Serazym*<sup>®</sup> Bovine Serum Albumin sensitive. No cross reactivity was observed.


### Lower detection limit

The lower detection limit was determined to be 0.5 ng/ml.


Typical calibration curve in the *Serazym*<sup>®</sup> Bovine Serum Albumin sensitive

### Dilution linearity and analytical sensitivity

### High-Dose-Hook-Effect

An analytical sensitivity of 5 ng BSA/ml for the reliable detection of a difference in concentration of BSA between samples was determined.

Development of signal intensity in dependence of the BSA concentration reveals a high-dose hook effect which starts by exceeding a concentration of about 1 µg/ml (see “limitations of the procedure”).

## Common Advices and Precautions

**This kit is for *in-vitro* use only.** Follow the working instructions carefully. The kit should be performed by trained technical staff only. The expiry date on the respective labels is appropriate to note. The same applies for the specified shelf life of the reconstituted reagents. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.** Do not use reagents from other manufacturers. Some of the reagents (2, 3, 4, 5, 6, 7) contain small amounts of Thimerosal (< 0.1% w/v) and Kathon (1.0% v/v) as preservative. They may neither be swallowed nor come into contact with skin or mucous membranes. All reagents should be kept at 2...8°C / 36...46°F before use. Precisely adhere to the prescribed incubation times and temperatures. Avoid time shift during dispensing of reagents. Prepare a data sheet that defines the positions of samples, standards and control before starting the assay.

Handle all samples as potentially infectious. Use a new pipette tip for each sample to avoid contaminations and absolutely BSA-free equipment to get correct results. Since the kit contains potentially hazardous materials, the following precautions must always be observed:

**Do not smoke, eat or drink while handling kit material!**

**Always use protective gloves!**

**Never pipette material by mouth!**

**Note safety precautions of the single test components!**

# Incubation Scheme *Serazym*<sup>®</sup> Bovine Serum Albumin sensitive (E-108)

1.  100 µl **CONJ HRP** (6)

+  
100 µl **diluted sample**

100 µl **STD S1 – S7** (4) alternatively S1 – S8

100 µl **CONTROL** (5)

 60 min incubation (20...25°C / 68...77°F)  
5 x wash with wash solution
  
2.  100 µl **SUBSTR TMB** (7)

15 min incubation (20...25°C / 68...77°F) protected from light
  
3.  100 µl **STOP** (8)

Read absorbances at 450 / ≥ 620 nm

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advisces

Use by

Storage temperature

Consult instructions for use