Serazym[®] Anti-Cardiolipin- β 2- GP I IgG/IgM

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IgA: Vaskulitis, Sicca Symptome, Thrombocytopenie.

Literatur:

1. Conrad K, Schöbner W, Hiepe F „Autoantibodies in Systemic Autoimmune Diseases “ A Diagnostic Reference, Vol. 2, 2002. Editors K. Conrad, U. Sack, Pabst Science Publishers Arthritis Rheum. 1999; 42:833-843.
2. McNeil HP, Simpson RJ, Chesterman CN, Krilis SA “Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding factor of coagulation: beta 2-glycoprotein I (apolipoprotein H)” Proc Natl Acad Sci USA 1990, 87:4120-4124.
3. Harris EN, Gharavi AE, Boey ML, Patel BM, Mackworth-Young GG, Loizou S and Hughes GRV “Anti-cardiolipin antibodies: detection by radioimmunoassay and association with thrombosis in systemic lupus erythematosus” Lancet 1983 11:1211

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Verfallsdatum

Lagertemperatur

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Serazym[®] Anti-Cardiolipin- β 2-GP I IgG/IgM

Enzyme immunoassay for detection of IgG and/or IgM antibodies to
cardiolipin- β 2-glycoprotein I -complex
in human serum or plasma

REF E-054  96  *In-vitro*-diagnostic device 

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Autoantibodies to cardiolipin are diagnostic markers for the antiphospholipid syndrome (APS), an autoimmune disease strongly associated with recurrent episodes of venous and arterial thrombosis, spontaneous abortion and thrombocytopenia. Anti-Cardiolipin antibodies are detectable in up to 60 % of patients with Systemic Lupus Erythematosus, in more than 50 % of patients with Juvenile Chronic Arthritis and Vasculitis, in 20-30 % of patients suffering from Rheumatoid Arthritis, in 18-46 % of Cerebral Ischemia in young adults, in 25 % of patients with Primary Sjögren's syndrome and in 21 % of young adults with Myocardial Infarction.

The determination of autoantibodies to cardiolipin requires the presence of the serum protein β 2-glycoprotein I (apolipoprotein H) as cofactor. The cardiolipin- β 2-glycoprotein I-complex is characterized by conformational changes thus building epitopes recognized by anti-cardiolipin autoantibodies. Antibodies to cardiolipin in the absence of β 2-glycoprotein-1 are often detected in patients with infectious diseases. Autoantibodies to cardiolipin can be of the IgG, IgM or IgA isotype with IgG being the most prevalent and clinically significant isotype.

In APS all cardiolipin isotypes may be detectable at all stages of the disease, but the individual isotypes are associated with specific manifestations:

IgG: Venous and/or arterial thrombosis, recurrent spontaneous abortion, neurological symptoms.

IgM: Venous thrombosis, thrombocytopenia, recurrent spontaneous abortion, neurological symptoms, livedo reticularis, endocardial lesions.

IgA: Vasculitis, sicca symptoms, thrombocytopenia.

References:

1. Conrad K, Schöbeler W, Hiepe F „Autoantibodies in Systemic Autoimmune Diseases “ A Diagnostic Reference, Vol. 2, 2002. Editors K. Conrad, U. Sack, Pabst Science Publishers Arthritis Rheum. 1999; 42:833-843.
2. McNeil HP, Simpson RJ, Chesterman CN, Krilis SA “Anti-phospholipid antibodies are directed against a complex antigen that includes a lipid-binding factor of coagulation: beta 2-glycoprotein I (apolipoprotein H)” Proc Natl Acad Sci USA 1990, 87:4120-4124.
3. Harris EN, Gharavi AE, Boey ML, Patel BM, Mackworth-Young GG, Loizou S and Hughes GRV “Anti-cardiolipin antibodies: detection by radioimmunoassay and association with thrombosis in systemic lupus erythematosus” Lancet 1983 11:1211

Intended Use

Serazym® Anti-Cardiolipin- β 2-GP I IgG/IgM is an *in-vitro* diagnostic device for quantitative determination of IgG and/or IgM antibodies to cardiolipin- β 2-glycoprotein I (β 2-GP I)-complex in human serum or plasma.

Principle of the test

In the first incubation step diluted samples, ready-to-use calibrators and positive control react with the solid-phase adsorbed cardiolipin- β 2-GP I-complex. Unbound components are removed after 60 minutes incubation at room temperature (RT, 18...25 °C) by aspirating the samples followed by a washing step.

In the second incubation step horseradish peroxidase (HRP)-labelled anti-human-IgG or -IgM antibodies bind to sample antibodies. Unbound conjugate is removed after 30 minutes incubation at RT by aspiration followed by a washing step.

In the third step the HRP converts the subsequently added colourless substrate solution of hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB) within a 15 min reaction time at RT into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow.

The absorbance of the solution read at 450/ \geq 620 nm (OD) is directly proportional to the amount of specifically bound antibodies. A calibration curve is established by plotting the concentrations of the antibodies of the calibrators (x-axis) and their corresponding absorbance values (y-axis). The antibody concentrations of the unknown specimens are directly read from the calibration curve.

Test components

		For 96 Wells
1	WELLS	Microtitration plate coated with cardiolipin-β2-GP I-complex
		12 single breakable 8-well strips vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold
		100 ml concentrate for 1000 ml solution white cap
3	DIL	Sample diluent
		100 ml · ready-to-use coloured red black cap
4	CAL 0 - 4	Calibrators 0-4 CAL 0 = 1 U/ml CAL 1 = 10 U/ml CAL 2 = 30 U/ml CAL 3 = 100 U/ml CAL 4 = 300 U/ml
		1.0 ml per calibrator ready-to-use coloured blue white cap
5	CONTROL +	Positive control Concentration see leaflet enclosed
		1.0 ml · ready-to-use coloured blue red cap
6	CONJ HRP IgG	HRP-conjugate HRP-labelled, anti-human-IgG antibodies
		15 ml · ready-to-use coloured red red cap
7	CONJ HRP IgM	HRP-conjugate HRP-labelled, anti-human-IgM antibodies
		15 ml · ready-to-use coloured green green cap
8	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide
		15 ml · ready-to-use blue cap
9	STOP	Stop solution 0.25 M sulphuric acid
		15 ml · ready-to-use yellow cap
10		Adhesive film
		2

Preparation and storage of samples

Collection and storage

Blood is taken by venipuncture. Serum or plasma is separated after clotting by centrifugation. Serum or plasma can be used in the *Serazym*[®] Anti-Cardiolipin-β2-GP I IgG/IgM. Contaminated samples should not be used. Repeated freezing and thawing should be avoided. The samples may be kept at 2...8 °C for up to two days. Long-term storage requires -20 °C.

Preparation before use

Allow samples to reach RT prior to assay. Take care to agitate samples gently in order to ensure homogeneity.

Note: Patient samples have to be diluted 1 : 101, e.g. 10 µl sample + 1000 µl sample diluent (3).

Materials required but not provided

Adjustable one-channel micropipette 0.100 - 1.000 ml and 0.010 - 0.100 ml · Adjustable 8-channel micropipette 0.050 - 0.200 ml · Pipette tips · Reagent containers for dispensing with 8-channel pipette · Graduated measuring flasks 10 ml and 100 ml · Test tubes 2.0 ml for sample dilution · Microtitration plate washer (automatic or hand wash head) · Microtitration plate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · Distilled or deionized water · Stop-watch

Preparation and storage of reagents

Kit size and expiry

The *Serazym*[®] Anti-Cardiolipin- β 2-GP I IgG/IgM has been designed for 12 x 8 determinations. The complete kit with unopened reagent bottles and microtitration strips is stable until the expiry date printed on the kit box in case of storage at 2...8 °C. After opening all kit components are stable for at least 2 months, provided proper storage. The ready-to-use wash solution can be used for at least one month when stored at 2...8 °C.

Reagent preparation

The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For example: 10 ml wash buffer (2) + 90 ml distilled or deionized water.

Besides wash buffer all components of the device are ready-to-use.

Assay procedure

Allow all components to reach room temperature prior to use in the assay.

Dilute patient samples with sample diluent (3) 1 : 101 (v/v), e.g. 10 μ l sample + 1000 μ l sample diluent (3).

The calibrators and the positive control are ready-to-use.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CAL 0** Calibrator 0 (4) if desired*
 100 µl **CAL 1 – 4** Calibrator 1, 2, 3, 4 (4)
 100 µl **CONTROL +** positive control (5)
 100 µl **diluted samples** resp.
 into the intended wells.
3. Cover plate (10) and incubate for 60 min at RT.
4. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
5. Dispense 100 µl **CONJ HRP IgG** HRP-conjugate (6) **or**
 100 µl **CONJ HRP IgM** HRP-conjugate (7) per well
6. Cover plate (10) and incubate for 30 min at RT.
7. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
8. Dispense 100 µl **SUBSTR TMB** substrate (8) per well.
9. Incubate for 15 min at RT, protected from light.
10. Dispense 100 µl **STOP** stop solution (9) per well, mix gently.
11. Read OD at 450 nm / ≥ 620 nm with a microplate reader within 30 min after reaction stop.

* Note:

Using CAL 0 (1 U/ml) enables quantification of concentrations below the value of CAL (10 U/ml), otherwise CAL 0 is not needed.

Result interpretation

Create a calibration curve by plotting the mean absorbances of the calibrators CAL 0-4 (Y-axis) to their corresponding antibody-concentrations (X-axis).
 Alternatively create a calibration curve by plotting the mean absorbances of the calibrators CAL 1-4 (Y-axis) to their corresponding antibody-concentrations (X-axis).

For the calibration curve use the 4-Parameter-regression-model.

Determine the antibody-concentrations of the unknown samples by referring their mean absorbances to the calibration curve. In case of using the prescribed sample dilution of 1 : 101, antibody-concentrations are not to be determined by multiplication with the dilution factor, since the calibrators are already prediluted.

If a sample after proper dilution of 1:101 produced a higher absorbance than calibrator CAL 4 please retest the sample at higher dilution and calculate its concentration in the following way:

Sample volume	Added volume of diluent	Final volume	Dilution of original sample	Multiplication factor to determine analyte concentration of the original sample
<i>10µl</i>	<i>1000µl</i>	<i>1010µl</i>	<i>1 : 101 = Sample 1</i>	<i>No multiplication required</i>
500µl Sample 1	500µl	1000µl	1 : 202 = Sample 2	factor 2
500µl Sample 2	500µl	1000µl	1 : 404 = Sample 3	factor 4
500µl Sample 3	500µl	1000µl	1 : 808 = Sample 4	factor 8

Please be aware that antibodies may not be stable in diluted samples resulting in activity loss. By that reason reinvestigation at higher dilution should be done at the same day.

Reference values

Serazym® Anti-Cardiolipin-β2-GP I IgG/IgM	
positive	≥ 10 U/ml
negative	< 10 U/ml

Test validity

A test run is valid if:

- the mean absorbance of calibrator CAL0 is < CAL 1
- the mean absorbance of calibrator CAL1 is ≤ 0.70
- the mean absorbance of calibrator CAL4 is ≥ 1.20

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation volumes, times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure

Any clinical diagnosis should not be based on the results of *in-vitro* diagnostic methods alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis. False results may be caused by cross contaminations of the kit reagents and samples, bacterial or fungal contaminations of reagents and/or samples, incorrect washing and incorrect incubation times. In case of background problems when using automatic microplate ELISA systems it may be recommendable to wash wells 5 times (instead of 3 times) in every wash cycle.

Sera of healthy individuals should be tested negative in the Serazym® Anti-Cardiolipin-β2-GP I IgG/IgM, but 3-5 % of healthy individuals exhibit anti-cardiolipin antibody-concentrations exceeding the cut-off value. Since anti-cardiolipin antibodies may also occur in infectious diseases, the test should be repeated after at least 6 weeks to ensure the diagnosis of APS. Anti-cardiolipin antibodies associated with infectious diseases do not require β2-GP I to recognize cardiolipin and they are not associated with the symptoms of APS.


Interference

Haemolytic, lipaemic and icteric samples did not interfere with the determination of IgG and IgM anti-cardiolipin-β2-I-antibodies up to a concentration of 500mg/dl (haemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors did not interfere with the determination of IgG anti- cardiolipin-β2-I-antibodies up to a concentration of 500 IU/ml, but can interfere with the determination of IgM anti- cardiolipin-β2-I-antibodies.

Performance characteristics

Frequency distribution

Frequency distribution of anti-Cardiolipin- β 2-GP I antibody concentrations in 400 blood donor serum samples measured in *Serazym*[®] Anti-Cardiolipin- β 2-GP I IgG/IgM:

Precision

Intra-assay coefficient of variation (CV [%]) in the *Serazym*[®] Anti-Cardiolipin- β 2-GP I IgG/IgM from 8-fold determinations of samples:


sample	IgG		IgM	
	concentration [U/ml]	CV [%]	Concentration [U/ml]	CV [%]
serum 1	26.9	2.3	27.4	4.1
serum 2	60.7	8.9	51.2	8.8
serum 3	143.4	6.0	109.0	4.0
serum 4	200.7	9.1	264.7	8.8

Inter-assay coefficient of variation (CV [%]) in the *Serazym*[®] Anti-Cardiolipin- β 2-GP I IgG/IgM from 8 different test runs with 3-fold determinations of samples:

sample	IgG		IgM	
	concentration [U/ml]	CV [%]	concentration [U/ml]	CV [%]
serum A	42.7	10.3	30.8	7.4
serum B	81.2	8.2	78.8	9.5
serum C	129.4	10.4	136.7	10.1
serum D	215.0	6.2	184.4	13.9

Typical reference curve

The typical reference curve in the *Serazym*[®] Anti-Cardiolipin- β 2-GP I IgG/IgM is illustrated:

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for diluted reagents. Do not use or mix reagents from different lots, damaged packages or bottles or reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8 °C before use.

Do not use or mix reagents from different lots except for sample diluent, wash buffer, substrate and stop solution.

Some of the reagents contain small amounts of Thimerosal (< 0.01 % w/v) and Kathon (1.0 % v/v) as preservative. They may neither be swallowed nor come into contact with skin or mucous membranes. In case of contact rinse with plenty of water. Handle all samples as potentially infectious. Since the kit contains potentially hazardous materials, the following precautions must always be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

Incubation scheme *Serazym*[®] Anti-Cardiolipin- β 2-GP I IgG/IgM

- | | | | |
|----|---|-------------|--|
| 1. |  | 100 μ l | CAL 0 – 4 (4) |
| | | 100 μ l | CONTROL + (5) |
| | | 100 μ l | diluted samples resp. |
| |  | 60 min | incubation (room temperature) |
| | | 3 x wash | with wash solution |
| 2. |  | 100 μ l | CONJ HRP IgG (6) or CONJ HRP IgM (7) |
| | | 30 min | incubation (room temperature) |
| | | 3 x wash | with wash solution |
| |  | | |
| 3. |  | 100 μ l | SUBSTR TMB (8) |
| | | 15 min | incubation at RT protected from light |
| 4. |  | 100 μ l | STOP (9) |

Read absorbances at 450 / \geq 620 nm

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use