## Serazym<sup>®</sup> H. pylori 2nd Gen.

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# Serazym<sup>®</sup> H. pylori 2nd Gen.

Enzyme immunoassay for the qualitative detection of *Helicobacter pylori*-specific antigen in stool samples of human origin

**REF** E-114-A  96  
**IVD** In-vitro-diagnostic medical device **CE**

Seramun Diagnostica GmbH • Spreehagener Str. 1 • 15754 Heidesee • Germany •  
T +49 33767 791-10 • [info@seramun.com](mailto:info@seramun.com) • [www.seramun.com](http://www.seramun.com)

- |   |   |  |
|---|---|--|
| <b>IVD</b> In-vitro diagnostic medical device   | <b>UDI</b> Unique device identifier   |  Manufacturer  |
|  Country of manufacture and date of manufacture | <b>REF</b> Article number   | <b>SN</b> Serial number  |
|  Keep away from sunlight                        |  Humidity limitation | <b>LOT</b> Batch code  |
|  Consult instructions for use                   |  Temperature limit   |  Do not reuse |
|  Sufficient for <i>n</i> tests                  |  Biohazard           |  Use-by date  |
|   |   |  Attention    |

## Intended Use

Serazym® *H. pylori* 2nd Gen. is an IVD test for the qualitative determination of *H. pylori* specific antigen in stool samples of human origin by a health professional in a laboratory environment.

The test can be used with the devices DS2® and DSX® (Dynerx Technologies) as well as ThunderBolt® (Gold Standard Diagnostics) for the semi-automatic performance.

It is intended to aid in the diagnosis of a *H. pylori* infection in samples from patients with symptoms of a *H. pylori* associated gastroenteritis, for the monitoring of an infection during eradication therapy and after eradication therapy.

The test must not be used with specimen materials other than stool samples of human origin, for diagnosis, screening, prediction, prognosis, as a companion diagnostic, in the near-patient setting and by lay persons.

## Principle of the Test

Serazym® *H. pylori* 2nd Gen. is a direct two-sided assay based on monoclonal antibodies against *Helicobacter pylori* antigen. Diluted, untreated stool samples as well as negative and positive control samples are dispensed into wells of the microtiter plate coated with monoclonal anti-*Helicobacter pylori* antibodies. After incubation unbound components are removed by a washing step, then peroxidase (HRP)-labelled monoclonal anti-*Helicobacter pylori* antibodies are dispensed into the wells. After incubation and a washing step, the colorless substrate solution in an enzymatic reaction is converted into a blue reaction product. The reaction is stopped by addition of the stop solution, resulting in a color change from blue to yellow. The optical density (OD) of the reaction product measured at 450 nm measuring and 620 nm reference filter, respectively, is directly proportional to the concentration of specifically bound *Helicobacter pylori* antigens.

## Test Components (Delivery Scope)

		<b>For 96 wells</b>	
1	<b>WELLS</b>	<b>Microtiter plate</b> coated with monoclonal anti- <i>H. pylori</i> antibodies (mouse)	12 single breakable 8-well strips, wine red color marking, vacuum-sealed with desiccant
2	<b>WASHBUF (10x)</b>	<b>Wash buffer (10x)</b> Seramun® Wash buffer A TRIS-based buffer	100 mL concentrate for 1000 mL solution, colorless, white cap
3	<b>DIL</b>	<b>Sample diluent</b> Seramun® Sample diluent A Phosphate-based buffer	100 mL, ready-to-use, colored yellow black cap
4	<b>CONTROL +</b>	<b>Positive control</b> <i>H. pylori</i> positive sample (inactivated)	2.0 mL, ready-to-use, colored blue, red cap
5	<b>CONTROL -</b>	<b>Negative control</b> TRIS-based buffer	2.0 mL, ready-to-use, colored blue, green cap
6	<b>CONJ HRP</b>	<b>HRP conjugate</b> HRP-labelled monoclonal anti- <i>H. pylori</i> antibodies (mouse)	15 mL, ready-to-use, colored green, violet cap
7	<b>SUBSTR</b>	<b>Substrate</b> SeramunBlau® automat fast < 0.1 % 3,3',5,5'-tetramethylbenzidine	15 mL, ready-to-use, colorless, blue cap
8	<b>STOP</b>	<b>Stop solution</b> SeramunBlau® stop 0.25 M sulphuric acid	15 mL, ready-to-use, colorless, yellow cap
9	<b>COVER</b>	<b>Covering film</b>	2 pieces

10 Certificate of Analysis 1 piece

11 Instructions for Use 1 piece

## Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel-micropipette • 8-channel pipette or multi-pipette with pipette tips • 8-channel wash comb with vacuum pump and waste bottle or microplate washer • microplate reader with 450 nm measuring filter and  $\geq 620$  nm reference filter • deionized water • measuring cylinder • tubes for sample preparation

## Important Information

**This device is for *in-vitro* diagnostic use only.** Follow the instructions carefully. The kit may be performed by laboratory professionals only.

Do not use reagents from damaged packages or bottles. Imprinted expiry dates must be observed. Do not mix components with reagents from other manufacturers.

**Mixing of test kit components of different lots is only allowed for sample diluent, wash buffer, substrate and stop solution.**

**Sample diluent, wash buffer, substrate and stop solution are universally applicable for Serazym<sup>®</sup> stool ELISA *Campylobacter* (E-093), *Clostridium difficile* GDH (E-107), *Clostridium difficile* Toxin A+B (E-040), *H. pylori* 2nd Gen. (E-114), Adenovirus (E-017), Astrovirus (E-045), Norovirus (E-061), Rotavirus (E-020), *Cryptosporidium parvum* (E-039), *Entamoeba histolytica* (E-018), and *Giardia* (E-106). In addition, wash buffer, substrate and stop solution of Serazym<sup>®</sup> Verotoxin (E-030) may be universally applied within the Serazym<sup>®</sup> stool ELISA.**

All serious incidents occurring in relation with Serazym<sup>®</sup> *H. pylori* 2nd Gen. must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

### Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all test components to room temperature before use. Reagents that appear contaminated should not be used.

Each well of a microtiter plate can be used once only.

Positive and negative control are ready-to-use.

For larger sample series, pipetting reagents from liquid reservoirs using a multichannel pipette is recommended to avoid time delays and contaminations. Follow the pipetting scheme and time schedules of the protocol.

Each sample and control should be pipetted with a new pipette tip.

Wash solution should be allowed a minimum reaction time of 5 s in the wells per wash cycle. Remove wash buffer residues by thoroughly aspirating or tapping out the cavities! The aspiration and washing steps can be performed manually or with the help of a microplate washer or waterjet pump.

Protect substrate from light!

### Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	EUH208 EUH210	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3 : 1). May produce an allergic reaction. Safety data sheet available on request.
DIL	-	Contains material of animal origin.
CONTROL +	-	Contains material of microbial and animal origin.
CONTROL -	-	Contains material of animal origin.
CONJ HRP	Hazard components H360D P202  P280  P308+P313 P501	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone  May damage the unborn child. Do not handle until all safety precautions have been read and understood. Wear protective gloves/protective clothing/eye protection/face protection. IF exposed or concerned: Get medical advice/attention. Dispose of contents/container in accordance with local regulations. Restricted to professional users. Contains material of animal origin.
SUBSTR	Hazard components H319 P264 P280 P305+P351 +P338  P337+P313	2-Pyrrolidone  Causes serious eye irritation. Wash hands thoroughly after handling. Wear protective gloves/protective clothing/eye protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
STOP	Hazard components H290	Sulphuric acid 2.5 %  May be corrosive to metals.

### Limitations of the Procedure

There is no correlation between absorbance measured and severity of the infection. Also, it is not allowed to correlate absorbances of samples with the absorbance of the positive control.

Cross contamination of reagents and samples may result in false positive results. Incorrect dilutions, insufficiently homogenized samples, and particles not sedimented by centrifugation may cause false negative as well as false positive test results. A negative test result obtained with Serazym® H. pylori 2nd Gen. does not exclude an infection. False negative tests may result from improper timing of sample collection or inhomogeneous antigen distribution in the sample. The overall interpretation of the ELISA test result should be made in the context of the clinical manifestations. Individual cases may require repeat testing at intervals of several weeks.

## Sample Treatment

### Sample Collection

Collect stool sample in suitable sampling container.

### Sample Shelf Life and Storage

Stool samples should be stored at 2...8 °C immediately after collection and examined within 72 h or stored frozen at -20 °C. Repeated freezing (> 3 x) and thawing of samples should be avoided due to the risk of incorrect results. Stool samples that have already been diluted in Seramun® Sample diluent

A according to the instructions for use can be stored at 2...8 °C for 72 h and subsequently analyzed by ELISA.

### Sample Preparation

Mix untreated stool samples well and dilute 1 : 6 with sample buffer.

Pipette 500 µL sample buffer into a reaction tube. For solid or semi-solid stool samples transfer 100 mg (diameter approx. 2 – 3 mm) with a disposable stick, for liquid stool samples transfer 100 µL into the sample buffer and mix thoroughly. If necessary, sediment suspended particles by centrifugation in a microcentrifuge for 1 min at maximum speed.

## Reagent Treatment

### Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips can be stored at 2...8 °C until the printed expiration date. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

### Reagent Preparation

Microtiter plate with breakable 8-well strips vacuum sealed with desiccant. Allow packaging to reach room temperature before opening. Protect unused wells from moisture and store refrigerated with desiccant in the original bag carefully resealed. Dilute wash buffer (10x) 1 : 10 with deionized water.

Example: 10 mL Seramun® Wash buffer A (10x) + 90 mL deionized water.

## Assay Procedure

1. Allow test reagents and required number of wells to reach room temperature (RT). Shake reagents gently before use. Avoid foaming.
2. Pipette 100 µL **CONTROL +** Positive control  
100 µL **CONTROL -** Negative control  
100 µL diluted stool specimen each.
3. Cover the plate and incubate for 60 min at RT.
4. Decant, then wash each well 5 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
5. Add 3 drops (or 100 µL) **CONJ HRP** HRP conjugate per well.
6. Cover plate and incubate for 30 min at RT.
7. Decant, then wash each well 5 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
8. Add 3 drops (or 100 µL) **SUBSTR** substrate per cavity.
9. Incubate for 10 min at RT **protected from light**.
10. Add 3 drops (or 100 µL) **STOP** stop solution per cavity, mix gently.
11. Read OD at 450 nm and ≥ 620 nm with a microplate reader within 30 min following reaction stop.

## Evaluation of Results

### Qualitative Evaluation:

Cut-off determination: OD negative control + 0.10

Samples showing OD values equal with or higher than the cut-off are considered positive, samples with OD values below cut-off are considered negative for *Helicobacter pylori* antigens.

The test run is valid if:

- mean OD of the negative control is  $\leq 0.20$  (manual performance)  
 $\leq 0.30$  (automatic performance)
- mean OD of the positive control is  $\geq 1.20$

If the above-mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the quality criteria contact the manufacturer.

## Interpretation of Results

Positive	$\geq$ cut-off
Negative	$<$ cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges.

## Performance Characteristics

### Precision

To determine precision, 3 stool samples were measured multiple times. For the determination of the intra-assay coefficient of variation samples were measured in a 24-fold determination in one test run. The determination of the inter-assay coefficient of variation was done by a 3-fold determination on 2 days in 10 different test runs. The lot-to-lot coefficient of variation was determined by 3-fold determination in 3 lots of the product.

Sample	Intra-assay-coefficient of variation		Inter-assay-coefficient of variation		Lot-to-Lot-coefficient of variation	
	$\bar{x}$ OD	CV (%)	$\bar{x}$ OD	CV (%)	$\bar{x}$ OD	CV (%)
1	2.377	5.3	2.483	6.3	2.626	10.4
2	0.871	7.4	1.044	13.2	1.265	17.4
3	0.124	11.0	0.142	10.3	0.162	18.5

### Detection Limit

The lower detection limit of the Serazym<sup>®</sup> H. pylori 2nd Gen. has been determined by titration of faecal samples spiked with purified *H. pylori* antigen. The lower detection limit is  $< 31.3$  ng/mL.

### Determination of the cut-off Value

The frequency distribution of the absorbances of 169 faecal samples has been investigated in order to fix the cut-off value for the Serazym<sup>®</sup> H. pylori 2nd Gen. the cut-off has been determined with absorbance negative control + 0.1.

## Sensitivity and Specificity

Sensitivity and specificity of Serazym® H. pylori 2nd Gen. have been determined in a retrospective study in comparison to two commercially available ELISAs.

n = 233	ELISA 1 positive	ELISA 1 negative
<b>Serazym® ELISA positive</b>	140	6**
<b>Serazym® ELISA negative</b>	2*	85

Sensitivity: 98.6 %

Specificity: 93.4 %

Samples labeled with \* and \*\* were re-tested in a commercially available lateral flow assay. This results in a corrected sensitivity of 99.3 % and a corrected specificity of 97.7 %. Compared to ELISA 1 the accuracy is 98.7 %.

n = 79	ELISA 2 positive	ELISA 2 negative
<b>Serazym® ELISA positive</b>	53	5***
<b>Serazym® ELISA negative</b>	0	21

Sensitivity: 100 %

Specificity: 80,8 %

Samples labeled with \*\*\* were confirmed positive by comparative ELISA 1. This results in a corrected specificity of 100 %. Compared to ELISA 2 the accuracy is 100 %.

## Cross Reactivity

The following microorganisms were diluted in sample diluent with  $\geq 10^8$  colony-forming units (CFU) per mL and tested negative with Serazym® H. pylori 2nd Gen. (OD 450 / 620 nm < cut-off):

<i>Aeromonas hydrophila</i>	ATCC 7966	<i>Campylobacter coli</i>	ATCC 33559
<i>Bacillus cereus</i>	ATCC 11778	<i>Campylobacter jejunii</i>	ATCC 32291
<i>Bacillus subtilis</i>	ATCC 6633	<i>Campylobacter fetus</i>	ATCC 27374
<i>Bacteroides fragilis</i>	ATCC 25285	<i>Campylobacter upsaliensis</i>	ATCC 43954
<i>Citrobacter freundii</i>	ATCC 8090	<i>Campylobacter lari</i>	ATCC 35221
<i>Clostridium sordelli</i>	ATCC 9714	<i>Vibrio cholerae</i>	clinical isolate
<i>Enterobacter aerogenes</i>	ATCC 13048	<i>Yersinia enterocolitica</i> 0:3	clinical isolate
<i>Enterobacter cloacae</i>	ATCC 13047	<i>Yersinia enterocolitica</i> 0:9	clinical isolate
<i>Enterococcus faecalis</i>	ATCC 29212	<i>Yersinia enterocolitica</i> Y11	clinical isolate
<i>Escherichia coli</i>	ATCC 25922	<i>Yersinia enterocolitica</i> RKI 0803733	clinical isolate
<i>Klebsiella pneumonia</i>	ATCC 13883	<i>Clostridium difficile</i>	VPI 10463
<i>Peptostreptococcus anaerobius</i>	ATCC 27337	<i>Salmonella infantis</i>	ATCC 51741
<i>Proteus vulgaris</i>	ATCC 8427	<i>Salmonella anatum</i>	ATCC 9270
<i>Pseudomonas aeruginosa</i>	ATCC 10145	<i>Salmonella paratyphi</i> A	ATCC 11511
<i>Salmonella enterica</i> Serovar <i>typhimurium</i>	ATCC 14028	<i>Salmonella paratyphi</i> B	ATCC 8759
<i>Salmonella enterica</i> ssp. <i>galolyticus/enteritidis</i>	ATCC 13076	<i>Salmonella paratyphi</i> C	Nr. 2 Pasteur
<i>Shigella flexneri</i>	ATCC 12022	<i>Lactococcus lactis</i>	DSM 20481
<i>Shigella sonnei</i>	ATCC 25931	<i>Proteus mirabilis</i>	ATCC 29906
<i>Staphylococcus aureus</i>	ATCC 25923	<i>Pseudomonas fluorescens</i>	ATCC 13525
<i>Staphylococcus epidermidis</i>	ATCC 12228	<i>Pseudomonas putida</i>	ATCC 49128
<i>Vibrio parahaemolyticus</i>	ATCC 17802	<i>Streptococcus agalactiae</i>	ATCC 13813
<i>Candida albicans</i>	ATCC 10231	<i>Morganella morganii</i>	ATCC 25830

## Interference

None of the following substances in the indicated concentrations added to *Helicobacter pylori* positive and negative stool samples did show a significant impact on the test result. The mentioned concentrations refer to 100 mg faeces:

(-)-Scopolamine N-butyl bromide (0.5 %, Buscopan®), barium sulfate (5 %), bismuth(III) subsalicylate (0.5 %, Pepto-Bismol), Cyclamat (5 %), Diclofenac (0.5 %), hemoglobin human (5 %), blood human (5 %), Hylak® N (5 %), Iberogast® (5 %), Immodium® akut duo (0.03/1.9 %), Loperamid hydrochloride (5 %, Loperamid-CT akut), Metronidazole (0.5 %), Mucin (5 %), Nexium® (0.03 %), Nifuroxazide (0.5 %, Pentofuryl®), palmitic acid (20 %), Perenterol forte (0.5 %), Rennie® (20 %), Simagel® (1 %), stearic acid (20 %), Vancomycin (0.5 %).

## Clinical performance

Seven clinical studies have been carried out with comparative ELISA 1 in Europe and the Arabian region between 2001 and 2005, which clearly demonstrated the clinical performance of comparative ELISA 1. The comparative tests with Serazym® H. pylori 2nd Gen. therefore result in the following data:

<b>Diagnostic sensitivity</b>	99.3 %
<b>Diagnostic specificity</b>	97.7 %
<b>Positive predictive value</b>	98.6 %
<b>Negative predictive value</b>	98.9 %
<b>Likelihood ratio</b>	43.2 (LR+) 0.007 (LR-)

Four clinical studies have been carried out with comparative ELISA 2 in Europe and Asia in the period from 2008 to 2020, which clearly confirmed the clinical performance of comparative ELISA 2. Therefore, the comparative tests with the Serazym® H. pylori 2nd Gen. yield in the following results:

<b>Diagnostic Sensitivity</b>	100 %
<b>Diagnostic Specificity</b>	100 %
<b>Positive Predictive Value</b>	100 %
<b>Negative Predictive Value</b>	100 %
<b>Likelihood ratio</b>	LR+ → +∞ LR- → 0


## Application

### Automatic Processing


Performing Serazym® H. pylori 2nd Gen. on fully automated microplate processors (e.g., DS2®, DSX®, Dynex Technologies or ThunderBolt®, Gold Standard Diagnostics) may cause elevated absorbance values in comparison to the manual procedure caused by differences in the wash procedures and technical specifications of the equipment. In these cases, a maximum value of OD = 0.3 is permissible for the negative control. It is recommended to program a wash protocol with 10 s soak time per strip and wash step. A final wash step with deionized water and a soak time of 10 s is recommended after each wash cycle. If necessary, the number of washing steps may be increased to 7 x or 8 x.

### Correlation: manual – automatic processing

A panel of 93 stool specimens was processed manually and automatically in parallel (DS2®, Dynex Technologies). The correlation coefficient was calculated at  $r = 0.928$ .

A panel of 94 stool specimens was processed manually and automatically in parallel (ThunderBolt®, Gold Standard Diagnostics). The correlation coefficient was calculated at  $r = 0.995$ .

## Change History

Version	Section	Modifications
2022-08_v01_DE_EN	Entire document	Updating of the intended use Conversion of subsections Insertion of safety instructions Updating of the performance characteristics Inclusion of ThunderBolt® data in chapter "Application" Removal of pipetting scheme

## References

1. Fischbach, W. et al. S2k-Leitlinie Helicobacter pylori und gastroduodenale Ulkuserkrankheiten. Z Gastro-enterol 2016; 54: 327–363
2. Malfertheiner P, et al. Management of Helicobacter pylori infection—the Maastricht V/Florence Consensus Report. Gut 2016; 0: 1–25