Serazym[®] Anti-dsDNA IgG

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Literatur:

1. Tan, E.M., Cohen, A.S., Fries, J.F. et al. "The 1982 revised criteria for the classification of systemic lupus erythematosus"
Arthritis Rheum. 25, 1271-1277, 1982
2. Schwartz, R.S. "Anti-DNA and the problem of autoimmunity"
Cell. Immunol. 99, 38-43, 1986
3. Isenberg, D. and Shoenfeld, Y. "The origin and significance of anti-DNA antibodies"
Immunol. Today 8, 279-282, 1987
4. Conrad, K., Schöblier, W., Hiepe, F., Fritzler, M.J. „Autoantibodies in systemic autoimmune diseases – a diagnostic reference“, volume 2, 2015, Pabst Science Publishers, Lengerich Germany
5. Feltkamp D.E.W. et al "The first international standard for antibodies to double stranded DNA"
Ann Rheum Dis 47, 740-746, 1988

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

Serazym[®] Anti-dsDNA IgG

Enzyme immunoassay for detection of IgG antibodies to dsDNA
in human serum or plasma

 E-058  96  *In-vitro*-diagnostic device 

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Systemic Lupus erythematosus (SLE) is a chronic inflammatory autoimmune disease affecting multiple organ systems. Beside clinical findings (Classification criteria of the American College of Rheumatology (ACR) for SLE) a large variety of autoantibodies may be supportive markers for the diagnosis of SLE. Among these autoantibodies directed to nuclear antigens (ANA), antibodies to double stranded DNA (anti-dsDNA) are highly pathogenic by formation of circulating immune complexes.

Two populations of anti-dsDNA antibodies exist:

- Antibodies to the phosphorylated ribose site-chain, which cross-react with single stranded DNA (ssDNA)
- Antibodies to the deoxyribose phosphate skeleton, which do not cross-react.

Anti-dsDNA antibodies are not only of superior diagnostic significance (ACR criterion), their titre is also correlated with disease activity. Moreover anti-dsDNA antibodies are considered as prognostic marker with regard to renal and cerebral manifestations of the disease.

Therefore the frequency of the determination of anti-dsDNA antibodies depends on the disease activity and the involvement of organs especially the kidney (lupus nephritis):

Inactive SLE	< 40 %
Active SLE without nephritis	> 50 – 70 %
Active SLE with nephritis	> 95 % (3)

Availability of an international standard for antibodies to dsDNA is a decisive factor for standardization of their determination by immunoassays.

Reference:

1. Tan, E.M., Cohen, A.S., Fries, J.F. et al. "The 1982 revised criteria for the classification of systemic lupus erythematosus" *Arthritis Rheum.* 25, 1271-1277, 1982
2. Schwartz, R.S. "Anti-DNA and the problem of autoimmunity" *Cell. Immunol.* 99, 38-43, 1986
3. Isenberg, D. and Shoenfeld, Y. "The origin and significance of anti-DNA antibodies" *Immunol. Today* 8, 279-282, 1987
4. Conrad, K., Schöblier, W., Hiepe, F., Fritzler, M.J. „Autoantibodies in systemic autoimmune diseases – a diagnostic reference”, volume 2, 2015, Pabst Science Publishers, Lengerich Germany
5. Feltkamp D.E.W. et al "The first international standard for antibodies to double stranded DNA" *Ann Rheum Dis* 47, 740-746, 1988

Intended use

Serazym® Anti-dsDNA IgG is an *in-vitro* diagnostic device for quantitative determination of IgG antibodies to double stranded DNA in human serum or plasma.

Principle of the test

In the first incubation step diluted samples, ready-to-use calibrators and positive control react with the solid-phase adsorbed double stranded DNA. Unbound components are removed after 60 minutes incubation at room temperature (RT, 18...25 °C) by aspirating the samples followed by a washing step.

In the second incubation step horseradish peroxidase (HRP)-labelled anti-human-IgG antibodies bind to sample antibodies. Unbound conjugate is removed after 30 minutes incubation at RT by aspiration followed by a washing step.

In the third step the HRP converts the subsequently added colourless substrate solution of hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB) within a 15 min reaction time at RT into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells turning the solution from blue to yellow.

The absorbance of the solution read at 450/≥620 nm (OD) is directly proportional to the amount of specifically bound antibodies. A calibration curve is established by plotting the concentrations of the antibodies of the calibrators (x-axis) and their corresponding absorbance values (y-axis). The antibody concentrations of the unknown specimens are directly read from the calibration curve.

Test components

			For 96 Wells
1	WELLS	Microtitration plate coated with antigen	12 single breakable 8-well strips vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready-to-use coloured red black cap
4	CAL 0 - 4	Calibrator 0 - 4 diluted serum samples CAL 0 = 1 IU/ml CAL 1 = 10 IU/ml CAL 2 = 30 IU/ml CAL 3 = 100 IU/ml CAL 4 = 300 IU/ml	1.0 ml per calibrator ready-to-use coloured blue white cap
5	CONTROL +	Positive control diluted serum sample concentration see certificate of analysis enclosed	1.0 ml · ready-to-use coloured blue red cap
6	CONJ HRP	HRP-conjugate HRP-labelled, anti-human-IgG antibodies	15 ml · ready-to-use coloured red red cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready-to-use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready-to-use yellow cap
9		Adhesive film	2

Preparation and storage of samples

Collection and storage

Blood is taken by venipuncture. Serum or plasma is separated after clotting by centrifugation. Serum or plasma can be used in the *Serazym*[®] Anti-dsDNA IgG. Contaminated samples should not be used. Repeated freezing and thawing should be avoided. The samples may be kept at 2...8 °C for up to two days. Long-term storage requires -20 °C.

Preparation before use

Allow samples to reach RT prior to assay. Take care to agitate samples gently in order to ensure homogeneity.

Note: Patient samples have to be diluted 1 : 101, e.g. 10 µl sample + 1000 µl sample diluent (3).

Materials required but not provided

Adjustable one-channel micropipette 0.100 - 1.000 ml and 0.010 - 0.100 ml · Adjustable 8-channel micropipette 0.050 - 0.200 ml · Pipette tips · Reagent containers for dispensing with 8-channel pipette · Graduated measuring flasks 10 ml and 100 ml · Test tubes 2.0 ml for sample dilution · Microtitration plate washer (automatic or hand wash head) · Microtitration plate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · Distilled or deionized water · Stop-watch

Preparation and storage of reagents

Kit size and expiry

The *Serazym*[®] Anti-dsDNA IgG has been designed for 12 x 8 determinations.

The complete kit with unopened reagent bottles and microtitration strips is stable until the expiry date printed on the kit box in case of storage at 2...8 °C. After opening all kit components are stable for at least 2 months, provided proper storage. The ready-to-use wash solution can be used for at least one month when stored at 2...8°C.

Reagent preparation

The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed.

Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For example: 10 ml wash buffer (2) + 90 ml distilled or deionized water.

Besides wash buffer all components of the device are ready-to-use.

Assay procedure

Allow all components to reach room temperature prior to use in the assay.

Dilute patient samples with sample diluent (3) 1 : 101 (v/v), e.g. 10 μ l sample + 1000 μ l sample diluent (3).

The calibrators and the positive control are ready-to-use.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash solution in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 100 µl **CAL 0** calibrator 0 (4) if desired*,
 100 µl **CAL 1 - 4** calibrators 1, 2, 3, 4 (4)
 100 µl **CONTROL +** positive control (5)
 100 µl **diluted samples** resp.
 into the intended wells.
3. Cover plate (9) and incubate for 60 min at RT.
4. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
5. Dispense 100 µl **CONJ HRP** HRP-conjugate (6) per well.
6. Cover plate (9) and incubate for 30 min at RT.
7. Decant, then wash each well 3x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper.
8. Dispense 100 µl **SUBSTR TMB** substrate (7) per well.
9. Incubate for 15 min at RT protected from light.
10. Dispense 100 µl **STOP** stop solution (8) per well, mix gently.
11. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction stop.

* Note:

Using CAL 0 (1 IU/ml) enables quantification of concentrations below the value of CAL 1 (10 IU/ml), otherwise CAL 0 is not needed.

Result interpretation

Create a calibration curve by plotting the mean absorbances of the calibrators CAL 0-4 (Y-axis) to their corresponding antibody-concentrations (X-axis).
 Alternatively create a calibration curve by plotting the mean absorbances of the calibrators CAL 1-4 (Y-axis) to their corresponding antibody-concentrations (X-axis).

For the calibration curve use the 4-Parameter-regression-model.

Determine the antibody-concentrations of the unknown samples by referring their mean absorbances to the calibration curve. In case of using the prescribed sample dilution of 1 : 101, antibody-concentrations are not to be determined by multiplication with the dilution factor, since the calibrators are already prediluted.

If a sample after proper dilution of 1:101 produced a higher absorbance than calibrator CAL 4 please retest the sample at higher dilution and calculate its concentration in the following way:

Sample volume	Added volume of diluent	Final volume	Dilution of original sample	Multiplication factor to determine analyte concentration of the original sample
<i>10µl</i>	<i>1000µl</i>	<i>1010µl</i>	<i>1 : 101 = Sample 1</i>	<i>No multiplication required</i>
500µl Sample 1	500µl	1000µl	1 : 202 = Sample 2	factor 2
500µl Sample 2	500µl	1000µl	1 : 404 = Sample 3	factor 4
500µl Sample 3	500µl	1000µl	1 : 808 = Sample 4	factor 8

Please be aware that antibodies may not be stable in diluted samples resulting in activity loss. By that reason reinvestigation at higher dilution should be done at the same day.

Reference values

Serazym® Anti-dsDNA IgG	
positive	> 35 IU/ml
Grey zone	30 – 35 IU/ml
negative	< 30 IU/ml

Samples giving borderline test results should be tested repeatedly.

Test validity

The test run is valid if:

- absorbance of calibrator CAL 0 is < CAL 1
- absorbance of calibrator CAL 1 is ≤ 0.50
- absorbance of calibrator CAL 4 is ≥ 1.20

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation volumes, times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the method

Any clinical diagnosis should not be based on the results of *in-vitro* diagnostic methods alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis. False results may be caused by cross contaminations of the kit reagents and samples, bacterial or fungal contaminations of reagents and/or samples, incorrect washing and incorrect incubation times. In case of background problems when using automatic microplate ELISA systems it may be recommendable to wash wells 5 times (instead of 3 times) in every wash cycle.

Interference

Haemolytic, lipaemic and icteric samples do not interfere with the test up to a concentration of 500 mg/dl (haemoglobin), 1000 mg/dl (lipids) and 20 mg/dl (bilirubin C and bilirubin F). Rheumatoid factors do not interfere up to a concentration of 500 IU/ml.


Performance characteristics

Calibration

Since the *Serazym*® Anti-dsDNA IgG is calibrated according to the WHO standard (Wo/80) the determined anti-dsDNA IgG antibody concentrations are expressed as IU/ml.

Frequency distribution

The frequency distribution of the anti-dsDNA antibody concentrations measured in 211 samples in the *Serazym*[®] Anti-dsDNA IgG is illustrated:

Precision

Intra-assay coefficient of variation (CV [%]) in the *Serazym*[®] Anti-dsDNA IgG from 8-fold determinations of samples:


sample	concentration [IU/ml]	CV [%]
1	180	3.34
2	92	3.96
3	67	4.48
4	41	2.36

Inter-assay coefficient of variation (CV [%]) in the *Serazym*[®] Anti-dsDNA IgG in 8 different test runs from twofold determination of samples:

sample	concentration [IU/ml]	CV [%]
A	157	7.60
B	130	8.42
C	92	5.02
D	41	10.6

Typical reference curve

The typical reference curve in the *Serazym*[®] Anti-dsDNA IgG is illustrated:

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. The expiration dates stated on the respective labels are to be observed. The same relates to the stability stated for diluted reagents. Do not use or mix reagents from different lots, damaged packages or bottles or reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8 °C before use.

Do not use or mix reagents from different lots except for sample diluent, wash buffer, substrate and stop solution.

Some of the reagents contain small amounts of Thimerosal (< 0.01 % w/v) and Kathon (1.0 % v/v) as preservative. They may neither be swallowed nor come into contact with skin or mucous membranes. In case of contact rinse with plenty of water. Handle all samples as potentially infectious. Since the kit contains potentially hazardous materials, the following precautions must always be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

Incubation scheme *Serazym*[®] Anti-dsDNA IgG

- | | | | |
|----|---|----------|--|
| 1. |  | 100 µl | CAL 0 - 4 (4) |
| | | 100 µl | CONTROL + (5) |
| | | 100 µl | diluted samples resp. |
| |  | 60 min | incubation (room temperature) |
| | | 3 x wash | with wash solution |
| 2. |  | 100 µl | CONJ HRP (6) |
| | | 30 min | incubation (room temperature) |
| |  | 3 x wash | with wash solution |
| 3. |  | 100 µl | SUBSTR TMB (7) |
| | | 15 min | incubation (room temperature) protected from light |
| 4. |  | 100 µl | STOP (8) |

Read OD at 450 / ≥ 620 nm

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use