
# Serazym<sup>®</sup> ANA screen

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# Serazym<sup>®</sup> ANA screen

Enzyme immunoassay for the qualitative or semiquantitative detection of antinuclear antibodies (ANA) of the IgG isotype against nuclear antigens and Jo-1 in serum or plasma of human origin

<b>REF</b>	E-065		96
<b>IVD</b>	In-vitro-diagnostic medical device		

 **Seramun Diagnostica GmbH** • Sprenhagener Str. 1 • 15754 Heidesee • Germany •  
T +49 33767 791-10 • [info@seramun.com](mailto:info@seramun.com) • [www.seramun.com](http://www.seramun.com)

<b>IVD</b> In-vitro diagnostic medical device	<b>UDI</b> Unique device identifier	 Manufacturer
 Country of manufacture and date of manufacture	<b>REF</b> Article number	<b>SN</b> Serial number
 Keep away from sunlight	 Humidity limitation	<b>LOT</b> Batch code
 Consult instructions for use	 Temperature limit	 Do not reuse
 Sufficient for <i>n</i> tests	 Biohazard	 Use-by date
		 Attention

## Intended Use

Serazym® ANA screen is an IVD test for the qualitative or semiquantitative determination of antinuclear antibodies (ANA) of the IgG isotype against nuclear antigens (dsDNA, RNP, Sm, Ro/SS-A (60kDa), La/SS-B, Scl-70, CENP-B) and Jo-1 in serum or plasma (citrate, EDTA, heparin) of human origin by a laboratory professional user.

It is intended to aid in the diagnosis of a collagenosis in specimen materials from patients with suspicion of a systemic autoimmune disease.

The test must not be used with specimen materials other than serum or plasma (citrate, EDTA, heparin) of human origin, for diagnosis, monitoring, screening, prediction, prognosis, as companion diagnostic, in the near-patient setting and by lay persons.

## Principle of the Test

Serazym® ANA screen is an enzyme immunoassay in which the diluted serum or plasma samples, the ready-to-use calibrator and the negative control react with antigens adsorbed to the solid phase. After a 60-minute incubation at room temperature (RT, 18...25 °C) unbound components are removed by a wash step. In the second incubation step horseradish peroxidase (HRP)- labeled anti-human IgG antibodies bind to the sample antibodies. After a 30-minute incubation at RT unbound conjugate is removed by a wash step. HRP converts the colorless substrate solution to a blue reaction product in the following 15-minute enzymatic reaction step. The reaction is stopped by addition of the stop solution, resulting in a color change from blue to yellow. The optical density (OD) of the reaction product measured at 450 nm measuring and  $\geq 620$  nm reference filter is directly proportional to the concentration of specifically bound antibodies. The cut-off value is calculated by multiplying the measured absorbances of the calibrator by the corresponding factor. The binding index (BI) is calculated by dividing the measured absorbance of the patient sample by the calculated cut-off value.

## Test Components (Delivery Scope)

			<b>For 96 wells</b>
1	<b>WELLS</b>	<b>Microtiter plate</b> coated with nuclear antigens (ANA) and Jo-1	12 single breakable 8-well strips, colorless, vacuum-sealed with desiccant
2	<b>WASHBUF (10x)</b>	<b>Wash buffer (10x)</b> Seramun® Wash buffer K TRIS based buffer	100 mL concentrate for 1000 mL solution, colorless, white cap
3	<b>DIL</b>	<b>Sample diluent</b> Seramun® Sample diluent M TRIS based buffer	100 mL, ready to use, colored red, black cap
4	<b>CAL</b>	<b>Calibrator</b> Factor see Certificate of Analysis	1.0 mL per calibrator ready to use, colored blue, red cap
5	<b>CONTROL -</b>	<b>Positive control</b> Concentration see Certificate of Analysis	1.0 mL, ready to use, colored blue, green cap
6	<b>CONJ HRP</b>	<b>HRP conjugate</b> HRP labeled anti-human IgG antibody (sheep)	15 mL, ready to use, colored red, red cap
7	<b>SUBSTR</b>	<b>Substrate</b> SeramunBlau® slow2 70 < 0.1 % 3,3',5,5'-tetramethylbenzidine	15 mL, ready to use, colorless, blue cap
8	<b>STOP</b>	<b>Stop solution</b> SeramunBlau® stop 0.25 M sulphuric acid	15 mL, ready to use, colorless, yellow cap

9	COVER	Covering film	2 pieces
10		Certificate of Analysis	1 piece
11		Instructions for Use	1 piece

## Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micropipette • 8-channel pipette or multi-pipette with pipette tips • reagent container for multi-channel micro-pipettes • 8-channel wash comb with vacuum pump and waste bottle or microplate washer • microplate reader with 450 nm measuring filter and  $\geq 620$  nm reference filter • deionized water • measuring cylinder • tubes for sample preparation

## Important Information

**This device is for *in-vitro* diagnostic use only.** Follow the instructions carefully. The kit may be performed by laboratory professionals only.

Do not use reagents from damaged packages or bottles. The shelf life specified must be observed. Do not mix components with reagents from other manufacturers.

**Mixing of test kit components of different lots is only allowed for wash buffer (10x), substrate and stop solution.**

All serious incidents occurring in relation with Serazym<sup>®</sup> ANA screen must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

### Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all test components to room temperature before use. For larger sample series, pipetting reagents from liquid reservoirs using a multichannel pipette is recommended to avoid time delays. Follow the pipetting scheme and time schedules of the protocol. The aspiration and washing steps can be performed manually or with the help of a microplate washer or waterjet pump. Wash solution should be allowed a minimum reaction time of 5s in the wells per wash cycle. Remove wash buffer residues by thoroughly aspirating or tapping out the cavities! Protect substrate from light! The test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

### Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
	EUH210	Safety data sheet available on request.
DIL	-	Contains material of animal origin.
CAL	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
	EUH210	Safety data sheet available on request.
CONTROL -	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
	EUH210	Safety data sheet available on request.
CONJ HRP	Hazard components H360D P202	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone  May damage the unborn child. Do not handle until all safety precautions have been read and understood.
	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P308+P313 P501	If exposed or concerned: Get medical advice/attention. Dispose of contents/container in accordance with local regulations.
	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
		Restricted to professional users. Contains material of animal origin.
SUBSTR	EUH210	Safety data sheet available on request
STOP	Hazard components H290	Sulphuric acid 2.5 %  May be corrosive to metals.

### Limitations of the Procedure

Contamination of reagents and samples by bacteria or fungi as well as cross-contamination of the test kit reagents and samples can lead to incorrect results. Incorrect washing to separate unbound components from the sample and test reagents as well as incorrect incubation times of samples, negative control and calibrator can cause incorrect results as well.

The overall interpretation of an ELISA test result should consider all clinical findings.

## Sample Treatment

### Sample Collection

Collect serum or plasma (citrate, EDTA, heparin) of human origin in suitable sampling container.

### Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 2 days. For longer periods samples have to be stored at -20 °C. Repeated (> 3x) freezing and thawing of samples should be avoided.

### Sample Preparation

Bring samples to room temperature before use. Homogeneity should be ensured by shaking components briefly. Samples have to be diluted 1 : 101 (v/v) with sample diluent.

Example: 10 µL sample and 1000 µL sample diluent.

## Reagent Treatment

### Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips can be stored at 2...8 °C until the printed expiration date. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

### Reagent Preparation

The microtiter plate with breakable 8-well strips is vacuum sealed with desiccant. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored at 2...8 °C and protected from moisture in the original cover carefully resealed.

Dilute wash buffer (10x) 1 : 10 with deionized water.

Example: 10 mL wash buffer (10x) + 90 mL deionized water.

Except for the wash buffer all components included in the test kit are ready to use.

The substrate must be protected from direct light. The test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

## Assay Procedure

1. Allow test reagents and required number of wells to reach room temperature (RT). Shake reagents gently before use. Avoid foaming.
2. Pipette 100 µL **CAL** Calibrator  
100 µL **CONTROL -** Negative control  
100 µL diluted sample each.
3. Cover the plate and incubate for 60 min at RT
4. Decant, then wash each cavity 3 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
5. 100 µL **CONJ HRP** HRP conjugate per cavity.
6. Cover the plate and incubate for 30 min at RT.
7. Decant, then wash each cavity 3 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
8. 100 µL **SUBSTR** substrate per cavity.
9. Incubate for 15 min at RT **protected from light**.
10. 100 µL **STOP** stop solution per cavity, mix gently.
11. Read OD at 450 nm and ≥ 620 nm with a microplate reader within 30 min following reaction stop.

## Evaluation of Results

The evaluation of the results is performed qualitatively by calculating a cut-off value or semi-quantitatively by determining a binding index (BI) for each sample based on the cut-off value using the following formula:

$$Absorbance_{cut-off} = absorbance_{calibrator} \times factor$$

The factor is lot-specific and shown on the Certificate of Analysis.

The binding index (ratio) is calculated according to the following formula:

$$BI = \frac{absorbance_{sample}}{absorbance_{cut-off}}$$

## Interpretation of Results

Serazym® ANA screen	Binding index
positive	≥ 1.0
negative	< 1.0

### Test validation

The test run is valid if

- absorbance of calibrator CAL is ≥ 0.80
- absorbance of negative control is ≤ 0.35

If the above-mentioned quality criteria are not met, the test should be repeated strictly following the test procedure (reagent preparation, incubation times and temperatures, wash steps, etc.). In case of repeated failure of the quality criteria contact the manufacturer.

## Performance Characteristics

For lack of international reference preparations results are interpreted by means of the binding index (BI).

### Precision

To determine precision 4 samples were measured multiple times. For the determination of the intra-assay coefficient of variation (CV) samples were measured in a 6-fold determination in one test run. The determination of the inter-assay coefficient of variation was done by a 2-fold determination in 16 different test runs:

Sample	Intra-assay coefficient of variation		Inter-assay coefficient of variation	
	BI	CV [%]	BI	CV [%]
1	4.87	3.26	4.49	5.21
2	3.76	2.01	3.66	4.07
3	2.00	1.80	2.22	3.48
4	1.31	2.71	1.74	5.16

## Change History

<b>Version</b>	<b>Section</b>	<b>Modifications</b>
2023-03_v01_de_en	Entire document	Updating of the Intended Use Conversion of subsections Insertion of Safety Instructions Updating Important Information

## References

1. Conrad, K, Schößler, W et al. (2015): "Autoantibodies in systemic autoimmune diseases: A diagnostic reference", Pabst Science Publishers, Vol. 3, Lengerich, Germany.
2. Peter, JB, Shoenfeld Y (1996): "Autoantibodies". Elsevier, Vol. 1, Amsterdam, The Netherlands.
3. Tan, EM (1989): "Antinuclear antibodies: diagnostic markers for autoimmune disease and probes for cell biology", Adv. Immunol., Vol. 44, p. 93-151.