Serazym[®] Nucleo ANAscreen plus

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Serazym[®] Nucleo ANAscreen plus

Enzyme immunoassay for the quantitative detection of autoantibodies of the IgG-, IgA- and IgM-isotype against nuclear antigens (antinuclear antibodies – ANA) in serum or plasma of human origin

REF	E-100		96
IVD	In-vitro-diagnostic medical device		

 **Seramun Diagnostica GmbH** • Sprenhagener Str. 1 • 15754 Heidesee • Germany •
T +49 33767 791-10 • info@seramun.com • www.seramun.com

IVD In-vitro diagnostic medical device	UDI Unique device identifier	 Manufacturer
 Country of manufacture and date of manufacture	REF Article number	SN Serial number
 Keep away from sunlight	 Humidity limitation	LOT Batch code
 Consult instructions for use	 Temperature limit	 Do not reuse
 Sufficient for <i>n</i> tests	 Biohazard	 Use-by date
		 Attention

Intended Use

Serazym® Nucleo ANAScreen plus is an IVD test for the quantitative determination of autoantibodies of the IgG-, IgA- and IgM-isotype against nuclear antigens (antinuclear antibodies – ANA) in serum or plasma (citrate, EDTA) of human origin by a laboratory professional user.

It is intended to aid in the diagnosis of a collagenosis in specimen materials from patients with suspicion of a systemic autoimmune disease.

The test must not be used with specimen materials other than serum or plasma (citrate, EDTA) of human origin, for diagnosis, monitoring, screening, prediction, prognosis, as companion diagnostic, in the near-patient setting and by lay persons.

Principle of the Test

Serazym® Nucleo ANAScreen plus is an enzyme immunoassay in which the diluted serum or plasma samples, the ready-to-use calibrators and the positive control react with the cell based and selective enriched antigens adsorbed to the solid phase. After a 60-minute incubation at room temperature (RT, 18...25 °C) unbound components are removed by a wash step. In the second incubation step horseradish peroxidase (HRP)- labeled anti-human IgG/IgA/IgM antibodies bind to the sample antibodies. After a 30-minute incubation at RT unbound conjugate is removed by a wash step. HRP converts the colorless substrate solution to a blue reaction product in the following 15-minute enzymatic reaction step. The reaction is stopped by addition of the stop solution, resulting in a color change from blue to yellow. The optical density (OD) of the reaction product measured at 450 nm measuring and ≥ 620 nm reference filter is directly proportional to the concentration of specifically bound antibodies. Using the measured absorbances of the calibrators and their corresponding antibody concentrations, a calibration curve is generated from which the concentrations of the patient samples can be read.

Test Components (Delivery Scope)

			For 96 wells
1	WELLS	Microtiter plate coated with nuclear antigens (ANA)	12 single breakable 8-well strips, colorless, vacuum-sealed with desiccant
2	WASHBUF (10x)	Wash buffer (10x) Seramun® Wash buffer G TRIS based buffer	100 mL concentrate for 1000 mL solution, colorless, white cap
3	DIL	Sample diluent Seramun® Sample diluent O TRIS based buffer	100 mL, ready to use, colored orange, black cap
4	CAL 0 - 4	Calibrators 0 – 4 CAL 0 = 1 U/mL CAL 1 = 10 U/mL CAL 2 = 30 U/mL CAL 3 = 100 U/mL CAL 4 = 300 U/mL	1.0 mL per calibrator ready-to-use colored blue white cap
5	CONTROL +	Positive control Concentration see certificate of analysis	1.0 mL, ready to use, colored blue, red cap
6	CONJ HRP	HRP conjugate HRP labeled anti-human IgG, anti-human IgA and anti-human IgM antibodies (sheep)	15 mL, ready to use, colored red, red cap
7	SUBSTR	Substrate SeramunBlau® slow2 70 < 0.1 % 3,3',5,5'-tetramethylbenzidine	15 mL, ready to use, colorless, blue cap

8	STOP	Stop solution SeramunBlau® stop 0.25 M sulphuric acid	15 mL, ready to use, colorless, yellow cap
9	COVER	Covering film	2 pieces
10		Certificate of Analysis	1 piece
11		Instructions for Use	1 piece

Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micropipette • 8-channel pipette or multi-pipette with pipette tips • reagent container for multi-channel micro-pipettes • 8-channel wash comb with vacuum pump and waste bottle or microplate washer • microplate reader with 450 nm measuring filter and ≥ 620 nm reference filter • deionized water • measuring cylinder • tubes for sample preparation

Important Information

This device is for *in-vitro* diagnostic use only. Follow the instructions carefully. The kit may be performed by health professionals only.

Do not use reagents from damaged packages or bottles. The shelf life specified must be observed. Do not mix components with reagents from other manufacturers.

Mixing of test kit components of different lots is only allowed for wash buffer (10x), substrate and stop solution.

All serious incidents occurring in relation with Serazym® Nucleo ANAscreen plus must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all test components to room temperature before use. For larger sample series, pipetting reagents from liquid reservoirs using a multichannel pipette is recommended to avoid time delays. Follow the pipetting scheme and time schedules of the protocol. The aspiration and washing steps can be performed manually or with the help of a microplate washer or waterjet pump. Wash solution should be allowed a minimum reaction time of 5 s in the wells per wash cycle. Remove wash buffer residues by thoroughly aspirating or tapping out the cavities! Protect substrate from light!

Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
	EUH210	Safety data sheet available on request.
DIL	-	Contains material of animal origin.
CAL 0 - 4	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
	EUH210	Safety data sheet available on request.
CONTROL +	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.
	EUH210	Safety data sheet available on request.
CONJ HRP	Hazard components H360D P202	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone May damage the unborn child. Do not handle until all safety precautions have been read and understood.
	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P308+P313 P501	If exposed or concerned: Get medical advice/attention. Dispose of contents/container in accordance with local regulations.
	EUH208	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Restricted to professional users.
		Contains material of animal origin.
SUBSTR	EUH210	Safety data sheet available on request
STOP	Hazard components H290	Sulphuric acid 2.5 % May be corrosive to metals.

Limitations of the Procedure

Contamination of reagents and samples by bacteria or fungi as well as cross-contamination of the test kit reagents and samples can lead to incorrect results. Incorrect washing to separate unbound components from the sample and test reagents as well as incorrect incubation times of the samples and calibrators can cause incorrect results as well.

The overall interpretation of an ELISA test result should consider all clinical findings.

Sample Treatment

Sample Collection

Collect serum or plasma (citrate, EDTA) of human origin in suitable sampling container.

Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 2 days. For longer periods samples have to be stored at -20 °C. Repeated (> 3x) freezing and thawing of samples should be avoided.

Sample Preparation

Bring samples to room temperature before use. Homogeneity should be ensured by shaking components briefly. Samples have to be diluted 1 : 101 (v/v) with sample diluent.

Example: 10 µL sample and 1000 µL sample diluent.

Reagent Treatment

Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips can be stored at 2...8 °C until the printed expiration date. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

Reagent Preparation

The microtiter plate with breakable 8-well strips is vacuum sealed with desiccant. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored at 2...8 °C and protected from moisture in the original cover carefully resealed.

Dilute wash buffer (10x) 1 : 10 with deionized water.

Example: 10 mL wash buffer (10x) + 90 mL deionized water.

Except for the wash buffer all components included in the test kit are ready to use.

The substrate must be protected from direct light. The test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

Assay Procedure

1. Allow test reagents and required number of wells to reach room temperature (RT). Shake reagents gently before use. Avoid foaming.
2. Pipette 100 µL

CAL	0 - 4
-----	-------

 Calibrators 0, 1, 2, 3 and 4
100 µL

CONTROL	+
---------	---

 Positive control
100 µL diluted sample each.
3. Cover the plate and incubate for 60 min at RT
4. Decant, then wash each well 3 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
5. 100 µL

CONJ HRP

 HRP conjugate per cavity.
6. Cover the plate and incubate for 30 min at RT.
7. Decant, then wash each well 3 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
8. 100 µL

SUBSTR

 substrate per cavity.
9. Incubate for 15 min at RT **protected from light**.
10. 100 µL

STOP

 stop solution per cavity, mix gently.
11. Read OD at 450 nm and ≥ 620 nm with a microplate reader within 30 min following reaction stop.

Evaluation of Results

A calibration curve is created by plotting the measured absorbances of the calibrators CAL 0 - 4 (y-axis) against the corresponding antibody concentrations (x-axis).

The 4-parameter regression model is recommended to create the calibration curve.

The absorbances of the diluted samples are converted to antibody concentrations U/mL by extrapolating from the calibration curve. The indicated concentrations of the calibrators already take into account the regular dilution factor of the serum samples of 1 : 101.

Measurements of samples with a higher absorbance than that of calibrator CAL 4 have to be repeated in a higher predilution. Results have to be multiplied with the factor of the additional dilution:

Sample volume	Added volume of sample diluent	Final volume	Dilution of original sample	Multiplication factor
10 µL	1000 µL	1010 µL	1 : 101 = Sample 1	-
500 µL Sample 1	500 µL	1000 µL	1 : 202 = Sample 2	2
500 µL Sample 2	500 µL	1000 µL	1 : 404 = Sample 3	4

Note that antibodies in dilutions may not be stable, which may lead to a loss of activity. For this reason, a reinvestigation should be performed in a higher dilution at the same day.

Interpretation of Results

Positive	> 10 U/mL
Gray zone	8 – 10 U/mL
Negative	< 8 U/mL

Test should be repeated if read-out lays within the gray zone.

Test validation

The test run is valid if

- adsorbance of calibrator CAL 0 < CAL 1
- adsorbance of calibrator CAL 1 ≤ 0.70
- adsorbance of calibrator CAL 4 ≥ 1.20

If the above-mentioned quality criteria are not met, the test should be repeated strictly following the test procedure (reagent preparation, incubation times and temperatures, wash steps, etc.). In case of repeated failure of the quality criteria contact the manufacturer.

Performance Characteristics

Sensitivity and Specificity

Sensitivity

A total of 62 serum samples from patients were tested to determine the sensitivity of Serazym® Nucleo ANAscreen plus compared to another commercially available ELISA. Each patient of this collective presented at least 4 ACR criteria for rheumatic diseases. This group of patients included 28 cases of systemic lupus erythematosus (SLE), 16 cases of primary Sjögren syndrome (PSS) and 18 different rheumatic diseases (myositis, CREST syndrome, MCTD).

The associated data are shown in the following table. Patient samples determined to be ranged in the grey zone were considered to be positive.

Diagnosis	n	Comparative ELISA		Serazym® Nucleo ANAscreen plus	
		positive	sensitivity	positive	sensitivity
SLE	28	27	96.4 %	27	96.4 %
PSS	16	16	100 %	16	100 %
Others	18	17	94.4 %	17	94.4 %
Total	62	60	96.8 %	60	96.8 %

Specificity

A total of 62 serum samples from blood donors were tested in Serazym® Nucleo ANAscreen plus in comparison to another commercially available ELISA.

	n	Comparative ELISA		Serazym® Nucleo ANAscreen plus	
		negative	specificity	negative	specificity
Blood donors	55	48	87.3 %	55	100 %

Control groups

A total of 80 serum samples were tested in Serazym® Nucleo ANAscreen plus. The 8 control groups consisted of 10 serum samples each with antibody titers for the following parameters: EBV, CMV, hypergammaglobulinemia, RF-IgA, RF-IgM, CCP, PR3 and MPO. 10 serum samples from this collective were detected false positive.

Comparative test

A total of 207 serum samples were tested in Serazym® Nucleo ANAscreen plus and in another commercially available ELISA. Results are shown in the following table:

		Comparative ELISA		
		positive	gray zone	negative
Serazym® Nucleo ANAscreen plus	positive	118	4	1
	gray zone	9	1	0
	negative	8	11	55


Precision

To determine precision 4 samples were measured multiple times. For the determination of the intra-assay coefficient of variation (CV) samples were measured in an 8-fold determination in one test run. The determination of the inter-assay coefficient of variation was done by a 2-fold determination in 6 different test runs:

Sample	Intra-assay coefficient of variation		Inter-assay coefficient of variation	
	concentration [U/mL]	CV [%]	concentration [U/mL]	CV [%]
1	210.5	5.3	100.5	10.8
2	121.9	7.1	23.6	9.1
3	48.7	4.1	16.2	10.9
4	4.2	7.5	5.1	4.5

Calibration Curve

A typical calibration curve in Serazym® Nucleo ANAscreen plus is shown below:

Interference

Hemolytic, lipemic and icteric samples do not interfere in a concentration of up to 500 mg/dL (hemoglobin), 1000 mg/dL (lipids) and 20 mg/dL (bilirubin C and bilirubin F). Rheumatoid factors do not interfere in a concentration of up to 500 IU/mL.

Change History

Version	Section	Modifications
2022-11_v01_DE_EN	Entire document	Updating of the Intended Use Conversion of subsections Insertion of Safety Instructions Updating Important Information

References

1. Conrad, K, Schößler, W et al. (2015): "Autoantibodies in systemic autoimmune diseases: A diagnostic reference", Pabst Science Publishers, Vol. 3, Lengerich, Germany.
2. Peter, JB, Shoenfeld Y (1996): "Autoantibodies". Elsevier, Vol. 1, Amsterdam, The Netherlands.
3. Tan, EM (1989): "Antinuclear antibodies: diagnostic markers for autoimmune disease and probes for cell biology", Adv. Immunol., Vol. 44, p. 93-151.