## SeraSpot® Vaskulitis-3 IgG

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## Literatur:

1. Jenette JC, Falk RJ, Andrassy K. et al, Nomenclature of systemic vasculitides: proposal of an international conference. *Arthritis Rheumatol* 1994, 37, 187-192.
2. Charles LA, Caldas MLR, Falk RJ et al, Antibodies against granule proteins activate neutrophils in vitro. *J Leukocyte Biol* 1991, 50, 539-546.
3. Churg J, Churg A, Idiopathic and secondary vasculitis: a review. *Mod. Pathol.* 1989, 2, 144-160.
4. Cohen Tervaert JW, van der Waude FJ, Fauci, AS et al, Association between active Wegener's granulomatosis and anticytoplasmic antibodies. *Arch Intern Med* 1989, 149, 2461-2465.
5. Weiss SJ, Tissue destruction by neutrophils. *New Engl. J Med* 1989, 320, 365-376.
6. Ozaki S, ANCA associated vasculitis: Diagnostic and therapeutic strategy. *Allergol. Int.* 2007, 56, 87-96.
7. Chen M, Yu F, Wang SX et al, Renal histology in Chinese patients with anti-myeloperoxidase autoantibody-positive Wegener's granulomatosis. *Nephrol. Dial. Transplant.* 2007, 22, 139-145.
8. Papiris SA, Manali ED, Kalomenidis I et al, Bench-to-bedside review: Pulmonary-renal syndromes- an update for the intensivist. *Critical Care* 2007, 213.
9. Salama AD, Levy J, Lightstone L et al, Goodpasture's disease. *Lancet* 2001, 358, 917-920.
10. Hudson BG, Tryggvason K, Sundaramoorthy M et al, Alport's syndrome, Goodpasture's syndrome, and type IV collagen. *N Engl J Med*, 2003, 348, 2543- 2556.

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## SeraSpot® Vaskulitis-3 IgG

Spotimmunoassay for detection of IgG antibodies in in systemic vasculitis  
in human serum or plasma

 SP-003-3 G-S6  48  SP-003-3 G-S12  96  SP-003-3 G-S24  2x 96  
 *In-vitro*- diagnostic device 

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · [info@seramun.com](mailto:info@seramun.com)

### Introduction

Primary systemic vasculitis of small and middle large vessels (1) of unknown etiopathogenesis is associated with anti-neutrophil cytoplasmic antibodies (ANCA). Because ANCA seem to activate directly neutrophil granulocytes which are accumulated in the lesions of small vessels and thereby initiate an inflammatory destructive process no immunoglobulins or complement components are seen in the vascular lesions (so-called “pauci-immune” vasculitis) (2, 3).

Major antigens of ANCA are:

**Proteinase 3 (PR3)** (synonym: azurophilic granule protein) a member of the serine proteases family, is the main target of cytoplasmic ANCA (cANCA) seen by immunofluorescence and is strongly associated with Wegener’s granulomatosis (85 % anti-PR3 antibody positive) (4).

**Myeloperoxidase (MPO)** is a strongly cationic protein of the azurophilic granules and is involved in the generation of reactive oxygen species (5). It is the main target of perinuclear ANCA (pANCA) which are serological markers for microscopic polyangiitis (45 %), Churg-Strauss syndrome (60 %) and idiopathic glomerulonephritis (65 %). MPO-ANCA associated vasculitis is more common in Japanese and Chinese and PR3-ANCA associated vasculitis more common in Europeans (6). In PR3-ANCA associated vasculitis glomerulonephritis is more common and more severe than in patients with MPO-associated vasculitis (7). Although MPO is the main target of pANCA a similar immuno-fluorescence picture can be produced by antibodies to elastase, cathepsin and lysozyme. Therefore methods (ELISA or LIA) with purified antigens as targets should be used to determine fine specificity of antibodies classified by immunofluorescence.

**Glomerular basement membrane (GBM)** is the target of autoantibodies which bind to the non-collagenous domain of the  $\alpha 3$  chain of collagen type IV (10). This structure is also present in other membranes. Anti-GBM antibodies detectable in about 90 % of patients activate complement after binding. Cellular proteases are attracted and cause membrane disruption and destruction of the filtration barrier. Due to knowledge of the target antigen highly purified or recombinant antigen in ELISAs or LIAs enables a rapid diagnosis and an early onset of adequate therapy which dramatically improves prognosis. The good pasture's syndrome is extremely rare (one case per 1.000.000 inhabitants per year, 3.8 % of patients with pulmonary-renal syndrome in intensive care units) (8) and predominantly affects Caucasians. However, due to rapidly progressive glomerulonephritis it is responsible for about 20 % of acute renal failure cases (9). Untreated, mortality is between 75 and 90 %. Genetic (> 80 % of patients are positive for HLA DR15 or DR4) and environmental factors (smoking, infections) have been implicated in the pathogenesis. Sometimes the good pasture's syndrome can be superimposed to ANCA positive systemic vasculitis.

#### **Literature:**

1. Jenette JC, Falk RJ, Andrassy K. et al, Nomenclature of systemic vasculitides: proposal of an international conference. *Arthritis Rheumatol* 1994, 37, 187-192.
2. Charles LA, Caldas MLR, Falk RJ et al, Antibodies against granule proteins activate neutrophils in vitro. *J Leukocyte Biol* 1991, 50, 539-546.
3. Churg J, Churg A, Idiopathic and secondary vasculitis: a review. *Mod. Pathol.* 1989, 2, 144-160.
4. Cohen Tervaert JW, van der Waude FJ, Fauci, AS et al, Association between active Wegener's granulomatosis and anticytoplasmic antibodies. *Arch Intern Med* 1989, 149, 2461-2465.
5. Weiss SJ, Tissue destruction by neutrophils. *New Engl. J Med* 1989, 320, 365-376.
6. Ozaki S, ANCA associated vasculitis: Diagnostic and therapeutic strategy. *Allergol. Int.* 2007, 56, 87-96.
7. Chen M, Yu F, Wang SX et al, Renal histology in Chinese patients with anti-myeloperoxidase autoantibody-positive Wegener's granulomatosis. *Nephrol. Dial. Transplant.* 2007, 22, 139-145.
8. Papisir SA, Manali ED, Kalomenidis I et al, Bench-to-bedside review: Pulmonary-renal syndromes- an update for the intensivist. *Critical Care* 2007, 213.
9. Salama AD, Levy J, Lightstone L et al, Goodpasture's disease. *Lancet* 2001, 358, 917-920.
10. Hudson BG, Tryggvason K, Sundaramoorthy M et al, Alport's syndrome, Goodpasture's syndrome, and type IV collagen. *N Engl J Med*, 2003, 348, 2543- 2556.

## Intended use

**SeraSpot® Vaskulitis-3 IgG is an *in vitro* diagnostic device (spotimmunoassay, SIA) for the detection of autoantibodies of the IgG isotype directed against the following antigens: Proteinase 3 (PR3), myeloperoxidase (MPO) and glomerular basement membrane (GBM) in human serum or plasma samples.**

## Principle of the test

SeraSpot® Vaskulitis-3 IgG test is a solid phase immunoassay based on the use of recombinant or purified native proteins as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96well-microtitration plates. The antigens serve as capture molecules for autoantibodies against proteinase 3, myeloperoxidase and glomerular basement membrane. Bound antibodies are detected by horseradish peroxidase-(HRP)-labeled antibodies against human antibodies of IgG-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed as precipitated product from colorless substrate solution. Color intensity is correlated to the antibody concentration. Pale blue to dark blue spots are visible by eye.

**The test procedure for antibody detection is performed in 3 steps:**

### Step 1

Incubation of 1 : 101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are sucked off and unbound components are removed by 3 wash cycles with wash solution.

### Step 2

Incubation of wells with the anti-human IgG-HRP-conjugate for 30 minutes at room temperature. After incubation the antibody is sucked off and unbound components are removed by 3 wash cycles with wash solution.

### Step 3

Incubation of wells with substrate solution SeramunBlau® spot dark for 30 minutes at room temperature. The reaction is stopped by sucking off the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the Seramun SpotSight® plate or Seramun SpotSight® strip scanner. The obtained images are interpreted by using the Software Seramun SpotSight® scan. Taking account of the test-specific evaluation criteria, the test result is determined.

## Used Autoantigens

Nomenclature	Clinical relevance / Specificity		Description
PR3	Wegener's granulomatosis	95 %	Proteinase 3
MPO	Mikroskopische polyangiitis (MPA); Churg-Strauss-syndrom (CSS); Idiopathic membranous glomerulonephritis	45 % 60 % 65 %	Myeloperoxidase
GBM	Goodpasture-Syndrom	90 %	"Glomerular basement membrane antigen" (GBM), $\alpha$ 3 chain of human collagen type IV

## Kit components

		For 48 determinations	For 96 determinations	For 2x 96 determinations
1	<b>WELLS</b> Wells with arrays	<b>6 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>12 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>2x 12- single breakable 8-well strips in frame</b> vacuum-sealed with desiccant
<b>Color coding: red</b>				
2	<b>WASHBUF</b> <b>CONC</b> <b>10X</b> Wash buffer	<b>100 ml concentrate</b> transparent bottle, white cup	<b>100 ml concentrate</b> transparent bottle, white cup	<b>2x 100 ml concentrate</b> transparent bottle, white cup
3	<b>DIL</b> Sample diluent	<b>55 ml ready-to-use solution</b> , colored red transparent bottle, black cup	<b>2x 55 ml ready-to-use solution</b> , colored red, transparent bottle, black cup	<b>4x 55 ml ready-to-use solution</b> , colored red, transparent bottle, black cup
4	<b>CONJ</b> <b>HRP</b> <b>IgG</b> Anti-human IgG-HRP-conjugate	<b>8 ml ready-to-use solution</b> , colored red, transparent bottle, red cup	<b>8 ml ready-to-use solution</b> , colored red, transparent bottle, red cup	<b>2x 8 ml ready-to-use solution</b> , colored red, transparent bottle, red cup
5	<b>SUBSTR</b> <b>TMB</b> Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	<b>8 ml ready-to-use solution</b> , black bottle, blue cup	<b>8 ml ready-to-use solution</b> , black bottle, blue cup	<b>2x 8 ml ready-to-use solution</b> , black bottle, blue cup
6	<b>COVER</b> Adhesive film	<b>2x</b>	<b>2x</b>	<b>4x</b>
7	<b>SWAB</b> Swab with 70 % (v/v) Isopropyl alcohol	<b>3x 2</b>	<b>6x 2</b>	<b>12x 2</b>

## Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 48 hours. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided!

**Serum or plasma** (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (10 µl sample and 1000 µl buffer) with the ready to use sample diluent.

## Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- stop watch
- Washer for 96-well microtiter plates (optional)
- Collecting devices for infectious material
- Seramun *SpotSight*<sup>®</sup> plate or Seramun *SpotSight*<sup>®</sup> strip scanner with evaluation software Seramun *SpotSight*<sup>®</sup> scan

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready to use wash solution can be used for at least 1 month when stored at 2...8 °C.

## Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtiter plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or de-ionized water. For Example: 10 ml **WASHBUF 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

## Workplace requirements

**Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.**

## Assay procedure

Execution at room temperature (18 ... 25 °C) and at the specified incubation times.

For the handling of SeraSpot® tests the use of the starter kit (order number: SP-000-1) is recommended. The kit is designed for multiple use.

### Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere to the underside of the bottom of the wells using the swab **SWAB** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

### Working steps


1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Add **100 µl** of 1 : 101 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 1000 µl sample dilution buffer.
3. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **three times** using **400 µl** wash solution per well (made of (2)). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **CONJ** **HRP** **IgG** to each well.
6. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard samples), then wash wells **three times** using **400 µl** wash solution per well (made of (2)). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **SUBSTR** **TMB** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent.
11. Clean bottom of wells with Swab **SWAB**.
12. Perform image acquisition using the Seramun *SpotSight*® plate or the Seramun *SpotSight*® strip scanner and evaluate the results with the Seramun *SpotSight*® scan software. If image acquisition is performed with the scanner Seramun *SpotSight*® strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

**Aspiration and washing steps can be performed manually or via washer for 96-well microtiter plates.**

**After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when stored in the dark.**

## Test evaluation

### Array layout

Parameter
05 GBM
07 MPO
09 PR3

Controls
01 Positive control (PC)
02 Negative control (NC), 0 rel. units
03 Cut-off control (CO), 30 rel. units
04 Reference 3 (R3), 60 rel. units
06 Reference 2 (R2), 100 rel. units
08 Reference 1 (R1), 300 rel. units
10 Serum control (SC)
P Well position marker

### Validity criteria for the test

The *SeraSpot*<sup>®</sup> Vaskulitis-3 IgG test includes the following control spots:

1. Positive control (PC). Intensively stained spot, darker than cut-off control. Always stained.
2. Cut-off control (CO), corresponding to 30 rel. units. Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC), corresponding to 0 rel. units. Pale spot with intensity lower than cut-off control.
4. Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.
5. Reference 1 (R1), corresponding to 300 rel. units. Intensively stained spot. Always stained.
6. Reference 2 (R2), corresponding to 100 rel. units. Intensity lower than R1. Always stained.
7. Reference 3 (R3), corresponding to 60 rel. units. Intensity lower than R2. Always stained.

The spots NC, CO, R1 ... R3 are used to create a reference curve (four parameter non-linear regression) to calculate the staining intensities of the spots in relative units (rel. units).

The test cannot be evaluated if one of the validity criteria listed in points 1 to 7 is not fulfilled.

### Result interpretation

Evaluation of the test has to be performed by using the Seramun SpotSight<sup>®</sup> scan software. A visual plausibility check is possible with the aid of the magnifying lens contained in the starter kit (order number: SP-000-1).

Judgment of spot staining is performed in accordance to the following classification:

Result Interpretation	IgG
<b>negative</b>	color intensity of antigen spots $\leq$ Cut-off control
<b>positive</b>	color intensity of one or more antigen spots $>$ Cut-off control

## **Limitations of the procedure**

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

### **Interference**

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) do not interfere with the test. Rheumatoid factors also do not interfere with the test.

### **Assay procedure**

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere underside the bottom of the wells or to the array area cause erroneous results.

## Performance characteristics

### Diagnostic sensitivity

Serologically predetermined samples (n = 109, reference test 1) of patients with suspected inflammatory vascular disease were tested in the *SeraSpot*<sup>®</sup> Vaskulitis-3 IgG test ("Initial sensitivity"). Samples with discrepant results were retested in a second assay (reference test 2) for detection of autoantibodies in systemic vasculitis. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results ("Amended sensitivity").

GBM	n = 105	Reference test 1	
		positive (n = 13)	negative (n = 92)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	10	1
	negative	3	91

Initial sensitivity: **76.9 %**

Amended sensitivity: **100.0 %**

MPO	n = 105	Reference test 1	
		positive (n = 17)	negative (n = 88)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	17	2
	negative	0	86

Sensitivity: **100.0 %**

PR3	n = 105	Reference test 1	
		positive (n = 75)	negative (n = 30)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	74	3
	negative	1	27

Initial sensitivity: **98.7 %**

Amended sensitivity: **100.0 %**

### Diagnostic specificity

A total of 189 blood donor samples were tested by the *SeraSpot*<sup>®</sup> Vaskulitis-3 IgG test ("Initial specificity"). Samples with positive results were retested in a second assay for detection of autoantibodies in systemic vasculitis. The results obtained were used for determination of the final status of the samples ("Amended specificity").

GBM	n = 189	positive (n = 0)	negative (n = 189)
		<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive
negative	0		186
negative	0		186

Initial specificity: **98.4 %**

Amended specificity: **99.5 %**

MPO	n = 189	positive (n = 0)	negative (n = 189)
		<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive
negative	0		188

Initial specificity: **99.5 %**

Amended specificity: **100.0 %**

PR3	n = 189	positive (n = 0)	negative (n = 189)
		<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive
negative	0		179

Initial specificity: **94.7 %**

Amended specificity: **94.7 %**

### Cross-reactivity

Potential cross-reacting samples were tested in the *SeraSpot*<sup>®</sup> Vaskulitis-3 IgG test ("Initial specificity") and followed up with a reference test ("Amended specificity").

### CCP-positive samplings (antibodies against cyclic citrullinated peptide)

GBM	n = 29	positive (n = 0)	negative (n = 29)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	0	0
	negative	0	29

Specificity: 100.0 %

MPO	n = 29	positive (n = 0)	negative (n = 29)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	0	0
	negative	0	29

Specificity: 100.0 %

PR3	n = 29	positiv (n = 0)	negativ (n = 29)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positiv	0	2
	negativ	0	27

Initial specificity: 93.1 %

Amended specificity: 93.1 %

### ANA-positive samplings (antinuclear antibodies)

GBM	n = 51	positive (n = 0)	negative (n = 51)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	0	3
	negative	0	48

Initial specificity: 94.1 %

Amended specificity: 96.0 %

MPO	n = 51	positive (n = 0)	negative (n = 51)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	0	1
	negative	0	50

Initial specificity: 98.0 %

Amended specificity: 98.0 %

PR3	n = 51	positive (n = 0)	negative (n = 51)
<i>SeraSpot</i> <sup>®</sup> Vaskulitis-3 IgG	positive	0	1
	negative	0	50

Initial specificity: 98.0 %

Amended specificity: 98.0 %

## Precision

Samples of known antibody titer were assayed by SeraSpot® Vaskulitis-3 IgG test. The staining intensities of spots were determined as relative units (rel. units) according to the reference curve and the variation coefficients (CV) were calculated in order to determine the within-run precision (intra-assay-CV), between-run precision (inter-assay-CV) and Lot-to-Lot precision (Interbatch-CV).

	Intra-assay CV		Inter-assay CV		Interbatch CV	
	rel. units [ $\bar{x}$ ] n = 8	CV [%]	rel. units [ $\bar{x}$ ] n = 48	CV [%]	rel. units [ $\bar{x}$ ] n = 144	CV [%]
<b>GBM</b>	91.6	<b>6.8</b>	90.9	<b>6.2</b>	92.8	<b>7.8</b>
<b>MPO</b>	62.5	<b>3.9</b>	65.3	<b>7.9</b>	61.0	<b>17.5</b>
<b>PR3</b>	71.8	<b>6.9</b>	74.8	<b>9.6</b>	71.7	<b>8.7</b>
<b>Procedure</b>	1 operator, 8x determination, 1 batch		1 operator, 8x determinations, 3 days, 1 batch		1 operator, 8x determinations, 3 days, 3 batches	

## Automation

Samples (n = 96) of known antibody titer were assayed manually and (in parallel) by use of the microplate processor DS2® (Dynex Technologies; manual sample dilution). Sample dilution was performed manually for both methods. Color intensity (OD) of spots was recorded and used for calculation of the coefficient of determination  $r^2$  for selected antibody reactivities and the different methods of operation.

	Dynex DS2® vs. manual processing
	$r^2$
<b>GBM</b>	<b>0.997</b>
<b>MPO</b>	<b>0.996</b>
<b>PR3</b>	<b>0.996</b>

## Common advices and precautions

**This kit is for *in-vitro* use only.** Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

**Do not use or mix reagents from different lots except for sample diluent, washing buffer and TMB substrate solution. These can be used across batches and parameters.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

## History of changes

Version	Section	Modifications
2018-03-01	Common advices and precautions	Update

## Incubation scheme *SeraSpot*<sup>®</sup> Vaskulitis-3 IgG

-  100 µl diluted sample (1 : 101)  
 30 min incubation (room temperature)
   
 3 x wash with washing solution, each 400 µl per well
-  50 µl conjugate **CONJ** **HRP** **IgG** (4)  
 30 min incubation (room temperature)
   
 3 x wash with washing solution, each 400 µl per well
-  50 µl substrate **SUBSTR** **TMB** (5)  
 30 min incubation (room temperature, protected from light)
   
 Aspiration
- imaging of wells and image analysis scanner Seramun *SpotSight*<sup>®</sup> plate / Seramun *SpotSight*<sup>®</sup> strip and software Seramun *SpotSight*<sup>®</sup> scan

SeramunBlau<sup>®</sup>, *SeraSpot*<sup>®</sup> und Seramun *SpotSight*<sup>®</sup> are registered trademarks of Seramun Diagnostica GmbH.

Manufacturer

Catalogue-No.

Lot-No.

Number of determinations

Biohazard

Notice advices

Use by

Storage temperature

Consult instructions for use