## **SeraSpot® HepAk-7 IgG**

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## SeraSpot® HepAk-7 IgG

Spotimmunoassay for detection of IgG antibodies in autoimmune liver diseases  
in human serum or plasma

REF SP-004-7 G-S6 ▽ 48 REF SP-004-7 G-S12 ▽ 96 REF SP-004-7 G-S24 ▽ 2x 96  
IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Spreehagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 33767 791-10 · fax +49 33767 791-99 · [info@seramun.com](mailto:info@seramun.com)

### Introduction

Primary autoimmune liver disease (PAL) can be subdivided into autoimmune hepatitis (AIH) type I and II, primary biliary cirrhosis (PBC) and primary sclerosing cholangitis (PSC). Beside these clearly defined diseases, variant forms described as overlap syndrome additionally exist. These forms include characteristics of AIH and PBC or PSC. Clinical symptoms of PAL resemble that of other chronic liver diseases. The progress of chronic liver diseases may result in the development of liver cirrhosis. Only an early diagnosis and treatment may prevent such an outcome. In early stages liver diseases are characterized by increased levels of serum transaminases. Unspecific symptoms as sickness, fatigue and arthralgia can additionally occur. About 15 % of chronic liver disease cases underlie an autoimmune pathogenesis. Therefore the determination of different autoantibodies is recommended after exclusion of a virus hepatitis. The early therapy as a result of a correct diagnosis correlates with a positive prognosis.

The *SeraSpot*® HepAk-7 IgG test enables the fast determination of autoantibodies with diagnostic relevance for autoimmune liver diseases.

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## Intended use

The *SeraSpot*<sup>®</sup> HepAk-7 IgG test is an *in vitro* diagnostic medical device (spotimmunoassay, SIA) for the detection of autoantibodies of the IgG isotype to the following antigens: AMA-M2, LKM1, LC1, SLA, F-Actin, gp210 and Sp100 in human serum or plasma.

## Principle of the test

The *SeraSpot*<sup>®</sup> HepAk-7 IgG test is a solid phase immunoassay based on the use of recombinant or purified native proteins / protein complexes as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96-well microtiter plates. The antigens serve as capture molecules for autoantibodies, which are formed in autoimmune liver diseases. Bound antibodies are detected by horseradish peroxidase-(HRP)-labeled antibodies against human antibodies of IgG-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed as precipitated product from colorless substrate solution. Pale blue to dark blue spots are visible by eye.

**The test procedure for antibody detection is performed in 3 steps:**

### Step 1

Incubation of 1 : 101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are sucked off and unbound components are removed by 3 wash cycles with wash solution.

### Step 2

Incubation of wells with the anti-human IgG-HRP-conjugate for 30 minutes at room temperature. After incubation the antibody is sucked off and unbound components are removed by 3 wash cycles with wash solution.

### Step 3

Incubation of wells with substrate solution SeramunBlau<sup>®</sup> spot dark for 30 minutes at room temperature. The reaction is stopped by sucking off the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the Seramun *SpotSight*<sup>®</sup> plate or Seramun *SpotSight*<sup>®</sup> strip scanner. The obtained images are interpreted by using the Software Seramun *SpotSight*<sup>®</sup> scan. Taking account of the test-specific evaluation criteria, the test result is determined.

## Used Autoantigens

Nomenclature	Clinical relevance / Specificity		Description
AMA-M2	PBC	95 %	M2 protein of the mitochondrial pyruvate dehydrogenase (PDH) complex
LKM1	AIH II	95 – 100 %	"Liver-kidney microsomal antigen 1", cytochrome P450 2D6 (LKM 1 hp)
LC1	AIH II	50 %	Liver cytosol antigen type 1, formiminotransferase cyclodeaminase (62 kDa)
SLA	AIH I		Soluble liver antigen, 50 kDa UGA-suppressor serine tRNA associated protein
F-Actin	AIH I	Less specific	Polymerized filaments of the Actin protein.
gp210	PBC	10 – 25 % Highly specific	210 kDa glycoprotein of nuclear pore complex
Sp100	PBC	30 %	„Speckled pattern 100“, protein of “nuclear dots”

## Kit components

		For 48 determinations	For 96 determinations	For 2x 96 determinations
1	<b>WELLS</b>  Wells with arrays  Color coding: green	<b>6 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>12 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>2x 12 single breakable 8-well strips in frame</b> separately vacuum-sealed with desiccant
2	<b>WASHBUF CONC 10X</b>  Wash buffer Seramun® Wash buffer A	<b>100 ml concentrate</b> transparent bottle, white cap	<b>100 ml concentrate</b> transparent bottle, white cap	<b>2x 100 ml concentrate</b> transparent bottle white cap
3	<b>DIL</b>  Sample diluent Seramun® Sample diluent B	<b>55 ml</b> ready-to-use solution, colored red transparent bottle black cap	<b>2x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cap	<b>4x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cap
4	<b>CONJ HRP IgG</b>  Anti-human IgG- HRP- conjugate	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap	<b>2x 8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap
5	<b>SUBSTR TMB</b>  Substrate SeramunBlau® spot dark 3,3',5,5'- Tetramethylbenzidine and hydrogen peroxide	<b>8 ml</b> ready-to-use solution, black bottle, blue cap	<b>8 ml</b> ready-to-use solution, black bottle, blue cap	<b>2x 8 ml</b> ready-to-use solution, black bottle, blue cap
6	<b>COVER</b>  Adhesive film	<b>2x</b>	<b>2x</b>	<b>4x</b>
7	<b>SWAB</b>  Swab with 70 % (v/v) Isopropyl alcohol	<b>3x2</b>	<b>6x2</b>	<b>12x2</b>

## Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 8 days. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided!

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (e.g. 10 µl sample and 1000 µl buffer) with the ready-to-use sample diluent.

## Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- stopwatch
- Washer for 96-well microtiter plates (optional)
- Collecting devices for infectious material
- Seramun *SpotSight*® plate or Seramun *SpotSight*® strip scanner with evaluation software Seramun *SpotSight*® scan

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used up to 4 weeks when stored at 2...8 °C and up to 2 weeks when stored at room temperature.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtiter plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or deionized water. For Example: 10 ml **WASHBUF CONC 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

### Workplace requirements

**Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.**

## Assay procedure

Execution at room temperature (18...25 °C) and at the specified incubation times.

For the handling of SeraSpot® tests the use of the starter kit (order number: SP-000-1) is recommended. The kit is designed for multiple use.

### Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere to the underside of the bottom of the wells using the swab **SWAB** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

### Working steps


1. Warm all reagents and required wells **WELL** to room temperature (RT) before use. Mix gently without causing foam.
2. Add **100 µl** of 1 : 101 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 1000 µl sample dilution buffer **DIL**.
3. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **CONJ HRP IgG** to each well.
6. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard conjugate), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **SUBSTR TMB** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent paper.
11. Clean bottom of wells with swab **SWAB**.
12. Perform image acquisition using the Seramun *SpotSight*® plate or the Seramun *SpotSight*® strip scanner and evaluate the results with the Seramun *SpotSight*® scan software. If image acquisition is performed with the scanner Seramun *SpotSight*® strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

**Aspiration and washing steps can be performed manually or via washer for 96-well microtiter plates.**

**After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when stored in the dark.**

## Test evaluation

### Array layout

Parameter	
05	LKM1
06	AMA-M2
07	LC1
09	SLA
10	Sp100
12	F-Actin
13	gp210

Controls	
01	Positive control (PC)
02	Negative control (NC), 0 rel. units
03	Cut-off control (CO), 30 rel. units
04	Reference 3 (R3), 60 rel. units
08	Reference 2 (R2), 100 rel. units
11	Reference 1 (R1), 300 rel. units
14	Serum control (SC)

**P** Well position marker

### Validity criteria for the test

The *SeraSpot*<sup>®</sup> HepAk-7 IgG test includes the following control spots:

1. Positive control (PC). Intensively stained spot, darker than cut-off control. Always stained.
2. Cut-off control (CO), corresponding to 30 rel. units. Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC), corresponding to 0 rel. units. Pale spot with intensity lower than cut-off control.
4. Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.
5. Reference 1 (R1), corresponding to 300 rel. units. Intensively stained spot. Always stained.
6. Reference 2 (R2), corresponding to 100 rel. units. Intensity lower than R1. Always stained.
7. Reference 3 (R3), corresponding to 60 rel. units. Intensity lower than R2. Always stained.

The spots NC, CO, R1 ... R3 are used to create a reference curve (four parameter non-linear regression) to calculate the staining intensities of the spots in relative units (rel. units).

The test cannot be evaluated if one of the validity criteria listed in points 1 to 7 is not fulfilled.

### Result interpretation

Evaluation of the test has to be performed by using the Seramun *SpotSight*<sup>®</sup> scan software. A visual plausibility check is possible with the aid of the magnifying lens contained in the starter kit (order number: SP-000-1).

Judgment of spot staining is performed in accordance to the following classification:

Result Interpretation	IgG
<b>negative</b>	color intensity of antigen spot ≤ Cut-off control
<b>positive</b>	color intensity of antigen spot > Cut-off control

When evaluating the *SeraSpot*<sup>®</sup> HepAk-7 IgG (SP-004-7 G) results it should be noted that the S2k guideline on autoimmune liver diseases from 2017 recommends the detection of anti-F-Actin autoantibodies of the IgG type by indirect immunofluorescence (5).

## Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

### Interference

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) do not interfere with the test. Rheumatoid factors also do not interfere with the test.

### Assay procedure

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere underside the bottom of the wells or to the array area cause erroneous results.

## Performance characteristics

### Diagnostic sensitivity

Serologically predetermined samples (reference test 1) of patients with suspected inflammatory vascular disease were tested in the *SeraSpot*<sup>®</sup> HepAk-7 IgG test ("Initial sensitivity"). Samples with discrepant results were retested in a second assay (reference test 2) for detection of autoantibodies in systemic vasculitis. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results ("Amended sensitivity").

LKM1	n = 8	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	8	0
	negative	0	0

Initial sensitivity: 100 %

Amended sensitivity: 100 %

AMA-M2	n = 29	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	29	0
	negative	0	0

Initial sensitivity: 100 %

Amended sensitivity: 100 %

LC1	n = 7	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	7	0
	negative	0	0

Initial sensitivity: 100 %

Amended sensitivity: 100 %

SLA	n = 7	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	7	0
	negative	0	0

Initial sensitivity: 100 %

Amended sensitivity: 100 %

Sp100	n = 18	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	18	0
	negative	0	0

Initial sensitivity: 100 %

Amended sensitivity: 100 %

F-Actin	n = 13	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	12	0
	negative	1	0

Initial sensitivity: 92.3 %

Amended sensitivity: 92.3 %

gp210	n = 19	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	18	0
	negative	1	0

Initial sensitivity: 94.7 %

Amended sensitivity: 100 %

## Diagnostic specificity

A collective of blood donor samples were tested by the *SeraSpot*<sup>®</sup> HepAk-7 IgG test ("Initial specificity"). Samples with positive results were retested in a second assay for detection of hepatitis autoantibodies. The results obtained were used for determination of the final status of the samples ("Amended specificity").

LKM1	n = 192	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	1
	negative	0	191

Initial specificity: 99.5 %

Amended specificity: 99.5 %

AMA-M2	n = 192	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	1
	negative	0	191

Initial specificity: 99.5 %

Amended specificity: 99.5 %

LC1	n = 192	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	1	2
	negative	0	189

Initial specificity: 99.0 %

Amended specificity: 99.0 %

SLA	n = 192	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	0
	negative	0	192

Initial specificity: 100 %

Amended specificity: 100 %

Sp100	n = 191	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	0
	negative	1	190

Initial specificity: 100 %

Amended specificity: 100 %

F-Actin	n = 182	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	0
	negative	0	182

Initial specificity: 100 %

Amended specificity: 100 %

gp210	n = 191	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	1
	negative	0	190

Initial specificity: 99.5 %

Amended specificity: 99.5 %

## Cross-reactivity

Potential cross-reacting samples were tested in the *SeraSpot*<sup>®</sup> HepAk-7 IgG test ("Initial specificity") and followed up with a reference test ("Amended specificity").

### HBV-positive samplings (hepatitis B virus)

HBV	n = 27	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	0
	negative	0	27

Initial specificity: 100 %

Amended specificity: 100 %

### HCV-positive samplings (hepatitis C virus)

HCV	n = 43	reference test 1	
		positive	negative
<i>SeraSpot</i> <sup>®</sup> HepAk-7 IgG	positive	0	0
	negative	0	43

Initial specificity 100 %

Amended specificity: 100 %

## Precision

Samples of known antibody titer were assayed by SeraSpot® HepAk-7 IgG test. The staining intensities of spots were determined as relative units (rel. units) according to the reference curve and the variation coefficients (CV) were calculated in order to determine the within-run precision (intra-assay-CV), between-run precision (inter-assay-CV) and Lot-to-Lot precision (Interbatch-CV).

	Intra-assay CV		Inter-assay CV		Interbatch CV	
	rel. units $[\bar{x}]$ n = 40	CV [%]	rel. units $[\bar{x}]$ n = 80	CV [%]	rel. units $[\bar{x}]$ n = 240	CV [%]
<b>LKM1</b>	58.2	<b>5.3</b>	84.8	<b>12.3</b>	86.2	<b>12.6</b>
<b>AMA-M2</b>	36.8	<b>6.8</b>	38.9	<b>14.9</b>	42.2	<b>16.3</b>
<b>LC1</b>	61.5	<b>10.1</b>	60.2	<b>16.7</b>	62.5	<b>18.1</b>
<b>SLA</b>	66.3	<b>6.2</b>	64.9	<b>8.6</b>	77.7	<b>18.1</b>
<b>Sp100</b>	88.6	<b>6.2</b>	84.7	<b>11.4</b>	92.4	<b>18.7</b>
<b>F-Actin</b>	53.8	<b>14.7</b>	59.0	<b>15.3</b>	51.2	<b>19.9</b>
<b>gp210</b>	52.9	<b>6.4</b>	51.2	<b>8.8</b>	49.9	<b>13.8</b>
<b>Procedure</b>	1 operator, 40x determinations, 1 batch		5 operators, 2x determinations, 2x testings per day, 20 days, 1 batch		5 operators, 2x determinations, 2x testings per day, 20 days, 3 batches	

## Automation

Samples (n = 96) of known antibody titer were assayed manually or by use of the microplate processor DS2® (Dynex Technologies; manual sample dilution). Color intensity of antigen spots was recorded and used for calculation of the coefficient of determination  $r^2$  for selected antibody reactivities and the different methods of operation.

	Dynex DS2® vs. manual processing
	$r^2$
<b>LKM1</b>	<b>0.984</b>
<b>AMA-M2</b>	<b>0.988</b>
<b>LC1</b>	<b>0.971</b>
<b>SLA</b>	<b>0.990</b>
<b>Sp100</b>	<b>0.954</b>
<b>F-Actin</b>	
<b>gp210</b>	<b>0.994</b>

$r^2$  was not determined if more than 90 % of the total number (n = 96) of the samples showed no antibody reactivity for the particular antigens.

## Common advices and precautions

**This kit is for *in-vitro* use only.** Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

**Do not use or mix reagents from different lots except for sample diluent, wash buffer and TMB substrate solution. These can be used across batches and parameters.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

## History of changes

Version	Section	Modifications
2020-01-20	General	General revision
	Introduction	Update
	Preparation and storage of samples	Stability of samples
	Preparation and storage of reagents	Storage of ready-to-use wash solution
	Kit components	Update
	Result interpretation	detection of anti-F-Actin autoantibodies
	Performance characteristics	Update of data

## Incubation scheme *SeraSpot*<sup>®</sup> HepAk-7 IgG

-  100 µl  
30 min  
 3 x wash  
diluted sample (1 : 101)  
incubation (room temperature)  
with wash solution, each 400 µl per well
-  50 µl  
30 min  
 3 x wash  
conjugate **CONJ HRP IgG**  
incubation (room temperature)  
with wash solution, each 400 µl per well
-  50 µl  
30 min  
 aspiration  
substrate **SUBSTR TMB**  
incubation (room temperature, protected from light)
4. imaging of wells and image analysis  
scanner Seramun *SpotSight*<sup>®</sup> plate / Seramun *SpotSight*<sup>®</sup> strip and software Seramun *SpotSight*<sup>®</sup> scan

 Manufacturer	 Date of manufacture	 Use by	<b>LOT</b> Batch code	<b>REF</b> Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	<b>IVD</b> <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

SeramunBlau®, SeraSpot® und Seramun SpotSight® are registered trademarks of Seramun Diagnostica GmbH.