SeraSpot[®] HepAk-7 IgG

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SeraSpot[®] HepAk-7 IgG

Spot immunoassay for the quantitative detection of IgG autoantibodies against 7 antigens of autoimmune liver diseases in serum or plasma of human origin

REF	SP-004-7 G-S6		48
REF	SP-004-7 G-S12		96
REF	SP-004-7 G-S24		2 x 96
IVD	In-vitro-diagnostic medical device		CE

Seramun Diagnostica GmbH • Sprenhagener Str. 1 • 15754 Heidesee • Germany •
T +49 33767 791-10 • info@seramun.com • www.seramun.com

In-vitro diagnostic medical device

Unique device identifier

Manufacturer

Country of manufacture and date of manufacture

Article number

Serial number

Keep away from sunlight

Humidity limitation

Batch code

Consult instructions for use

Temperature limit

Do not reuse

Sufficient for n tests

Biohazard

Use-by date

Attention

Intended Use

SeraSpot® HepAk-7 IgG is an IVD test for the quantitative determination of autoantibodies of the IgG isotype against the antigens AMA-M2, LKM1, LC1, SLA, F-Actin, gp210 and Sp100 in serum or plasma (citrate, EDTA, heparin) of human origin by a laboratory professional user.

The test is used in combination with the Seramun SpotSight® plate mono / Seramun SpotSight® strip device for image acquisition and the software Seramun SpotSight® scan for image analysis.

It is intended to aid in the diagnosis of autoimmune liver disease in specimen materials from patients with suspicion of autoimmune liver disease.

The test must not be used with specimen materials other than serum or plasma (citrate, EDTA, heparin) of human origin, for diagnosis, monitoring, screening, prediction, prognosis, as companion diagnostic, in the near-patient setting and by lay persons.

Principle of the Test

SeraSpot® HepAk-7 IgG is a solid phase immunoassay (spot immunoassay) using recombinant or native purified antigens printed in array format (spot array) on the bottom of wells of 96-well microtiter plates. Antibodies will bind to the immobilized antigens AMA-M2, LKM1, LC1, SLA, F-Actin, gp210 and Sp100. After incubation and wash steps bound antibodies are detected by horseradish peroxidase (HRP)-labelled antibodies against human antibodies of IgG-type by substrate reaction with 3,3',5,5'-tetramethylbenzidine (TMB) and hydrogen peroxide. Blue spots indicate the formation of immune complexes. Spots range from pale to dark blue and are visible by eye.

Detection of Specific Antibodies is Performed in 3 Steps:

Step 1

Incubation of 1 : 101 diluted samples for 30 min at room temperature in selected wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with diluted wash buffer.

Step 2

Incubation of wells with the HRP-labeled anti-human IgG for 30 min at room temperature. Removal of unbound conjugate antibodies by aspiration and 3 wash cycles with diluted wash buffer.

Step 3

Incubation of wells with substrate SeramunBlau® spot dark for 30 min at room temperature. The reaction is stopped by aspirating the substrate, followed by tapping the plate / strip dry onto lint-free absorbent paper. The developed arrays are to be stored light-protected until image analysis.

Test Components (Delivery Scope)

			For 48 determinations	For 96 determinations	For 2 x 96 determinations
1	WELLS	Microtiter plate (wells with arrays) specific antigens and controls immobilized as spots in array layout	6 single breakable 8-well strips in frame color coding: green vacuum sealed with desiccant	12 single breakable 8-well strips in frame color coding: green vacuum sealed with desiccant	2 x 12 single breakable 8-well strips in frame color coding: green vacuum sealed with desiccant
2	WASHBUF (10x)	Wash buffer (10x) Seramun® Wash buffer A TRIS-based buffer	100 mL concentrate for 1000 mL buffer colorless white cap	100 mL concentrate for 1000 mL buffer colorless white cap	2 x 100 mL concentrate for each 1000 mL buffer colorless white cap
3	DIL	Sample diluent Seramun® Sample diluent B	55 mL ready-to-use solution colored red black cap	2 x 55 mL ready-to-use solution colored red black cap	4 x 55 mL ready-to-use solution colored red black cap
4	CONJ HRP IgG	Conjugate anti-Human IgG-HRP conjugate (sheep)	8.0 mL ready-to-use solution colored red red cap	8.0 mL ready-to-use solution colored red red cap	2 x 8.0 mL ready-to-use solution colored red red cap
5	SUBSTR	Substrate SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine	8.0 mL ready-to-use solution colorless blue cap	8.0 mL ready-to-use solution colorless blue cap	2 x 8.0 mL ready-to-use solution colorless blue cap
6	COVER	Covering film	2 pieces	2 pieces	4 pieces
7	SWAB	Swab 70 % (v/v) Isopropyl alcohol	3 x 2 pieces	6 x 2 pieces	12 x 2 pieces
8		Certificate of Analysis	1 piece	1 piece	1 piece
9		Instructions for Use	1 piece	1 piece	1 piece

Antigens Used

Name	Description	Clinical relevance
AMA-M2	M2-protein of the mitochondrial pyruvate dehydrogenase (PDH)-complex	primary biliary cholangitis (PBC)
LKM1	Liver-kidney microsomal antigen type 1, cytochrom P450 2D6 (LKM 1 hp)	Autoimmune Hepatitis type II (AIH II)
LC1	Liver cytosol antigen type 1, formiminotransferase cyclodeaminase (62 kDa)	Autoimmune Hepatitis type II (AIH II)
SLA	Soluble liver antigen, 50 kDa UGA-Suppressor Serin-tRNA associated Protein	Autoimmune Hepatitis type III/I (AIH III/I)
F-Actin	Polymerized filaments of the Actin Protein	Autoimmunhepatitis Typ I (AIH I), less specific
gp210	210 kDa glycoprotein of nuclear pore complex	primary biliary cholangitis (PBC), high specific
Sp100	„Speckled pattern 100“, protein of „nuclear dots“	primary biliary cholangitis (PBC)

Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micro-pipettes • adjustable 8-channel micro-pipette or multi-pipette • reagent container for multi-channel micro-pipettes • measuring cylinder • test tubes for sample dilution • test tube rack • washer for 96-well microtiter plates • laboratory system: Seramun SpotSight® plate mono / strip scanner with evaluation software Seramun SpotSight® scan • deionized water • lint-free absorbent paper • stop watch • collecting devices for infectious material • opaque cover (substrate reaction)

Important Information

This device is for *in-vitro* diagnostic use only. The kit may be performed by laboratory professional user only.

Follow the instructions carefully. The shelf life specified must be observed. Do not use reagents from damaged packages or bottles. Do not use reagents from other manufacturers.

Mixing of test kit components of different lots is only allowed for sample diluent, wash buffer and substrate.

All serious incidents occurring in relation with SeraSpot® HepAk-7 IgG must be reported to the manufacturer and the competent authority of the EU member state in which the user and/or the patient are located.

Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all components to room temperature before use.

For larger sample series, pipetting reagents from liquid reservoirs using a multi-channel pipette is recommended to avoid time delays. Reagents that appear contaminated should not be used.

Avoid time delays when dispensing reagents.

Substrate must be protected from direct light!

Incorrect dilution of samples, incorrect washing of wells and incorrect timing may lead to erroneous results.

Air bubbles generated by forceful pipetting of samples and/or reagents may cause uneven signal formation of arrays. These arrays cannot be evaluated.

Damaged arrays, e.g. by **scratching** the bottom of the well with pipet tips or washer needles, are not suitable for evaluation.

Before taking images of the wells make sure to remove particles or fibers which may adhere to the bottom side of the wells using the swab provided in the kit.

Images of dried spots may appear more intense, which can lead to slight deviations of the measured values when scanned repeatedly. Assessment of the individual parameters is not changed.

Workplace Requirements

Processing spot immunoassays requires a clean workplace. Adhesion of fibers to the plastic surface of the wells should be avoided. Fibers may cause interferences when taking images of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	EUH208 EUH210	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request.
DIL	EUH208 EUH210 -	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request. Contains material of animal origin.
CONJ HRP IgG	Hazard components EUH208 H360D P202 P280 P308+P313 P501 -	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Restricted to professional users. May damage the unborn child. Do not handle until all safety precautions have been read and understood. Wear protective gloves / protective clothing / eye protection / face protection. IF exposed or concerned: Get medical advice / attention. Dispose of contents / container in accordance with local regulations. Contains material of animal origin.
SWAB	Hazard components H225 H319 H336 P210 P280 P304+P340 P305+P351+P338	Propan-2-ol; isopropyl alcohol; isopropanol Highly flammable liquid and vapour. Causes serious eye irritation. May cause drowsiness or dizziness. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Wear protective gloves / protective clothing / eye protection / face protection / hearing protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Limitations of the Procedure

Result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks.

Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

Interference

In rare cases, specimens may contain antibodies to BSA (bovine serum albumin) and / or AGE (advanced glycation end products), which may cause nonspecific reactions, so that the test result is read as "not available" ("n.a.") by the Seramun SpotSight® scan software.

Sample Treatment

Sample Collection

Collect serum or plasma (citrate, EDTA, heparin) of human origin in suitable sampling container.

Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 7 days. For longer periods samples have to be stored at < -15 °C. Repeated freezing and thawing of samples should be avoided!

Sample Preparation

Bring samples to room temperature before use. Homogeneity should be ensured by shaking components briefly.

Samples have to be diluted 1 : 101 (v/v) with sample diluent.

Example: 10 µL sample and 1000 µL sample diluent.

Reagent Treatment

Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips is stable until the printed expiration date when stored at 2...8 °C. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

Reagent Preparation

Allow all components to reach **room temperature** prior to use. The **microtiter plate** is vacuum sealed with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored at 2...8 °C and protected from moisture in the original cover carefully resealed.

Dilute **wash buffer (10x)** 1 : 10 with deionized water.

Example: 10 mL wash buffer (10x) + 90 mL deionized water. The prepared wash buffer must be thoroughly mixed before use.

The **substrate** must be protected from direct light. The test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

Assay Procedure

Performance at room temperature (RT, 18...25 °C). Bring **all** ready-to-use test reagents, prepared wash buffer and microplate to **room temperature**. Mix all reagents before use by gently shaking, avoid foaming.

The aspiration and washing steps can be carried out manually or with the aid of a microtiter plate washer.

Important Notes on Test Procedure:

Avoid mechanical contact (**scratching**) on the bottom of the wells with pipette tips or washer needles. This will irreparably damage the array!

All liquid reagents (sample, conjugate and substrate) have to be pipetted **without** causing **air bubbles** into the wells.


1. Provide the required number of wells **WELLS** according to the number of samples.
2. Preparation of working dilution of the sample: 1 : 101 dilution
e.g. 10 µL sample to 1000 µL sample diluent **DIL**
3. Pipette **100 µL** diluted sample per well.
4. Cover plate with covering film **COVER** and incubate for **30 min** at RT.
5. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,
Tap plate dry onto lint-free absorbent paper.
6. **Add 50 µL** of the conjugate **CONJ HRP IgG** (anti-Human IgG-HRP) to each well.
7. Cover plate with covering film **COVER** and incubate for **30 min** at RT.
8. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,
Tap plate dry onto lint-free absorbent paper.
9. **Add 50 µL** of the substrate **SUBSTR** to each well.
10. Cover plate and protect from light. Incubate for **30 min** at RT.
11. Remove liquid. Tap plate dry onto lint-free absorbent paper.
12. Clean bottom of wells with swab **SWAB** shortly before scanning the images.
13. **Take images** using the Seramun SpotSight® plate mono / strip und evaluate the results with the Seramun SpotSight® scan software.

If image acquisition is performed with the scanner Seramun SpotSight® strip, transfer single 8-well microtiter plate strips to the scanner's strip holder.

After aspiration of the substrate the color of developed spots is stable for 24 h when the plate is stored protected from light.

Evaluation of Results

Arraylayout

Antigens

05	LKM1
06	AMA-M2
07	LC1
09	SLA
10	Sp100
12	F-Actin
13	gp210

Controls

01	Positive control (PC)
02	Negative control (NC), 0 rel.Units
03	Cut-off control (CO), 30 rel.Units
04	Reference 3 (R3), 60 rel.Units
08	Reference 2 (R2), 100 rel.Units
11	Reference 1 (R1), 300 rel.Units
14	Serum control (SC)
P	Well position marker

Quantitative evaluation

Validity Criteria for the Test

The spots NC, CO, R1 ... R3 are used to create a reference curve (4-parameter non-linear regression) to calculate the staining intensities of the antigen spots in relative units (rel.Units).

SeraSpot® HepAk-7 IgG includes the following control spots:

1. Positive control (PC). Intense spot, stained darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than cut-off control
4. Serum control (SC). Intense stained spot, always stained in presence of sample. Absence of spot indicates absence of sample.
5. Reference 1 (R1) corresponds to 300 rel.Units. Intense spot, always stained.
6. Reference 2 (R2) corresponds to 100 rel.Units. Weaker staining than R1, always stained.
7. Reference 3 (R3) corresponds to 60 rel.Units. Weaker staining than R2, always stained.

The test cannot be evaluated if one of the quality criteria listed in points 1. to 7. is not met.

If the above mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the validity criteria contact the manufacturer.

Interpretation of Results

Evaluation of the test has to be performed by using the Seramun SpotSight® plate mono / strip and in combination with the Seramun SeraSpot® scan software.

The results are interpreted as follows:

Evaluation	Conditions
Positive	Color intensity of one or more antigen spots* > cut-off control
Negative	Color intensity of the antigen spots ≤ cut-off control

*When evaluating the result it should be noted that the S2k guideline on autoimmune liver diseases (2017) recommends the detection of anti-F-Actin autoantibodies of the IgG type by indirect immunofluorescence (5).

Performance Characteristics

Precision

Samples with known antibody reactivity were examined in the SeraSpot® HepAk-7 IgG and the rel.Units were determined. These values were then used to determine the coefficients of variation (CV) as a measure of the precision within a test run (intra-assay CV), between different test runs (inter-assay CV) and between different test lots (lot-to-lot CV).

Antigen	Intra-assay coefficient of variation		Inter-assay coefficient of variation		Lot-to-lot coefficient of variation	
	\bar{x} rel.Units n=40	CV [%]	\bar{x} rel.Units n=80	CV [%]	\bar{x} rel.Units n=240	CV [%]
AMA-M2	36.8	6.8	38.9	14.9	42.2	16.3
LKM1	58.2	5.3	84.8	12.3	86.2	12.6
LC1	61.5	10.1	60.2	16.7	62.5	18.1
SLA	66.3	6.2	64.9	8.6	77.7	18.1
F-Actin	53.8	14.7	59.0	15.3	51.2	19.9
gp210	52.9	6.4	51.2	8.8	49.9	13.8
Sp100	88.6	6.2	84.7	11.4	92.4	18.7
Procedure	1 operator 40 determinations 1 batch		5 operators 2 determinations 2 testings per day 20 days 1 batch		5 operators 2 determinations 2 testings per day 20 days 3 batches	

Determination of the Cut-off Value

Data from sensitivity and specificity studies compared to reference tests were used to determine the value for the cut-off control by ROC analysis (Receiver Operating Characteristics). The cut-off range is determined specifically for each test.

Interfering Substances

All tested substances do not have a significant impact on the test results even if present in elevated concentrations in the serum: 20 mg/dL bilirubin C and F each (simulation of icteric samples). 500 mg/dL hemoglobin (simulation of hemolytic samples). 1000 mg/dL lipids (triglycerides) (simulation of lipemic samples) and 500 IU/mL rheumatoid factor.

Sensitivity

Pre-characterized samples were analyzed to determine the sensitivity and specificity of SeraSpot® HepAk-7 IgG in comparison with IgG reference tests. Samples with discrepant results were retested in a second reference test.

AMA-M2		Reference test	
		positive	negative
SeraSpot®	positive	29	0
	negative	0	0
Samples total		29	0
Sensitivity:		100 %	

LKM1		Reference test	
		positive	negative
SeraSpot®	positive	8	0
	negative	0	0
Samples total		8	0
Sensitivity:		100 %	

LC1		Reference test	
		positive	negative
SeraSpot®	positive	7	0
	negative	0	0
Samples total		7	0
Sensitivity:		100 %	

SLA		Reference test	
		positive	negative
SeraSpot®	positive	7	0
	negative	0	0
Samples total		7	0
Sensitivity:		100 %	

F-Actin		Reference test	
		positive	negative
SeraSpot®	positive	12	0
	negative	1	0
Samples total		13	0
Sensitivity:		92.3 %	

gp210		Reference test	
		positive	negative
SeraSpot®	positive	19	0
	negative	0	0
Samples total		19	0
Sensitivity:		100 %	

Sp100		Reference test	
		positive	negative
SeraSpot®	positive	18	0
	negative	0	0
Samples total		18	0
Sensitivity:		100 %	

Specificity

Blood donor samples were examined for the specificity in comparison to reference tests in IgG detection of SeraSpot® HepAk-7 IgG. Samples with discrepant results were retested in a second reference test.

AMA-M2		Reference test	
		positive	positive
SeraSpot®	positive	0	1
	negative	0	191
Samples total		0	192
Specificity:		99.5 %	

LKM1		Reference test	
		positive	negative
SeraSpot®	positive	0	1
	negative	0	191
Samples total		0	192
Specificity:		99.5 %	

LC1		Reference test	
		positive	negative
SeraSpot®	positive	1	2
	negative	0	189
Samples total		1	191

Specificity: 99.0 %

SLA		Reference test	
		positive	negative
SeraSpot®	positive	0	0
	negative	0	192
Samples total		0	192

Specificity: 100 %

F-Actin		Reference test	
		positive	negative
SeraSpot®	positive	0	0
	negative	0	182
Samples total		0	182

Specificity: 100 %

gp210		Reference test	
		positive	negative
SeraSpot®	positive	0	1
	negative	0	190
Samples total		0	191

Specificity: 99.5 %

Sp100		Reference test	
		positive	negative
SeraSpot®	positive	0	0
	negative	0	191
Samples total		0	191

Specificity: 100 %

Cross Reactivity

Potential cross-reacting samples were tested in the SeraSpot® HepAk-7 IgG and re-tested in a reference test. Samples with discrepant results were re-tested in a second reference test. Results obtained with the second reference test were used for the final status determination of the discrepant samples.

Hepatitis-B-Virus (HBV)- positive samples (n=27)

HBV		Reference test	
		positive	negative
SeraSpot®	positive	0	0
	negative	0	27
Samples total		0	27

Specificity: 100 %

Hepatitis-C-Virus (HCV)- positive samples (n=43)

HCV		Reference test	
		positive	negative
SeraSpot®	positive	0	0
	negative	0	43
Samples total		0	43

Specificity: 100 %

Application

Automated Workflow

Manually diluted positive samples (n=96) were assayed side by side in SeraSpot® HepAk-7 IgG. by hand and by use of an automatic microplate processor. Data were generated with the device DS2® (Dynex Technologies). The correlation was calculated with $R^2 > 0.9$ for the antigens AMA-M2, LKM1, LC1, SLA, gp210 and Sp100. Assay procedure using an automatic microplate processor such as DS2® is possible. R^2 was not determined if more than 90 % of the samples tested showed no antibody reactivity for the respective antigens.

Change History

Version	Section	Modifications
2022-10_v01_DE_EN	Entire document	Use of new IFU template Updating of the intended use Conversion of the subsections Insertion of the safety instructions Removal of pipetting scheme
	Important Information	Insertion note on dealing with serious incidents
	Entire document	Updating of wording and correction of spelling mistakes

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