SeraSpot® Anti-Yersinia-6 IgG SeraSpot® Anti-Yersinia-6 IgA

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Literatur:

1. Cremer, J., Putzker, M., Faulde, M., Zöller, L., Immunoblotting of Yersinia plasmid-encoded released proteins: A tool for serodiagnosis, *Electrophoresis* 1993, 14: 952-959
2. Iriarte, M., Cornelis, G.R, Identification of SycN, YscX, and YscY, Three new elements of the Yersinia Yop Virulon, *J. Bacteriol.* 1999, 181/2: 675-680
3. Heesemann, J., Gross, U., Schmidt, N., Laufs, R., Immunochemical analysis of plasmid-encoded proteins released by enteropathogenic Yersinia sp. grown in calcium-deficient media *Infect. Immun.* 1986, 54/2: 561-567
4. Heesemann, J., Grüter, L., Genetic evidence that the outer membrane protein Yop1 of *Yersinia enterocolitica* mediates adherence and phagocytosis resistance to human epithelial cells *FEMS Microbiol. Lett.* 1987, 40, 37-41
5. Kendrick, C.J., Baker, B., Morris, A.J., O'Toole, P.W., Identification of Yersinia-infected blood donors by anti-Yop IgA immunoassay, *Transfusion*, 2001, 41: 1365-1372
6. Brandis, H., Köhler, W., Eggers, H.J., Pulverer, G., *Lehrbuch der Medizinischen Mikrobiologie* Gustav Fischer Verlag Stuttgart Jena New York, 1994, 7. Auflage
7. Cornelis, G.R., Boland, A., Boyd, A.P., Geuijen, C. Iriarte, M., Neyt, C., Sory, M.P., Stainier, I. The virulence plasmid of Yersinia, an antihost genome *Microbiol. Mol. Biol. Rev.* 1998, 62/4, 1315-1352
8. Paerregaard, Interaction between *Yersinia enterocolitica* and the host with special reference to virulence plasmid encoded adhesion and humoral immunity *Dan. Med. Bull.* 1992, 39/2, 155-172
9. Rastawicki, W., Humoral response to selected antigens of Yersinia enterocolitica and Yersinia pseudotuberculosis in the course of yersiniosis in humans. I. Occurrence of antibodies to Yersinia lipopolysacharydes and Yop proteins by ELISA *Med Dosw Mikrobiol* 2006; 58(4): 303-319
10. Autenrieth, S.E. and Autenrieth, I.B., Yersinia enterocolitica: subversion of adaptive immunity and implications for vaccine development. *Int J Med Microbiol* 2008; 298(1-2): 69-77
11. Thomas, L. (Hrsg.), *Labor und Diagnose*. TH Books Verlagsgesellschaft FFM, 5. Auflage 1998
12. Rosner, B.M. et al., Clinical aspects and self-reported symptoms of sequelae of Yersinia enterocolitica infections in a population-based study, Germany 2009-2013. *BMC Infectious Diseases* 2013, 13, 236
13. Atkinson, S. and Williams, P., Yersinia virulence factors-a sophisticated arsenal for combating host defences. *F1000 Research* 2016, 5 (F1000 Faculty Rev.): 1370, Last updated 14 JUN 2016
14. Cambronne, E. D. and Schneewind, O., Yersinia enterocolitica Type III Secretion: yscM1 and yscM2 Regulate yop Gene Expression by a Posttranscriptional Mechanism That Targets the S'Untranslated Region of yop mRNA. *J. Bacteriol.* 2002, Vol. 184(21), 5880-5893
15. Rosqvist, R. et al., Intracellular Targeting of the Yersinia YopE Cytotoxin in Mammalian Cells Induces Actin Microfilament Disruption. *Infection and Immunity* 1991, Vol. 59(12), 4562-4569
16. Black, D.S. and Bliska, J.B., The Rho GAP activity of the Yersinia pseudotuberculosis cytotoxin YopE is required for antiphagocytic function and virulence. *Mol. Microbiol.* 2000, 37, 515-527
17. Guan, K.L. and Dixon, J.E., Protein tyrosine phosphatase activity of an essential virulence determinant in Yersinia. *Science* 1990, 249, 553-556
18. Ruckdeschel, K. et al., Differential Contribution of Yersinia enterocolitica Virulence Factors to Evasion of Microbicidal Action of Neutrophils. *Infection and Immunity* 1996, Vol.64(3), 724-733
19. Yao, T. et al., Suppression of T and B Lymphocyte Activation by a Yersinia pseudotuberculosis Virulence Factor, YopH. *J. Exp. Med.* 1999, Vol. 190(9), 1343-1350
20. Leung, K.Y. et al., YopM inhibits platelet aggregation and is necessary for virulence of Yersinia pestis in mice. *Infection and Immunity* 1990, 58, 3262-3271
21. Pawel-Ramminger, U.v. et al., GAP activity of the Yersinia YopE Cytotoxin Specifically targets the Rho pathway: a mechanism for disruption of actin microfilaments. *Mol. Microbiol.* 2000, 36(3), 737-748

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SeraSpot® Anti-Yersinia-6 IgG

SeraSpot® Anti-Yersinia-6 IgA

Spotimmunoassay for detection of IgG or IgA antibodies to *Yersinia enterocolitica*
in human serum or plasma

REF SP-005-6 G-S6 ▽ 48	REF SP-005-6 G-S12 ▽ 96	REF SP-005-6 G-S24 ▽ 2x 96
IVD <i>In-vitro</i> -diagnostic device CE		
REF SP-005-6 A-S6 ▽ 48	REF SP-005-6 A-S12 ▽ 96	REF SP-005-6 A-S24 ▽ 2x 96
IVD <i>In-vitro</i> -diagnostic device CE		

Seramun Diagnostica GmbH · Sprenghagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 33767 791-10 · fax +49 33767 791-99 · info@seramun.com

Introduction

The genus *Yersinia* comprises gram-negative bacteria belonging to the family Enterobacteriaceae. There are at least three *Yersinia* species which are pathogenic to human: *Y. pestis*, *Y. pseudotuberculosis* and *Y. enterocolitica*. More than 60 serovariants are known from *Y. enterocolitica* (6, 11). The infection of humans is transmitted by the consumption of contaminated meat or sausages, unpasteurized milk or contaminated drinking water. The transmission by direct contact from person to person is rather seldom. The pathogens enter the Peyer's patches of the small intestine, disseminate to lymphoid tissues and initiate a humoral immune response (10). In the first 2-3 weeks after the infection IgA and IgM class antibodies are dominating, then IgG antibodies are mainly detectable, which usually persist for several years (9,11). In cases of immunopathological complications and chronic forms of *Yersinia* infections IgA antibodies are detectable for a longer period of time (1). In patients with reactive arthritis IgA antibodies to specific *Yersinia* antigens may persist for up to 1 year. In these patients the corresponding antibody class is also detectable in the synovial fluid (9).

The human pathogenic *Yersinia* species are characterized by a 70 kb plasmid encoding for cytosolic, membraneous and releasing proteins. The main virulence markers are the *Yersinia* outer proteins (Yop) which are released by the bacteria (2, 8). These proteins are highly immunogenic and induce the generation of specific antibodies of the IgA, IgG and IgM isotype (3, 4).

Although no precise data are available the assumed incubation period for *Y. enterocolitica* ranges from 2 to 10 days. Many mammalian species are known to serve as reservoir hosts, but rodents, domestic and farm animals are of major epidemiological importance (6).

The clinical symptoms of an infection with *Y. enterocolitica* and *Y. pseudotuberculosis* are very similar. The clinical spectrum of *Yersinia enterocolitica* infections comprises acute intestinal symptoms (enterocolitis, watery diarrhea) as well as extra intestinal complications like reactive arthritis, Erythema nodosum and Reiter's syndrome as late sequelae of *Yersinia* enteritis. In infants an infection induces a self-limiting gastroenteritis frequently accompanied by fever and vomiting. Juveniles often react with a lymphadenitis characterized by abdominal pain often mistaken as appendicitis (8,11,12).

Literature:

1. Cremer, J., Putzker, M., Faulde, M., Zöller, L., Immunoblotting of Yersinia plasmid-encoded released proteins: A tool for serodiagnosis, *Electrophoresis* 1993, 14: 952-959
2. Iriarte, M., Cornelis, G.R, Identification of SycN, YscX, and YscY, Three new elements of the Yersinia Yop Virulon, *J. Bacteriol.* 1999, 181/2: 675-680
3. Heesemann, J., Gross, U., Schmidt, N., Laufs, R., Immunochemical analysis of plasmid-encoded proteins released by enteropathogenic Yersinia sp. grown in calcium-deficient media *Infect. Immun.* 1986, 54/2: 561-567
4. Heesemann, J., Grüter, L., Genetic evidence that the outer membrane protein Yop1 of *Yersinia enterocolitica* mediates adherence and phagocytosis resistance to human epithelial cells *FEMS Microbiol. Lett.* 1987, 40, 37-41
5. Kendrick, C.J., Baker, B., Morris, A.J., O'Toole, P.W., Identification of Yersinia-infected blood donors by anti-Yop IgA immunoassay, *Transfusion*, 2001, 41: 1365-1372
6. Brandis, H., Köhler, W., Eggers, H.J., Pulverer, G., *Lehrbuch der Medizinischen Mikrobiologie* Gustav Fischer Verlag Stuttgart Jena New York, 1994, 7. Auflage
7. Cornelis, G.R., Boland, A., Boyd, A.P., Geuijen, C. Iriarte, M., Neyt, C., Sory, M.P., Stainier, I. The virulence plasmid of Yersinia, an antihost genome *Microbiol. Mol. Biol. Rev.* 1998, 62/4, 1315-1352
8. Paerregaard, Interaction between *Yersinia enterocolitica* and the host with special reference to virulence plasmid encoded adhesion and humoral immunity *Dan. Med. Bull.* 1992, 39/2, 155-172
9. Rastawicki, W., Humoral response to selected antigens of Yersinia enterocolitica and Yersinia pseudotuberculosis in the course of yersiniosis in humans. I. Occurrence of antibodies to Yersinia lipopolysacharydes and Yop proteins by ELISA *Med Dosw Mikrobiol* 2006; 58(4): 303-319
10. Autenrieth, S.E. and Autenrieth, I.B., Yersinia enterocolitica: subversion of adaptive immunity and implications for vaccine development. *Int J Med Microbiol* 2008; 298(1-2): 69-77
11. Thomas, L. (Hrsg.), *Labor und Diagnose*. TH Books Verlagsgesellschaft FFM, 5. Auflage 1998
12. Rosner, B.M. et al., Clinical aspects and self-reported symptoms of sequelae of Yersinia enterocolitica infections in a population-based study, Germany 2009-2013. *BMC Infectious Diseases* 2013, 13, 236
13. Atkinson, S. and Williams, P., Yersinia virulence factors-a sophisticated arsenal for combating host defences. *F1000 Research* 2016, 5 (F1000 Faculty Rev.): 1370, Last updated 14 JUN 2016
14. Cambronne, E. D. and Schneewind, O., Yersinia enterocolitica Type III Secretion: yscM1 and yscM2 Regulate yop Gene Expression by a Posttranscriptional Mechanism That Targets the 5' Untranslated Region of yop mRNA. *J. Bacteriol.* 2002, Vol. 184(21), 5880-5893
15. Rosqvist, R. et al., Intracellular Targeting of the Yersinia YopE Cytotoxin in Mammalian Cells Induces Actin Microfilament Disruption. *Infection and Immunity* 1991, Vol. 59(12), 4562-4569
16. Black, D.S. and Bliska, J.B., The Rho GAP activity of the Yersinia pseudotuberculosis cytotoxin YopE is required for antiphagocytic function and virulence. *Mol. Microbiol.* 2000, 37, 515-527
17. Guan, K.L. and Dixon, J.E., Protein tyrosine phosphatase activity of an essential virulence determinant in Yersinia. *Science* 1990, 249, 553-556
18. Ruckdeschel, K. et al., Differential Contribution of Yersinia enterocolitica Virulence Factors to Evasion of Microbicidal Action of Neutrophils. *Infection and Immunity* 1996, Vol.64(3), 724-733
19. Yao, T. et al., Suppression of T and B Lymphocyte Activation by a Yersinia pseudotuberculosis Virulence Factor, YopH. *J. Exp. Med.* 1999, Vol. 190(9), 1343-1350
20. Leung, K.Y. et al., YopM inhibits platelet aggregation and is necessary for virulence of Yersinia pestis in mice. *Infection and Immunity* 1990, 58, 3262-3271
21. Pawel-Ramminger, U.v. et al., GAP activity of the Yersinia YopE Cytotoxin Specifically targets the Rho pathway: a mechanism for disruption of actin microfilaments. *Mol. Microbiol.* 2000, 36(3), 737-748

Intended use

The **SeraSpot® Anti-Yersinia-6 IgG / SeraSpot® Anti-Yersinia-6 IgA test** is an *in vitro* diagnostic device (spotimmunoassay, SIA) for the detection of *Yersinia*-specific IgG or IgA antibodies in human serum, plasma or synovial fluid samples. It can be used as screening assay as well as confirmatory assay after a positive ELISA result.

Principle of the test

The **SeraSpot® Anti-Yersinia-6 IgG / SeraSpot® Anti-Yersinia-6 IgA test** is a solid phase immunoassay based on the use of recombinant *Yersinia* proteins as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96-well microtiter plates. The recombinant antigens serve as capture molecules for antibodies against *Yersinia* antigens. Bound antibodies are detected by horseradish peroxidase-(HRP)-labeled antibodies against human antibodies of IgG- or IgA-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed by precipitated product from colorless substrate solution. Color intensity is correlated to the antibody concentration. Pale blue to dark blue spots are visible by eye.

The test procedure for antibody detection is performed in 3 steps:

Step 1

Incubation of 1 : 51 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are sucked off and unbound components are removed by 3 wash cycles with wash solution.

Step 2

Incubation of wells with the HRP-labeled anti-human IgG- or IgA for 30 minutes at room temperature. After incubation the antibody is sucked off and unbound components are removed by 3 wash cycles with wash solution.

Step 3

Incubation of wells with substrate solution **SeramunBlau®** spot dark for 30 minutes at room temperature. The reaction is stopped by sucking off the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the **Seramun SpotSight®** plate or **Seramun SpotSight®** strip scanner; the obtained images are interpreted by using the Software **Seramun SpotSight®** scan. Taking account of the test-specific evaluation criteria, the test result is determined.

Used Yersinia antigens

Nomenclature	Characterization / Function		Relevance
YopB – Yersinia outer protein B (12)	42 kD	Part of the Translocon	Yersinia outer proteins (Yop) are structure- and effector proteins of the Type III secretion systems of human pathogenic Yersinia. All listed Yop are high specific for human pathogenic Yersinia.
YopD – Yersinia outer protein D (12, 13)	36 kD	Part of the Translocon, Regulation of Yop expression	
YopE – Yersinia outer protein E (14, 15, 20)	27 kD	Inhibition of phagocytosis	
YopH – Yersinia outer protein H (16, 17, 18)	48 kD	Inhibition of phagocytosis	
YopM – Yersinia outer protein M (19)	45 kD	Activation of eukaryotic kinases	
YopN – Yersinia outer protein M (7)	34 kD	Regulation of Yop release	

*The molecular weight of the antigens differs in the literature.

Kit components

		For 48 determinations	For 96 determinations	For 2x 96 determinations	
1	WELLS	Wells with arrays	6 single breakable 8-well strips in frame vacuum-sealed with desiccant	12 single breakable 8-well strips in frame vacuum-sealed with desiccant	2x 12- single breakable 8-well strips in frame separately vacuum-sealed with desiccant
IgG detection: color coding dark brown; IgA detection: color coding light brown					
2	WASHBUF CONC 10X	Wash buffer Seramun® Wash buffer A	100 ml concentrate transparent bottle, white cap	100 ml concentrate transparent bottle, white cap	2x 100 ml concentrate transparent bottle white cap
3	DIL	Sample diluent Seramun® Sample diluent B	55 ml ready-to-use solution, colored red transparent bottle black cap	2x 55 ml ready-to-use solution, colored red, transparent bottle, black cap	4x 55 ml ready-to-use solution, colored red, transparent bottle, black cap
4	CONJ HRP IgG	Anti-human IgG-HRP-conjugate	8 ml ready-to-use solution, colored red, transparent bottle, red cap	8 ml ready-to-use solution, colored red, transparent bottle, red cap	2x 8 ml ready-to-use solution, colored red, transparent bottle, red cap
	CONJ HRP IgA	Anti-human IgA-HRP-conjugate	8 ml ready-to-use solution, colored purple, transparent bottle, purple cap	8 ml ready-to-use solution, colored purple, transparent bottle, purple cap	2x 8 ml ready-to-use solution, colored purple, transparent bottle, purple cap
5	SUBSTR TMB	Substrate SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	8 ml ready-to-use solution, black bottle, blue cap	8 ml ready-to-use solution, black bottle, blue cap	2x 8 ml ready-to-use solution, black bottle, blue cap
6	COVER	Adhesive film	2x	2x	4x
7	SWAB	Swab with 70 % (v/v) Isopropyl alcohol	3x2	6x2	12x2

Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 48 hours. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided!

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 51 (e.g. 10 µl sample and 500 µl sample diluent) with the ready-to-use sample diluent.

For individual cases the investigation of synovial fluid is possible in principle. A sample dilution of 1:51 is recommended (e.g. 10 µl sample and 500 µl sample diluent). Due to the limited number of tested samples a common recommendation for result interpretation cannot be given. The responsibility lies with the laboratory's own authority.

Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- stopwatch
- Washer for 96-well microtiter plates (optional)
- Collecting devices for infectious material
- Seramun *SpotSight*® plate or Seramun *SpotSight*® strip scanner with evaluation software Seramun *SpotSight*® scan

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used up to 4 weeks when stored at 2...8 °C and up to 2 weeks when stored at room temperature.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtiter plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or de-ionized water. For Example: 10 ml **WASHBUF CONC 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

Workplace requirements

Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

Assay procedure

Execution at room temperature (18...25 °C) and at the specified incubation times.

For the handling of *SeraSpot*[®] tests the use of the starter kit (order number: SP-000-1) is recommended. The kit is designed for multiple use.

Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere to the underside of the bottom of the wells using the swab **[SWAB]** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

Working steps


1. Warm all reagents to room temperature before use. Mix gently without causing foam.
2. Add **100 µl** of 1 : 51 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 500 µl sample dilution buffer.
3. Cover plate with adhesive film **[COVER]** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of (2)). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **[CONJ HRP IgG]** or **[CONJ HRP IgA]** to each well.
6. Cover plate with adhesive film **[COVER]** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard conjugate), then wash wells **3 x** using **400 µl** wash solution per well (made of (2)). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **[SUBSTR TMB]** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent paper.
11. Clean bottom of wells with Swab **[SWAB]**.
12. Perform image acquisition using the Seramun *SpotSight*[®] plate or the Seramun *SpotSight*[®] strip scanner and evaluate the results with the Seramun *SpotSight*[®] scan software. If image acquisition is performed with the scanner Seramun *SpotSight*[®] strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

Aspiration and washing steps can be performed manually or via washer for 96-well microtiter plates.

After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when stored the plate in the dark.

Test evaluation

Array layout

Parameter	Controls
04 YOP B	01 Positive control (PC)
05 YOP D	02 Cut-off control (CO)
06 YOP E	03 Negative control (NC)
07 YOP H	10 IgG conjugate control (GC)
08 YOP M	11 IgA conjugate control (AC)
09 YOP N	12 Serum control (SC)
	P Well position marker

Validity criteria for the test

The *SeraSpot*[®] Anti-Yersinia-6 IgG / *SeraSpot*[®] Anti-Yersinia-6 IgA tests include the following control spots:

1. Positive control (PC). Intensively stained spot, darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than cut-off control
4. IgG, IgA conjugate control (GC, AC). Intensively stained spots with different position for IgG and IgA antibody detection. Serving as antibody isotype control.
5. Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.

The test cannot be evaluated if one of the validity criteria listed in points 1 to 5 is not fulfilled.

Result interpretation

Evaluation of the test has to be performed by using the *Seramun SpotSight*[®] scan software. A visual plausibility check is possible with the aid of the magnifying lens contained in the starter kit (order number: SP-000-1).

Judgment of spot staining is performed according to the following classification:

Result Interpretation	IgG	IgA
negative	all antigen spots < Cut-off control	all antigen spots < Cut-off control
borderline	YOP D-Spot = Cut-off control or one antigen-Spot except YOP D > Cut-off control	one antigen-Spot except YOP D ≥ Cut-off control
positive	YOP D-Spot > Cut-off control or at least two different antigen spots except YOP D > Cut-off control	YOP D-Spot > Cut-off control or at least two different antigen spots except YOP D > Cut-off control

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

Interference

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) and rheumatoid factors (500 IU/ml) do not interfere with the test.

Investigation of synovial fluid

In individual cases the detection of anti-*Yersinia enterocolitica* antibodies in synovial fluid with the SeraSpot® Anti-*Yersinia*-6 IgG/IgA test is possible in principle. In the course of performance evaluation studies a total of 24 joint punctate samples were tested in parallel by SeraSpot® Anti-*Yersinia*-6 IgG/IgA and a commercially available line immunoassay. The overall agreement for IgG and IgA antibody detection was 91.7 % and 95.8 % resp. Due to this limited number of samples there is no common recommendation for a final result interpretation. Result interpretation should be conducted by the laboratory in consideration of further laboratory data and clinical findings.

Assay procedure

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere to the underside of the bottom of the wells or to the array area cause erroneous results.

Performance characteristics

Diagnostic sensitivity

Serologically pre-determined samples (n = 217 (IgG determination), n = 218 (IgA determination), reference test 1) of patients with suspect of Yersiniosis were tested in the *SeraSpot*[®] Anti-Yersinia-6 IgG / *SeraSpot*[®] Anti-Yersinia-6 IgA test ("Initial sensitivity"). Samples with discrepant results were retested in a second assay for detection of anti-Yersinia antibodies. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results ("Amended sensitivity").

IgG	n = 217	Reference test 1		
		positive (n = 209)	borderline (n = 0)	negative (n = 8)
<i>SeraSpot</i> [®] Anti-Yersinia-6 IgG	positive	200	0	2
	borderline	1	0	0
	negative	8	0	6

Borderline samples were valued positive.

Initial sensitivity: **96.2 %**

Amended sensitivity: **97.6 %**

IgA	n = 218	Reference test 1		
		positive (n = 112)	borderline (n = 0)	negative (n = 106)
<i>SeraSpot</i> [®] Anti-Yersinia-6 IgA	positive	95	0	28
	borderline	3	0	14
	negative	14	0	64

Borderline samples were valued positive.

Initial sensitivity: **87.5 %**

Amended sensitivity: **95.2 %**

Diagnostic specificity

A total of 205 predetermined blood donors (Reference test 1) were tested in the *SeraSpot*[®] Anti-Yersinia-6 IgG / *SeraSpot*[®] Anti-Yersinia-6 IgA test. Samples with discrepant results were retested and assessed as described in "Diagnostic sensitivity".

IgG	n = 205	Reference test 1		
		positive (n = 33)	borderline (n = 13)	negative (n = 159)
<i>SeraSpot</i> [®] Anti-Yersinia-6 IgG	positive	30	9	20
	borderline	3	4	11
	negative	0	0	128

Borderline samples were valued positive.

Initial specificity: **80.8 %**

Amended specificity: **87.1 %**

IgA	n = 205	Reference test 1		
		positive (n = 13)	borderline (n = 6)	negative (n = 186)
<i>SeraSpot</i> [®] Anti-Yersinia-6 IgA	positive	13	2	5
	borderline	0	4	6
	negative	0	0	175

Borderline samples were valued positive.

Initial specificity: **94.1 %**

Amended specificity: **95.6 %**

Precision

Samples of known antibody titer were assayed by *SeraSpot*[®] Anti-Yersinia-6 IgG / *SeraSpot*[®] Anti-Yersinia-6 IgA test. Color intensity (OD) of spots was recorded and the coefficients of variation (CV) were calculated.

IgG detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	OD [\bar{x}] n = 8	CV* [%]	OD [\bar{x}] n = 48	CV* [%]	OD [\bar{x}] n = 144	CV* [%]
YOP B	58.2	1.4	58.2	6.1	57.4	5.7
YOP D	70.4	2.2	72.4	8.7	72.5	9.0
YOP E	59.4	1.4	56.4	4.5	54.3	4.5
YOP H	37.1	4.2	32.2	14.1	32.3	11.1
YOP M	39.7	1.9	38.9	5.4	38.5	5.0
YOP N	35.7	2.1	32.3	10.6	31.6	11.0
Procedure	1 operator, 8x determinations, 1 batch		1 operator, 8x determinations, 2x testings per day, 3 days, 1 batch		1 operator, 8x determinations, 2x testings per day, 3 days, 3 batches	

CV*: The CV for OD-values < 20 has not been determined.

IgA detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	OD [\bar{x}] n = 8	CV* [%]	OD [\bar{x}] n = 48	CV* [%]	OD [\bar{x}] n = 144	CV* [%]
YOP B	22.4	2.6	22.4	9.5	23.6	8.6
YOP D	121.3	2.6	130.5	12.1	133.4	11.7
YOP M	26.8	6.0	28.3	6.2	28.4	4.9
Procedure	1 operator, 8x determinations, 1 batch		1 operator, 8x determinations, 2x testings per day, 3 days, 1 batch		1 operator, 8x determinations, 2x testings per day, 3 days, 3 batches	

CV*: The CV for OD-values < 20 has not been determined.

Automation

Samples (n = 96) of known antibody titer (IgG and IgA) were assayed manually or by use of the microplate processor DS2[®] (Dynerx Technologies; manual sample dilution). Color intensity (OD) of spots was recorded and used for calculation of the coefficient of determination r^2 for selected antibody reactivities and the different methods of operation.

Ig isotype detection	Dynerx DS2 [®] vs. manual processing	
	r^2 IgG	r^2 IgA
YOP B	0.979	0.960
YOP D	0.993	0.978
YOP E	0.984	
YOP H	0.960	
YOP M	0.978	0.966
YOP N	0.966	

r^2 was not determined if more than 90 % of the total number (n = 96) of the samples showed no antibody reactivity for the particular antigens.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

Do not use or mix reagents from different lots except for sample diluent, washing buffer and TMB substrate solution. These can be used across batches and parameters.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

History of changes

Version	Section	Modifications
2017-10-20	General	Adding of synovial fluid as sample matrix
	Assay procedure	Update
	Test evaluation	Update
	Performance characteristics	Update of precision and automation
	Common advices and precautions	Update
2018-03-01	Common advices and precautions	Update
2019-02-15	Kit components	Update
	Preparation and storage of reagents	Addition

Incubation scheme *SeraSpot*[®] Anti-Yersinia-6 IgG / IgA

1.  100 µl diluted sample (1 : 51)
30 min incubation (room temperature)

 3 x wash with wash solution, each 400 µl per well
2.  50 µl conjugate **CONJ HRP IgG** or **CONJ HRP IgA** (4)
30 min incubation (room temperature)

 3 x wash with wash solution, each 400 µl per well
3.  50 µl substrate **SUBSTR TMB** (5)
30 min incubation (room temperature, protected from light)

 aspiration
4. imaging of wells and image analysis scanner Seramun *SpotSight*[®] plate / Seramun *SpotSight*[®] strip and software Seramun *SpotSight*[®] scan

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

SeramunBlau[®], *SeraSpot*[®] und Seramun *SpotSight*[®] are registered trademarks of Seramun Diagnostica GmbH.