SeraSpot[®] Anti-Borrelia-10 IgG

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SeraSpot[®] Anti-Borrelia-10 IgG

Spot immunoassay for the qualitative detection of IgG antibodies against 10 *Borrelia burgdorferi* sensu lato specific antigens in serum or plasma of human origin.

REF	SP-006-10 G-S6		48
REF	SP-006-10 G-S12		96
REF	SP-006-10 G-S24		2 x 96
REF	SP-006-10 G-S120		10 x 96
IVD	In-vitro-diagnostic medical device		

Seramun Diagnostica GmbH • Sprenhagener Str. 1 • 15754 Heidesee • Germany •
T +49 33767 791-10 • info@seramun.com • www.seramun.com

In-vitro diagnostic medical device

Unique device identifier

Manufacturer

Country of manufacture and date of manufacture

Article number

Serial number

Keep away from sunlight

Humidity limitation

Batch code

Consult instructions for use

Temperature limit

Do not reuse

Sufficient for *n* tests

Biohazard

Use-by date

Attention

Intended Use

SeraSpot® Anti-Borrelia-10 IgG is an IVD test for the qualitative determination of antibodies of the IgG isotype against *Borrelia burgdorferi* sensu lato specific antigens VlsE Bafz, p39 Bafz, p58 Bgar, p100 Bafz, OspC Bafz, OspC Bgar, OspC Bbur, DbpA Bafz, DbpA Bgar, DbpA Bbur in serum or plasma (citrate, EDTA, heparin) of human origin by a laboratory professional user.

The test is used in combination with the Seramun SpotSight® plate mono / Seramun SpotSight® strip device for image acquisition and the software Seramun SpotSight® scan for image analysis.

It is intended to aid in the diagnosis of Lyme disease in specimen materials from patients with suspicion of borrelia infection and as confirmatory test following a positive result.

The test must not be used with specimen materials other than serum or plasma (citrate, EDTA, heparin) of human origin, for diagnosis, monitoring, screening, prediction, prognosis, as companion diagnostic, in the near patient setting and by lay persons.

Principle of the Test

SeraSpot® Anti-Borrelia-10 IgG is a solid phase immunoassay (spot immunoassay) using recombinant or native purified antigens printed in array format (spot array) on the bottom of wells of 96-well microtiter plates. The antibodies will bind to the immobilized borrelia antigens. After incubation and wash steps bound antibodies are detected by horseradish peroxidase (HRP)-labeled antibodies against human antibodies of IgG-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). Blue spots indicate the formation of immune complexes. Spots range from pale to dark blue and are visible by eye.

Detection of Specific Antibodies is Performed in 3 Steps:

Step 1

Incubation of 1 : 101 diluted samples for 30 min at room temperature in selected wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with diluted wash buffer.

Step 2

Incubation of wells with the HRP-labeled anti-Human IgG for 30 min at room temperature. Removal of unbound conjugate antibodies by aspiration and 3 wash cycles with diluted wash buffer.

Step 3

Incubation of wells with substrate SeramunBlau® spot dark for 30 min at room temperature. The reaction is stopped by aspirating the substrate, followed by tapping the plate / strip dry onto lint-free absorbent paper. The developed arrays are to be stored light-protected until image analysis.

Test Components (Delivery Scope)

	For 48 determinations	For 96 determinations	For 2 x 96 determinations	For 10 x 96 determinations
1	<p>Microtiter plate (wells with arrays) borrelia-specific antigens and controls immobilized as spots in array layout</p> <p>6 single breakable 8-well strips in frame color coding: dark blue vacuum sealed with desiccant</p>	<p>12 single breakable 8-well strips in frame color coding: dark blue vacuum sealed with desiccant</p>	<p>2 x 12 single breakable 8-well strips in frame color coding: dark blue vacuum sealed with desiccant</p>	<p>10 x 12 single breakable 8-well strips in frame color coding: dark blue vacuum sealed with desiccant</p>
2	<p>Wash buffer Seramun® Wash buffer A 10x TRIS-based buffer</p>	<p>100 mL each concentrate for 1000 mL buffer colorless white cap</p>	<p>2 x 100 mL concentrate for each 1000 mL buffer colorless white cap</p>	<p>10 x 100 mL concentrate for each 1000 mL buffer colorless white cap</p>
3	<p>Sample diluent Seramun® Sample diluent B</p>	<p>55 mL ready-to-use solution colored red black cap</p>	<p>4 x 55 mL ready-to-use solution colored red black cap</p>	<p>20 x 55 mL ready-to-use solution colored red black cap</p>
4	<p>Conjugate anti-Human IgG-HRP conjugate (Sheep)</p>	<p>8.0 mL ready-to-use solution colored red red cap</p>	<p>2 x 8.0 mL ready-to-use solution colored red red cap</p>	<p>10 x 8.0 mL ready-to-use solution colored red red cap</p>
5	<p>Substrate SeramunBlau® spot dark 3,3',5,5' - Tetramethylbenzidine</p>	<p>8.0 mL ready-to-use solution colorless blue cap</p>	<p>2 x 8.0 mL ready-to-use solution colorless blue cap</p>	<p>10 x 8.0 mL ready-to-use solution colorless blue cap</p>
6	<p>Sealing tapes</p>	<p>2 pieces</p>	<p>4 pieces</p>	<p>20 pieces</p>
7	<p>Swab 70 % (v/v) Isopropyl alcohol</p>	<p>3 x 2 pieces</p>	<p>12 x 2 pieces</p>	<p>60 x 2 pieces</p>
8	<p>Certificate of analysis</p>	<p>1 piece</p>	<p>1 piece</p>	<p>1 piece</p>
9	<p>Instructions for use</p>	<p>1 piece</p>	<p>1 piece</p>	<p>1 piece</p>

Used Specific Antigens

Name*	Description / Characterization		Clinical relevance	Origin
VlsE Bafz	Variable major protein-like sequence Expressed	specific	Frequently detected in an early stage of infection, often IgG-answer	<i>B. afzelii</i>
p39 Bafz	Flagella complex Borrelia membrane protein A	highly specific	Frequently already detected in an early stage of infection	<i>B. afzelii</i>
p58 Bgar	not characterized	specific	Frequently detected in late stage of infection	<i>B. garinii</i>
p100 Bafz	Protein of surface membrane-vesicles	highly specific	Marker for late infection stage	<i>B. afzelii</i>
OspC (p23)**	Outer surface protein C	highly specific	Marker of an early stage of infection, often IgM-answer	<i>B. afzelii</i> , <i>B. garinii</i> , <i>B. burgdorferi sensu stricto</i>
DbpA (p18)**	Outer surface protein, Decorin binding protein A	specific	Specific marker for neuroborreliosis and lyme arthritis, IgM- and IgG-answer	<i>B. afzelii</i> , <i>B. garinii</i> , <i>B. burgdorferi sensu stricto</i>

* The nomenclature of some antigens regarding molecular mass and specificity is not consistently used in scientific literature.

** Each antigen of the used *Borrelia* species is printed separately.

Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micro-pipettes • adjustable 8-channel micro-pipette or multi-pipette • reagent container for multi-channel micro-pipettes • measuring cylinder • test tubes for sample dilution • test tube rack • washer for 96-well microtiter plates • laboratory system: Seramun SpotSight® plate mono / strip scanner with evaluation software Seramun SpotSight® scan • deionized water • lint-free absorbent paper • stop watch • collecting devices for infectious material • opaque cover (substrate reaction)

Important Information

This device is for *in-vitro* diagnostic use only. The kit may be performed by laboratory professionals only.

Follow the instructions carefully. The shelf life specified must be observed. Do not use reagents from damaged packages or bottles. Do not use reagents from other manufacturers.

Mixing of test kit components of different lots is only allowed for sample diluent, wash buffer and substrate.

All serious incidents occurring in relation with SeraSpot® Anti-Borrelia-10 IgG must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all components to room temperature before use.

For larger sample series, pipetting reagents from liquid reservoirs using a multichannel pipette is recommended to avoid time delays.

Avoid time delays when dispensing reagents.

Substrate must be protected from direct light!

Incorrect dilution of samples, incorrect washing of wells and incorrect timing may lead to erroneous results.

Air bubbles generated by forceful pipetting of samples and/or reagents may cause uneven signal formation of arrays. These arrays cannot be evaluated.

Damaged arrays, e.g. by **scratching** the bottom of the well with pipet tips or washer needles, are not suitable for evaluation.

Before taking images of the wells make sure to remove particles or fibers which may adhere to the bottom side of the wells using the swab provided in the kit!

Images of dried spots may appear more intense, which can lead to slight deviations of the measured values when scanned repeatedly. Assessment of the individual parameters is not changed.

Workplace Requirements

Processing spot immunoassays requires a clean workplace. Adhesion of fibers to the plastic surface of the wells must be avoided. Fibers may cause interferences when taking images of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	EUH208 EUH210	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request.
DIL	EUH208 EUH210 -	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request. Contains material of animal origin.
CONJ HRP IgG	Hazard components EUH208 H360D P202 P280 P308+P313 P501 -	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Restricted to professional users. May damage the unborn child. Do not handle until all safety precautions have been read and understood. Wear protective gloves / protective clothing / eye protection / face protection. IF exposed or concerned: Get medical advice / attention. Dispose of contents / container in accordance with local regulations. Contains material of animal origin.
SWAB	Hazard components H225 H319 H336 P210 P280 P304+P340 P305+P351+P338	Propan-2-ol; isopropyl alcohol; isopropanol Highly flammable liquid and vapour. Causes serious eye irritation. May cause drowsiness or dizziness. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Wear protective gloves / protective clothing / eye protection / face protection / hearing protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Limitations of the Procedure

Result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. In an early stage of infection antibodies may not be present yet or below detection limit.

Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

Interference

In rare cases, specimens may contain antibodies to BSA (bovine serum albumin) and / or AGE (advanced glycation end products), which may cause nonspecific reactions, so that the test result is read as "not available" ("n.a.") by the Seramun SpotSight® scan software.

Sample Treatment

Sample Collection

Collect serum or plasma (citrate, EDTA, heparin) of human origin in suitable sampling container.

Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 7 days. For longer periods samples have to be stored at < -15 °C. Repeated freezing and thawing of samples should be avoided!

Sample Preparation

Bring samples to room temperature before use. Homogeneity should be ensured by shaking components briefly.

Samples have to be diluted 1 : 101 (v/v) with sample diluent, e.g.: 10 µL sample and 1000 µL sample diluent.

Reagent Treatment

Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips is stable until the printed expiration date when stored at 2...8 °C. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

Reagent Preparation

Allow all components to reach room temperature prior to use. The microtiter plate is vacuum sealed with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored at 2...8 °C and protected from moisture in the original cover carefully resealed.

Dilute **wash buffer (10x)** 1 : 10 with deionized water.

Example: 10 mL Wash buffer (10x) + 90 mL deionized water. The diluted wash buffer must be thoroughly mixed before use.

The **substrate** must be protected from direct light. The test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

Assay Procedure

Performance at room temperature (RT, 18...25 °C). Bring **all** ready-to-use test reagents, diluted wash buffer and microplate to **room temperature**. Mix all reagents before use by gently shaking, avoid foaming.

The aspiration and washing steps can be carried out manually or with the aid of a microtiter plate washer.

Important notes on test procedure:

Avoid mechanical contact (**scratching**) on the bottom of the wells with pipette tips or washer needles. This will irreparably damage the array!

All liquid reagents (sample, conjugate, and substrate) have to be pipetted **without** causing **air bubbles** into the wells!


1. Provide the required number of **WELLS** according to the number of samples.
2. Preparation of working dilution of the sample: 1 : 101 dilution
e.g. 10 µL sample to 1000 µL sample diluent **DIL**
3. Pipette **100 µL** diluted sample per well.
4. Cover plate with adhesive film **COVER** and incubate for **30 min** at RT.
5. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well, Tap plate dry onto lint-free absorbent paper.
6. **Add 50 µL** of the conjugate **CONJ HRP IgG** (anti-Human IgG-HRP) to each well.
7. Cover plate with adhesive film **COVER** and incubate for **30 min** at RT.
8. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well, Tap plate dry onto lint-free absorbent paper.
9. **Add 50 µL** of the substrate **SUBSTR** to each well.
10. Cover plate and protect from light. Incubate for **30 min** at RT.
11. Remove liquid. Tap plate dry onto lint-free absorbent paper.
12. Clean bottom of wells with swab **SWAB** shortly before scanning the images.
13. **Take images** using the Seramun SpotSight® plate mono / strip und evaluate the results with the Seramun SpotSight® scan software.

If image acquisition is performed with the scanner Seramun SpotSight® strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

After aspiration of the substrate the color of developed spots is stable for 24 h when the plate is stored protected from light.

Evaluation of Results

Arraylayout

Antigens

- 04 VlsE Bafz
- 05 p39 Bafz
- 06 p58 Bgar
- 07 p100 Bafz
- 08 OspC Bafz
- 09 OspC Bgar
- 10 OspC Bbur
- 11 DbpA Bafz
- 12 DbpA Bgar
- 13 DbpA Bbur

Controls

- 01 Positive control (PC)
 - 02 Cut-off control (CO)
 - 03 Negative control (NC)
 - 14 IgG conjugate control (GC)
 - 15 Serum control (SC)
- P Well position marker

Qualitative Evaluation

Validity criteria for the test

SeraSpot® Anti-Borrelia-10 IgG test includes the following control spots:

1. Positive control (PC). Intense spot, stained darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than cut-off control
4. Serum control (SC). Intense spot, always stained in presence of sample. Absence of spot indicates absence of sample.
5. IgG conjugate control (GC). Intense spot. Serves as antibody isotype control.

The test cannot be evaluated if one of the validity criteria listed in points 1 to 5 is not met.

If the above-mentioned validity criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the validity criteria contact the manufacturer.

Interpretation of Results

Evaluation of the test has to be performed by using the Seramun SpotSight® plate mono / strip and in combination with the Seramun SeraSpot® scan software.

The results are interpreted as follows:

Evaluation	Condition
Positive	Color intensity of 2 or more antigen spots* > cut-off control
Borderline	Color intensity of VlsE Bafz spot > cut-off control
Negative	Color intensity of max. 1 antigen spot* (except VlsE Bafz spot) > cut-off control

* Several OspC spots count as 1 spot / several DbpA spots count as 1 spot

Performance Characteristics

Precision

Samples with known antibody reactivity were examined in the SeraSpot® Anti-Borrelia-10 IgG and RAWs were determined. These values were then used to determine the coefficients of variation (CV) as a measure of the precision within a test run (intra-assay CV), between different test runs (inter-assay CV) and between different test lots (lot-to-lot CV).

Antigen	Intra-assay coefficient of variation		Inter-assay coefficient of variation		Lot-to-lot coefficient of variation	
	RAW n=8	CV* [%]	RAW n=48	CV* [%]	RAW n=144	CV* [%]
VlsE Bafz	153.3	1.5	151.8	1.7	154.5	2.0
p39 Bafz	114.2	2.7	111.6	4.9	117.4	6.8
p58 Bgar	143.0	2.8	145.5	3.0	154.4	6.0
p100 Bafz	53.3	5.4	50.5	8.2	55.8	8.9
OspC Bafz	34.9	4.9	34.2	5.4	35.5	4.0
DbpA Bafz	159.9	1.2	165.3	2.4	167.4	3.0
Procedure	1 operator 8 determinations 1 batch		1 operator 8 determinations 2 testings per day 3 days 1 batch		1 operator 8 determinations 2 testings per day 3 days 3 batches	

Determination of the Cut-off Value

The cut-off range is determined specifically for each test.

Interfering Substances

All tested substances do not have a significant impact on the test results even if present in elevated concentrations in the serum: 20 mg/dL bilirubin C and F each (simulation of icteric samples), 500 mg/dL hemoglobin (simulation of hemolytic samples), 1000 mg/dL lipids (triglycerides) (simulation of lipemic samples).

Rheumatoid factors may affect test results in concentrations of 500 IU/mL or more.

Sensitivity

Pre-characterized positive patient samples (n=370) were analyzed to determine the sensitivity in comparison to reference test. Samples with discrepant results were re-tested in a second reference test.

IgG (n = 370)		Reference test	
		positive	negative
SeraSpot®	positive	222	3
	negative	4	141
Samples total		226	144

Sensitivity*: 98.2 %

*Borderline samples were considered positive.

Specificity

Blood donor samples (n=150) were examined for the specificity in comparison to reference test. Samples with discrepant results were re-tested in a second reference test.

IgG (n = 150)		Reference test	
		positive	negative
SeraSpot®	positive	4	2
	negative	1	143
Samples total		5	145

Specificity*: 98.6 %

*Borderline samples were considered positive.

Cross reactivity

Potential cross-reacting samples were tested in the SeraSpot® Anti-Borrelia-10 IgG and re-tested in a reference test. Samples with discrepant results were re-tested in a second reference test. Results obtained with the second reference test were used for the final status determination of the discrepant samples.

EBV (n = 23)		Reference test	
		positive	negative
SeraSpot®	positive	1	0
	negative	0	22
Samples total		1	22

Specificity*: 100 %

Rheumatoid factor (n = 30)		Reference test	
		positive	negative
SeraSpot®	positive	0	0
	negative	0	30
Samples total		0	30

Specificity*: 100 %

ANA (n = 51)		Reference test	
		positive	negative
SeraSpot®	positive	4	0
	negative	0	47
Samples total		4	47

Specificity*: 100 %

<i>Treponema pallidum</i> (n = 48)		Reference test	
		positive	negative
SeraSpot®	positive	2	0
	negative	0	46
Samples total		2	46

Specificity*: 100 %

Pregnancy (n = 50)		Reference test	
		positive	negative
SeraSpot®	positive	0	0
	negative	3	47
Samples total		3	47
Specificity*:		100 %	

*Borderline samples were considered positive.

Application

Automated Workflow

Manually diluted positive samples were assayed side by side in SeraSpot® Anti-Borrelia-10 IgG, by hand and by use of an automatic microplate processor. Data were generated with the devices DS2® (Dyner Technologies) and Thunderbolt® (GSD). The correlation was calculated with $R^2 > 0.9$ for the antigens (VlsE Bafz, p58 Bgar, p100 Bafz, OspC Bafz, DbpA Bafz). Assay procedure using an automatic microplate processor such as DS2® is possible. R^2 was not determined if more than 90 % of the total number (n=96) of samples tested showed no antibody reactivity for the respective antigens.

Change History

Version	Section	Modifications
2022-10_v01_DE_EN	Entire document	Use of new IFU template Updating of the intended use Conversion of the subsections Insertion of the safety instructions Removal of pipetting scheme

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