SeraSpot[®] Anti-Yersinia-6 IgA

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SeraSpot[®] Anti-Yersinia-6 IgA

Spot immunoassay for the qualitative detection of IgA antibodies against 6 *Yersinia enterocolitica* specific antigens in serum or plasma of human origin

REF	SP-005-6 A-S6		48
REF	SP-005-6 A-S12		96
REF	SP-005-6 A-S24		2 x 96
IVD	In-vitro-diagnostic medical device		CE

Seramun Diagnostica GmbH • Sprenhagener Str. 1 • 15754 Heidesee • Germany •
T +49 33767 791-10 • info@seramun.com • www.seramun.com

In-vitro diagnostic medical device

Unique device identifier

Manufacturer

Country of manufacture and date of manufacture

Article number

Serial number

Keep away from sunlight

Humidity limitation

Batch code

Consult instructions for use

Temperature limit

Do not reuse

Sufficient for *n* tests

Biohazard

Use-by date

Attention

Intended Use

SeraSpot® Anti-Yersinia-6 IgA is an IVD test for the qualitative determination of antibodies of the IgA isotype against *Yersinia enterocolitica* specific antigens Yop (Yersinia outer proteins) B, Yop D, Yop E, Yop H, Yop M and Yop N in serum or plasma (citrate, EDTA, heparin) of human by a laboratory professional user. Detection of anti-*Yersinia enterocolitica* IgA antibodies from synovial fluid of human origin is possible.

The test is used in combination with the Seramun SpotSight® plate mono / Seramun SpotSight® strip device for image acquisition and the software Seramun SpotSight® scan for image analysis.

It is intended to be used for screening, as a confirmatory assay following a positive initial test result, and aid in the diagnosis of yersiniosis in specimen materials from patients with suspicion of a *Yersinia enterocolitica* infection.

The test must not be used with specimen materials other than serum, plasma (citrate, EDTA, heparin) or synovial fluid of human origin, for diagnosis, monitoring, prediction, prognosis, as companion diagnostic, in the near patient setting and by lay persons.

Principle of the Test

SeraSpot® Anti-Yersinia-6 IgA is a solid phase immunoassay (spot immunoassay) using recombinant *Yersinia* antigens printed in array format (spot array) on the bottom of wells of 96-well microtiter plates. Antibodies will bind to the immobilized Yop-antigens. After incubation and wash steps bound antibodies are detected by horseradish peroxidase (HRP)-labeled antibodies against human antibodies of the IgA-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). Blue spots indicate the formation of immune complexes. Spots range from pale to dark blue and are visible by eye.

Detection of Specific Antibodies is Performed in 3 Steps:

Step 1

Incubation of 1 : 51 diluted samples for 30 min at room temperature in selected wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with diluted wash buffer.

Step 2

Incubation of wells with the HRP-labeled anti-Human IgA for 30 min at room temperature. Removal of unbound conjugate antibodies by aspiration and 3 wash cycles with diluted wash buffer.

Step 3

Incubation of wells with substrate SeramunBlau® spot dark for 30 min at room temperature. The reaction is stopped by aspirating the substrate, followed by tapping the plate / strip dry onto lint-free absorbent paper. The developed arrays are to be stored light-protected until image analysis.

Test Components (Delivery Scope)

			For 48 determinations	For 96 determinations	For 2 x 96 determinations
1	WELLS	Microtiter plate (wells with arrays) Yersinia-specific antigens and controls immobilized as spots in array layout	6 single breakable 8-well strips in frame color coding: light-brown vacuum sealed with desiccant	12 single breakable 8-well strips in frame color coding: light-brown vacuum sealed with desiccant	2 x 12 single breakable 8-well strips in frame color coding: light-brown vacuum sealed with desiccant
2	WASHBUF (10x)	Wash buffer (10x) Seramun® Wash buffer A 10x TRIS-based buffer	100 mL concentrated for 1000 mL buffer colorless white cap	100 mL concentrated for 1000 mL buffer colorless white cap	2 x 100 mL concentrated for each 1000 mL buffer colorless white cap
3	DIL	Sample diluent Seramun® Sample diluent B	55 mL ready-to-use solution colored red black cap	2 x 55 mL ready-to-use solution colored red black cap	4 x 55 mL ready-to-use solution colored red black cap
4	CONJ HRP IgA	Conjugate anti-Human IgA-HRP conjugate (Sheep)	8.0 mL ready-to-use solution colored purple purple cap	8.0 mL ready-to-use solution colored purple purple cap	2 x 8.0 mL ready-to-use solution colored purple purple cap
5	SUBSTR	Substrate SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine and Hydrogen peroxide	8.0 mL ready-to-use solution colorless blue cap	8.0 mL ready-to-use solution colorless blue cap	2 x 8.0 mL ready-to-use solution colorless blue cap
6	COVER	Covering film	2 pieces	2 pieces	4 pieces
7	SWAB	Swab 70 % (v/v) Isopropyl alcohol	3 x 2 pieces	6 x 2 pieces	12 x 2 pieces
8		Certificate of Analysis	1 piece	1 piece	1 piece
9		Instructions for Use	1 piece	1 piece	1 piece

Antigens Used

Name	Characterization / Function	Clinical relevance
Yop B	Part of the translocon	Yersinia outer proteins (Yop) are structure and effector proteins of the type III secretion system of human pathogenic <i>Yersinia</i> species. All Yops listed are highly specific for human pathogenic <i>Yersinia</i> species.
Yop D	Part of the translocon, Regulation of Yop expression	
Yop E	Inhibition of phagocytosis	
Yop H	Inhibition of phagocytosis	
Yop M	Activation of eukaryotic kinases	
Yop N	Regulation of Yop release	

Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micro-pipettes • adjustable 8-channel micro-pipette or multi-pipette • reagent container for multi-channel micro-pipettes • measuring cylinder • test tubes for sample dilution • test tube rack • washer for 96-well microtiter plates • laboratory system: Seramun SpotSight® plate mono / strip scanner with evaluation software Seramun SpotSight® scan • deionized water • lint-free absorbent paper • stop watch • collecting devices for infectious material • opaque cover (substrate reaction)

Important Information

This device is for *in-vitro* diagnostic use only. The kit may be performed by health professionals only.

Follow the instructions carefully. The shelf life specified must be observed. Do not use reagents from damaged packages or bottles. Do not use reagents from other manufacturers.

Mixing of test kit components of different lots is only allowed for sample diluent, wash buffer and substrate.

All serious incidents occurring in relation with SeraSpot® Anti-Yersinia-6 IgA must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all components to room temperature before use.

For larger sample series, pipetting reagents from liquid reservoirs using a multichannel pipette is recommended to avoid time delays.

Avoid time delays when dispensing reagents.

Substrate must be protected from direct light!

Incorrect dilution of samples, incorrect washing of wells and incorrect timing may lead to erroneous results.

Air bubbles generated by forceful pipetting of samples and/or reagents may cause uneven signal formation of arrays. These arrays cannot be evaluated.

Damaged arrays, e.g., by **scratching** the bottom of the well with pipet tips or washer needles, are not suitable for evaluation.

Before taking images of the wells make sure to remove particles or fibers which may adhere to the bottom side of the wells using the swab provided in the kit.

Images of dried spots may appear more intense, which can lead to slight deviations of the measured values when scanned repeatedly. Assessment of the individual parameters is not changed.

Workplace Requirements

Processing spot immunoassays requires a clean workplace. Adhesion of fibers to the plastic surface of the wells must be avoided. Fibers may cause interferences when taking images of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	EUH208 EUH210	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request.
DIL	EUH208 EUH210 -	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request. Contains material of animal origin.
CONJ HRP IgA	Hazard components EUH208 H360D P202 P280 P308+P313 P501 -	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Restricted to professional users. May damage the unborn child. Do not handle until all safety precautions have been read and understood. Wear protective gloves / protective clothing / eye protection / face protection. IF exposed or concerned: Get medical advice / attention. Dispose of contents / container in accordance with local regulations. Contains material of animal origin.
SWAB	Hazard components H225 H319 H336 P210 P280 P304+P340 P305+P351+P338	Propan-2-ol; isopropyl alcohol; isopropanol Highly flammable liquid and vapour. Causes serious eye irritation. May cause drowsiness or dizziness. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Wear protective gloves / protective clothing / eye protection / face protection / hearing protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Limitations of the Procedure

Result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. In an early stage of infection antibodies may not be present yet or below detection limit.

Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

Interference

In rare cases, specimens may contain antibodies to BSA (bovine serum albumin) and / or AGE (advanced glycation end products), which may cause nonspecific reactions, so that the test result is read as "not available" ("n.a.") by the Seramun SpotSight® scan software.

Sample Treatment

Sample Collection

Collect serum or plasma (citrate, EDTA, heparin) of human origin in suitable sampling container.

Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 7 days. For longer periods samples have to be stored at < -15 °C. Repeated freezing and thawing of samples should be avoided!

Sample Preparation

Bring samples to room temperature before use. Homogeneity should be ensured by shaking components briefly.

Samples have to be diluted 1 : 51 (v/v) with sample diluent.

Example: 10 µL sample and 500 µL sample diluent.

Reagent Treatment

Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips is stable until the printed expiration date when stored at 2...8 °C. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

Reagent Preparation

Allow all components to reach **room temperature** prior to use. The **microtiter plate** is vacuum sealed with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored at 2...8 °C and protected from moisture in the original cover carefully resealed.

Dilute **wash buffer (10x)** 1 : 10 with deionized water.

For example: 10 mL wash buffer (10x) + 90 mL deionized water. The prepared wash buffer must be thoroughly mixed before use.

The **substrate** must be protected from direct light. Test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

Assay Procedure

Performance at room temperature (RT, 18...25 °C). Bring **all** ready-to-use test reagents, prepared wash buffer and microplate to **room temperature**. Mix all reagents before use by gently shaking, avoiding foaming.

The aspiration and washing steps can be carried out manually or with the aid of a microtiter plate washer.

Important Notes on Test Procedure:

Avoid mechanical contact (**scratching**) on the bottom of the wells with pipette tips or washer needles. This will irreparably damage the array!

All liquid reagents (sample, conjugate and substrate) have to be pipetted **without** causing **air bubbles** into the wells!


1. Provide the required number of wells **WELLS** according to the number of samples.
2. Preparation of working dilution of the sample: 1 : 51 dilution
e.g. 10 µL sample to 500 µL sample diluent **DIL**
3. Pipette **100 µL** diluted sample per well.
4. Cover plate with covering film **COVER** and incubate for **30 min** at RT.
5. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,
Tap plate dry onto lint-free absorbent paper.
6. **Add 50 µL** of the conjugate **CONJ HRP IgA** (anti-Human IgA-HRP) to each well.
7. Cover plate with covering film **COVER** and incubate for **30 min** at RT.
8. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,
Tap plate dry onto lint-free absorbent paper.
9. **Add 50 µL** of the substrate **SUBSTR** to each well.
10. Cover plate and protect from light. Incubate for **30 min** at RT.
11. Remove liquid. Tap plate dry onto lint-free absorbent paper.
12. Clean bottom of wells with swab **SWAB** shortly before scanning the images.
13. **Take images** using the Seramun SpotSight® plate mono / strip und evaluate the results with the Seramun SpotSight® scan software.

If image acquisition is performed with the scanner Seramun SpotSight® strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

After aspiration of the substrate the color of developed spots is stable for 24 h when the plate is stored protected from light.

Evaluation of Results

Arraylayout

Antigens

- 04 Yop B
- 05 Yop D
- 06 Yop E
- 07 Yop H
- 08 Yop M
- 09 Yop N

Controls

- 01 Positive control (PC)
 - 02 Cut-off control (CO)
 - 03 Negative control (NC)
 - 10 IgA conjugate control (AC)
 - 11 Serum control (SC)
- P Well position marker

Qualitative Evaluation

Validity Criteria for the Test

SeraSpot® Anti-Yersinia-6 IgA includes the following control spots:

1. Positive control (PC). Intense spot, stained darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than cut-off control
4. Serum control (SC). Intense stained spot, always stained in presence of sample. Absence of spot indicates absence of sample.
5. IgA conjugate control (AC). Intensively stained spot. Serving as antibody isotype control.

The test cannot be evaluated if one of the validity criteria listed in points 1. to 5. is not met.

If the above mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the validity criteria contact the manufacturer.

Interpretation of Results

Evaluation of the test has to be performed by using the Seramun SpotSight® plate mono / strip and in combination with the Seramun SeraSpot® scan software.

Results are interpreted as follows:

Evaluation	Conditions
Positive	Color intensity of Yop D > cut-off control or Color intensity of at least two antigen spots > cut-off control
Borderline	Color intensity of one antigen spot except Yop D \geq cut-off control
Negative	Color intensity of the antigen spots except < cut-off control and Color intensity of Yop D \leq cut-off control

Performance Characteristics

Precision

Samples with known antibody reactivity were examined in the SeraSpot® Anti-Yersinia-6 IgA and the color intensities of spots (RAW) were determined. These values were then used to determine the coefficients of variation (CV) as a measure of the precision within a test run (intra-assay CV), between different test runs (inter-assay CV) and between different test lots (lot-to-lot CV).

Antigen	Intra-assay coefficient of variation		Inter-assay coefficient of variation		Lot-to-lot coefficient of variation	
	\bar{x} RAW n = 8	CV [%]	\bar{x} RAW n = 48	CV [%]	\bar{x} RAW n = 144	CV [%]
Yop B	22.4	2.6	22.4	9.5	23.6	8.6
Yop D	121.3	2.6	130.5	12.1	133.4	11.7
Yop M	26.8	6.0	28.3	6.2	28.4	4.9
Procedure	1 operator 8 determination 1 batch		1 operator 8 determinations 2 testings per day 3 days 1 batch		1 operator 8 determinations 2 testings per day 3 days 3 batches	

Determination of the Cut-off Value

The cut-off range is determined specifically for each test.

Interfering Substances

All tested substances do not have a significant impact on the test results up to present concentrations in the serum of 20 mg/dL bilirubin C and F each (simulation of icteric samples), 500 mg/dL hemoglobin (simulation of hemolytic samples), 1000 mg/dL lipids (triglycerides) (simulation of lipemic samples) and 500 IU/mL rheumatoid factor.

Sensitivity

Pre-characterized samples (n=218) were analyzed to determine the sensitivity in comparison to reference tests of SeraSpot® Anti-Yersinia-6 IgA. Samples with discrepant results were retested in a second reference test.

IgA		Reference test	
		positive	negative
SeraSpot®	positive	119	21
	negative	6	72
Samples total		125	93

Sensitivity*: 95.2 %

*Borderline samples were considered positive.

Specificity

Blood donor samples (n=205) were examined for the specificity in comparison to reference tests of SeraSpot® Anti-Yersinia-6 IgA. Samples with discrepant results were retested in a second reference test.

IgA		Reference test	
		positive	negative
SeraSpot®	positive	22	8
	negative	0	175
Proben gesamt		22	183
Specificity*:		95.6 %	

*Borderline samples were considered positive.

Application

Automated Workflow

Manually diluted positive samples were assayed side by side in SeraSpot® Anti-Yersinia-6 IgA, by hand and by use of an automatic microplate processor. Data were generated with the device DS2® (Dynex Technologies). The correlation was calculated with $R^2 > 0.9$ for the antigens Yop B, Yop D, and Yop M. Assay procedure using an automatic microplate processor such as DS2® is possible. R^2 was not determined if more than 90 % of the total number (n=96) of samples tested showed no antibody reactivity for the respective antigens.

Synovial Fluid

Detection of anti-*Yersinia enterocolitica* antibodies in synovial fluid using SeraSpot® Anti-Yersinia-6 IgA is possible. In the performance evaluation test, 24 joint punctate samples were examined in parallel, in the SeraSpot® test and a commercial line immunoassay. The agreement was 95.8 % for the detection of IgA. Due to the small number of samples, no general recommendation for result interpretation can be given. This should be worked out by the laboratory, taking into account additional data and clinical findings.

Changes History

Version	Section	Modifications
2022-10_v02_DE_EN	Important Information	Insertion note on dealing with serious incidents
	Interpretation of Results	Adaption of evaluation conditions

References

1. Cambronne E. D. and Schneewind O., Yersinia enterocolitica Type III Secretion: yscM1 and yscM2 Regulate yop Gene Expression by a Posttranscriptional Mechanism That Targets the S'Untranslated Region of yop mRNA. J. Bacteriol. 2002, 184 (21), 5880-5893
2. Cremer J., Putzker M., Faulde M., Zöller L., Immunoblotting of Yersinia plasmid-encoded released proteins: A tool for serodiagnosis. Electrophoresis 1993, 14, 952-959
3. Kendrick C.J., Baker B., Morris A.J., O'Toole P.W., Identification of Yersinia-infected blood donors by anti-Yop IgG immunoassay. Transfusion, 2001, 41, 1365-1372
4. Kist M., Ackermann A., Autenrieth I.B., von Eichel-Streiber Chr., Frick J., Fruth A., et al., MIQ 09: Gastrointestinale Infektionen. 2nd Edition, Qualitätsstandards in der mikrobiologisch-infektiologischen Diagnostik, 2013
5. Rastawicki W., Humoral response to selected antigens of Yersinia enterocolitica and Yersinia pseudotuberculosis in the course of yersiniosis in humans. I. Occurrence of antibodies to Yersinia lipopolysacharydes and Yop proteins by ELISA, Med Dosw Mikrobiol 2006; 58 (4), 303-319