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SeraSpot® Anti-Borrelia-10 IgG

SeraSpot® Anti-Borrelia-10 IgM

Spotimmunoassay for detection of IgG or IgM antibodies to
Borrelia burgdorferi sensu lato in human serum or plasma

REF SP-006-10 G-S6 ▽ 48 REF SP-006-10 G-S12 ▽ 96 REF SP-006-10 G-S24 ▽ 2x 96
REF SP-006-10 G-S120 ▽ 10x 96 IVD *In-vitro*-diagnostic medical device CE
REF SP-006-10 M-S6 ▽ 48 REF SP-006-10 M-S12 ▽ 96 REF SP-006-10 M-S24 ▽ 2x 96
REF SP-006-10 M-S120 ▽ 10x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 33767 791-10 · fax +49 33767 791-99 · info@seramun.com

Introduction

Lyme borreliosis is a systemic infectious disease with various clinical manifestations. The disease can be divided into 3 stages involving different organ systems. Cardinal symptoms of the first stage (4 to 8 weeks p. i.) are the erythema chronicum migrans and local or general lymphadenopathies (lymphadenosis cutis benigna). The clinical manifestations of the second stage (1 to 12 months p. i.) range from meningitis, meningopolyneuritis, encephalitis up to hemiparesis, myoneuralgia and arthralgia, especially knee joint arthritis. Cardiac manifestations like life-threatening myocarditis / pancarditis are rarely detected. The third stage (months to years p. i.) is characterized by chronic infection of the nervous system (neuroborreliosis, progressive encephalomyelitis), the skin (acrodermatitis chronica atrophicans) and the joints (chronic erosive arthritis). Especially late manifestations of Lyme borreliosis can reduce quality of life and antibiotic treatment is difficult. Therefore an early diagnosis of borrelia infections is of great importance.

The infectious agent of Lyme disease was isolated from ticks in 1982 by BURGENDORFER. It belongs to spirochetes and represents the independent genus *Borrelia*. The tick species *Ixodes ricinus* is the main vector of *Borrelia burgdorferi* in Europe. The occurrence of ticks within densely wooded rural areas correlates with the frequency of Lyme disease cases with a peak in the summer and autumn months. While the tick is taking blood from the host the *Borrelia* are transmitted from the contents of the tick's intestine. Persons at risk for *Borrelia* infections are people who work in forestry or agriculture, hunters and campers, soldiers and all people who occasionally spend time in wooded areas. The incidence of Lyme disease reaches from 2 to 40 / 100 000 in middle Germany to 300 / 100 000 in Austria. The estimated number of unrecognized cases is presumably higher, because a tick bite is not always recognized and the typical erythema chronicum migrans only develops in about 50 % of *Borrelia* infections.

The direct detection of the infectious agent via culture from blood, cerebral spinal fluid (CSF) or skin biopsies undoubtedly represents the safest method for diagnosis of Lyme disease. Unfortunately this method is hampered by the long reproduction time of the *Borrelia*, the complexity of the culture medium

and the relatively low sensitivity in case of culture from patient material. Thus culture as well as Borrelia detection from tissue sections via immunohistological methods are not suitable for routine diagnosis. Methods of choice are detection of IgG and/or IgM antibodies via immunofluorescence or enzyme immunoassay (ELISA). The ELISAs use Borrelia sonicates, partially purified or recombinant antigens for plate coating. Because of the regionally differing occurrence of subtypes (genospecies) and the known variability of the cell surface proteins of *Borrelia burgdorferi*, antigen mixtures are often preferred. Confirmatory tests like western blot, dot blot or line assay are essential and the serological test results have to be interpreted together with the clinical picture. Negative test results do not exclude Lyme disease (early stage of infection, seronegative cases).

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Intended use

The *SeraSpot*[®] Anti-*Borrelia*-10 IgG / *SeraSpot*[®] Anti-*Borrelia*-10 IgM test is an *in vitro* diagnostic medical device (spotimmunoassay, SIA) for the detection of *Borrelia*-specific IgG or IgM antibodies in human serum or plasma samples.

Principle of the test

The *SeraSpot*[®] Anti-*Borrelia*-10 IgG / *SeraSpot*[®] Anti-*Borrelia*-10 IgM test is a solid phase immunoassay based on the use of recombinant *Borrelia* proteins as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96well-microtiter plates. The recombinant antigens serve as capture molecules for antibodies against *Borrelia* antigens. Bound antibodies are detected by horseradish peroxidase-(HRP)-labeled antibodies against human antibodies of IgG- or IgM-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed by precipitated product from colorless substrate solution. Pale blue to dark blue spots are visible by eye.

The test procedure for antibody detection is performed in 3 steps:

Step 1

Incubation of 1 : 101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with wash solution.

Step 2

Incubation of wells with the anti-human IgG- or IgM-HRP-conjugate for 30 minutes at room temperature. After incubation the antibody is aspirated and unbound components are removed by 3 wash cycles with wash solution.

Step 3

Incubation of wells with substrate solution SeramunBlau[®] spot dark for 30 minutes at room temperature. The reaction is stopped by aspirating the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the Seramun *SpotSight*[®] plate or Seramun *SpotSight*[®] strip scanner; the obtained images are interpreted by using the Software Seramun *SpotSight*[®] scan. Taking account of the test-specific evaluation criteria, the test result is determined.

Used *Borrelia* antigens

Nomenclature*	Properties / Specificity	Relevance	Origin	
VisE (7, 8, 9, 10, 13)	variable major protein-like sequence Expressed	specific	frequently already detectable during early stage, frequently IgG- answer	<i>B. afzelii</i>
p39 (BmpA) (1, 11, 14, 16)	Flagella complex <i>Borrelia</i> membrane protein A	high specific	frequently already detectable during early stage	<i>B. afzelii</i>
p58 (1, 11)	not characterized	specific	frequently in stage 3	<i>B. garinii</i>
p100 (1, 11, 15)	Protein of surface membrane- vesicles	high specific	marker for late infection stage	<i>B. afzelii</i>
p23 (OspC)** (1, 6, 11, 12, 17)	Outer surface protein C	high specific	marker for early stage, frequently IgM-answer	<i>B. afzelii</i> , <i>B. garinii</i> , <i>B.</i> <i>gurgdorferi sensu stricto</i>
p18 (DbpA)** (10, 11, 13)	Outer surface protein, Decorin binding protein A	specific	marker especially for neuroborreliosis and lyme arthritis, IgM- and IgG-answer	<i>B. afzelii</i> , <i>B. garinii</i> , <i>B.</i> <i>burgdorferi sensu stricto</i>

* The nomenclature for some antigens is not consistently used in the literature concerning molecular mass and specificity.

** Each antigen species is printed separately.

Kit components

	For 48 determinations	For 96 determinations	For 2x 96 determinations	For 10x 96 determinations
1	WELLS Wells with arrays 6 single breakable 8-well strips in frame vacuum-sealed with desiccant	12 single breakable 8-well strips in frame vacuum-sealed with desiccant	2x 12 single breakable 8-well strips in frame vacuum-sealed with desiccant	10x 12 single breakable 8-well strips in frame vacuum-sealed with desiccant
InG detection: color coding light blue				
2	WASHBUF CONC 10X Wash buffer Seramun® wash buffer A 10-fold	100 ml concentrate transparent bottle, white cap	2x 100 ml concentrate transparent bottle white cap	10x 100 ml concentrate transparent bottle white cap
3	DIL Sample diluent Seramun® sample diluent B	55 ml ready-to-use solution, colored red transparent bottle black cap	4x 55 ml ready-to-use solution, colored red, transparent bottle, black cap	20x 55 ml ready-to-use solution, colored red, transparent bottle, black cap
4	CONJ HRP IgG Anti-human IgG- HRP-conjugate or Anti-human IgM- HRP-conjugate	8 ml ready-to-use solution, colored red, transparent bottle, red cap	2x 8 ml ready-to-use solution, colored red, transparent bottle, red cap	10x 8 ml ready-to-use solution, colored red, transparent bottle, red cap
5	SUBSTR TMB Substrate SeramunBlau® spot dark 3,3'-5,5'- Tetramethylbenzidine and hydrogen peroxide	8 ml ready-to-use solution, black bottle, blue cap	2x 8 ml ready-to-use solution, colored green, transparent bottle, green cap	10x 8 ml ready-to-use solution, colored green, transparent bottle, green cap
6	COVER Adhesive film	2x	4x	20x
7	SWAB Swab with 70 % (v/v) Isopropyl alcohol	3x2	12x2	60x2

Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 7 days. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided! Bring the samples to room temperature and mix well.

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (10 µl sample and 1000 µl diluent) with the ready-to-use sample diluent.

Materials required but not provided

• Adjustable single-channel micro-pipettes incl. pipette tips • Adjustable multi-channel micro-pipettes incl. pipette tips • Reagent container for multi-channel micro-pipettes • Distilled or deionized water • lint-free absorbent paper • Test tubes (1 ml) for sample dilution • Beaker and measuring cylinder • stop watch • Washer for 96-well microtiter plates (optional) • Collecting devices for infectious material • Seramun *SpotSight*[®] plate or Seramun *SpotSight*[®] strip scanner with evaluation software Seramun *SpotSight*[®] scan

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 1x 48, 1x 96, 2x 96 or 10x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used up to 4 weeks when stored at 2...8 °C and up to 2 weeks when stored at room temperature.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash buffer by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or deionized water. For Example: 10 ml **WASHBUF CONC 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

Workplace requirements

Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

Assay procedure

Performance at room temperature (18...25 °C) and at the specified incubation times.

Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere to the underside of the bottom of the wells using the swab **SWAB** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

Working steps


1. Warm all reagents and required wells **WELL** to room temperature (RT) before use. Mix gently without causing foam.
2. Add **100 µl** of 1:101 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 1000 µl sample diluent **DIL**.
3. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **CONJ HRP IgG** or **CONJ HRP IgM** to each well.
6. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **SUBSTR TMB** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent.
11. Clean bottom of wells with swab **SWAB**.
12. Perform image acquisition using the Seramun *SpotSight*[®] plate or the Seramun *SpotSight*[®] strip scanner and evaluate the results with the Seramun *SpotSight*[®] scan software. If image acquisition is performed with the scanner Seramun *SpotSight*[®] strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

Aspiration and washing steps can be performed manually or via washer for 96-well microtiter plates.

After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when stored the plate in the dark.

Test evaluation

Array layout

Parameter	Controls
04 VlsE (AFZ)	01 Positive control (PC)
05 p39 (AFZ)	02 Cut-off control (CO)
06 p58 (GAR)	03 Negative control (NC)
07 p100 (AFZ)	14 IgG conjugate control (GC)
08 OspC (AFZ)	15 IgM conjugate control (MC)
09 OspC (GAR)	16 Serum control (SC)
10 OspC (BUR)	
11 DbpA (AFZ)	
12 DbpA (GAR)	
13 DbpA (BUR)	
	P Well position marker

Validity criteria for the test

The *SeraSpot*[®] Anti-Borrelia-10 IgG / *SeraSpot*[®] Anti-Borrelia-10 IgM test includes the following control spots:

1. Positive control (PC). Intensively stained spot, darker than Cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than Cut-off control
4. IgG, IgM conjugate control (GC, MC). Intensively stained spots with different position for IgG and IgM antibody detection. Serving as antibody isotype control.
5. Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.

The test cannot be evaluated if one of the validity criteria listed in points 1 to 5 is not fulfilled.

Result interpretation

Evaluation of the test has to be performed by using the Seramun *SpotSight*[®] scan software.

Judgment of spot staining is performed in accordance with the MIQ guidelines (12) as follows:

Result Interpretation	IgG	IgM
negative	max. 1 antigen spot * (except VlsE-Spot) > Cut-off control	no antigen spot > Cut-off control
borderline	VlsE-Spot > Cut-off control	1 antigen spot* (except OspC*) > Cut-off control
positive	2 or more antigen spots * > Cut-off control	OspC* spot(s) or 2 other antigen spots* > Cut-off control

* Several OspC spots count as 1 spot / several DbpA spots count as 1 spot

Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

Interference

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) do not interfere with the test. Rheumatoid factors may influence the results.

In rare cases, specimens may contain antibodies to BSA (bovine serum albumin) and/or AGE (advanced glycation endproducts), which may cause nonspecific reactions, so that the test result is evaluated as "not available" ("n.a.") by the Seramun *SpotSight*[®] scan software.

Assay procedure

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere to the underside of the bottom of the wells or to the array area cause erroneous results.

Performance characteristics

Diagnostic sensitivity

1. Investigation of serologically pre-determined samples

Serologically pre-determined samples (n = 370 (IgG determination), n = 370 (IgM determination), reference test 1) of patients with suspect of Borreliosis were tested in the *SeraSpot*[®] Anti-Borreliia-10 IgG / *SeraSpot*[®] Anti-Borreliia-10 IgM test ("Initial sensitivity"). Samples with discrepant results were retested in a second assay for detection of anti-*Borreliia* antibodies. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results ("Amended sensitivity").

IgG	n = 370	Reference test 1		
		positive (n = 205)	borderline (n = 7)	negative (n = 158)
<i>SeraSpot</i> [®] Anti-Borreliia-10 IgG	positive	193	4	11
	borderline	7	1	9
	negative	5	2	138

Borderline samples were valued positive.

Initial sensitivity: **96.7 %**

Amended sensitivity: **98.2 %**

IgM	n = 370	Reference test 1		
		positive (n = 160)	borderline (n = 2)	negative (n = 208)
<i>SeraSpot</i> [®] Anti-Borreliia-10 IgM	positive	151	0	22
	borderline	0	1	8
	negative	9	1	178

Borderline samples were valued positive.

Initial sensitivity: **93.8 %**

Amended sensitivity: **97.2 %**

Diagnostic specificity

A total of 150 predetermined blood donors (Reference test 1) were tested in the *SeraSpot*[®] Anti-Borrelia-10 IgG / *SeraSpot*[®] Anti-Borrelia-10 IgM test. Samples with discrepant results were retested and assessed as described in "Diagnostic sensitivity".

IgG	n = 150	Reference test 1		
		positive (n = 4)	borderline (n = 0)	negative (n = 146)
<i>SeraSpot</i> [®] Anti-Borrelia-10 IgG	positive	3	0	2
	borderline	0	0	1
	negative	1	0	143

Borderline samples were valued positive.

Initial specificity: **97.9 %**

Amended specificity: **98.6 %**

IgM	n = 150	Reference test 1		
		positive (n = 4)	borderline (n = 0)	negative (n = 146)
<i>SeraSpot</i> [®] Anti-Borrelia-10 IgM	positive	4	0	0
	borderline	0	0	0
	negative	0	0	146

Borderline samples were valued positive.

Specificity: **100.0 %**

Cross-reactivity

Potential cross-reacting samples were tested in the *SeraSpot*[®] Anti-*Borrelia*-10 IgG / *SeraSpot*[®] Anti-*Borrelia*-10 IgM test, and retested in a reference test for detection of anti-*Borrelia* antibodies. Samples with discrepant results were retested and assessed as described in "Diagnostic sensitivity".

SeraSpot[®]	EBV (n = 23)	IgG, Reference test 1			IgM, Reference test 1		
		positive (n = 1)	borderline (n = 0)	negative (n = 22)	positive (n = 2)	borderline (n = 0)	negative (n = 21)
	positive	1	0	0	1	0	0
	borderline	0	0	0	0	0	1
	negative	0	0	22	1	0	20
	Borderline samples were valued positive.						
Specificity: 100.0 %				Initial Specificity: 95.2 %			
				Amended Specificity: 100.0 %			
SeraSpot[®]	Rheumatic factor (n = 30)	IgG, Reference test 1			IgM, Reference test 1		
		positive (n = 0)	borderline (n = 0)	negative (n = 30)	positive (n = 0)	borderline (n = 0)	negative (n = 30)
	positive	0	0	0	0	0	0
	borderline	0	0	0	0	0	0
	negative	0	0	30	0	0	30
	Borderline samples were valued positive.						
Specificity: 100.0 %				Specificity: 100.0 %			
SeraSpot[®]	Autoimmune (n = 51)	IgG, Reference test 1			IgM, Reference test 1		
		positive (n = 3)	borderline (n = 0)	negative (n = 48)	positive (n = 2)	borderline (n = 0)	negative (n = 49)
	positive	3	0	0	2	0	1
	borderline	0	0	1	0	0	2
	negative	0	0	47	0	0	46
	Borderline samples were valued positive.						
Initial Specificity: 97.9 %				Initial Specificity: 93.9 %			
Amended Specificity: 100.0 %				Amended Specificity: 95.8 %			
SeraSpot[®]	Syphilis (n = 48)	IgG, Reference test 1			IgM, Reference test 1		
		positive (n = 2)	borderline (n = 0)	negative (n = 46)	positive (n = 2)	borderline (n = 0)	negative (n = 46)
	positive	2	0	0	2	0	1
	borderline	0	0	0	0	0	1
	negative	0	0	46	0	0	44
	Borderline samples were valued positive.						
Specificity: 100.0 %				Initial Specificity: 95.7 %			
				Amended Specificity: 95.7 %			

Investigation of sera from pregnant women

SeraSpot®	Pregnancy (n = 50)	IgG, Reference test 1			IgM, Reference test 1		
		positive (n = 3)	borderline (n = 0)	negative (n = 47)	positive (n = 2)	borderline (n = 3)	negative (n = 45)
	positive	0	0	0	2	0	1
	borderline	0	0	0	0	0	4
	negative	3	0	47	0	3	40
Borderline samples were valued positive.							
Specificity: 100.0 %			Specificity: 88.9 %				
Amended Specificity: 88.9 %							

Precision

Samples of known antibody titer were assayed by SeraSpot® Anti-Borrelia-10 IgG / SeraSpot® Anti-Borrelia-10 IgM test. Color intensity (RAW) of spots was recorded and the coefficients of variation (CV) were calculated.

IgG detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	RAW [\bar{x}] n = 8	CV* [%]	RAW [\bar{x}] n = 48	CV* [%]	RAW [\bar{x}] n = 144	CV* [%]
VisE (AFZ)	153.3	1.5	151.8	1.7	154.5	2.0
p39 (AFZ)	114.2	2.7	111.6	4.9	117.4	6.8
p58 (GAR)	143.0	2.8	145.5	3.0	154.4	6.0
p100 (AFZ)	53.3	5.4	50.5	8.2	55.8	8.9
OspC (AFZ)	34.9	4.9	34.2	5.4	35.5	4.0
DbpA (AFZ)	159.9	1.2	165.3	2.4	167.4	3.0
Procedure	1 operator, 8x determination, 1 batch		1 operator, 8x repeat determinations, 2x testing per day, 3 days, 1 batch		1 operator, 8x repeat determinations, 2x testing per day, 3 days, 3 batches	

CV*: The CV for RAW-values < 20 has not been determined.

AFZ: *B. afzelii*; GAR: *B. garinii*

IgM detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	RAW [\bar{x}] n = 8	CV* [%]	RAW [\bar{x}] n = 48	CV* [%]	RAW [\bar{x}] n = 144	CV* [%]
OspC (AFZ)	103.5	5.8	108.5	2.1	107.1	2.6
OspC (GAR)	105.4	3.5	107.8	3.8	108.1	2.2
OspC (BUR)	109.2	3.8	106.7	4.8	108.5	2.5
DbpA (AFZ)	16.4	7.4	16.2	8.6	16.0	7.4
Procedure	1 operator, 8x determination, 1 batch		1 operator, 8x repeat determinations, 2x testing per day, 3 days, 1 batch		1 operator, 8x repeat determinations, 2x testing per day, 3 days, 3 batches	

CV*: The CV for RAW-values < 20 has not been determined.

AFZ: *B. afzelii*; GAR: *B. garinii*; BUR: *B. burgdorferi* s.s.

Automation

Samples (n = 96) of known antibody titer (IgG and IgM) were assayed manually or by use of the microplate processors DS2® (Dynex Technologies; manual sample dilution) and Thunderbolt® (GSD; processor-performed sample dilution). Color intensity (RAW) of spots was recorded and used for calculation of the coefficient of determination r^2 for selected antibody reactivities and the different methods of operation.

Ig-isotype detection	Dynex DS2® vs. manual processing		GSD Thunderbolt® vs. manual processing		Dynex DS2® vs. GSD Thunderbolt® processing	
	r^2 IgG	r^2 IgM	r^2 IgG	r^2 IgM	r^2 IgG	r^2 IgM
VlsE (AFZ)	0.973	0.952	0.963	0.923	0.968	0.910
p58 (GAR)	0.971		0.947		0.971	
p100 (AFZ)	0.965		0.948		0.954	
OspC (AFZ)	0.955	0.974	0.950	0.908	0.975	0.927
OspC (GAR)		0.971		0.907		0.740
OspC (BUR)		0.966		0.934		0.923
DbpA (GAR)	0.988		0.922		0.967	

r^2 was not determined if more than 90 % of the total number (n = 96) of the samples showed no antibody reactivity for the particular antigens.

AFZ: *B. afzelii*; GAR: *B. garinii*; BUR: *B. burgdorferi* s.s.

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

Do not use or mix reagents from different lots except for sample diluent, wash buffer and TMB substrate solution. These can be used across batches and parameters.

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

History of changes

Version	Section	Modifications
2019-11-11	General	General revision
	Kit components	Adding of S120 version
	Preparation and storage of samples	Change of stability of samples
	Preparation and storage of reagents	Expiry of ready-to-use wash solution
2022-02-01	Result interpretation	„Starter kit“ SP-000-1 removed
	Limitations of the procedure	Interference supplemented

Incubation scheme *SeraSpot*[®] Anti-Borrelia-10 IgG / IgM

- 

100 µl
30 min

3 x wash

diluted sample (1 : 101)
incubation (room temperature)
with wash solution, each 400 µl per well
- 

50 µl
30 min

3 x wash

conjugate **CONJ HRP IgG** or **CONJ HRP IgM**
incubation (room temperature)
with wash solution, each 400 µl per well
- 

50 µl
30 min

aspiration

substrate **SUBSTR TMB**
incubation (room temperature, protected from light)
- imaging of wells and image analysis

scanner Seramun *SpotSight*[®] plate / Seramun *SpotSight*[®] strip and software Seramun *SpotSight*[®] scan

 Manufacturer	 Date of manufacture	 Use by	LOT Batch code	REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

SeramunBlau®, SeraSpot® und Seramun SpotSight® are registered trademarks of Seramun Diagnostica GmbH.