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# SeraSpot<sup>®</sup> Anti-Helicobacter-6 IgG

# SeraSpot<sup>®</sup> Anti-Helicobacter-6 IgA

Spotimmunoassay for detection of IgG or IgA antibodies to *Helicobacter pylori* in human serum or plasma

<b>REF</b> SP-007-6 G-S6 ▾ 48	<b>REF</b> SP-007-6 G-S12 ▾ 96	<b>REF</b> SP-007-6 G-S24 ▾ 2x 96
<b>IVD</b> <i>In-vitro</i> -diagnostic device <b>CE</b>		
<b>REF</b> SP-007-6 A-S6 ▾ 48	<b>REF</b> SP-007-6 A-S12 ▾ 96	<b>REF</b> SP-007-6 A-S24 ▾ 2x 96
<b>IVD</b> <i>In-vitro</i> -diagnostic device <b>CE</b>		

Seramun Diagnostica GmbH · Sprenghagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
 phone +49 33767 791-10 · fax +49 33767 791-99 · [info@seramun.com](mailto:info@seramun.com)

## Introduction

*Helicobacter pylori* (*H. pylori*) is a gram-negative, spiral shaped, human pathogenic bacterium which has been found in the stomachs of humans in all parts of the world. In developing countries, 70 to 90% of the population carries *H. pylori* mostly already acquired before the age of 10 years. In developed countries the prevalence of infection ranges from 25 to 50%. *H. pylori* infections are transmitted via the oral-oral, the faecal-oral route or iatrogenic. An animal reservoir for *H. pylori* aside from the human stomach has not been confirmed so far (1).

*H. pylori* is characterized by a strong urease activity and some strains (type I) additionally produce toxins like VacA (vacuolating toxin) and CagA (cytotoxin-associated antigen A) (2, 3). Further studies have identified HcpC (Helicobacter cysteine-rich protein C) and GroEL (Chaperonin) as new independent virulence factors (4, 5), whereas GroEL is more important than HcpC.

These extracellular products contribute to the pathogenesis by direct damage of the gastric epithelium accompanied by a chronic inflammation with enhanced levels of inflammation mediators like interleukin 8, interleukin 1 or tumor necrosis factor alpha. Patients infected with *H. pylori* often develop autoantibodies to parietal cells. The key-role of these autoantibodies in the development of atrophic gastritis as precursor of gastric cancer has been reviewed (1, 6).

About 10% of the infected persons develop *H. pylori* associated gastritis and secondary diseases like chronic active gastritis, Ulcus ventriculi, Ulcus duodeni, gastric cancer and MALT lymphoma. Patients suffering from gastric malignomas and Ulcus duodeni are infected with *H. pylori* in nearly 100% of the cases. About 80% of patients with chronic gastritis, 70% of patients with Ulcus ventriculi and 60% of patients with gastric cancer are infected with *H. pylori* (6).

Diagnosis is usually based on gastroendoscopy combined with the detection of the pathogen in biopsy material by culture, histology, molecular assays or rapid urease test. Culture from biopsy material is difficult and not always successful. The <sup>13</sup>C breath test which is often performed as follow up test, detects CO<sub>2</sub> released by the bacterial urease from radioactive labeled urea in patient's breath. This

method is not invasive but the need for special equipment and the uptake of radioactive urea by the patients are disadvantageous.

Immunological tests on the basis of specific anti-*H. pylori* antibodies (mono- or polyclonal) are available. They enable the direct detection of *H. pylori* antigens from faecal specimens and may be used for therapeutic surveillance.

Serology is the third method commonly used as a non-invasive method to diagnose *H. pylori* infection and is concerning sensitivity and specificity comparable to other methods. Furthermore, serology is the only test that is not affected by proton pump inhibitor treatment that could lead to false-negative results in other tests.

Invasive as well as non-invasive tests for detection of *H. pylori* infections contain methodological shortcomings. Therefore, the S2k-guideline (8) recommend to include at least two positive test results to ensure a safe diagnosis of an infection with *H. pylori*.

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## Intended use

The **SeraSpot® Anti-Helicobacter-6 IgG / SeraSpot® Anti-Helicobacter-6 IgA test** is an *in vitro* diagnostic device (spotimmunoassay, SIA) for the detection of *Helicobacter*-specific IgG or IgA antibodies in human serum or plasma. It can be used as screening assay as well as confirmatory assay after a positive ELISA result.

## Principle of the test

The **SeraSpot® Anti-Helicobacter-6 IgG / SeraSpot® Anti-Helicobacter-6 IgA test** is a solid phase immunoassay based on the use of recombinant *Helicobacter* proteins as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96-well microtiter plates. The recombinant antigens serve as capture molecules for antibodies against *Helicobacter* antigens. Bound antibodies are detected by horseradish peroxidase-(HRP)-labeled antibodies against human antibodies of IgG- or IgA-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed by precipitated product from colorless substrate solution. Color intensity is correlated to the antibody concentration. Pale blue to dark blue spots are visible by eye.

**The test procedure for antibody detection is performed in 3 steps:**

### Step 1

Incubation of 1 : 101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with wash solution.

### Step 2

Incubation of wells with the HRP-labeled anti-human IgG- or IgA for 30 minutes at room temperature. After incubation the antibody is aspirated and unbound components are removed by 3 wash cycles with wash solution.

### Step 3

Incubation of wells with substrate solution SeramunBlau® spot dark for 30 minutes at room temperature. The reaction is stopped by aspirating the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the Seramun *SpotSight*® plate or Seramun *SpotSight*® strip scanner; the obtained images are interpreted by using the Software Seramun *SpotSight*® scan. Taking account of the test-specific evaluation criteria, the test result is determined.

## Used Helicobacter antigens

Nomenclature	Characterization	Incidence*	
CagA	120-145 kD	cytotoxin-associated antigen A	type I (high virulent)
gGT	59-61 kD	γ-glutamyl transpeptidase	type II (virulent)
GroEL	54-62 kD	chaperone	type II (virulent)
UreA	26 kD	urease (subunit A)	type II (virulent)
VacA	87-97 kD	vacuolating cytotoxin	type I (high virulent)
HcpC	28-32 kD	Helicobacter cysteine-rich protein C	type II (virulent)

\*strain typing only in IgG-detection possible

## Kit components

		For 48 determinations	For 96 determinations	For 2x 96 determinations	
1	<b>WELLS</b>	<b>Wells with arrays</b>	<b>6 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>12 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>2x 12- single breakable 8-well strips in frame</b> separately vacuum-sealed with desiccant
<b>IgG detection: color coding green; IgA detection: color coding light green</b>					
2	<b>WASHBUF CONC 10X</b>	<b>Wash buffer</b> Seramun® Wash buffer A	<b>100 ml concentrate</b> transparent bottle, white cap	<b>100 ml concentrate</b> transparent bottle, white cap	<b>2x 100 ml concentrate</b> transparent bottle, white cap
3	<b>DIL</b>	<b>Sample diluent</b> Seramun® Sample diluent B	<b>55 ml</b> ready-to-use solution, colored red transparent bottle, black cap	<b>2x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cap	<b>4x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cap
4	<b>CONJ HRP IgG</b>	<b>Anti-human IgG-HRP-conjugate</b>	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap	<b>2x 8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap
	<b>CONJ HRP IgA</b>	<b>Anti-human IgA-HRP-conjugate</b>	<b>8 ml</b> ready-to-use solution, colored purple, transparent bottle, purple cap	<b>8 ml</b> ready-to-use solution, colored purple, transparent bottle, purple cap	<b>2x 8 ml</b> ready-to-use solution, colored purple, transparent bottle, purple cap
5	<b>SUBSTR TMB</b>	<b>Substrate</b> SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	<b>8 ml</b> ready-to-use solution, black bottle, blue cap	<b>8 ml</b> ready-to-use solution, black bottle, blue cap	<b>2x 8 ml</b> ready-to-use solution, black bottle, blue cap
6	<b>COVER</b>	<b>Adhesive film</b>	<b>2x</b>	<b>2x</b>	<b>4x</b>
7	<b>SWAB</b>	<b>Swab</b> with 70 % (v/v) Isopropyl alcohol	<b>3x2</b>	<b>6x2</b>	<b>12x2</b>

## Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 48 hours. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided!

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (e.g. 10 µl sample and 1000 µl sample diluent) with the ready-to-use sample diluent.

## Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- stop watch
- Washer for 96-well microtiter plates (optional)
- Collecting devices for infectious material
- Seramun *SpotSight*® plate or Seramun *SpotSight*® strip scanner with evaluation software Seramun *SpotSight*® scan

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used up to 4 weeks when stored at 2...8 °C and up to 2 weeks when stored at room temperature.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtiter plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or de-ionized water. For Example: 10 ml **WASHBUF CONC 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

## Workplace requirements

**Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.**

## Assay procedure

Performance at room temperature (18...25 °C) and at the specified incubation times.

For the handling of *SeraSpot*<sup>®</sup> tests the use of the starter kit (order number: SP-000-1) is recommended. The kit is designed for multiple use.

### Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere to the underside of the bottom of the wells using the swab **[SWAB]** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

### Working steps


1. Warm all reagents to room temperature before use. Mix gently without causing foam.
2. Add **100 µl** of 1 : 101 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 1000 µl sample dilution buffer.
3. Cover plate with adhesive film **[COVER]** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **[WASHBUF CONC 10X]**). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **[CONJ HRP IgG]** or **[CONJ HRP IgA]** to each well.
6. Cover plate with adhesive film **[COVER]** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard conjugate), then wash wells **3 x** using **400 µl** wash solution per well (made of **[WASHBUF CONC 10X]**). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **[SUBSTR TMB]** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent paper.
11. Clean bottom of wells with Swab **[SWAB]**.
12. Perform image acquisition using the Seramun *SpotSight*<sup>®</sup> plate or the Seramun *SpotSight*<sup>®</sup> strip scanner and evaluate the results with the Seramun *SpotSight*<sup>®</sup> scan software. If image acquisition is performed with the scanner Seramun *SpotSight*<sup>®</sup> strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

**Aspiration and washing steps can be performed manually or via washer for 96-well microtiter plates.**

**After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when the plate is stored in the dark.**

## Test evaluation

### Array layout

Parameter	Controls
04 CagA	01 Positive control (PC)
05 gGT	02 Cut-off control (CO)
06 GroEL	03 Negative control (NC)
07 UreA	10 IgG conjugate control (GC)
08 VacA	11 IgA conjugate control (AC)
09 HcpC	12 Serum control (SC)
	P Well position marker

### Validity criteria for the test

The *SeraSpot*<sup>®</sup> Anti-Helicobacter-6 IgG / *SeraSpot*<sup>®</sup> Anti-Helicobacter-6 IgA test includes the following control spots:

1. Positive control (PC). Intensively stained spot, darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than Cut-off control
4. IgG, IgA conjugate control (GC, AC). Intensively stained spots with different position for IgG and IgA antibody detection. Serving as antibody isotype control.
5. Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.

The test cannot be evaluated if one of the validity criteria listed in points 1 to 5 is not fulfilled.

### Result interpretation

Evaluation of the test has to be performed by using the Seramun *SpotSight*<sup>®</sup> scan software. A visual plausibility check is possible with the aid of the magnifying lens contained in the starter kit (order number: SP-000-1).

Judgment of spot staining is performed according to the following classification:

Result Interpretation	IgG	IgA
<b>negative</b>	no antigen spot $\geq$ Cut-off control	no antigen spot $\geq$ Cut-off control
<b>borderline</b>	UreA- or HcpC- or gGT spot $\geq$ Cut-off control	UreA- or HcpC- or gGT spot $\geq$ Cut-off control
<b>positive</b>	CagA- or VacA- or GroEL spot or two or more antigen spots $\geq$ Cut-off control	CagA- or VacA- or GroEL spot or two or more antigen spots $\geq$ Cut-off control

## **Limitations of the procedure**

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

### **Interference**

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) and rheumatoid factors (500 IU/ml) do not interfere with the test.

### **Assay procedure**

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere to the underside of the bottom of the wells or to the array area cause erroneous results.

## Performance characteristics

### Diagnostic sensitivity and specificity

Serologically pre-determined samples (reference test 1) were tested in the *SeraSpot*<sup>®</sup> Anti-Helicobacter-6 IgG / *SeraSpot*<sup>®</sup> Anti-Helicobacter-6 IgA test ("Initial sensitivity" / "Initial specificity"). Samples with discrepant results were retested in a second assay for detection of *Helicobacter pylori* antibodies. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results ("Amended sensitivity" / "Amended specificity").

IgG		Reference test 1		
		positive (n = 65)	borderline (n = 7)	negative (n = 48)
<i>SeraSpot</i> <sup>®</sup> Anti-Helicobacter-6 IgG	positive	65	2	1
	borderline	0	0	0
	negative	0	5	47

Borderline samples were valued positive.

Initial sensitivity: **93.1** %

Amended sensitivity: **100.0** %

Initial specificity: **97.9** %

Amended specificity: **100.0** %

IgA		Reference test 1		
		positive (n = 30)	borderline (n = 5)	negative (n = 85)
<i>SeraSpot</i> <sup>®</sup> Anti-Helicobacter-6 IgA	positive	28	2	14
	borderline	0	2	4
	negative	2	1	65

Borderline samples were valued positive.

Initial sensitivity: **91.4** %

Amended sensitivity: **95.9** %

Initial specificity: **78.3** %

Amended specificity: **94.2** %

## Precision

Samples of known antibody titer were assayed by SeraSpot® Anti-Helicobacter-6 IgG / SeraSpot® Anti-Helicobacter-6 IgA test. Color intensity of antigen spots divided by the color intensity of the Cut-off spot (Ratio) was recorded and the coefficients of variation (CV) were calculated.

IgG detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	Ratio [ $\bar{x}$ ] n = 40	CV [%]	Ratio [ $\bar{x}$ ] n = 80	CV [%]	Ratio [ $\bar{x}$ ] n = 240	CV [%]
<b>CagA</b>	3.67	<b>8.09</b>	3.90	<b>7.91</b>	4.03	<b>10.80</b>
<b>gGT</b>	4.04	<b>8.19</b>	2.70	<b>11.07</b>	2.75	<b>12.60</b>
<b>GroEL</b>	3.43	<b>8.04</b>	3.13	<b>11.83</b>	3.25	<b>11.13</b>
<b>UreA</b>	1.45	<b>7.67</b>	1.15	<b>10.05</b>	1.18	<b>11.70</b>
<b>VacA</b>	1.89	<b>7.41</b>	1.55	<b>10.54</b>	1.54	<b>9.71</b>
<b>HcpC</b>	3.91	<b>5.12</b>	3.09	<b>7.28</b>	3.19	<b>10.74</b>
<b>Procedure</b>	1 operator, 40x determination, 1 batch		5 operator, 2 laboratories 2x determination, 2x testing per day, 20 days,1 batch		5 operator, 2 laboratories 2x determination, 2x testing per day, 20 days, 3 batches	

IgA detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	Ratio [ $\bar{x}$ ] n = 40	CV [%]	Ratio [ $\bar{x}$ ] n = 80	CV [%]	Ratio [ $\bar{x}$ ] n = 240	CV [%]
<b>CagA</b>	2.92	<b>8.91</b>	2.52	<b>8.86</b>	2.61	<b>11.68</b>
<b>gGT</b>	1.00	<b>9.34</b>	1.28	<b>14.98</b>	1.31	<b>14.42</b>
<b>GroEL</b>	2.88	<b>9.78</b>	2.54	<b>9.22</b>	2.66	<b>13.96</b>
<b>UreA</b>	1.86	<b>9.59</b>	2.31	<b>10.49</b>	2.34	<b>14.11</b>
<b>VacA</b>	3.60	<b>8.48</b>	3.29	<b>13.74</b>	3.40	<b>16.70</b>
<b>HcpC</b>	1.54	<b>8.55</b>	1.63	<b>6.84</b>	1.67	<b>9.76</b>
<b>Procedure</b>	1 operator, 8x determination, 1 batch		1 operator, 2 laboratories 8x determination, 2x testing per day, 3 days,1 batch		1 operator, 2 laboratories 8x determination, 2x testing per day, 3 days, 3 batches	

## Automation

Samples (n = 96) of known antibody titer (IgG and IgA) were assayed manually or by use of the microplate processor DS2® (Dy nex Technologies; manual sample dilution). Color intensity of antigen spots divided by the color intensity of the Cut-off spot (Ratio) was recorded and used for calculation of the coefficient of determination  $r^2$  for selected antibody reactivities and the different methods of operation.

Ig isotype detection	Dy nex DS2® vs. manual processing	
	$r^2$ IgG	$r^2$ IgA
<b>CagA</b>	<b>0.958</b>	<b>0.956</b>
<b>gGT</b>	<b>0.930</b>	
<b>GroEL</b>	<b>0.924</b>	<b>0.954</b>
<b>UreA</b>	<b>0.965</b>	<b>0.929</b>
<b>VacA</b>	<b>0.967</b>	<b>0.961</b>
<b>HcpC</b>	<b>0.941</b>	<b>0.958</b>

$r^2$  was not determined if more than 90 % of the total number (n = 96) of the samples showed no antibody reactivity for the particular antigens.

## Common advices and precautions

**This kit is for *in-vitro* use only.** Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

**Do not use or mix reagents from different lots except for sample diluent, washing buffer and TMB substrate solution. These can be used across batches and parameters.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. The reagents must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

## History of changes

Version	Section	Modifications
2019-09-09	General Preparation and storage of reagents	General revision Storage of ready-to-use wash solution

## Incubation scheme *SeraSpot*<sup>®</sup> Anti-Helicobacter-6 IgG / IgA

- |    |   |   |
|----|---|---|
| 1. | <br>100 µl<br>30 min | diluted sample (1 : 101)<br>incubation (room temperature)   |
|    | <br>3 x wash         | with wash solution, each 400 µl per well  |
| 2. | <br>50 µl<br>30 min  | conjugate <b>CONJ HRP IgG</b> or <b>CONJ HRP IgA</b><br>incubation (room temperature)   |
|    | <br>3 x wash         | with wash solution, each 400 µl per well  |
| 3. | <br>50 µl<br>30 min  | substrate <b>SUBSTR TMB</b><br>incubation (room temperature, protected from light)  |
|    | <br>aspiration       |   |
| 4. | imaging of wells and image analysis   | scanner Seramun <i>SpotSight</i> <sup>®</sup> plate / Seramun <i>SpotSight</i> <sup>®</sup> strip and software Seramun <i>SpotSight</i> <sup>®</sup> scan |

 Manufacturer	 Date of manufacture	 Use by	<b>LOT</b> Batch code	<b>REF</b> Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	<b>IVD</b> <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

SeramunBlau<sup>®</sup>, *SeraSpot*<sup>®</sup> und Seramun *SpotSight*<sup>®</sup> are registered trademarks of Seramun Diagnostica GmbH.