SeraSpot[®] Anti-Helicobacter-6 IgG

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SeraSpot[®] Anti-Helicobacter-6 IgG

Spot immunoassay for the qualitative detection of IgG antibodies against 6 *Helicobacter pylori* specific antigens in serum or plasma of human origin

REF	SP-007-6 G-S6		48
REF	SP-007-6 G-S12		96
REF	SP-007-6 G-S24		2 x 96
IVD	In-vitro-diagnostic medical device		CE

Seramun Diagnostica GmbH • Sprenhagener Str. 1 • 15754 Heidesee • Germany •
T +49 33767 791-10 • info@seramun.com • www.seramun.com

In-vitro diagnostic medical device

Unique device identifier

Manufacturer

Country of manufacture and date of manufacture

Article number

Serial number

Keep away from sunlight

Humidity limitation

Batch code

Consult instructions for use

Temperature limit

Do not reuse

Sufficient for n tests

Biohazard

Use-by date

Attention

Intended Use

SeraSpot® Anti-*Helicobacter*-6 IgG is an IVD test for the qualitative determination of antibodies of the IgG isotype against *Helicobacter pylori* specific antigens cytotoxin-associated antigen A (CagA), γ -glutamyl-transpeptidase (gGT), GroEL, urease subunit A (UreA), vacuolating cytotoxin (VacA) and *Helicobacter* cysteine-rich protein C (HcpC) in serum or plasma (citrate, EDTA, heparin) of human origin by a laboratory professional user.

The test is used in combination with the Seramun SpotSight® plate mono / Seramun SpotSight® strip device for image acquisition and the software Seramun SpotSight® scan for image analysis.

It is intended to aid in the diagnosis of a past or persistent *Helicobacter pylori* infection in specimen materials from patients with symptoms of *H. pylori* associated gastritis.

The test must not be used with specimen materials other than serum or plasma (citrate, EDTA, heparin) of human origin, for diagnosis, monitoring, screening (especially for gastric carcinomas), prediction, prognosis, as companion diagnostic, in the near patient setting and by lay persons.

Principle of the Test

SeraSpot® Anti-*Helicobacter*-6 IgG is a solid phase immunoassay (spot immunoassay) using recombinant *Helicobacter pylori* antigens printed in array format (spot array) on the bottom of wells of 96-well microtiter plates. Antibodies will bind to the immobilized *Helicobacter pylori* antigens. After incubation and wash steps bound antibodies are detected by horseradish peroxidase (HRP)-labeled antibodies against human antibodies of the IgG-isotype by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). Blue spots indicate the formation of immune complexes. Spots range from pale to dark blue and are visible by eye.

Detection of Specific Antibodies is performed in 3 Steps:

Step 1

Incubation of 1 : 101 diluted samples for 30 min at room temperature in selected wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with diluted wash buffer.

Step 2

Incubation of wells with the HRP-labeled anti-Human IgG for 30 min at room temperature. Removal of unbound conjugate antibodies by aspiration and 3 wash cycles with diluted wash buffer.

Step 3

Incubation of wells with substrate SeramunBlau® spot dark for 30 min at room temperature. The reaction is stopped by aspirating the substrate, followed by tapping the plate / strip dry onto lint-free absorbent paper. The developed arrays are to be stored light-protected until image analysis.

Test Components (Delivery Scope)

		For 48 determinations	For 96 determinations	For 2 x 96 determinations	
1	WELLS	Microtiter plate (wells with arrays) <i>Helicobacter pylori</i> specific antigens and controls immobilized as spots in array layout	6 single breakable 8-well strips in frame color coding: green vacuum sealed with desiccant	12 single breakable 8-well strips in frame color coding: green vacuum sealed with desiccant	2 x 12 single breakable 8-well strips in frame color coding: green vacuum sealed with desiccant
2	WASHBUF (10x)	Wash buffer (10x) Seramun® Wash buffer A 10x TRIS-based buffer	100 mL concentrated for 1000 mL buffer colorless white cap	100 mL concentrated for 1000 mL buffer colorless white cap	2 x 100 mL concentrated for 1000 mL buffer each colorless white cap
3	DIL	Sample diluent Seramun® Sample diluent B	55 mL ready-to-use solution colored red black cap	2 x 55 mL ready-to-use solution colored red black cap	4 x 55 mL ready-to-use solution colored red black cap
4	CONJ HRP IgG	Conjugate anti-Human IgG-HRP conjugate (sheep)	8.0 mL ready-to-use solution colored red red cap	8.0 mL ready-to-use solution colored red red cap	2 x 8.0 mL ready-to-use solution colored red red cap
5	SUBSTR	Substrate SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine	8.0 mL ready-to-use solution colorless blue cap	8.0 mL ready-to-use solution colorless blue cap	2 x 8.0 mL ready-to-use solution colorless blue cap
6	COVER	Covering film	2 pieces	2 pieces	4 pieces
7	SWAB	Swab 70 % (v/v) Isopropyl alcohol	3 x 2 pieces	6 x 2 pieces	12 x 2 pieces
8		Certificate of Analysis	1 piece	1 piece	1 piece
9		Instructions for Use	1 piece	1 piece	1 piece

Antigens Used

Name	Description		Occurrence*
CagA	120-145 kDa	cytotoxin-associated antigen A	type I (highly virulent)
gGT	59-61 kDa	γ-glutamyl transpeptidase	type II (virulent)
GroEL	54-62 kDa	chaperone	type II (virulent)
UreA	26 kDa	urease (subunit A)	type II (virulent)
VacA	87-97 kDa	vacuolating cytotoxin	type I (highly virulent)
HcpC	28-32 kDa	<i>Helicobacter</i> cysteine-rich protein C	type II (virulent)

*Strain-typing possible when detecting IgG only

Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel micro-pipettes • adjustable 8-channel micro-pipette or multi-pipette • reagent container for multi-channel micro-pipettes • measuring cylinder • test tubes for sample dilution • test tube rack • washer for 96-well microtiter plates • laboratory system: Seramun SpotSight® plate mono / strip scanner with evaluation software Seramun SpotSight® scan • deionized water • lint-free absorbent paper • stop watch • collecting devices for infectious material • opaque cover (substrate reaction)

Important Information

This device is for *in-vitro* diagnostic use only. The kit may be performed by laboratory professional user only.

Follow the instructions carefully. The shelf life specified must be observed. Do not use reagents from damaged packages or bottles. Do not use reagents from other manufacturers.

Mixing of test kit components of different lots is only allowed for sample diluent, wash buffer and substrate.

All serious incidents occurring in relation with SeraSpot® Anti-Helicobacter-6 IgG must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all components to room temperature before use.

For larger sample series, pipetting reagents from liquid reservoirs using a multichannel pipette is recommended to avoid time delays.

Avoid time delays when dispensing reagents.

Substrate must be protected from direct light!

Incorrect dilution of samples, incorrect washing of wells and incorrect timing may lead to erroneous results.

Air bubbles generated by forceful pipetting of samples and/or reagents may cause uneven signal formation of arrays. These arrays cannot be evaluated.

Damaged arrays, e.g. by **scratching** the bottom of the well with pipet tips or washer needles, are not suitable for evaluation.

Before taking images of the wells make sure to remove particles or fibers which may adhere to the bottom side of the wells using the swab provided in the kit.

Images of dried spots may appear more intense, which can lead to slight deviations of the measured values when scanned repeatedly. Assessment of the individual parameters is not changed.

Workplace Requirements

Processing spot immunoassays requires a clean workplace. Adhesion of fibers to the plastic surface of the wells must be avoided. Fibers may cause interferences when taking images of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.

Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	EUH208 EUH210	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request.
DIL	EUH208 EUH210 -	Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Safety data sheet available on request. Contains material of animal origin.
CONJ HRP IgG	Hazard components EUH208 H360D P202 P280 P308+P313 P501 -	N-methyl-2-pyrrolidone; 1-methyl-2-pyrrolidone Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. Restricted to professional users. May damage the unborn child. Do not handle until all safety precautions have been read and understood. Wear protective gloves / protective clothing / eye protection / face protection. IF exposed or concerned: Get medical advice / attention. Dispose of contents / container in accordance with local regulations. Contains material of animal origin.
SWAB	Hazard components H225 H319 H336 P210 P280 P304+P340 P305+P351+P338	Propan-2-ol; isopropyl alcohol; isopropanol Highly flammable liquid and vapour. Causes serious eye irritation. May cause drowsiness or dizziness. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Wear protective gloves / protective clothing / eye protection / face protection / hearing protection. IF INHALED: Remove person to fresh air and keep comfortable for breathing. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Limitations of the Procedure

Result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. In an early stage of infection antibodies may not be present yet or below detection limit.

Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

Interference

In rare cases specimens may contain antibodies to BSA (bovine serum albumin) and / or AGE (advanced glycation end products), which may cause nonspecific reactions, so that the test result is read as "not available" ("n.a.") by the Seramun SpotSight® scan software.

Sample Treatment

Sample Collection

Collect serum or plasma (citrate, EDTA, heparin) of human origin in suitable sampling container.

Sample Shelf Life and Storage

Samples can be stored at 2...8 °C for a maximum of 7 days. For longer periods samples have to be stored at < -15 °C. Repeated freezing and thawing of samples should be avoided!

Sample Preparation

Bring samples to room temperature before use. Homogeneity should be ensured by shaking components briefly.

Samples have to be diluted 1 : 101 (v/v) with sample diluent.

Example: 10 µL sample and 1000 µL sample diluent.

Reagent Treatment

Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtitration strips is stable until the printed expiration date when stored at 2...8 °C. All opened test kit components are stable for up to 2 months when stored properly at 2...8 °C. The diluted wash buffer can be stored at 2...8 °C for up to 1 month.

Reagent Preparation

Allow all components to reach **room temperature** prior to use. The **microtiter plate** is vacuum sealed with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored at 2...8 °C and protected from moisture in the original cover carefully resealed.

Dilute **wash buffer (10x)** 1 : 10 with deionized water.

Example: 10 mL wash buffer (10x) + 90 mL deionized water. The prepared wash buffer must be thoroughly mixed before use.

The **substrate** must be protected from direct light. Test should not be performed with a substrate solution that is colored dark or contains colored precipitates.

Assay Procedure

Performance at room temperature (RT, 18...25 °C). Bring **all** ready-to-use test reagents, prepared wash buffer and microplate to **room temperature**. Mix all reagents before use by gently shaking, avoid foaming.

The aspiration and washing steps can be carried out manually or by means of a microtiter plate washer.

Important Notes on Test Procedure:

Avoid mechanical contact (**scratching**) on the bottom of the wells with pipette tips or washer needles. This will irreparably damage the array!

All liquid reagents (sample, conjugate and substrate) have to be pipetted **without** causing **air bubbles** in the wells.


1. Provide the required number of wells **WELLS** according to the number of samples.
2. Preparation of working dilution of the sample: 1 : 101 dilution
e.g. 10 µL sample to 1000 µL sample diluent **DIL**
3. Pipette **100 µL** diluted sample per well.
4. Cover plate with sealing tape **COVER** and incubate for **30 min** at RT.
5. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,
Tap plate dry onto lint-free absorbent paper.
6. **Add 50 µL** of the conjugate **CONJ HRP IgG** (anti-Human IgG-HRP) to each well.
7. Cover plate with sealing tape **COVER** and incubate for **30 min** at RT.
8. Remove liquid, then wash wells **3x** using **400 µL** diluted wash buffer per well,
Tap plate dry onto lint-free absorbent paper.
9. **Add 50 µL** of the substrate **SUBSTR** to each well.
10. Cover plate and protect from light. Incubate for **30 min** at RT.
11. Remove liquid. Tap plate dry onto lint-free absorbent paper.
12. Clean bottom of wells with swab **SWAB** shortly before scanning the images.
13. **Take images** using the Seramun SpotSight® plate mono / strip und evaluate the results with the Seramun SpotSight® scan software.

If image acquisition is performed with the scanner Seramun SpotSight® strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

After aspiration of the substrate the color of developed spots is stable for 24 h when the plate is stored protected from light.

Evaluation of Results

Arraylayout

Antigens

- 04 CagA
- 05 gGT
- 06 GroEL
- 07 UreA
- 08 VacA
- 09 HcpC

Controls

- 01 Positive control (PC)
- 02 Cut-off control (CO)
- 03 Negative control (NC)
- 10 IgG conjugate control (GC)
- 11 Serum control (SC)

P Well position marker

Qualitative Evaluation

Validity Criteria for the Test

SeraSpot® Anti-Helicobacter-6 IgG includes the following control spots:

1. Positive control (PC). Intense spot, stained darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot, intensity lower than cut-off control.
4. Serum control (SC). Intense spot, always stained in presence of a sample. Absence of spot indicates absence of sample.
5. IgG conjugate control (GC). Intense spot. Serves as antibody isotype control.

The test cannot be evaluated if one of the validity criteria listed in points 1. to 5. is not met.

If the above mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the validity criteria contact the manufacturer.

Interpretation of Results

Evaluation of the test has to be performed by using the Seramun SpotSight® plate mono / strip in combination with the Seramun SeraSpot® scan software.

Results are interpreted as follows:

Evaluation	Conditions
Positive	Color intensity of the antigen spots CagA or VacA or GroEL \geq cut-off control or Color intensity of two or more antigen spots \geq cut-off control
Borderline	Color intensity of the antigen spots UreA or HcpC or gGT \geq cut-off control
Negative	Color intensity of the antigen spots $<$ cut-off control

Performance Characteristics

Precision

Samples with known antibody reactivity were examined in the SeraSpot® Anti-Helicobacter-6 IgG and ratios were determined. These values were then used to determine the coefficients of variation (CV) as a measure of the precision within a test run (intra-assay CV), between different test runs (inter-assay CV) and between different test lots (lot-to-lot CV).

Antigen	Intra-assay coefficient of variation		Inter-assay coefficient of variation		Lot-to-lot coefficient of variation	
	\bar{x} Ratio n=40	CV [%]	\bar{x} Ratio n=80	CV [%]	\bar{x} Ratio n=240	CV [%]
CagA	3.67	8.09	3.90	7.91	4.03	10.80
gGT	4.04	8.19	2.70	11.07	2.75	12.60
GroEL	3.43	8.04	3.13	11.83	3.25	11.13
UreA	1.45	7.67	1.15	10.05	1.18	11.70
VacA	1.89	7.41	1.55	10.54	1.54	9.71
HcpC	3.91	5.12	3.09	7.28	3.19	10.74
Procedure	1 operator 40 determinations 1 batch		5 operators 2 determinations 2 testings per day 20 days 1 batch		5 operators 2 determinations 2 testings per day 20 days 3 batches	

Determination of the Cut-off Value

The cut-off range is determined specifically for each test.

Interfering Substances

All tested substances do not have a significant impact on the test results even if present in elevated concentrations in the serum: 20 mg/dL bilirubin C and F each (simulation of icteric samples), 500 mg/dL hemoglobin (simulation of hemolytic samples), 1000 mg/dL lipids (triglycerides) (simulation of lipemic samples) and 500 IU/mL rheumatoid factor.

Sensitivity and Specificity

Pre-characterized samples (n=120) were analyzed to determine the sensitivity and specificity of SeraSpot® Anti-Helicobacter-6 IgG in comparison with IgG reference tests. Samples with discrepant results were retested in a second reference test.

IgG		Reference test	
		positive	negative
SeraSpot®	positive	68	0
	negative	0	52
Samples total		68	52

Sensitivity*: 100 %

Specificity*: 100 %

*Borderline samples were considered positive.

Application

Automated Workflow

Manually diluted positive samples were assayed side by side in SeraSpot® Anti-Helicobacter-6 IgG, by hand and by use of an automatic microplate processor, e.g. DS2® (Dynex Technologies). For determination of the linear regression the coefficient of determination R^2 was used. The correlation was calculated with $R^2 > 0.9$ for the antigens CagA, GroEL, UreA, VacA and HcpC. Assay procedure using an automatic microplate processor such as DS2® is possible. R^2 was not determined if more than 90 % of the total number ($n=96$) of samples tested showed no antibody reactivity for the respective antigens.

Change History

Version	Section	Modifications
2022-12_v01_DE_EN	Entire document	Updating of the intended use Conversion of the subsections Insertion of safety instructions Removal of pipetting scheme
	Sample treatment	Insertion of sample shelf life and storage
	Test procedure	Starter-Kit SP-000-1 removed
	Limitations of the procedure	Update of interference
	Characteristic performance	Update of analytical and diagnostic performance

References

1. Dunn B E, Cohen H, Blaser M J, Helicobacter pylori. Clinical Microbiology Reviews, Oct 1997 p 720-741.
2. Fischbach W, Bornschein J, Hoffmann J C et al, S2k-Leitlinie "Helicobacter pylori und gastroduodenale Ulkuskrankheit" der Deutschen Gesellschaft für Gastroenterologie, Verdauungs- und Stoffwechselkrankheiten (DGVS). Z Gastroenterol 2022
3. Gao L, Weck M N, Michel A et al, Association between chronic atrophic gastritis and serum antibodies to 15 Helicobacter pylori proteins measured by multiplex serology. Cancer Res 2009; 69:2973-80.
4. Gao L, Michel A, Weck M N et al, Helicobacter pylori Infection and Gastric Cancer Risk: Evaluation of 15 H. pylori proteins determined by novel multiplex serology. Cancer Res 2009; 69:2973-80.
5. Hahn H, Falke D, Ullmann U (Hrsg.), Medizinische Mikrobiologie und Infektiologie. Springer Verlag Berlin, Heidelberg, New York, 3. Auflage 1999
6. Hocker M, Hohenberger P, Helicobacter pylori virulence factors-one part of a big picture. Lancet 2003; 362:1231-3.
7. Huang J O, Zheng G F, Sumanac K et al, Meta-analysis of the relationship between cagA seropositivity and gastric cancer. Gastroenterology 2003; 125:1636-44.
8. Malfertheiner P, Megraud F, O'Morain C A et al, Management of Helicobacter pylori infection – the Maastricht V/ Florence Consensus Report. Gut 2017; 66:6-30.