# SeraSpot® Anti-Parvovirus-6 IgG

# SeraSpot® Anti-Parvovirus-5 IgM

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#### **Literatur:**

1. Bundesministerium für Gesundheit (2010). Parvovirus B19 - Stellungnahmen des Arbeitskreises Blut des Bundesministeriums für Gesundheit, Bundesgesundheitsblatt, Springer Verlag, 53: 944-956.
2. Cossart, Y. E., Field, A. M., Cant, B., Widdows, D. (1975). Parvovirus-like particles in human sera. *Lancet*, 1(7898):72-73.
3. Anderson, M. J., Jones, S. E., Fisher-Hoch S. P., Lewis, E., Hall, S. M., Bartlett, C. L. R. Cohen, B. J. Mortimer, P. P., Pereira, M. S. (1983). Human parvovirus, the cause of erythema infectiosum (fifth disease)?, *Lancet*, 1(8338):1378.
4. Arbeitsgemeinschaft der Wissenschaftlichen Medizinischen Fachgesellschaften (AWMF) online (2014). Labordiagnostik schwangerschaftsrelevanter Virusinfektionen – Sk2-Leitlinie, AWMF Registernummer 0093/001, Kapitel C.II.7 Ringelröteln, S.161-173.
5. Schleunig (1996). Parvovirus-B19-Infektionen – Sind es nur harmlose Ringelröteln?, *Deutsches Ärzteblatt*, 93(43): A-2781-284.
6. Modrow, S., Gärtner, B. (2006). Parvovirus-B-19-Infektion in der Schwangerschaft. *Deutsches Ärzteblatt*, 103(43): A-2869-2876.
7. Servey, J. T., Reamy, B. V., Hodge, J. (2007). Clinical presentations of parvovirus B19 infection. *Am Fam Physician*, 75(3):373-376.
8. Giorgio, E., de Oronzo, M. A., Iozza, I., di Natale, A., Cianci, S., Garofalo, G., Giacobbe, A. M., Politi, S. (2010). Parvovirus B19 during pregnancy: a review.
9. Staroselsky, A., Klieger-Grossmann, C., Garcia-Bournissen, F., Koren, G. (2009). Exposure to fifth disease in pregnancy. *Can Fam Physician*. 55(12):1195-8.
10. Kaufmann, B., Simpson, A.A., Rossmann, M.G. (2004). The structure of human parvovirus B19. *PNAS*, 101(32):11628-11633
11. Zhi, N., Mills, I. P., Lu J., Wong, S., Filippone, Brown, K. E. (2006). Molecular and functional analyses of a human parvovirus B19 infectious clone demonstrates essential roles for NS1, VP1, and the 11- Kilodalton Protein in Virus Replication and Infectivity. *J Virol*, 80(12):5941-5950.

12. Heegaard, E. D.; Rasksen, C. J. & Christensen, J. (2002). Detection of parvovirus B19 NS1-specific antibodies by ELISA and western blotting employing recombinant NS1 protein as antigen. *J Med Virol*, 67(3):375-383.
13. Von Poblotski, A., Gigler, A., Lang, B., Wolf, H., Modrow, S. (1995). Antibodies to parvovirus B19 NS-1 protein in infected individuals. *J Gen Virol*, 76:519-527.
14. Hemauer A., Gigler A., Searle K., Beckenlehner K., Raab U., Broliden K., Wolf H., Enders G., Modrow S. (2000). Seroprevalence of parvovirus B19 NS1-specific IgG in B19-infected and uninfected individuals and in infected pregnant women. *J Med Virol*, 60(1):48-55.
15. Dembinski J., Eis-Hübinger A.M., Maar J., Schild R., Bartmann P. (2003). Long term follow up of serostatus after maternofetal parvovirus B19 infection. *Arch Dis Child*, 88:219-221.

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## SeraSpot® Anti-Parvovirus-6 IgG SeraSpot® Anti-Parvovirus-5 IgM

Spotimmunoassay for detection of IgG or IgM antibodies to Parvovirus B19  
in human serum or plasma

<b>REF</b> SP-012-6 G-S6 ▽ 48	<b>REF</b> SP-012-6 G-S12 ▽ 96	<b>REF</b> SP-012-6 G-S24 ▽ 2x 96
	<b>IVD</b> <i>In-vitro</i> -diagnostic device	<b>CE</b>
<b>REF</b> SP-012-5 M-S6 ▽ 48	<b>REF</b> SP-012-5 M-S12 ▽ 96	<b>REF</b> SP-012-5 M-S24 ▽ 2x 96
	<b>IVD</b> <i>In-vitro</i> -diagnostic device	<b>CE</b>

Seramun Diagnostics GmbH · Spreenhagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 33767 791-10 · fax +49 33767 791-99 · [info@seramun.com](mailto:info@seramun.com)

### Introduction

In 1975, the human parvovirus B19 was discovered by chance in a blood donor test. The designation B19 is a result of the nomenclature used in this screening: the corresponding sample was the 19th serum in row B (1, 2). In 1983, the virus was detected in patients who had flu-like symptoms or aplastic crises. Since then, the human parvovirus B19 virus has been known to cause ringworm (Erythema infectiosum, also called Quintus disease or Fifth disease) (1, 3). It has now been proven that numerous other conditions may be associated with parvovirus B19 infection. These include arthritis and arthralgia (in 8% of infected children and 80% in adults), chronic anemia (in patients with leukemia and immunosuppressed patients), congenital anemia and rare forms of myocarditis, vasculitis and glomerulonephritis (1).

Parvovirus B19 infection during pregnancy can cause complications in women without immune protection. Fetal deaths or hydrops fetalis (accumulation of fluid in compartments of the fetus) may be the result (1, 4, 5, 6, 7, 8, 9, 10). Parvovirus infections in children are 30% asymptomatic (6). Even in healthy adults, infection often has no apparent symptoms (5).

The human parvovirus B19 belongs to the family Parvoviridae. It is the only known human pathogenic parvovirus (10). The virus proliferates primarily in erythroid cells. It is a virus that has no lipid envelope and a diameter of 18 to 26 nm. The virus particles consist of 60 copies of the capsid protein and single-stranded DNA (5,596 nucleotides). The capsid protein is composed of the structural proteins VP-1 and VP-2. Portions of these structural proteins are referred to as VP-C and VP-N. In addition, the viral DNA encodes non-structural proteins. The major non-structural protein is designated NS-1 (1). This plays a role in viral replication and replication (12, 13).

Parvovirus transmission usually takes place via droplet infections, but infection routes via blood or blood products are also possible (1, 5). In the course of an infection, antibodies against the structural proteins VP-1 and VP-2 partly also form against the nonstructural protein NS-1 (11, 12, 13). IgM antibodies are present only at an early stage of infection, while parvovirus IgG antibodies are maintained throughout their lifetime (5). The prevalence rate in Germany is very high and is for example

already 35 % in children aged 4 to 6 years, 65 % in young adults (18-25 years) and 80 % in the elderly (65-75 years) (4).

Laboratory diagnostic tests for the detection of parvovirus B19 infection take place among other things in the course of pregnancy examinations in order to determine the immune status of the women. Especially pregnant women, who are in contact with children due to their job or who live in the household with children under the age of six, should have their immune status determined at an early stage of pregnancy. Seronegative pregnant women should follow strict hygiene management to avoid infection during pregnancy (possible consequences of infection: abortion, hydrodrops fetalis).

The S2k guideline published by the Society for Virology (GfV) for the laboratory diagnosis of pregnancy-related viral infections (AWMF registry number 0093/001; C.II.7 ring rubella) (4) recommends molecular biology methods (PCR) or serological methods or a combination of these for the determination of acute and / or recent parvovirus B19 infections. IgM antibodies can be detected approximately 10 days after infection. After 30 to 60 days, this immunoglobulin class is no longer present. The time of infection can be narrowed down by application of Western or immunoblots by testing antibody reactivity against linear and denatured epitopes of VP2 proteins. Furthermore, an avidity determination can give an indication of the infection period. Low avidity antibodies indicate a recent infection.

To determine immunity, the detection of parvovirus IgG antibodies is recommended. Negative parvovirus IgM detection confirms a recent infection (4).

Several test methods, e.g. ELISA tests, based on the use of recombinant VP-1 and VP-2 polypeptides, are available. The infection status can then be interpreted with a line immunoassay. These use different recombinant polypeptides derived from sequence sections of the viral proteins VP-1 and VP-2 (VP-C, VP-N, VP-2r - linear epitope, VP-2p - conformational epitope), and the protein NS-1 in order to characterize the course of the infection by detecting IgG or IgM antibodies.

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## Intended use

The *SeraSpot*<sup>®</sup> Anti-Parvovirus-6 IgG / *SeraSpot*<sup>®</sup> Anti-Parvovirus-5 IgM test is an *in vitro* diagnostic device (spotimmunoassay, SIA) for the detection of *Parvovirus*-specific IgG or IgM antibodies in human serum or plasma samples.

## Principle of the test

The *SeraSpot*<sup>®</sup> Anti-Parvovirus-6 IgG / *SeraSpot*<sup>®</sup> Anti-Parvovirus-5 IgM test is a solid phase immunoassay based on the use of recombinant Parvovirus B19 proteins as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96well-microtiter plates. The recombinant antigens serve as capture molecules for antibodies against Parvovirus B19 antigens. Bound antibodies are detected by horseradish peroxidase-(HRP)-labeled antibodies against human antibodies of IgG- or IgM-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed by precipitated product from colorless substrate solution. Color intensity correlates to the antibody concentration. Pale blue to dark blue spots are visible by eye.

**The test procedure for antibody detection is performed in 3 steps:**

### Step 1

Incubation of 1 : 101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with wash solution.

### Step 2

Incubation of wells with the anti-human IgG- or IgM-HRP-conjugate for 30 minutes at room temperature. After incubation the antibody is aspirated and unbound components are removed by 3 wash cycles with wash solution.

### Step 3

Incubation of wells with substrate solution SeramunBlau<sup>®</sup> spot dark for 30 minutes at room temperature. The reaction is stopped by aspirating the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the Seramun *SpotSight*<sup>®</sup> plate or Seramun *SpotSight*<sup>®</sup> strip scanner; the obtained images are interpreted by using the Software Seramun *SpotSight*<sup>®</sup> scan. Taking account of the test-specific evaluation criteria, the test result is determined.

## Used Parvovirus antigens

Nomenclature		Epitope structure	Origin
NS-1	C-terminus of non-capsid protein NS-1 (aa 403 - 671)	linear	recombinant
VP-2p	capsid protein VP2 (aa 1 - 554)	conformation	recombinant
VP-2r	capsid protein VP2 (aa 1 - 554)	linear	recombinant
VP-1S	N-terminus of capsid protein VP1 (aa 1 - 227)	linear	recombinant
VP-N	N-terminus of capsid protein VP1 (aa 1 - 486)	linear	recombinant
VP-C	C-terminus of capsid protein VP2 (aa 150 - 554)	linear	recombinant

## Kit components

		For 48 determinations	For 96 determinations	For 2x 96 determinations	
1	<b>WELLS</b>	<b>Wells with arrays</b>	<b>6 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>12 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>2x 12- single breakable 8-well strips in frame</b> separately vacuum-sealed with desiccant
<b>IqG detection: color coding dark purple; IgM detection: color coding light purple</b>					
2	<b>WASHBUF CONC 10X</b>	<b>Wash buffer</b> Seramun® wash buffer A	<b>100 ml concentrate</b> transparent bottle, white cap	<b>100 ml concentrate</b> transparent bottle, white cap	<b>2x 100 ml concentrate</b> transparent bottle white cap
3	<b>DIL</b>	<b>Sample diluent</b> Seramun® sample diluent B	<b>55 ml</b> ready-to-use solution, colored red transparent bottle black cap	<b>2x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cap	<b>4x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cap
4	<b>CONJ HRP IgG</b>	<b>Anti-human IgG-HRP-conjugate</b>	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap	<b>2x 8 ml</b> ready-to-use solution, colored red, transparent bottle, red cap
	<b>CONJ HRP IgM</b>	<b>Anti-human IgM-HRP-conjugate</b>	<b>8 ml</b> ready-to-use solution, colored green, transparent bottle, green cap	<b>8 ml</b> ready-to-use solution, colored green, transparent bottle, green cap	<b>2x 8 ml</b> ready-to-use solution, colored green, transparent bottle, green cap
5	<b>SUBSTR TMB</b>	<b>Substrate</b> SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	<b>8 ml</b> ready-to-use solution, black bottle, blue cap	<b>8 ml</b> ready-to-use solution, black bottle, blue cap	<b>2x 8 ml</b> ready-to-use solution, black bottle, blue cap
6	<b>COVER</b>	<b>Adhesive film</b>	<b>2x</b>	<b>2x</b>	<b>4x</b>
7	<b>SWAB</b>	<b>Swab</b> with 70 % (v/v) Isopropyl alcohol	<b>3x2</b>	<b>6x2</b>	<b>12x2</b>

## Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 48 hours. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided!

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (10 µl sample and 1000 µl sample dilution buffer) with the ready to use sample diluent.

## Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- stop watch
- Washer for 96-well microtiter plates (optional)
- Collecting devices for infectious material
- Seramun *SpotSight*® plate or Seramun *SpotSight*® strip scanner with evaluation software Seramun *SpotSight*® scan

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used up to 4 weeks when stored at 2...8 °C and up to 2 weeks when stored at room temperature.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash buffer by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or de-ionized water. For Example: 10 ml **WASHBUF CONC 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

## Workplace requirements

**Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.**

## Assay procedure

Performance at room temperature (18...25 ° C) and at the specified incubation times.

For the handling of *SeraSpot*<sup>®</sup> tests the use of the starter kit (order number: SP-000-1) is recommended. The kit is designed for multiple use.

### Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere to the underside of the bottom of the wells using the swab **SWAB** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

### Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Add **100 µl** of 1 : 101 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 1000 µl sample dilution buffer **DIL**.
3. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **CONJ HRP IgG** or **CONJ HRP IgM** to each well.
6. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **SUBSTR TMB** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent.
11. Clean bottom of wells with Swab **SWAB**.
12. Perform image acquisition using the Seramun *SpotSight*<sup>®</sup> plate or the Seramun *SpotSight*<sup>®</sup> strip scanner and evaluate the results with the Seramun *SpotSight*<sup>®</sup> scan software. If image acquisition is performed with the scanner Seramun *SpotSight*<sup>®</sup> strip, transfer single 8-well microtiter plate strips to the scanners strip holder.


**Aspiration and washing steps can be performed manually or via washer for 96-well microtiter plates.**

**After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when the plate is stored in the dark.**

## Test evaluation

### Arraylayout IgG detection

SP-012-6 G

#### Parameter

04	NS-1
05	VP-2p
06	VP-2r
07	VP-1S
08	VP-C
09	VP-N


#### Controls

01	Positive control (PC)
02	Cut-off control (CO)
03	Negative control (NC)
10	IgG conjugate control (GC)
11	Serum control (SC)

P Well position marker

### Arraylayout IgM detection

SP-012-5 M

#### Parameter

04	VP-2p
05	VP-2r
06	VP-1S
07	VP-C
08	VP-N

#### Controls

01	Positive control (PC)
02	Cut-off control (CO)
03	Negative control (NC)
09	IgM conjugate control (MC)
10	Serum control (SC)

P Well position marker

### Validity criteria for the test

The *SeraSpot*<sup>®</sup> Anti-Parvovirus-6 IgG / *SeraSpot*<sup>®</sup> Anti-Parvovirus-5 IgM test includes the following control spots:

1. Positive control (PC). Intensively stained spot, darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than Cut-off control
4. IgG, IgM conjugate control (GC, MC). Intensively stained spots with different position for IgG and IgM antibody detection. Serving as antibody isotype control.
5. Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.

The test cannot be evaluated if one of the validity criteria listed in points 1 to 5 is not fulfilled.

## Result interpretation

Evaluation of the test has to be performed by using the Seramun *SpotSight*<sup>®</sup> scan software. A visual plausibility check is possible with the aid of the magnifying lens contained in the starter kit (order number: SP-000-1).

Judgment of spot staining is performed according to the following classification:

Result Interpretation	IgG*	IgM
<b>negative</b>	no antigen spot > Cut-off control <b>or</b> VP-1S isolated > Cut-off control	no antigen spot > Cut-off control <b>or</b> VP-2p isolated > Cut-off control
<b>borderline</b>	VP-N <b>or</b> VP-2r <b>or</b> VP-C > Cut-off control	VP-N <b>or</b> VP-1S <b>or</b> VP-2r <b>or</b> VP-C > Cut-off control
<b>positive</b>	VP-2p isolated <b>or</b> two or more antigen spots > Cut-off control	two or more antigen spots > Cut-off control

\* The parameter NS-1 is not evaluated. An existing NS-1 reactivity in IgG detection indicates a past infection (6 weeks). However, this only applies to about 20% of the cases. In conjunction with the corresponding symptoms and clinical findings, an NS-1 titer may indicate a possible persistence of this virus (12-15). However, the diagnosis "persistent parvovirus B19 infection" must never be made solely on the basis of the serological findings.

## Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

### Interference

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) and rheumatoid factors (500 IU/ml) do not interfere with the test.

### Assay procedure

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere to the underside of the bottom of the wells or to the array area cause erroneous results.

## Performance characteristics

### Sensitivity and specificity

#### 1. Comparative testing of serologically pre-characterized samples

Serologically pre-determined samples of patients with suspect of Parvovirus B19 infection (n = 111 (IgG detection), n = 91 (IgM detection) were tested in the *SeraSpot*<sup>®</sup> Anti-Parvovirus-6 IgG / *SeraSpot*<sup>®</sup> Anti-Parvovirus-5 IgM test.

IgG	n = 111	Reference test		
		positive (n = 62)	borderline (n = 0)	negative (n = 49)
<i>SeraSpot</i> <sup>®</sup> Anti-Parvovirus-6 IgG	positive	60	0	0
	borderline	0	0	0
	negative	2	0	49

Borderline samples were valued positive.

Sensitivity: **97 %**

Specificity: **100 %**

IgM	n = 91	Reference test 1		
		positive (n = 24)	borderline (n = 0)	negative (n = 67)
<i>SeraSpot</i> <sup>®</sup> Anti-Parvovirus-5 IgM	positive	23	0	1
	borderline	0	0	0
	negative	1	0	66

Borderline samples were valued positive.

Sensitivity: **96 %**

Specificity: **99 %**

#### 2. Determination of sensitivity by means of the WHO standard

The *SeraSpot*<sup>®</sup> Anti-Parvovirus-6 IgG assay was compared to a reference test (IgG detection) by use of the WHO International Standard (2nd International Standard for Anti-Parvovirus B 19 plasma, human, NIBSC 01/602, 77 IU/ml). The dilution of the WHO standard to 19.25 IU / ml in both tests showed clearly positive signals for IgG antibodies against VP-1s and VP-2p.

#### 3. Examination of sera of women in pregnancy

IgG-positive (n = 53) and IgG-negative samples (n = 75) pre-determined in a comparative ELISA, as well as IgM-positive (n = 2) and IgM-negative (n = 100) samples were tested with the *SeraSpot*<sup>®</sup> Anti-Parvovirus-6 IgG / *SeraSpot*<sup>®</sup> Anti-Parvovirus-5 IgM assay. This resulted in a sensitivity of 100% in IgG and IgM detection, as well as a specificity of 99% in IgG detection and 98% in IgM detection.

## Precision

Samples of known antibody titer were assayed by *SeraSpot*<sup>®</sup> Anti-Parvovirus-6 IgG / *SeraSpot*<sup>®</sup> Anti-Parvovirus-5 IgM test. Color intensity of antigen spots divided by the color intensity of the cut-off control (Ratio) was recorded and the coefficients of variation (CV) were calculated.

IgG detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	Ratio [ $\bar{x}$ ] n = 48	CV [%]	Ratio [ $\bar{x}$ ] n = 80	CV [%]	Ratio [ $\bar{x}$ ] n = 240	CV [%]
<b>NS-1</b>	1.67	<b>8.00</b>	1.70	<b>7.38</b>	1.71	<b>10.82</b>
<b>VP-2p</b>	1.37	<b>4.45</b>	1.30	<b>7.54</b>	1.31	<b>9.80</b>
<b>VP-2r</b>	1.52	<b>4.88</b>	1.53	<b>8.72</b>	1.54	<b>7.22</b>
<b>VP-1S</b>	1.23	<b>6.31</b>	1.33	<b>8.92</b>	1.29	<b>9.70</b>
<b>VP-C</b>	1.36	<b>6.14</b>	1.37	<b>7.37</b>	1.37	<b>7.11</b>
<b>VP-N</b>	1.53	<b>8.05</b>	1.58	<b>8.68</b>	1.55	<b>8.10</b>
<b>Procedure</b>	1 operator, 48x determination, 1 batch		2 operator, 2x determination, 2x testing per day, 20 days,1 batch		2 operator, 2x determination, 2x testing per day, 20 days,3 batches	

IgM detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	Ratio [ $\bar{x}$ ] n = 48	CV [%]	Ratio [ $\bar{x}$ ] n = 80	CV [%]	Ratio [ $\bar{x}$ ] n = 240	CV [%]
<b>VP-2p</b>	1.83	<b>6.59</b>	1.67	<b>11.47</b>	1.64	<b>8.93</b>
<b>VP-2r</b>	1.82	<b>7.20</b>	1.88	<b>12.79</b>	1.74	<b>13.53</b>
<b>VP-1S</b>	1.47	<b>8.66</b>	1.59	<b>11.70</b>	1.49	<b>12.89</b>
<b>VP-C</b>	1.43	<b>5.87</b>	1.52	<b>11.85</b>	1.67	<b>14.96</b>
<b>VP-N</b>	1.41	<b>8.08</b>	1.28	<b>11.22</b>	1.20	<b>13.91</b>
<b>Procedure</b>	1 operator, 48x determination, 1 batch		5 operator, 2x determination, 2x testing per day, 20 days,1 batch		5 operator, 2x determination, 2x testing per day, 20 days,3 batches	

## Automation

Samples of known antibody titer were assayed manually or by use of the microplate processors DS2® (Dynex Technologies; manual sample dilution). Color intensity of antigen spots divided by the color intensity of the cut-off spot (Ratio) was recorded and used for calculation of the coefficient of determination  $r$  for selected antibody reactivities and the different methods of operation.

Ig-isotype detection	Dynex DS2® vs. manual processing	
	r IgG [n = 96]	r IgM [n = 88]
NS-1	0.92	./.
VP-2p	0.98	0.97
VP-2r	0.98	0.97
VP-1S	0.97	0.93
VP-C	0.99	0.98
VP-N	0.96	0.95

## Common advices and precautions

**This kit is for *in-vitro* use only.** Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

**Do not use or mix reagents from different lots except for sample diluent, washing buffer and TMB substrate solution. These can be used across batches and parameters.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**


## History of changes

Version	Section	Modifications
2019-09-09	General	General revision

## Incubation scheme **SeraSpot®** Anti-Parvovirus-6 IgG / Anti-Parvovirus 5 IgM

1.  100 µl diluted sample (1 : 101)  
30 min incubation (room temperature)  
 3 x wash with washing solution, each 400 µl per well
2.  50 µl conjugate **CONJ HRP IgG** or **CONJ HRP IgM**  
30 min incubation (room temperature)  
 3 x wash with washing solution, each 400 µl per well
3.  50 µl substrate **SUBSTR TMB**  
30 min incubation (room temperature, protected from light)  
 aspiration
4. imaging of wells and image analysis scanner Seramun *SpotSight*® plate / Seramun *SpotSight*® strip and software Seramun *SpotSight*® scan

 Manufacturer	 Date of manufacture	 Use by	<b>LOT</b> Batch code	<b>REF</b> Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	<b>IVD</b> <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

SeramunBlau®, SeraSpot® und Seramun SpotSight® are registered trademarks of Seramun Diagnostica GmbH.