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#### Literatur:

1. Andersson, A., Vetter, V., Kreutzer, L., Bauer, G.: Avidities of IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr virus serology, *J. Med. Virol.* 43: 238-244, 1994
2. Bauer, G.: Simplicity through complexity: Immunoblot with recombinant antigens as the new gold standard in Epstein-Barr virus serology, *Clin. Lab.* 47: 223-230, 2001
3. Bauer, G.: The rational basis for efficient Epstein-Barr virus (EBV) serology, *Clin. Lab.* 41: 623-634, 1995
4. Buisson, M., Fleurent, B., Mak, M., Morand, P., Chan, L., Ng, A., Guan, M., Chin, D., Seignerin, J. M.: Novel immunoblot assay using four recombinant antigens for diagnosis of Epstein-Barr virus primary infection and reactivation, *J. Clin. Microbiol.* 37(8): 2709-2714, 1999
5. Callan, M. F., Tan, L., Annel, N., Ogg, G. S., Wilson, J. D., O'Callaghan, C. A., Steven, N., McMichael, A. J., Rickinson, A. B.: Direct visualization of antigen-specific CD8<sup>+</sup> T cells during the primary immune response to Epstein-Barr virus *In vivo*, *J. Exp. Med.* 187: 1395-1402, 1998
6. Coyle, P. V., Wyatt, D., Connolly, J. H., Lynch, G. A.: Antibodies to Epstein-Barr virus in patients with nasopharyngeal carcinoma in Northern Ireland, *Ir. J. Med.* 156: 182-184, 1987
7. De Paschale, M., Agrappi, C., Manco, M. T., Mirri, P., Viganò, E. F., Clerici, P.: Seroepidemiology of EBV and interpretation of the „isolated VCA IgG“ pattern, *J. Med. Virol.* 81: 325-331, 2009
8. De Paschale, M., Clerici, P.: Serological diagnosis of Epstein-Barr virus infection: Problems and solutions, *World J. Virol.* 1(1): 31-43, 2012
9. Evans, A. S., Niedermann, J. C., Cenabre, L. C., West, B., Richards, V. A.: A prospective evaluation of heterophile and EBV specific Ig; antibody tests in clinical and subclinical infectious mononucleosis, *J. Infect. Dis.* 132: 546-554, 1975
10. Fleisher, G., Henle, W., Henle, G., Lennette, E. T., Biggar, R. J.: Primary infection with Epstein-Barr virus in infants in the United States: clinical and serologic observations, *J. Infect. Dis.* 139: 553-558, 1975
11. Gärtner, B. C., Hess, R. D., Bandt, D., Kruse, A., Rethwilm, A., Römer, K., Müller-Lantzsch, N.: Evaluation of four commercially available Epstein-Barr virus enzyme immunoassays with an immunofluorescence assay as the reference method, *Clin. Diagn. Lab. Immunol.* 10: 78-82, 2003
12. Gärtner, B. C., Kortmann, K., Schäfer, M., Müller-Lantzsch, N., Sester, U., Kaul, H., Pees, H.: No correlation in Epstein-Barr virus reactivation between serological parameters and viral load, *J. Clin. Microbiol.* 38: 2458, 2000
13. Gärtner, B. C., Müller-Lantzsch, N.: Epstein-Barr-Virus, Mikrobiologische Diagnostik, 2009
14. Gulley, M. L., Tang, W.: Laboratory assays for Epstein-Barr virus related disease, *J. Mol. Diagn.* 10: 279-292, 2008
15. Henle, G., Henle, W.: Immunofluorescence in cells derived from Burkitt's lymphoma, *J. Bacteriol.* 91: 1248-1256, 1966
16. Henle, W., Henle, G.: Epstein-Barr virus-specific serology in immunologically compromised individuals, *Cancer Res.* 41, 4222-4225, 1981
17. Henle, W., Henle, G., Andersson, J., Ernberg, I., Klein, G., Horwitz, C. A., Marklund, G., Rymo, L., Wellinder, C., Straus, S. E.: Antibody response to Epstein-Barr virus-determined nuclear antigen (EBNA)-1 and EBNA-2 in acute and chronic Epstein-Barr virus infection, *Proc. Natl. Acad. Sci. USA* 84: 570-574, 1987
18. Hess, R. D.: Routine Epstein-Barr virus diagnostics from the laboratory perspective: still challenging after 35 years, *J. Clin. Microbiol.* 42(8): 3381-3387, 2004
19. Hislop, A. D., Taylor, G. S., Sauce, D., Rickinson, A. B.: Cellular response to viral infection in humans: lessons from Epstein-Barr virus, *Annu. Rev. Immunol.* 25: 587-617, 2007
20. Holmes, R. D., Sokol, R. J.: Epstein-Barr virus and post-transplant lymphoproliferative disease. *Pediatr. Transplant.* 6: 456-464, 2002
21. Inoue, N. et al.: Use of enzyme-linked immunosorbent assays with chimeric fusion proteins to titrate antibodies against Epstein-Barr virus nuclear antigen 1, *Clin. Micro.* 30(6): 1442-1448, 1992
22. Jenson, H. B.: Virologic Diagnosis, Viral Monitoring, and Treatment of Epstein-Barr Virus Infectious Mononucleosis, *Curr. Infect. Dis. Rep.* 6: 200-207, 2004
23. Klemola, E., von Essen, R., Henle, G., Henle, W.: IM-like disease with negative heterophile agglutination test, *J. Infect. Dis.* 121: 608-614, 1970
24. Lamy, M. E., Favart, A. M., Cornu, C., Mendez, M., Segas, M., Burtonboy, G.: Study of Epstein-Barr virus (EBV) antibodies: IgG and IgM anti-VCA, IgG anti-EA and Ig anti-EBNA obtained with an original microtiter technique: -serological criteria of primary and recurrent EBV infections and follow-up of infectious mononucleosis—seroepidemiology of EBV in Belgium based on 5178 sera from patients, *Acta. Clin. Belg.* 37: 281-298, 1982

25. Leight, E. R., Sugden, B.: EBNA-1: a protein pivotal to latent infection by Epstein-Barr virus, *Rev. Med. Virol.* 10: 83-100, 2000
26. Lennette, E. T., Henle, W.: Epstein-Barr virus infections: clinical and serological features, *Lab. Manage.* 25: 23-28, 1987
27. Linde, A.: Diagnosis and pathogenesis of infectious mononucleosis and other Epstein-Barr virus associated diseases, *Rev. Med. Microbiol.* 3: 43-51, 1992
28. Linde, A.: Diagnosis of Epstein-Barr virus-related diseases, *Scand. J. Infect. Dis. Suppl.* 100: 83-88, 1996
29. Linde, A., Kallin, B., Dillner, J., Andersson, J., Jägdahl, L., Lindvall, A., Wahren, B.: Evaluation of enzyme-linked immunosorbent assays with two synthetic peptides of Epstein-Barr virus for infectious mononucleosis, *J. Infect. Dis.* 161: 903-909, 1990
30. Kampmann, M., Henninger, K., Bauer, G.: Determination of antibodies directed specifically against Epstein-Barr virus nuclear antigen-1 (EBNA-1) by anticomplementary immunofluorescence (ACIF), *Med. Microbiol.* 2: 1-8, 1993
31. Niedobitek, G., Hamilton-Dutoit, S., Herbst, H., Finn, T., Vetner, M., Pallesen, G., Stein, H.: Identification of Epstein-Barr virus-infected cells in tonsils of acute infectious mononucleosis by *in situ* hybridization, *Hum. Pathol.* 20: 796-799, 1989
32. Nikoskelainen, J., Lelkola, J., Lelkola E.: IgM antibodies specific for EBV in IM without heterophile agglutination test, *Br. Med. J.* 4: 72-75, 1974
33. Nikoskelainen, J., Neel, E. U., Stevens, D. A.: Epstein-Barr virus-specific serum immunoglobulin A as an acute-phase antibody in infectious mononucleosis, *J. Clin. Microbiol.* 10: 75-79, 1979
34. Reedman, B. M., Klein, G.: Cellular localization of an Epstein-Barr virus (EBV)-associated complement-fixing antigen in producer and non-producer lymphoblastoid cell lines, *Int. J. Cancer.* 11: 499-520, 1973
35. Rickinson, A. B., Kieff, E.: Epstein-Barr virus, In: Knipe, D. M., Howley, P. M., Griffin, D. E., Martin, M. A., Lamb, R. A., Roizman, B., editors, *Fields virology*. Philadelphia: Lippincott Williams and Wilkins: 2575-2627, 2001
36. Schillinger, M., Kampmann, M., Henninger, K., Murray, G., Hanselmann, I., Bauer, G.: Variability of humoral immune response to acute Epstein-Barr virus (EBV) infection: evaluation of the significance of serological marker, *Med. Microbiol. Letters* 2: 296-303, 1993
37. Schmitz, H., Volz, D., Krainlek-Richert, C. H., Schere, M.: Acute EBV infections in Children, *Med. Microbiol. Immunol.* 158: 58-63; 1972
38. Schubert, J., Zens, W., Weissbrich, B.: Comparative evaluation of the use of immunoblots and of IgG avidity assays as confirmatory tests for the diagnosis of acute EBV infections, *J. Clin. Virol. Methods* 11: 161-172, 1998
39. Straus, S. E., Cohen, J. I., Tosato, G., Meier, J., NIH conference: Epstein-Barr virus infections: biology, pathogenesis, and management, *Ann. Intern. Med.* 118: 45-58, 1993
40. Svahn, A., Magnusson, M., Jägdahl, L., Schloss, L., Kahlmeter, G., Linde, A.: Evaluation of three commercial enzyme-linked immunosorbent assays and two latex agglutination assays for diagnosis of primary Epstein-Barr virus infection, *J. Clin. Microbiol.* 35: 2728-2832, 1997
41. Sumya, C. V., Ench, Y.: Epstein-Barr virus infectious mononucleosis in children, II, Heterophil antibody and viral-specific responses, *Paediatrics.* 75:1011-1019, 1985
42. Thorley-Lawson, D. A.: Epstein-Barr virus: exploiting the immune system, *Nat. Rev. Immunol.* 1: 75-82, 2001
43. Tierney, R. J., Steven, N., Young, L. S., Rickinson, A. B.: Epstein-Barr virus latency in blood mononuclear cells: analysis of viral gene transcription during primary infection and in the carrier state, *J. Virol.* 68: 7374-7385, 1994
44. Yao, Q. Y., Rickinson, A. B., Epstein, M. A.: A re-examination of the Epstein-Barr virus carrier state in healthy seropositive individuals, *Int. J. Cancer* 35: 35-42, 1985
45. Van Gunsven, W. M. J., A. Nabbe, and J.M. Middeldorp: Identification and molecular characterization of two diagnostically relevant marker proteins of the Epstein-Barr virus capsid antigen, *J. Med. Virol.* 40:161-169, 1993
46. Van Grunsven, W. M. J., W. J. M. Spaan, and Middeldorp: Localization and diagnostic application of immunodominant domains of the BFRF3-encoded Epstein-Barr virus capsid protein, *J. Infect. Dis.* 170:13-19, 1994
47. Tranchand-Bunel, D., C. Auriault, E. Diesis, and H. Gras-Masse: Detection of human antibodies using "convergent" combinatorial peptide libraries or "mixotopes" designed from non-variable antigen: application to the EBV viral capsid antigen p18, *J. Pept. Res.* 52:495-508, 1998
48. Tranchand-Bunel, D., H. Gras-Masse, B. Bourez, L. Dedecker, and C. Auriault: Evaluation of an Epstein-Barr Virus (EBV) Immunoglobulin M Enzyme-Linked Immunosorbent Assay Using a Synthetic Convergent Peptide Library, or Mixotope, for Diagnosis of Primary EBV Infection, *J. Clin. Microbiol.* 37 (1), 2366-2368, 1999

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## SeraSpot<sup>®</sup> Anti-EBV-4 IgG SeraSpot<sup>®</sup> Anti-EBV-3 IgM

Spotimmunoassay for detection of IgG or IgM antibodies to *Epstein-Barr-Virus* (EBV)  
in human serum or plasma

<b>REF</b> SP-013-4 G-S6 ▽ 48	<b>REF</b> SP-013-4 G-S12 ▽ 96	<b>REF</b> SP-013-4 G-S24 ▽ 2x 96
	<b>IVD</b> <i>In-vitro</i> - diagnostic device	<b>CE</b>
<b>REF</b> SP-013-3 M-S6 ▽ 48	<b>REF</b> SP-013-3 M-S12 ▽ 96	<b>REF</b> SP-013-3 M-S24 ▽ 2x 96
	<b>IVD</b> <i>In-vitro</i> - diagnostic device	<b>CE</b>

**Seramun Diagnostica GmbH** · Sprenhagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 33767 791-10 · fax +49 33767 791-99 · [info@seramun.com](mailto:info@seramun.com)

### Introduction

The *Epstein-Barr-Virus* (EBV, also known as *Human gammaherpesvirus 4*) is a small, enveloped, double-stranded DNA virus in the family of *Herpesviridae*. The DNA contains internal and terminal repetitive sequences, whereby the terminal sequences enable a ring formation for the DNA replication.

The EBV is the aetiologic agent of the infectious mononucleosis (glandular fever). After a symptomatic or asymptomatic EBV infection the virus could establish a lifelong persistent infection. The virus is mainly transmitted through saliva (contact infection), a transmission through transplantations or blood transfusions is rare. At the age of 4, the infection rate is about 50%, from the age of 20 around 90% and from the age of 50 up to 99%.

An infection during childhood is often asymptomatic, but in case of young adults and adults the infection could lead in 30 to 60 % of all cases to the onset of an infectious mononucleosis. After an incubation time of 10 to 50 days (the average incubation time is 10 to 15 days) the first symptoms could occur. The common symptoms are fever, lymphadenopathy (lymphoma), hepatosplenomegaly (enlargement of the liver and spleen), pharyngitis (sore throat) and tonsillitis. More respiratory symptoms could be dominant in children. In immunocompetent persons the disease is self-limiting and last no longer than 3 weeks.

Especially EBV is associated in immunosuppressed patients with different lymphoproliferative diseases, such as the Burkitt's lymphoma (lymph node cancer) and nasopharyngeal carcinoma (nose and throat cancer).

During primary infection with EBV, the epithelial cells and B-lymphocytes in the oropharyngeal tract are infected. The virus is spread through the B-lymphocytes in the entire body. The infection could be lytic or latent.

The lytic infection is mainly in epithelial cells and leads to the lysis of the infected cells and the release of virus particles. During the lytic phase regulatory proteins, which belong to the group of "immediately early antigen" (IEA) and "early antigen" (EA), are synthesized to enable the production of the viral DNA,

structural components of the virion ("viral capsid antigen" (VCA)) and membrane proteins (MA). The gene products possess antigenic characteristics and initiate the formation of antibodies.

B-lymphocytes are usually latent infected. The viral genome is in the cell nucleus of the infected cell in an extrachromosomal form (episome). This leads to an immortalization of the B-cells and to the transformation of continuously growing blasts. In immunocompetent persons the growth of the transformed cells is controlled by cytotoxic T-lymphocytes and NK cells. Especially the CD8+ T-cells are directed against antigens of the lytic and latent phase. The immune response is insufficient to disrupt all infected cells. Therefore, the virus can persist in the organism.

During the latent phase, „EBV nuclear antigens" (EBNAs) and 3 latent membrane proteins (LMPs) are expressed in the infected cells. The EBNAs are a complex of at least 6 proteins (ENBA 1-6). EBNA-1 is responsible for the maintenance of the episomal status of the EBV-DNA.

The EBV reactivation, which could occur during the latent phase, is enabled through the recirculation of the infected B-cells in the lymphatic tissue, which are differentiated to plasma cells.

The humoral immune response includes antibodies against antigens of the lytic and latent phase:

Antibodies against EA showing 2 different patterns: a diffuse pattern (EA-D) und restricted pattern (EA-R). EA-D-specific IgG-antibodies are developed during the first 3-4 weeks and are undetectable after 3-4 months. Approximately 85% of the patients with a acute infection are positive for EA-D (up to 3 months after onset of symptoms). In some cases, the antibodies could be detectable over years after a primary infection. During an EBV reactivation (in immunosuppression) and in patients with a nasopharyngeal carcinoma high EA-D-specific antibody titer could be detected.

IgG antibodies against the capsid antigen VCA are developed with the first clinical symptoms of an acute infection remain in the body throughout the life. IgM antibodies are developed in the same time as the IgG antibodies for VCA but disappear within a few weeks.

EBNA-1 IgG antibodies are not detectable during the first 3-4 weeks after onset of the symptoms and are therefore a marker for a past infection. However, 5% of the patients are not producing anti-EBNA-1. A negative Anti-ENBA-1 IgG result in combination with a positive VCA-IgG, could lead to a false positive diagnosis of an acute infection. In 7% of samples out of the laboratory routine, VCA IgG is isolated reactive. p18 (VCA) is another late marker, which is in contrast with EBNA-1 always produced.

Using 3 parameters (VCA IgG, VCA IgM, EBNA IgG) a differentiation between acute and past infection is in immunocompetent patients usually possible. In immunosuppressed patients is the interpretation of the antibody reactivities is difficult. Due to a failure of the immune system and the fact that the production of the antibodies could differ (depending on the dynamic of the infection), atypical antibody patterns could occur.

In diagnostics commonly used screening tests are Immunofluorescence assays (IFAs), Enzyme Immunoassays (EIAs) and Chemiluminescence Immunoassays (CLIAs), which test for EA IgG, EBNA-1 IgG, VCA IgG and VCA IgM. In IFAs normally EBV-transformed B-cells from patients with Burkitt's lymphoma are used. In contrast to that, in EIAs and CLIAs purified, native or recombinant proteins, synthetic peptides or fusion proteins are used.

As confirmatory tests Immunoblot assays or Line-Immunoassays are applied, which are using viral lysates from EBV-transformed cells or recombinant antigens. Especially for the differentiation between acute and past infection (in the case, that VCA IgG is positive but EBNA-1 IgG and VCA IgM are negative) Immunoblot assays or Line-Immunoassays are well suitable.

#### **Literatur:**

50. Andersson, A., Vetter, V., Kreutzer, L., Bauer, G.: Avidities if IgG directed against viral capsid antigen or early antigen: useful markers for significant Epstein-Barr virus serology, *J. Med. Virol.* 43: 238-244, 1994
51. Bauer, G.: Simplicity through complexity: Immunoblot with recombinant antigens as the new gold standard in Epstein-Barr virus serology, *Clin. Lab.* 47: 223-230, 2001
52. Bauer, G.: The rational basis for efficient Epstein-Barr virus (EBV) serology, *Clin. Lab.* 41: 623-634, 1995
53. Buisson, M., Fleurent, B., Mak, M., Morand, P., Chan, L., Ng, A., Guan, M., Chin, D., Seignerin, J. M.: Novel immunoblot assay using four recombinant antigens for diagnosis of Epstein-Barr virus primary infection and reactivation, *J. Clin. Microbiol.* 37(8): 2709-2714, 1999
54. Callan, M. F., Tan, L., Annels, N., Ogg, G. S., Wilson, J. D., O'Callaghan, C. A., Steven, N., McMichael, A. J., Rickinson, A. B.: Direct visualization of antigen-specific CD8+ T cells during the primary immune response to Epstein-Barr virus *In vivo*, *J. Exp. Med.* 187: 1395-1402, 1998

55. Coyle, P. V., Wyatt, D., Connolly, J. H., Lynch, G. A.: Antibodies to Epstein-Barr virus in patients with nasopharyngeal carcinoma in Northern Ireland, *Ir. J. Med.* 156: 182-184, 1987
56. De Paschale, M., Agrappi, C., Manco, M. T., Mirri, P., Viganò, E. F., Clerici, P.: Seroepidemiology of EBV and interpretation of the „isolated VCA IgG“ pattern, *J. Med. Virol.* 81: 325-331, 2009
57. De Paschale, M., Clerici, P.: Serological diagnosis of Epstein-Barr virus infection: Problems and solutions, *World J. Virol.* 1(1): 31-43, 2012
58. Evans, A. S., Niedermann, J. C., Cenabre, L. C., West, B., Richards, V. A.: A prospective evaluation of heterophile and EBV specific Ig; antibody tests in clinical and subclinical infectious mononucleosis, *J. Infect. Dis.* 132: 546-554, 1975
59. Fleisher, G., Henle, W., Henle, G., Lennette, E. T., Biggar, R. J.: Primary infection with Epstein-Barr virus in infants in the United States: clinical and serologic observations, *J. Infect. Dis.* 139: 553-558, 1975
60. Gärtner, B. C., Hess, R. D., Bandt, D., Kruse, A., Rethwilm, A., Römer, K., Müller-Lantzsch, N.: Evaluation of four commercially available Epstein-Barr virus enzyme immunoassays with an immunofluorescence assay as the reference method, *Clin. Diagn. Immunol.* 10: 78-82, 2003
61. Gärtner, B. C., Kortmann, K., Schäfer, M., Müller-Lantzsch, N., Sester, U., Kaul, H., Pees, H.: No correlation in Epstein-Barr virus reactivation between serological parameters and viral load, *J. Clin. Microbiol.* 38: 2458, 2000
62. Gärtner, B. C., Müller-Lantzsch, N.: Epstein-Barr-Virus, *Mikrobiologische Diagnostik*, 2009
63. Gulley, M. L., Tang, W.: Laboratory assays for Epstein-Barr virus related disease, *J. Mol. Diagn.* 10: 279-292, 2008
64. Henle, G., Henle, W.: Immunofluorescence in cells derived from Burkitt's lymphoma, *J. Bacteriol.* 91: 1248-1256, 1966
65. Henle, W., Henle, G.: Epstein-Barr virus-specific serology in immunologically compromised individuals, *Cancer Res.* 41, 4222-4225, 1981
66. Henle, W., Henle, G., Andersson, J., Ernberg, I., Klein, G., Horwitz, C. A., Marklund, G., Rymo, L., Wellinder, C., Straus, S. E.: Antibody response to Epstein-Barr virus-determined nuclear antigen (EBNA)-1 and EBNA-2 in acute and chronic Epstein-Barr virus infection, *Proc. Natl. Acad. Sci. USA* 84: 570-574, 1987
67. Hess, R. D.: Routine Epstein-Barr virus diagnostics from the laboratory perspective: still challenging after 35 years, *J. Clin. Microbiol.* 42(8): 3381-3387, 2004
68. Hislop, A. D., Taylor, G. S., Sauce, D., Rickinson, A. B.: Cellular response to viral infection in humans: lessons from Epstein-Barr virus, *Annu. Rev. Immunol.* 25: 587-617, 2007
69. Holmes, R. D., Sokol, R. J.: Epstein-Barr virus and post-transplant lymphoproliferative disease, *Pediatr. Transplant.* 6: 456-464, 2002
70. Inoue, N. et al.: Use of enzyme-linked immunosorbent assays with chimeric fusion proteins to titrate antibodies against Epstein-Barr virus nuclear antigen 1, *Clin. Micro.* 30(6): 1442-1448, 1992
71. Jenson, H. B.: Virologic Diagnosis, Viral Monitoring, and Treatment of Epstein-Barr Virus Infectious Mononucleosis, *Curr. Infect. Dis. Rep.* 6: 200-207, 2004
72. Klemola, E., von Essen, R., Henle, G., Henle, W.: IM-like disease with negative heterophile agglutination test, *J. Infect. Dis.* 121: 608-614, 1970
73. Lamy, M. E., Favart, A. M., Cornu, C., Mendez, M., Segas, M., Burtonboy, G.: Study of Epstein-Barr virus (EBV) antibodies: IgG and IgM anti-VCA, IgG anti-EA and Ig anti-EBNA obtained with an original microtiter technique: -serological criteria of primary and recurrent EBV infections and follow-up of infectious mononucleosis—seroepidemiology of EBV in Belgium based on 5178 sera from patients, *Acta. Clin. Belg.* 37: 281-298, 1982
74. Leight, E. R., Sugden, B.: EBNA-1: a protein pivotal to latent infection by Epstein-Barr virus, *Rev. Med. Virol.* 10: 83-100, 2000
75. Lennette, E. T., Henle, W.: Epstein-Barr virus infections: clinical and serological features, *Lab. Manage.* 25: 23-28, 1987
76. Linde, A.: Diagnosis and pathogenesis of infectious mononucleosis and other Epstein-Barr virus associated diseases, *Rev. Med. Microbiol.* 3: 43-51, 1992
77. Linde, A.: Diagnosis of Epstein-Barr virus-related diseases, *Scand. J. Infect. Dis. Suppl.* 100: 83-88, 1996
78. Linde, A., Kallin, B., Dillner, J., Andersson, J., Jägdahl, L., Lindvall, A., Wahren, B.: Evaluation of enzyme-linked immunosorbent assays with two synthetic peptides of Epstein-Barr virus for infectious mononucleosis, *J. Infect. Dis.* 161: 903-909, 1990
79. Kampmann, M., Henninger, K., Bauer, G.: Determination of antibodies directed specifically against Epstein-Barr virus nuclear antigen-1 (EBNA-1) by anticomplementary immunofluorescence (ACIF), *Med. Microbiol.* 2: 1-8, 1993
80. Niedobitek, G., Hamilton-Dutoit, S., Herbst, H., Finn, T., Vetner, M., Pallesen, G., Stein, H.: Identification of Epstein-Barr virus-infected cells in tonsils of acute infectious mononucleosis by *in situ* hybridization, *Hum. Pathol.* 20: 796-799, 1989
81. Nikoskelainen, J., Lelkola, J., Lelkola E.: IgM antibodies specific for EBV in IM without heterophile agglutination test, *Br. Med. J.* 4: 72-75, 1974
82. Nikoskelainen, J., Neel, E. U., Stevens, D. A.: Epstein-Barr virus-specific serum immunoglobulin A as an acute-phase antibody in infectious mononucleosis, *J. Clin. Microbiol.* 10: 75-79, 1979
83. Reedman, B. M., Klein, G.: Cellular localization of an Epstein-Barr virus (EBV)-associated complement-fixing antigen in producer and non-producer lymphoblastoid cell lines, *Int. J. Cancer.* 11: 499-520, 1973
84. Rickinson, A. B., Kieff, E.: Epstein-Barr virus, In: Knipe, D. M., Howley, P. M., Griffin, D. E., Martin, M. A., Lamb, R. A., Roizman, B., editors, *Fields virology*. Philadelphia: Lippincott Williams and Wilkins: 2575-2627, 2001

85. Schillinger, M., Kampmann, M., Henninger, K., Murray, G., Hanselmann, I., Bauer, G.: Variability of humoral immune response to acute Epstein-Barr virus (EBV) infection: evaluation of the significance of serological marker, *Med. Microbiol. Letters* 2: 296-303, 1993
86. Schmitz, H., Volz, D., Krainlek-Richert, C. H., Schere, M.: Acute EBV infections in Children, *Med. Microbiol. Immunol.* 158: 58-63; 1972
87. Schubert, J., Zens, W., Weissbrich, B.: Comparative evaluation of the use of immunoblots and of IgG avidity assays as confirmatory tests for the diagnosis of acute EBV infections, *J. Clin. Virol. Methods* 11: 161-172, 1998
88. Straus, S. E., Cohen, J. I., Tosato, G., Meier, J., NIH conference: Epstein-Barr virus infections: biology, pathogenesis, and management, *Ann. Intern. Med.* 118: 45-58, 1993
89. Svahn, A., Magnusson, M., Jägdahl, L., Schloss, L., Kahlmeter, G., Linde, A.: Evaluation of three commercial enzyme-linked immunosorbent assays and two latex agglutination assays for diagnosis of primary Epstein-Barr virus infection, *J. Clin. Microbiol.* 35: 2728-2832, 1997
90. Sumya, C. V., Ench, Y.: Epstein-Barr virus infectious mononucleosis in children, II, Heterophil antibody and viral-specific responses, *Paediatrics.* 75:1011-1019, 1985
91. Thorley-Lawson, D. A.: Epstein-Barr virus: exploiting the immune system, *Nat. Rev. Immunol.* 1: 75-82, 2001
92. Tierney, R. J., Steven, N., Young, L. S., Rickinson, A. B.: Epstein-Barr virus latency in blood mononuclear cells: analysis of viral gene transcription during primary infection and in the carrier state, *J. Virol.* 68: 7374-7385, 1994
93. Yao, Q. Y., Rickinson, A. B., Epstein, M. A.: A re-examination of the Epstein-Barr virus carrier state in healthy seropositive individuals, *Int. J. Cancer* 35: 35-42, 1985
94. Van Gunsven, W. M. J., A. Nabbe, and J.M. Middeldorp: Identification and molecular characterization of two diagnostically relevant marker proteins of the Epstein-Barr virus capsid antigen. *J. Med. Virol.* 40:161-169, 1993
95. Van Grunsven, W. M. J., W. J. M. Spaan, and Middeldorp: Localization and diagnostic application of immunodominant domains of the BFRF3-encoded Epstein-Barr virus capsid protein. *J. Infect. Dis.* 170:13-19, 1994
96. Tranchand-Bunel, D., C. Auriault, E. Diesis, and H. Gras-Masse: Detection of human antibodies using "convergent" combinatorial peptide libraries or "mixotopes" designed from non-variable antigen: application to the EBV viral capsid antigen p18. *J. Pept. Res.* 52:495-508, 1998
97. Tranchand-Bunel, D., H. Gras-Masse, B. Bourez, L. Dedecker, and C. Auriault: Evaluation of an Epstein-Barr Virus (EBV) Immunoglobulin M Enzyme-Linked Immunosorbent Assay Using a Synthetic Convergent Peptide Library, or Mixotope, for Diagnosis of Primary EBV Infection. *J. Clin. Microbiol.* 37 (1), 2366-2368, 1999

## Intended use

The *SeraSpot*<sup>®</sup> Anti-EBV-4 IgG / *SeraSpot*<sup>®</sup> Anti-EBV-3 IgM test is an *in vitro* diagnostic device (spotimmunoassay, SIA) for the detection of EBV-specific IgG or IgM antibodies in human serum or plasma samples.

## Principle of the test

The *SeraSpot*<sup>®</sup> Anti-EBV-4 IgG / *SeraSpot*<sup>®</sup> Anti-EBV-3 IgM test is a solid phase immunoassay based on the use of recombinant EBV antigens as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96well-microtiter plates. The recombinant antigens serve as capture molecules for antibodies against EBV antigens. Bound antibodies are detected by horseradish peroxidase-(HRP)-labeled antibodies against human antibodies of IgG- or IgM-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed by precipitated product from colorless substrate solution. Color intensity correlates to the antibody concentration. Pale blue to dark blue spots are visible by eye.

**The test procedure for antibody detection is performed in 3 steps:**

### Step 1

Incubation of 1 : 101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation samples are aspirated and unbound components are removed by 3 wash cycles with wash solution.

### Step 2

Incubation of wells with the anti-human IgG- or IgM-HRP-conjugate for 30 minutes at room temperature. After incubation the antibody is aspirated and unbound components are removed by 3 wash cycles with wash solution.

### Step 3

Incubation of wells with substrate solution SeramunBlau<sup>®</sup> spot dark for 30 minutes at room temperature. The reaction is stopped by aspirating the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper. The coloring of the developed arrays is stable when stored protected from light.

After finishing step 3, developed assays can be scanned within 24 hours either by using the Seramun *SpotSight*<sup>®</sup> plate or Seramun *SpotSight*<sup>®</sup> strip scanner; the obtained images are interpreted by using the Software Seramun *SpotSight*<sup>®</sup> scan. Taking account of the test-specific evaluation criteria, the test result is determined.

## Used EBV antigens

Nomenclature	Properties	Relevance	Origin
EBNA-1	Epstein-Barr Nuclear Antigen-1	EBNA-1-specific IgG antibodies are specific marker for a past EBV infection. In some patients the EBNA-1-specific IgG antibody titer could be decreased (secondary EBNA-1 loss).	recombinant
VCA p18	Virus Capsid Antigen, p18	VCA p18-specific IgM and IgG antibodies are specific marker for a EBV infection. IgM antibodies are early detectable, later declining and replaced by the IgG antibody response. VCA p18-specific IgG antibodies are persistent and detectable lifelong.	recombinant
EA-D p54	Early Antigen-Diffuse, p54	EA-D is described as an early antigen. It is involved in the replication of the viral DNA. EA-D-specific antibodies could occur during primary infection. EA-D antibodies are often detected in patients with nasopharyngeal carcinoma.	recombinant
VCA gp125	Virus Capsid Antigen, Glycoprotein125	The glycoprotein 125 is one of the "virus capsid antigens". The proteins gp125 and p18 are immunodominant. VCA gp125-specific IgM antibodies are often produced very early during a primary infection. VCA gp125-specific IgG antibodies are persistent and detectable lifelong.	recombinant

## Kit components

		For 48 determinations	For 96 determinations	For 2x 96 determinations	
1	<b>WELLS</b>	<b>Wells with arrays</b>	<b>6 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>12 single breakable 8-well strips in frame</b> vacuum-sealed with desiccant	<b>2x 12- single breakable 8-well strips in frame</b> separately vacuum-sealed with desiccant
<b>IgG detection: color coding red; IgM detection: color coding orange</b>					
2	<b>WASHBUF CONC 10X</b>	<b>Wash buffer</b> Seramun® Wash buffer A	<b>100 ml concentrate</b> transparent bottle, white cup	<b>100 ml concentrate</b> transparent bottle, white cup	<b>2x 100 ml concentrate</b> transparent bottle white cup
3	<b>DIL</b>	<b>Sample diluent</b> Seramun® Sample diluent B	<b>55 ml</b> ready-to-use solution, colored red transparent bottle black cup	<b>2x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cup	<b>4x 55 ml</b> ready-to-use solution, colored red, transparent bottle, black cup
4	<b>CONJ HRP IgG</b>	<b>Anti-human IgG-HRP-conjugate</b>	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cup	<b>8 ml</b> ready-to-use solution, colored red, transparent bottle, red cup	<b>2x 8 ml</b> ready-to-use solution, colored red, transparent bottle, red cup
	<b>CONJ HRP IgM</b>	<b>Anti-human IgM-HRP-conjugate</b>	<b>8 ml</b> ready-to-use solution, colored green, transparent bottle, green cup	<b>8 ml</b> ready-to-use solution, colored green, transparent bottle, green cup	<b>2x 8 ml</b> ready-to-use solution, colored green, transparent bottle, green cup
5	<b>SUBSTR TMB</b>	<b>Substrate</b> SeramunBlau® spot dark 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	<b>8 ml</b> ready-to-use solution, black bottle, blue cup	<b>8 ml</b> ready-to-use solution, black bottle, blue cup	<b>2x 8 ml</b> ready-to-use solution, black bottle, blue cup
6	<b>COVER</b>	<b>Adhesive film</b>	<b>2x</b>	<b>2x</b>	<b>4x</b>
7	<b>SWAB</b>	<b>Swab</b> with 70 % (v/v) Isopropyl alcohol	<b>3x 2</b>	<b>6x 2</b>	<b>12x 2</b>

## Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 48 hours. For longer storage times samples have to be stored at -20 °C. Repeated freezing and thawing of samples should be avoided!

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (10 µl sample and 1000 µl sample dilution buffer) with the ready to use sample diluent.

## Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- stop watch
- Washer for 96-well microtiter plates (optional)
- Collecting devices for infectious material
- Seramun *SpotSight*® plate or Seramun *SpotSight*® strip scanner with evaluation software Seramun *SpotSight*® scan

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used up to 4 weeks when stored at 2...8 °C and up to 2 weeks when stored at room temperature.

### Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash buffer by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or de-ionized water. For Example: 10 ml **WASHBUF CONC 10X** + 90 ml distilled water. The substrate **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

## Workplace requirements

**Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.**

## Assay procedure

Performance at room temperature (18...25 ° C) and at the specified incubation times.

For the handling of *SeraSpot*<sup>®</sup> tests the use of the starter kit (order number: SP-000-1) is recommended. The kit is designed for multiple use.

### Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipet tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere to the underside of the bottom of the wells using the swab **SWAB** provided in the kit.**
- **Due to the drying of residual liquid, the images may appear more contrast-intensive, which can lead to slight deviations of the measured values during repeated scanning. Evaluation of the individual parameters is not changed!**

### Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Add **100 µl** of 1:101 diluted sample to each well. Preparation of working dilution of the sample e.g. 10 µl serum and 1000 µl sample dilution buffer.
3. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **CONJ HRP IgG** or **CONJ HRP IgM** to each well.
6. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of **WASHBUF CONC 10X**). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of the ready-to-use substrate solution **SUBSTR TMB** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent.
11. Clean bottom of wells with Swab **SWAB**.
12. Perform image acquisition using the Seramun *SpotSight*<sup>®</sup> plate or the Seramun *SpotSight*<sup>®</sup> strip scanner and evaluate the results with the Seramun *SpotSight*<sup>®</sup> scan software. If image acquisition is performed with the scanner Seramun *SpotSight*<sup>®</sup> strip, transfer single 8-well microtiter plate strips to the scanners strip holder.


**Aspiration and washing steps can be performed manually or via washer for 96-well microtiter plates.**

**After aspiration of the substrate solution, the color of developed spots is stable for 24 hours when the plate is stored in the dark.**

## Test evaluation

### Arraylayout IgG detection

SP-013-4 G

#### Parameter

04	EBNA-1
05	VCA p18
06	EA-D p54
07	VCA gp125


#### Kontrollen

01	Positive control (PC)
02	Cut-off control (CO)
03	Negative control (NC)
08	IgG conjugate control (GC)
09	Serum control (SC)

P Well position marker

### Arraylayout IgM detection

SP-013-3 M

#### Parameter

04	VCA p18
05	EA-D p54
06	VCA gp125

#### Kontrollen

01	Positive control (PC)
02	Cut-off control (CO)
03	Negative control (NC)
07	IgM conjugate control (MC)
08	Serum control (SC)

P Well position marker

### Validity criteria for the test

The *SeraSpot*<sup>®</sup> Anti-EBV-4 IgG / *SeraSpot*<sup>®</sup> Anti-EBV-3 IgM test includes the following control spots:

1. Positive control (PC). Intensely stained spot, darker than cut-off control. Always stained.
2. Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
3. Negative control (NC). Pale spot with intensity lower than cut-off control
4. IgG, IgM conjugate control (GC, MC). Intensely stained spots with different position for IgG and IgM antibody detection. Serving as antibody isotype control.
5. Serum control (SC). Intensely stained spot always stained in presence of sample. Absence of spot indicates absence of sample.

The test cannot be evaluated if one of the validity criteria listed in points 1 to 5 is not fulfilled.

## Result interpretation

Evaluation of the test has to be performed by using the Seramun *SpotSight*<sup>®</sup> scan software. A visual plausibility check is possible with the aid of the magnifying lens contained in the starter kit (order number: SP-000-1).

For a correct interpretation of the test results, the results of the IgG test should be considered in combination with the IgM test results.

Result Interpretation	IgG	IgM
<b>negative</b>	no antigen spot > Cut-off control	no antigen spot > Cut-off control or EA-D p54* (isolated positive) > Cut-off control
<b>positive</b>	one or more antigen spots > Cut-off control	one ( <b>except</b> EA-D p54* isolated positive) or more antigen spots > Cut-off control

\*In general, a screening for EA-D p54 IgG/IgM alone does not identify the stage of disease, but combining these antibodies with other parameters can be useful in laboratory diagnoses of EBV.

## Suggested interpretation of the test results.

IgG	IgM	Interpretation
EBNA-1 negative	positive	indication of an acute or primary EBV infection
EBNA-1 positive	negative	indication of a past EBV infection
EBNA-1 negative and VCA p18 + VCA gp125 positive	negative	indication of a past EBV infection
negative	negative	no antibodies against EBV antigens are detectable
all other reactivities*	all other reactivities*	results not clear (control recommended).

\*During EBV reactivation different antibody reactivities could occur. The EBV reactivation cannot be diagnosed with a specific antibody reactivity profile.

## Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

## Interference

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids and 20 mg/dl bilirubin C / bilirubin F) and rheumatoid factors (500 IU/ml) do not interfere with the test.

## Assay procedure

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may adhere to the underside of the bottom of the wells or to the array area cause erroneous results.

## Performance characteristics

### Diagnostic sensitivity and specificity

Serologically pre-determined samples (reference test 1) were tested in the *SeraSpot*® Anti-EBV-4 IgG / *SeraSpot*® Anti-EBV-3 IgM test (“Initial sensitivity” / “Initial specificity”). Samples with discrepant results were retested in a second assay for detection of EBV antibodies. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results (“Amended sensitivity” / “Amended specificity”).

IgG	n = 130	Reference test 1	
		positive (n = 84)	negative (n = 46)
<b>SeraSpot</b> ® Anti-EBV-4 IgG	positive	82	2
	negative	2	44

Initial sensitivity: **97.6 %**

Amended sensitivity: **97.7 %**

Initial specificity: **95.7 %**

Amended specificity: **97.8 %**

IgM	n = 132	Reference test 1	
		positive (n = 22)	negative (n = 110)
<b>SeraSpot</b> ® Anti-EBV-3 IgM	positive	20	3
	negative	2	107

Initial sensitivity: **90.9 %**

Amended sensitivity: **92.0 %**

Initial specificity: **97.3 %**

Amended specificity: **100 %**

## Precision

Samples of known antibody titer were assayed by *SeraSpot*<sup>®</sup> Anti-EBV-4 IgG / *SeraSpot*<sup>®</sup> Anti-EBV-3 IgM test. Color intensity of antigen spots divided by the color intensity of the Cut-off control (Ratio) was recorded and the coefficients of variation (CV) were calculated.

IgG detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	Ratio [ $\bar{x}$ ] n = 48	CV [%]	Ratio [ $\bar{x}$ ] n = 80	CV [%]	Ratio [ $\bar{x}$ ] n = 240	CV [%]
<b>EBNA-1</b>	1.23	<b>8.94</b>	1.30	<b>9.27</b>	1.39	<b>11.58</b>
<b>VCA p18</b>	1.83	<b>8.16</b>	1.63	<b>12.98</b>	1.63	<b>14.61</b>
<b>EA-D p54</b>	1.56	<b>9.03</b>	1.24	<b>19.18</b>	1.10	<b>24.20</b>
<b>VCA gp125</b>	1.04	<b>8.71</b>	1.34	<b>16.92</b>	1.28	<b>18.86</b>
<b>Procedure</b>	1 operator, 48x determination, 1 batch		2 operator, 2x determination, 2x testing per day, 20 days, 1 batch		2 operator, 2x determination, 2x testing per day, 20 days, 3 batches	

IgM detection						
	Intra-assay CV		Inter-assay CV		Interbatch CV	
	Ratio [ $\bar{x}$ ] n = 48	CV [%]	Ratio [ $\bar{x}$ ] n = 80	CV [%]	Ratio [ $\bar{x}$ ] n = 240	CV [%]
<b>EBNA-1</b>	1.77	<b>8.52</b>	1.20	<b>11.60</b>	1.19	<b>19.34</b>
<b>EA-D p54</b>	1.24	<b>7.49</b>	1.11	<b>14.32</b>	1.24	<b>14.50</b>
<b>VCA gp125</b>	1.02	<b>9.34</b>	1.03	<b>13.31</b>	1.29	<b>21.88</b>
<b>Procedure</b>	1 operator, 48x determination, 1 batch		2 operator, 2x determination, 2x testing per day, 20 days, 1 batch		2 operator, 2x determination, 2x testing per day, 20 days, 3 batches	

## Automation

Samples of known antibody titer were assayed manually or by use of the microplate processors DS2<sup>®</sup> (Dynex Technologies; manual sample dilution). Color intensity of antigen spots divided by the color intensity of the Cut-off spot (Ratio) was recorded and used for calculation of the coefficient of determination  $r^2$  for selected antibody reactivities and the different methods of operation.

Ig-isotype detection	Dynex DS2 <sup>®</sup> vs. manual processing	
	$r^2$ IgG [n = 141]	$r^2$ IgM [n = 168]
<b>EBNA-1</b>	<b>0.96</b>	<b>./.</b>
<b>VCA p18</b>	<b>0.95</b>	<b>0.80</b>
<b>EA-D p54</b>	<b>0.86</b>	<b>0.86</b>
<b>VCA gp125</b>	<b>0.96</b>	<b>0.88</b>

## Common advices and precautions

**This kit is for *in-vitro* use only.** Follow the working instruction carefully. The kit should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiration dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

**Do not use or mix reagents from different lots except for sample diluent, washing buffer and TMB substrate solution. These can be used across batches and parameters.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8 °C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**


## History of changes

Version	Section	Modifications
2019-06-01	-	Initial release

## Incubation scheme SeraSpot® Anti-EBV-4 IgG / Anti-EBV-3 IgM

-  100 µl diluted sample (1 : 101)  
30 min incubation (room temperature)  
 3 x wash with washing solution, each 400 µl per well
-  50 µl conjugate **CONJ HRP IgG** or **CONJ HRP IgM**  
30 min incubation (room temperature)  
 3 x wash with washing solution, each 400 µl per well
-  50 µl substrate **SUBSTR TMB**  
30 min incubation (room temperature, protected from light)  
 aspiration
4. imaging of wells and image analysis scanner Seramun SpotSight® plate / Seramun SpotSight® strip and software Seramun SpotSight® scan

 Manufacturer	 Date of manufacture	 Use by	<b>LOT</b> Batch code	<b>REF</b> Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	<b>IVD</b> <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	

SeramunBlau®, SeraSpot® und Seramun SpotSight® are registered trademarks of Seramun Diagnostica GmbH.