Serazym[®] Entamoeba histolytica

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Literatur:

1. Lunzen, J.v., Tannich, E., Burchard, G.-D. (1996): „Amöbenruhr und Amöbenleberabszess“. Deutsches Ärzteblatt-Ärztliche Mitteilungen. 93. Jahrgang/Heft 51-52, C: Seite 2390-2397
2. Gonzalez-Ruiz, A. et al. (1994): „Diagnosis of Amebic Dysentery by Detection of Entamoeba histolytica Fecal Antigen by an Invasive Strain-Specific, Monoclonal Antibody-Based Enzyme-Linked- Immunosorbent Assay“. Journal Of Clinical Microbiology, Vol. 32, No. 4, p. 964-970
3. Blakely, P., Sargeant, P.G. and Reed, S.L. (1990): „An immunogenic 30-kDa Surface Antigen of Pathogenic Clinical Isolates of Entamoeba histolytica“. The Journal of Infectious Diseases 162: 949-954
4. Stanley S.L. et al. (1995): „The Serine-rich Entamoeba histolytica Protein is a Phosphorylated Membrane Protein Containing O-Linked Terminal N-Acetylglucosamine Residues“. The Journal of Biological Chemistry Vol. 270, No. 8, p. 4121-4126
5. Tachibana H. et al. (1991): „Differences in genomic DNA Sequences between Pathogenic and Nonpathogenic Isolates of Entamoeba histolytica Identified by Polymerase Chain Reaction“. Journal Of Clinical Microbiology Vol. 29, No. 10, p. 2234-2239
6. Mirelman D., Nuchamowitz Y. and Stolarsky T. (1997): „Comparison of Use of Enzyme-Linked Immunosorbent Assay-Based Kits and PCR Amplification of rRNA Genes for Simultaneous Detection of Entamoeba histolytica and E. dispar“. Journal Of Clinical Microbiology Vol. 35, No. 9, p. 2405-2407

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In-vitro- Diagnostikum

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Serazym[®] Entamoeba histolytica

Enzyme immunoassay for detection of *Entamoeba histolytica* in faecal samples

REF E-018-A  96  *In-vitro*-diagnostic medical device 

Seramun Diagnostica GmbH · Sprenhagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Entamoeba histolytica is the causative agent of amoebiasis (amoebic dysentery, amoebic liver abscess). Only trophozoites of virulent strains invade the intestinal wall and cause ulcers releasing blood and mucous into the faeces. Invasive Amoebae reach the liver via the enteral vascular system and cause extensive abscesses (1). The infectious cysts are spread via faecal-oral transmission. The vegetative trophozoites are released by excystation and multiply by bisection.

Meanwhile molecular biological (PCR) and biochemical (isoenzyme analysis) methods confirm the classification of the pathogenic, invasive strains as separate species *Entamoeba histolytica*. The non-invasive strains are classified as *Entamoeba dispar*. The two species share an identical morphology and therefore a microscopic differentiation is impossible (1 - 6).

Entamoeba histolytica releases specific antigens into the intestine during its life cycle. These antigens are excreted with the faeces of infected persons. The antigen detection by enzyme immunoassay can serve as specific and easy to perform alternative to microscopy.

References:

1. Lunzen, J.v., Tannich, E., Burchard, G.-D. (1996): „Amöbenruhr und Amöbenleberabszess“. Deutsches Ärzteblatt-Ärztliche Mitteilungen. 93. Jahrgang/Heft 51-52, C: Seite 2390-2397
2. Gonzalez-Ruiz, A. et al. (1994): “Diagnosis of Amebic Dysentery by Detection of *Entamoeba histolytica* Fecal Antigen by an Invasive Strain-Specific, Monoclonal Antibody-Based Enzyme-Linked- Immunosorbent Assay“. Journal Of Clinical Microbiology, Vol. 32, No. 4, p. 964-970
3. Blakely, P., Sargeant, P.G. and Reed, S.L. (1990): „An immunogenic 30-kDa Surface Antigen of Pathogenic Clinical Isolates of *Entamoeba histolytica*“. The Journal of Infectious Diseases 162: 949-954
4. Stanley S.L. et al. (1995): „The Serine-rich *Entamoeba histolytica* Protein is a Phosphorylated Membrane Protein Containing O-Linked Terminal N-Acetylglucosamine Residues“. The Journal of Biological Chemistry Vol. 270, No. 8, p. 4121-4126
5. Tachibana H. et al. (1991): „Differences in genomic DNA Sequences between Pathogenic and Nonpathogenic Isolates of *Entamoeba histolytica* Identified by Polymerase Chain Reaction“. Journal Of Clinical Microbiology Vol. 29, No. 10, p. 2234-2239
6. Mirelman D., Nuchamowitz Y. and Stolarsky T. (1997): „Comparison of Use of Enzyme-Linked Immunosorbent Assay-Based Kits and PCR Amplification of rRNA Genes for Simultaneous Detection of *Entamoeba histolytica* and *E. dispar*“. Journal Of Clinical Microbiology Vol. 35, No. 9, p. 2405-2407

Intended use

The **Serazym®** *Entamoeba histolytica* is an *in-vitro*-diagnostic medical device for direct detection of *Entamoeba histolytica* specific antigen in faecal samples.

Principle of the test

Serazym® *Entamoeba histolytica* is a fast enzymometric two-step immunoassay based on polyclonal peptide antibodies recognizing two different epitopes of the serine-rich 30 kD membrane protein (SREHP) of *Entamoeba histolytica*.

Diluted stool specimens as well as positive and negative controls react in the first incubation step of 60 minutes at room temperature with the solid phase bound antibodies. In the following washing step unbound components are removed from the wells.

In the next step horseradish-peroxidase (HRP) labelled antibodies react with solid phase bound antibody-antigen-complexes within a reaction time of 30 minutes at room temperature. Non-bound material is separated from the solid-phase immune complexes by a subsequent washing step.

HRP converts the subsequently added colourless substrate solution of 3,3', 5,5'-tetramethylbenzidine (TMB) into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells after 10 min incubation at room temperature turning the solution from blue to yellow.

The optical density (OD) of the solution read at 450 / \geq 620 nm is directly proportional to the specifically bound amount of *Entamoeba histolytica* antigen.

Test components

			For 96 Wells
1	WELLS	Microtitration plate coated with polyclonal anti-SREHP- antibodies (rabbit)	12 single breakable 8-well strips colour coding light brown vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap
3	DIL	Sample diluent	100 ml · ready to use coloured yellow black cap
4	CONTROL +	Positive control SREHP peptide	2.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control <i>Entamoeba histolytica</i> negative sample	2.0 ml · ready to use coloured blue green cap
6	CONJ HRP	HRP-conjugate HRP-labelled, polyclonal anti-SREHP- antibodies (Sheep)	15 ml · ready to use coloured green violet cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap

Preparation and storage of samples

Collection and storage

Stool samples should be stored at 2...8°C immediately after collection and processed within 72 h. Longer storage is possible at -20°C. Repeated (> 3x) freezing and thawing of samples should be avoided. Stool samples already diluted with the *Serazym*[®] sample diluent can be stored for up to 72 h at 2...8°C before testing. Faecal samples that are already diluted in transportation media should not be used.

Preparation

Warm samples to room temperature and mix well. Pipette 1000 µl of sample diluent into a clean tube. Using a disposable stirring rod transfer about 200 mg (diameter about 4 - 6 mm) of faeces if solid or 200 µl if liquid into the tube and suspend thoroughly. If necessary, sediment floating particles by a centrifugation step with a micro centrifuge for one min at maximum speed.

Materials required but not provided

Micropipettes · multi-channel pipette or multi-pipette · reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with optical filters of 450 nm for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation

Preparation and storage of reagents

Kit size and expiry

One kit is designed for 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C, preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least one month when stored at 2...8°C.

Reagent preparation

Allow all components to reach room temperature prior to use in the assay.

The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay procedure

Dilute samples with sample diluent (3) 1 : 6, e.g. 200 mg or 200 µl faeces + 1.0 ml sample diluent (3). Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that residual fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette : 100 µl **CONTROL +** positive control (4)
100 µl **CONTROL -** negative control (5)
100 µl **diluted sample**
3. Cover plate and incubate for 60 min at RT.
4. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
5. Dispense 3 drops (or 100 µl) **CONJ HRP** HRP-conjugate (6) per well.
6. Cover plate and incubate for 30 min at RT.
7. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
8. Dispense 3 drops (or 100 µl) **SUBSTR TMB** substrate (7) per well.
9. Incubate for 10 min at RT protected from light.
10. Dispense 3 drops (or 100 µl) **STOP** stop solution (8) per well, mix gently.
11. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction stop.

Result interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.20

Samples with absorbances equal to or higher than the cut-off value are considered positive, samples with absorbances below the cut-off value are considered negative for *Entamoeba histolytica* antigen.

Reference values

Serazym® Entamoeba histolytica	
Positive	\geq Cut-off
Negative	$<$ Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters, too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

- the mean OD of the negative control is \leq 0.20 (manual performance)
 \leq 0.30 (automatic performance)
- the mean OD of the positive control is \geq 0.80

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure


There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions and not sufficiently homogenized samples can cause false results. A negative ELISA result does not exclude an *Entamoeba histolytica* infection, because the number of excreted cysts can decrease below the detection limit of the assay in invasive amoebiasis. Thus additional investigations (e. g. detection of specific antibodies or ultrasound) should be performed in case of a negative ELISA result but clinical suspect. Clinical findings have to be considered for a final result interpretation.

Automatic Processing

Performing the *Serazym*[®] *Entamoeba histolytica* on fully automated microplate processors (e.g. DS2, DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x - 8x.

Correlation: manual vs. automatic processing

A panel of 140 stool specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.998$.

Performance characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym*[®] Entamoeba histolytica calculated from 40-fold determinations of samples:

sample	mean OD	standard deviation	CV (%)
1	2.022	0.055	2.7
2	1.106	0.024	2.1
3	0.413	0.020	4.9

Inter-assay coefficient of variation (CV) in the *Serazym*[®] Entamoeba histolytica calculated from 2-fold determinations of samples in a total of 40 test runs, spread over 20 days and 2 operators:


sample	mean OD	standard deviation	CV (%)
1	1.448	0.133	9.2
2	0.883	0.067	7.6
3	0.374	0.035	9.4

Lower detection limit

The lower detection limit of *Serazym*[®] Entamoeba histolytica has been determined by titration of faecal samples spiked with *E. histolytica* HM-1:IMSS (ATCC 30459) trophozoites. The lower detection limit is < 150 trophozoites/ml.

Specificity

A panel of in all n = 188 stool samples (microbiological routine laboratory) were tested in the *Serazym*[®] Entamoeba histolytica. Of these, n = 178 samples were determined as negative and n = 10 samples as positive. The n = 10 samples were determined as negative with a reference test and a specificity of 94.7% is obtained for this panel of samples.

Specificity and sensitivity

Specificity and sensitivity of the *Serazym*[®] Entamoeba histolytica have been determined in a retrospective study in comparison to another commercially available ELISA.

	comparative ELISA positive	comparative ELISA negative
<i>Serazym</i> [®] ELISA positive	14	5**
<i>Serazym</i> [®] ELISA negative	2*	18

Specificity: 78.3% Sensitivity: 85.7%

The samples labeled with * and ** were re-tested at an external laboratory using PCR antigen detection. This results in a corrected specificity of 81.8% and a corrected sensitivity of 100%.

Cross reactivity

The following microorganisms were diluted in sample diluent with $\geq 10^8$ colony forming units and tested negative with the *Serazym*[®] *Entamoeba histolytica*. (OD 450 / 620 nm < Cut-Off):

<i>Aeromonas hydrophila</i>	ATCC (7966)
<i>Bacillus cereus</i>	ATCC (11778)
<i>Bacillus subtilis</i>	ATCC (6633)
<i>Bacteroides fragilis</i>	ATCC (25285)
<i>Citrobacter freundii</i>	ATCC (8090)
<i>Clostridium sordelli</i>	ATCC (9714)
<i>Enterobacter aerogenes</i>	ATCC (13048)
<i>Enterobacter cloacae</i>	ATCC (13047)
<i>Enterococcus faecalis</i>	ATCC (29212)
<i>Escherichia coli</i>	ATCC (25922)
<i>Klebsiella pneumonia</i>	ATCC (13883)
<i>Peptostreptococcus anaerobius</i>	ATCC (27337)
<i>Proteus vulgaris</i>	ATCC (8427)
<i>Pseudomonas aeruginosa</i>	ATCC (10145)
<i>Salmonella enterica</i> Serovar <i>thyphimurium</i>	ATCC (14028)
<i>Salmonella enterica</i> ssp. <i>galolyticus/enteritidis</i>	ATCC (13076)
<i>Shigella flexneri</i>	ATCC (12022)
<i>Shigella sonnei</i>	ATCC (25931)
<i>Staphylococcus aureus</i>	ATCC (25923)
<i>Staphylococcus epidermidis</i>	ATCC (12228)
<i>Vibrio parahaemolyticus</i>	ATCC (17802)
<i>Candida albicans</i>	ATCC (10231)

<i>Campylobacter coli</i>	ATCC (33559)
<i>Campylobacter jejuni</i>	ATCC (32291)
<i>Campylobacter fetus</i>	ATCC (27374)
<i>campylobacter upsaliensis</i>	ATCC (43954)
<i>Campylobacter lari</i>	ATCC (35221)
<i>Vibrio cholerae</i>	clinical isolate
<i>Yersinia enterocolitica</i> O:3	clinical isolate
<i>Yersinia enterocolitica</i> O:9	clinical isolate
<i>Yersinia enterocolitica</i> Y11	clinical isolate
<i>Yersinia enterocolitica</i> RKI 0803733	clinical isolate
<i>Clostridium difficile</i>	VPI 10463
<i>Salmonella infantis</i>	ATCC (51741)
<i>Salmonella anatum</i>	ATCC (9270)
<i>Salmonella paratyphi</i> A	ATCC (11511)
<i>Salmonella paratyphi</i> B	ATCC (8759)
<i>Salmonella paratyphi</i> C	Nr. 2 Pasteur
<i>Lactococcus lactis</i>	DSM (20481)
<i>Proteus mirabilis</i>	ATCC (29906)
<i>Pseudomonas fluorescens</i>	ATCC (13525)
<i>Pseudomonas putida</i>	ATCC (49128)
<i>Streptococcus agalactiae</i>	ATCC (13813)
<i>Morganella morganii</i>	ATCC (25830)

Faecal samples containing *Entamoeba dispar* in high concentrations (at least $3.9 \cdot 10^5$ cfu/ml sample) may occur a positiv signal in the *Serazym*[®] *Entamoeba histolytica* ELISA.

Interference

None of the following substances added to *E. histolytica* positive and negative stool samples showed a significant impact on the test result: The mentioned concentrations refer to 100 mg faeces.

(-)-Scopolamine N-butyl bromid (0.5%), Buscopan[®], Barium sulfate (5%), Bismuth(III) subsalicylate (0.5%, Pepto-Bismol), Cyclamat (5%), Diclofenac (0.5%), Hemoglobin human (5%), Blood human (5%), Hylak[®] N (5%), Iberogast[®] (5%), Immodium[®] akut duo (0.06/3.8%), Loperamid hydrochloride (5%, Loperamid-CT akut), Metronidazole (0.5%), Mucin (5%), Nexium[®] (0.06%), Nifuroxazide (0.5%, Pentofuryl[®]), Palmitic acid (20%), Perenterol forte (0.5%), Rennie[®] (20%), Simagel[®] (1%), Stearic acid (20%), Vancomycin (0.5%)

Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent, wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*[®] stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8°C before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

Do not smoke, eat or drink while handling kit material!

Always use protective gloves!

Never pipette material by mouth!

Note safety precautions of the single test components!

History of Changes

Version	Section	Modifications
2021-02-02	Assay procedure	Update
	Intended use	Correction

Incubation scheme *Serazym*[®] *Entamoeba histolytica* (E-018-A)

- | | | | |
|----|---|---------------------------------------|--|
| 1. |  | pipette
100 µl
100 µl
100 µl | <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-bottom: 5px;">CONTROL +</div> (4)
<div style="border: 1px solid black; padding: 2px; display: inline-block; margin-bottom: 5px;">CONTROL -</div> (5)
diluted stool sample
60 min incubation (room temperature)
 5 x wash with wash solution |
| 2. |  | 3 drops (or 100 µl) | <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-bottom: 5px;">CONJ HRP</div> (6)
30 min incubation (room temperature)
 5 x wash with wash solution |
| 3. |  | 3 drops (or 100 µl) | <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-bottom: 5px;">SUBSTR TMB</div> (7)
10 min incubation (room temperature) protected from light |
| 4. |  | 3 drops (or 100 µl) | <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-bottom: 5px;">STOP</div> (8) |

Read OD at 450 / ≥ 620 nm

	Manufacturer		Date of manufacture		Use by	<div style="border: 1px solid black; padding: 2px; display: inline-block;">LOT</div>	Batch code	<div style="border: 1px solid black; padding: 2px; display: inline-block;">REF</div>	Catalog number
	Keep away from sunlight		Temperature limits		Biological risks		Do not reuse		
	Consult instructions for use		Caution	<div style="border: 1px solid black; padding: 2px; display: inline-block;">IVD</div>	<i>In-vitro</i> -diagnostic medical device			Contains sufficient for <n> tests	