Serazym[®] Verotoxin 1+2

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Serazym[®] Verotoxin 1+2

Enzyme immunoassay for detection of Verotoxin 1+2 (Shiga-Toxin 1+2) in stool specimens

REF E-030 ▽ 96 REF E-030-A2 ▽ 2x 96 IVD *In-vitro*-diagnostic medical device CE

Seramun Diagnostica GmbH · Sprehagener Straße 1 · 15754 Heidesee · Germany · www.seramun.com
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · info@seramun.com

Introduction

Invasive and toxigenic *Escherichia coli* strains cause diarrhoea in infants and adults. Among pathogenic *E. coli* strains the group of enterohaemorrhagic *E. coli* (EHEC) can cause lifethreatening haemorrhagic colitis and haemolytic uraemic syndrome (HUS) leading to acute renal failure and haemolytic anaemia with thrombocytopenia (1, 2, 3). Strains like *E. coli* O:157; O:26; O:111 and other serovars are characterized by the production of cytotoxins (verotoxin 1 and 2 or shiga-toxin 1 and 2, shiga-toxin variants). The diagnosis of an EHEC infection is initially done by detection of the shiga-toxins. Diagnostic methods can be direct toxin detection by cytotoxicity test on vero-cells and subsequent neutralization test or the detection of the encoding genes with probes or polymerase chain reaction (PCR). These methods are time-consuming and not suited for a routine diagnostic laboratory. Immunological methods like enzyme immunoassay enable a fast and specific shiga-toxin detection in stool specimens. It is commonly recommended to enrich the EHEC bacteria in selective broth media prior to the test run to enhance the sensitivity of the method (4, 5, 6).

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Intended Use

The **Serazym® Verotoxin 1+2** is an *in-vitro*-diagnostic medical device for direct detection of Verotoxin 1 and 2 (shiga-toxin 1 and 2) in faecal specimens and stool culture supernatants.

Principle Of The Test

Serazym® Verotoxin 1+2 is an indirect two-site-immunoassay for the qualitative determination of verotoxin 1 and 2 based on polyclonal and monoclonal antibodies. Verotoxin 1 and/or 2 of specimens and the positive control react with polyclonal anti-verotoxin 1 and 2 antibodies coated on the solid phase of the microplate. After incubation for 60 minutes at room temperature (RT) non-bound material is removed by a washing step. Subsequently bound toxins specifically react with biotinylated monoclonal anti-verotoxin 1 and anti-verotoxin 2 antibodies during a second incubation period of 30 min at RT. Non-bound material is separated from the solid-phase immune complexes by a subsequent washing step. During the next incubation period of 30 min at RT horseradish peroxidase (HRP) conjugated streptavidin reacts with the bound biotinylated antibodies. Unbound conjugate is removed by a washing step. HRP converts the subsequently added colourless substrate solution of 3,3',5,5'-tetramethylbenzidine (TMB) into a blue product. The enzyme reaction is terminated by sulphuric acid dispensed into the wells after 15 min incubation at RT turning the solution from blue to yellow. The optical density (OD) of the solution read at 450 / \geq 620 nm is directly proportional to the specifically bound amount of verotoxin 1 and/or verotoxin 2. Considering the cut-off value results are interpreted as positive or negative.

Test Components

			For 96 Wells	For 2x 96 Wells
1	WELLS	Microtitration plate coated with polyclonal anti-Verotoxin 1+2 antibodies (sheep)	12 single breakable 8-well strips colour coding orange vacuum-sealed with desiccant	2x 12 single breakable 8-well strips colour coding orange vacuum-sealed with desiccant
2	WASHBUF CONC 10x	Wash buffer 10-fold	100 ml concentrate for 1000 ml solution white cap	2x 100 ml concentrate for 2x 1000 ml solution white cap
3	DIL	Sample diluent for <i>Serazym</i> [®] Verotoxin 1+2	100 ml · ready to use coloured orange black cap	2x 100 ml · ready to use coloured orange black cap
4	CONTROL +	Positive control Inactivated Verotoxin positive culture supernatant	2.0 ml · ready to use coloured blue red cap	4.0 ml · ready to use coloured blue red cap
5	CONTROL -	Negative control Verotoxin negative sample	2.0 ml · ready to use coloured blue green cap	4.0 ml · ready to use coloured blue green cap
6/1	CONJ BIOTIN	Biotin-conjugate Biotinylated, monoclonal anti-Verotoxin 1 and 2 antibodies (mouse)	15 ml · ready to use coloured green white cap	30 ml · ready to use coloured green white cap
6/2	CONJ STREPT	Streptavidin-poly-HRP-conjugate	15 ml · ready to use coloured red brown cap	30 ml · ready to use coloured red brown cap
7	SUBSTR TMB	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	15 ml · ready to use blue cap	28 ml · ready to use blue cap
8	STOP	Stop solution 0.25 M sulphuric acid	15 ml · ready to use yellow cap	28 ml · ready to use yellow cap

Preparation And Storage Of Samples

Collection and storage

Stool samples should be stored at 2...8°C (36...46°F) immediately after collection and processed within 72 hours. Longer storage is possible at least at -20°C (-4°F). In case of direct verotoxin detection from stool suspension in sample diluent, testing immediately after the sample has arrived in the laboratory should be preferred. Stool specimens for enrichment culture should be transferred to the enrichment broth within 1 - 2 hours after the sample has arrived. Samples from broth cultures already diluted with the *Serazym*[®] sample diluent can be stored for up to 24 h at 2...8°C before testing in the ELISA.

Preparation

Sample preparation for direct testing from diluted stool specimens

Quickly thaw frozen stool specimens and mix them well. Samples treated with transport media should also be mixed before testing. Pipette 500 µl of *Serazym*[®] Verotoxin 1+2 sample diluent into a clean tube. Transfer 200 mg (diameter about 4 - 6 mm) or 200 µl stool sample into the tube with the 500 µl of *Serazym*[®] Verotoxin 1+2 sample diluent and mix thoroughly.

Caution: the direct investigation of stool specimens with ELISA without prior enrichment should only serve as screening method for a fast preliminary result. A subsequent additional investigation of the

concerning sample after enrichment is absolutely necessary to reach a sufficient sensitivity. **A negative ELISA result does not necessarily exclude an infection with EHEC when stool samples are tested without enrichment culture.**

Sample preparation for testing from enrichment culture

Transfer about 200 mg or 200 µl of stool sample into a tube with 4 ml enrichment broth, e.g. EHEC-Direktmedium (Haipha) or mTSB (Mast) containing 50 ng/ml Mitomycin C and incubate for 18 to 20 hours at 37°C (98°F). If possible, use a shaker during sample incubation. Subsequently allow floating particles to sediment or if necessary, sediment floating particles by a centrifugation step. Dilute the culture supernatant 1 : 2 with *Serazym*[®] Verotoxin 1+2 sample diluent and dispense 100 µl diluted culture sample per well.

Materials Required But Not Provided

Adjustable one-channel micropipettes and pipette tips · adjustable multi-channel pipette or multi-pipette and pipette tips · Reagent container for multi-channel pipette · 8-channel wash comb with vacuum pump and waste bottle or microplate washer · microplate reader with 450 nm filter for measurement and ≥ 620 nm for reference · distilled or deionized water · glassware · tubes (2 ml) for sample preparation · enrichment broth, e. g. EHEC-Direktmedium (Haipha) or mTSB (Mast) containing 50 ng/ml Mitomycin C

Preparation And Storage Of Reagents

Kit size and expiry

One kit is designed for 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label, that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8°C (36...46°F), preferably in the original kit box. After opening all kit components are stable for at least 2 months, provided proper storage. The ready to use wash solution can be used for at least 1 month when stored at 2...8°C (36...46°F).

Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting the 10-fold concentrated wash buffer 1 + 9 with distilled or deionized water.

For Example: 10 ml wash buffer concentrate (2) + 90 ml distilled or deionized water.

Assay Procedure

Dilute stool samples with *Serazym*[®] Verotoxin 1+2 sample diluent (3) 1 + 2.5, e.g. 200 mg or 200 µl stool + 0.5 ml *Serazym*[®] Verotoxin 1+2 sample diluent (3).

For enrichment culture: Transfer 200 mg or 200 µl stool sample to 4 ml enrichment broth (e.g. EHEC-Direktmedium or mTSB + 50 ng/ml Mitomycin C) and incubate for 18 – 20 hours at 37°C (98°F) if possible on a shaker. Dilute culture supernatant 1:2 with *Serazym*[®] Verotoxin 1+2 sample diluent (3), mix thoroughly and use 100 µl / well for ELISA testing.

Avoid any time shift during dispensing of reagents and samples.

Make sure the soak time of the wash buffer in the wells is at least 5 seconds per wash cycle and that the remaining fluid is completely drained in every single wash cycle!

Avoid light exposure of the TMB substrate solution!

Working steps

1. Warm all reagents to room temperature (RT) before use. Mix gently without causing foam.
2. Pipette: 120 µl **CONTROL +** positive control (4)
120 µl **CONTROL -** negative control (5)
100 µl **diluted stool specimen** or **diluted culture supernatant**.
3. Cover plate and incubate for 60 min at RT.
4. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
5. Dispense 3 drops (or 120 µl) **CONJ BIOTIN** biotin-conjugate (6/1) per well.
6. Cover plate and incubate for 30 min at RT.
7. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
8. Dispense 3 drops (or 120 µl) **CONJ STREPT** streptavidin-poly-HRP-conjugate (6/2) per well.
9. Cover plate and incubate for 30 min at RT.
10. Decant, then wash each well 5x with 300 µl wash solution (diluted from (2)) and tap dry onto absorbent paper if necessary.
11. Dispense 3 drops (or 120 µl) **SUBSTR TMB** substrate (7) per well.
12. Incubate for 15 min at RT protected from light.
13. Dispense 3 drops (or 120 µl) **STOP** stop solution (8) per well, mix gently.
14. Read OD at 450 nm / \geq 620 nm with a microplate reader within 30 min after reaction stop.

Result Interpretation

Qualitative evaluation

Cut-off determination: OD negative control + 0.20

Samples with OD values equal with or higher than the cut-off are considered positive, samples with OD values below the cut-off are considered negative in the *Serazym*[®] Verotoxin 1+2.

Reference Values

<i>Serazym</i>[®] Verotoxin 1+2	
Positive	\geq Cut-off
Negative	$<$ Cut-off

It is recommended that each laboratory establishes its own normal and pathological reference ranges as usually done for other diagnostic parameters too. Therefore, the mentioned reference values provide a guide only to values which might be expected.

Test validity

The test run is valid if:

- the mean OD of the negative control is ≤ 0.20 (manual performance)
 ≤ 0.30 (automatic performance)
- the mean OD of the positive control is ≥ 1.00

If the above mentioned quality criteria are not met, repeat the test and make sure that the test procedure is followed correctly (incubation times and temperatures, sample and wash buffer dilution, wash steps etc.). In case of repeated failure of the quality criteria contact your supplier.

Limitations of the procedure


There is no correlation between measured absorbance and seriousness of the infection. It is also not allowed to correlate absorbances of the samples with that of the positive control. Cross contamination of reagents and samples can produce false positive results. Incorrect dilutions, not sufficiently homogenized samples or solid particles after centrifugation of the suspension can cause false results. A negative test result in the *Serazym*[®] Verotoxin 1+2 does not necessarily exclude an EHEC infection. Due to the usually low toxin concentrations in stool samples, the selective enrichment of EHEC bacteria in special culture media is of decisive influence on the sensitivity of the toxin detection with ELISA (6). The overall interpretation of the ELISA results should always consider the microbiological examination as well as clinical findings.

Automatic Processing

Performing the *Serazym*[®] Verotoxin 1+2 on fully automated microplate processors (e.g. DS2, DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment. In these cases a maximum value of 0.3 absorbance units is permissible for the negative control. It is recommended to use a wash procedure including 10 seconds soak time per strip and wash step followed by a wash step with distilled or deionized water with 10 seconds of soak time after the final wash step of each wash cycle. If necessary the number of washing steps can be enhanced from 5x to 7x-8x.

Correlation: manual – automatic processing

A panel of 188 specimens was investigated in parallel by manual and automatic processing method (DS2, Dynex Technologies) resp. The correlation was calculated with $r = 0.98$.

Performance Characteristics

Precision

Intra-assay coefficient of variation (CV) in the *Serazym*[®] Verotoxin 1+2 calculated from 12-fold determinations of samples:

verotoxin 2 (pg/ml)	mean OD	standard deviation	CV (%)
3125	2.268	0.051	2.3
800	0.799	0.037	4.7
200	0.262	0.013	5.1
0	0.056	0.006	11.3

Inter-assay coefficient of variation (CV) in the *Serazym*[®] Verotoxin 1+2 in 11 different test runs from 3-fold determinations of samples:

verotoxin 2 (pg/ml)	mean OD	standard deviation	CV (%)
3125	2.026	0.057	2.8
800	0.752	0.055	7.3
200	0.241	0.022	9.3
0	0.048	0.007	15.1

Lower detection limit

The lower detection limit of the *Serazym*[®] Verotoxin 1+2 was determined with < 100 pg / ml (<10 pg / well) by separate titration of verotoxin 1 and 2.

Specificity and sensitivity

A total of 825 stool specimens were tested in parallel with the vero-cell cytotoxicity assay and the *Serazym*[®] Verotoxin 1+2. The sample material consisted of 795 specimens sent for TPE-group diagnosis and 30 specimens that were already characterized by shiga-toxin gen PCR and culture and stored at -20°C (-4°F) until testing. Enrichment culture (mTSB; 18 - 20 hours at 37°C (98°F) of all samples was carried out prior to ELISA testing.

		Vero-cell assay	
		positive	negative
<i>Serazym</i> [®] ELISA	positive	0	2
	negative	0	793

Specificity: 98.8%

		Vero-cell assay	
		positive	negative
<i>Serazym</i> [®] ELISA	positive	11	0
	negative	2	17

Sensitivity: 84.6%

Cross reactivity

The routine diagnosis for TPE-group bacteria in this study revealed a positive pathogen detection in 141 samples. These pathogens belonged to 12 different bacterial species:

Staphylococcus aureus, enterotoxin negative; *Staphylococcus aureus*, enterotoxin positive; EHEC; *Pseudomonas aeruginosa*; *Salmonella typhimurium*; *Salmonella enteritidis*; *Salmonella spec.* *Aeromonas hydrophila*; *Aeromonas caviae*; *Campylobacter spec.*; *Hafnia alvei*; *Yersinia enterocolitica* O:3. None of them caused false positive reactions in the *Serazym*[®] Verotoxin 1+2.

Common Advices and Precautions

This kit is for *in-vitro* use only. Follow the working instructions carefully. The kit should be performed by trained technical staff only. Do not use reagents from damaged packages or bottles. The expiration dates stated on the respective labels are to be observed. **Do not use or mix reagents from different lots except for sample diluent, wash buffer, TMB/substrate solution and stop solution.**

The sample diluent IS NOT universally applicable!

The wash buffer, TMB/substrate solution and stop solution is universally applicable for the *Serazym*[®] stool ELISA Adenovirus (E-017), Rotavirus (E-020), Astrovirus (E-045), Norovirus (E-061), Clostridium difficile Toxin A+B (E-040), Clostridium difficile GDH (E-107), Campylobacter (E-093), H. pylori 2nd Gen. (E-114), Entamoeba histolytica (E-018), Cryptosporidium parvum (E-039), Giardia lamblia (E-038) and Giardia (E-106).

Do not use reagents from other manufacturers. Avoid time shift during dispensing of reagents. All reagents should be kept at 2...8 °C (36...46°F) before use. Some of the reagents may contain biocides as preservative. Further information can be found in the safety data sheet. They must not be swallowed or allowed to come into contact with skin or mucous membranes. Handle all components and all

patient samples as if potentially hazardous. Since the kit contains potentially hazardous materials, the following precautions should generally be observed:

- Do not smoke, eat or drink while handling kit material!**
- Always use protective gloves!**
- Never pipette material by mouth!**
- Note safety precautions of the single test components!**

History of Changes

Version	Section	Modifications
2021-02-02	Preparation And Storage Of Samples, Assay Procedure	Update

Incubation Scheme *Serazym*[®] Verotoxin 1+2 (E-030)

1.  pipette:
 120 µl **CONTROL +** (4)
 120 µl **CONTROL -** (5)
 100 µl **diluted stool specimen** or **culture supernatant**
 60 min incubation (room temperature)
 5 x wash with wash solution
2.  3 drops (or 120 µl) **CONJ BIOTIN** (6/1)
 30 min incubation (room temperature)
 5 x wash with wash solution
3.  3 drops (or 120 µl) **CONJ STREPT** (6/2)
 30 min incubation (room temperature)
 5 x wash with wash solution
4.  3 drops (or 120 µl) **SUBSTR TMB** (7)
 15 min incubation (room temperature) protected from light
5.  3 drops (or 120 µl) **STOP** (8)
Read OD at 450 / ≥ 620 nm

 Manufacturer	 Date of manufacture	 Use by	 LOT Batch code	 REF Catalog number
 Keep away from sunlight	 Temperature limits	 Biological risks	 Do not reuse	
 Consult instructions for use	 Caution	 IVD <i>In-vitro</i> -diagnostic medical device	 Contains sufficient for <n> tests	