

REF: AGX-5-000-0204

LOT: -

Product Identity: Full length Bacterial extracellular solute-binding protein, family 5
(aa 1-527) *Borrelia garinii*

Source: *E.coli*

Size: 86.7 kDa

Description: Protein is expressed as Glutathione-S-Transferase fusion and purified
with affinity chromatographic techniques

Purity: ≥ 90% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: PBS (0.15 M NaCl), pH 7.2 – 7.4

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only
Avoid repeated freeze-thaw cycles