

REF: AGX-5-000-0231

LOT: -

Product Identity: Full length Outer surface protein C (aa 1-207)

Borrelia garinii (strain PBi)

Source: *E.coli*

Size: 48.6 kDa

Description: Protein is expressed as Glutathione-S-Transferase fusion and purified
with affinity chromatographic techniques

Purity: ≥ 90% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: PBS (0.15 M NaCl), pH 7.2 – 7.4

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only

Avoid repeated freeze-thaw cycles