

REF: AGX-5-000-0136

LOT: -

Product Identity: Full length Flagellar filament 41 kDa core protein (aa 1-336)

Borrelia burgdorferi sensu stricto

Source: *E.coli*

Size: 62.0 kDa

Description: Protein is expressed as Glutathione-S-Transferase fusion and purified
with affinity chromatographic techniques

Purity: ≥ 95% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: PBS (0.15 M NaCl), pH 7.2 – 7.4

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only

Avoid repeated freeze-thaw cycles