

**REF:** AGX-5-000-0139

**LOT:** -

**Product Identity:** Full length Decorin-binding protein A (aa 1-170)

*Borrelia afzelii*

**Source:** *E.coli*

**Size:** 44.8 kDa

**Description:** Protein is expressed as Glutathione-S-Transferase fusion and purified with affinity chromatographic techniques

**Purity:**  $\geq 70\%$  by 12.5% SDS-PAGE

**Appearance:** clear solution

**Quantification:** BCA Protein Assay

**Concentration:** Lot specific

**Buffer:** PBS (0.15 M NaCl), pH 7.2 – 7.4

**Shelf Life:** 3 years at -70°C or below

**Note:** This product is intended for research use only

Avoid repeated freeze-thaw cycles