

**REF:**      AGX-5-000-0140

**LOT:**      -

**Product Identity:**      Full length Decorin-binding protein A (aa 1-191)

*Borrelia burgdorferi sensu stricto*

**Source:**                      *E.coli*

**Size:**                              46.8 kDa

**Description:**                  Protein is expressed as Glutathione-S-Transferase fusion and purified with affinity chromatographic techniques

**Purity:**                              ≥ 80% by 12.5% SDS-PAGE

**Appearance:**                      clear solution

**Quantification:**                      BCA Protein Assay

**Concentration:**                      Lot specific

**Buffer:**                              PBS (0.15 M NaCl), pH 7.2 – 7.4

**Shelf Life:**                              3 years at -70°C or below

**Note:**                                      This product is intended for research use only

Avoid repeated freeze-thaw cycles