

**REF:**     AGX-5-000-0141

**LOT:**     -

**Product Name:**     Full length Decorin-binding protein A (aa 1-175)

*Borrelia garinii*

**Source:**            *E.coli*

**Size:**               45.6 kDa

**Description:**     Protein is expressed as Glutathione-S-Transferase fusion and purified  
with affinity chromatographic techniques

**Purity:**             ≥ 80% by 12.5% SDS-PAGE

**Appearance:**     clear solution

**Quantification:**    BCA Protein Assay

**Concentration:**    Lot specific

**Buffer:**            PBS (0.15 M NaCl), pH 7.2 – 7.4

**Shelf Life:**        3 years at -70°C or below

**Note:**              This product is intended for research use only

Avoid repeated freeze-thaw cycles