

REF: AGX-5-000-0195

LOT: -

Product Identity: Outer surface protein C (aa 33-212)

Borrelia burgdorferi sensu stricto

Source: *E.coli*

Size: 20.0 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified
with affinity chromatographic techniques

Purity: ≥ 98% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: Acetate (pH 4.5), 150 mM NaCl

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only

Avoid repeated freeze-thaw cycles