

**REF:**      AGX-5-000-0199

**LOT:**      -

**Product Identity:**    Full length Outer surface protein C (aa 1-210)

*Borrelia spielmanii*

**Source:**                *E.coli*

**Size:**                    48.4 kDa

**Description:**        Protein is expressed as Glutathione-S-Transferase fusion and purified with affinity chromatographic techniques

**Purity:**                 ≥ 90% by 12.5% SDS-PAGE

**Appearance:**        clear solution

**Quantification:**     BCA Protein Assay

**Concentration:**      Lot specific

**Buffer:**                PBS (0.15 M NaCl), pH 7.2 – 7.4

**Shelf Life:**            3 years at -70°C or below

**Note:**                    This product is intended for research use only

Avoid repeated freeze-thaw cycles