

**REF:**      AGX-5-000-0376

**LOT:**      -

**Product Identity:**    Decorin-binding protein A (aa 22-170)

*Borrelia afzelii*

**Source:**                *E.coli*

**Size:**                    19.7 kDa

**Description:**        Protein is expressed with His-tag (N-terminally) and purified  
with affinity chromatographic techniques

**Purity:**                 ≥ 98% by 12.5% SDS-PAGE

**Appearance:**        clear solution

**Quantification:**     BCA Protein Assay

**Concentration:**     Lot specific

**Buffer:**                20 mM Tris (pH 8.0), 300 mM NaCl

**Shelf Life:**            3 years at -70°C or below

**Note:**                    This product is intended for research use only

Avoid repeated freeze-thaw cycles