

REF: AGX-5-000-0379

LOT: -

Product Identity: Basic membrane protein A (aa 19-339)

Borrelia afzelii

Source: *E.coli*

Size: 38.4 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified
with affinity chromatographic techniques

Purity: ≥ 98% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: 20 mM Tris (pH 8.0), 300 mM NaCl

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only

Avoid repeated freeze-thaw cycles