

REF: AGX-5-000-0389

LOT: -

Product Identity: Bacterial extracellular solute-binding protein, family 5 (aa 20-527)
Borrelia garinii

Source: *E.coli*

Size: 61.6 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified
with affinity chromatographic techniques

Purity: ≥ 95% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: 20 mM Tris (pH 8.0), 300 mM NaCl

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only
Avoid repeated freeze-thaw cycles