

REF: AGX-5-000-0477

LOT: -

Product Identity: Decorin-binding protein A (aa 26-191)

Borrelia burgdorferi sensu stricto

Source: *E.coli*

Size: 21.9 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified
with affinity chromatographic techniques

Purity: ≥ 95% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: 20 mM Tris (pH 8.0), 300 mM NaCl

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only

Avoid repeated freeze-thaw cycles