

**REF:** AGX-5-001-0710

**LOT:** -

**Product Identity:** Outer surface protein C (aa 19-209)

*Borrelia afzelii*

**Source:** *E.coli*

**Size:** 21.2 kDa

**Description:** Protein is expressed with His-tag (N-terminally) and purified with affinity chromatographic techniques

**Purity:**  $\geq 95\%$  by 12.5% SDS-PAGE

**Appearance:** clear solution

**Quantification:** BCA Protein Assay

**Concentration:** Lot specific

**Buffer:** 20 mM Tris (pH 8.0), 300 mM NaCl

**Shelf Life:** 3 years at -70°C or below

**Note:** This product is intended for research use only

Avoid repeated freeze-thaw cycles