

REF: AGX-5-000-0256

LOT: -

Product Identity: Full length NAD-specific glutatmate dehydrogenase (aa 1-421)
Clostridium difficile

Source: *E.coli*

Size: 47.4 kDa

Description: Protein is expressed with His-tag (N-terminally) and purified
with affinity chromatographic techniques

Purity: ≥ 98% by 12.5% SDS-PAGE

Appearance: clear solution

Quantification: BCA Protein Assay

Concentration: Lot specific

Buffer: PBS (0.15 M NaCl), pH 7.2 – 7.4

Shelf Life: 3 years at -70°C or below

Note: This product is intended for research use only
Avoid repeated freeze-thaw cycles